

11 - BioCell gold conjugates



Colloidal gold was first introduced as a cytochemical marker for TEM in 1971 by Faulk and Taylor*, and has become the universally accepted label of choice for electron microscope immunocytochemistry applications and detection of cellular and subcellular macromolecules. It is now a well established technique. Gold labels combine the benefit of unequivocal identification of antigens together with high resolution localisation. The wide choice of gold particle sizes available for labelling at different working magnifications, together with its quantitative facility, also make gold labelling a very flexible tool. One of the major advantages of colloidal gold for electron microscopy is its high electron density, offering easy detection. Procedures for using gold reagents are now well established, and these techniques have been proven to be simple, efficient and economical, as well as non-hazardous.

Colloidal gold has also been established as an excellent marker for light microscopic detection, and silver enhancing methods have allowed gold labels to provide a highly specific and sensitive method for visualising antigens in light microscopy. Gold labels have also become recognised as very important tools for detection and quantitation of proteins, antigens and nucleic acids when used with other techniques such as blotting.

British BioCell International (BBI) is a leading manufacturer of high quality immunogold and silver reagents for research and diagnostic applications, with over 30 years of experience in the manufacture of superior quality nanoparticles and secondary gold probes for the life science industry.

All products are manufactured to the highest specifications, and a quality assurance certificate is provided with each product indicating the mean particle size and coefficient of variation. This section includes many new additions to the BBI product range, such as high quality silver colloids, 2 nm gold conjugates, specialist gold and silver nanoparticles, BIOBOND™ adhesive and BIOMOUNT™ mounting medium.

BBI standard gold conjugates are supplied in a buffer containing 20 % glycerol, and are stable at 4 °C for at least 12 months from delivery, and for many years at -25 °C or below. Antibody gold conjugates are absorbed against certain species of serum proteins, minimising cross reactions to those species.

Detailed instructions for the use of all gold conjugates and silver reagents are given in the technical information booklet provided with each product.

* Faulk, W.P. & Taylor, G.M. *Immunocytochem* 8, 1081-1083, 1991.

Unconjugated gold and silver colloids

Nanometre sized gold and silver particles of uniform shape and size are invaluable tools in nanotechnology, eg. assembled arrays, light scattering etc.

The colloids can be conjugated to antibodies, proteins or many types of macromolecules for binding and reaction labelling studies. The colloids are manufactured with a very narrow size distribution.

Colloidal gold is available in sizes ranging from 2 - 250 nm: 2 - 20 nm optimised for electron microscopy; 30 - 80 nm optimised for lateral flow *in vitro* diagnostics and 100 - 250 nm for optical microscopy.

The silver colloids are available in size range 20 - 80 nm. These are not as well characterised as the gold but give alternative optical properties.

Gold colloids		Quantity		
Particle size (nm)	Particles per ml	100 ml	500 ml	20 ml
2	1.5 x 10 ¹⁴	R14076	R14076-1	R14076-2
5	5.0 x 10 ¹³	R14077	R14077-1	R14077-2
10	5.7 x 10 ¹²	R14078	R14078-1	R14078-2
15	1.4 x 10 ¹²	R14079	R14079-1	R14079-2
20	7.0 x 10 ¹¹	R14080	R14080-1	R14080-2
30	2.0 x 10 ¹¹	R14081	R14081-1	R14081-2
40	9.0 x 10 ¹⁰	R14082	R14082-1	R14082-2
50	4.5 x 10 ¹⁰	R14150	R14150-1	R14150-2
60	2.6 x 10 ¹⁰	R14151	R14151-1	R14151-2
80	1.1 x 10 ¹⁰	R14152	R14152-1	R14152-2
100	5.6 x 10 ⁹	R14153	R14153-1	R14153-2
150	1.7 x 10 ⁹	R14154	R14154-1	R14154-2
200	7.0 x 10 ⁸	R14155	R14155-1	R14155-2
250	3.6 x 10 ⁸	R14156	R14156-1	R14156-2

Silver colloids		Quantity		
Particle size (nm)	Particles per ml	100 ml	500 ml	20 ml
20	7.0 x 10 ¹¹	R14270	R14270-1	R14270-2
40	9.0 x 10 ¹⁰	R14271	R14271-1	R14271-2
60	2.6 x 10 ¹⁰	R14272	R14272-1	R14272-2
80	1.1 x 10 ¹⁰	R14273	R14273-1	R14273-2

All gold colloids are supplied at optical density 1.0 measured at 520 nm. They contain 0.01 % concentration of HAuCl₄, except 2 nm size which has a concentration of 0.002 %.

Gold conjugates for EM

Gold labels are the most satisfactory method for labelling antigens for visualisation in the electron microscope. Since their introduction, there has been a vast growth in the number of applications in both animal and plant biology. The main advantages of gold labelling lie in the high contrast, the unequivocal nature of the label, the range of particles sizes for different magnifications and the sensitivity and stability of the signal. Minimal particle clustering is also an important feature, especially where quantitative results are required.

EM grade conjugates are available in a range of sizes. All sizes are completely non-overlapping, with the narrowest available coefficient of variation to allow multi-labelling experiments. For low magnification work, larger particles (15 - 30 nm) are more easily seen. Small gold particles (5 - 10 nm) are well suited to high magnification studies and high intensity labelling. 10 nm gold conjugates are recommended for those just beginning immunogold labelling with EM and performing studies over a range of magnifications, whereas 2 nm gold particles offer ultra high sensitivity, but are usually visualised only after silver enhancing on the section.

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Goat anti-rabbit IgG (H+L) (human abs) <i>Used in immunocytochemistry of rabbit IgG primary</i>	2	R14265	R14265-1
	5	R14001	R14001-1
	10	R14002	R14002-1
	15	R14003	R14003-1
	20	R14004	R14004-1
	30	R14005	R14005-1
Goat anti-mouse IgG (H) (human abs) <i>Used in immunocytochemistry of mouse IgG primary</i>	2	R14266	R14266-1
	5	R14008	R14008-1
	10	R14009	R14009-1
	15	R14010	R14010-1
	20	R14011	R14011-1
Goat anti-mouse IgG (H+L) (human abs) <i>For use with mouse IgG primary</i>	5	R14162	R14162-1
	10	R14163	R14163-1
	15	R14164	R14164-1
	20	R14165	R14165-1
	30	R14166	R14166-1
Goat anti-mouse IgM (mu chain specific) <i>Used in immunocytochemistry of mouse IgM primary</i>	5	R14013	R14013-1
	10	R14014	R14014-1
	15	R14015	R14015-1
	20	R14016	R14016-1
	30	R14017	R14017-1
Goat anti-mouse IgG+IgM (H+L) (human abs) <i>Used in immunocytochemistry of combined mouse IgG and IgM primary</i>	5	R14018	R14018-1
	10	R14019	R14019-1
	15	R14020	R14020-1
	20	R14021	R14021-1
	30	R14172	R14172-1

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Goat anti-mouse IgA	10	R14173	R14173-1
Goat anti-mouse IgA+IgM+IgG (polyvalent) <i>Used in immunocytochemistry of combined mouse IgA, IgG and IgM primary</i>	10	R14174	R14174-1
	5	R14176	R14176-1
Goat anti-mouse IgG (rat abs) <i>Controls cross reaction. Used in immunocytochemistry of mouse IgG primary</i>	5	R14176	R14176-1
Goat anti-rat IgG (H+L) (human abs) <i>Used in immunocytochemistry of rat IgG primary</i>	5	R14024	R14024-1
	10	R14025	R14025-1
	15	R14026	R14026-1
	20	R14027	R14027-1
Goat anti-rat IgG (H+L) (mouse abs) <i>Controls cross reaction used in immunocytochemistry of rat IgG primary</i>	5	R14182	R14182-1
	10	R14183	R14183-1
Goat anti-human IgG (gamma chain specific) (mouse abs) <i>Controls cross reaction. Used in immunocytochemistry of human IgG primary</i>	10	R14029	R14029-1
	15	R14030	R14030-1
	20	R14031	R14031-1
Goat anti-human IgG (H+L) <i>Used in immunocytochemistry of human IgG primary</i>	5	R14191	R14191-1
	10	R14192	R14192-1
	15	R14193	R14193-1
	20	R14194	R14194-1

Gold conjugates for EM

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml	Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Goat anti-human IgM (mu chain specific) <i>Used in immunocytochemistry of human IgM primary</i>	5	R14198	R14198-1	Protein AG <i>Mimics a secondary antibody by binding to the FC part of the primary</i>	5	R14225	R14225-1
	10	R14199	R14199-1		10	R14226	R14226-1
					15	R14227	R14227-1
			20		R14228	R14228-1	
Goat anti-guinea pig IgG (H+L) <i>Used in immunocytochemistry of guinea pig IgG primary</i>	5	R14032	R14032-1	Cationic colloidal gold (poly-L-lysine) <i>Allows staining of negatively charged sites in cells and tissue</i>	5	R14060	R14060-1
	10	R14033	R14033-1		10	R14061	R14061-1
	15	R14034	R14034-1		15	R14062	R14062-1
	20	R14035	R14035-1		20	R14063	R14063-1
Rabbit anti-chicken IgG (H+L) <i>Used in immunocytochemistry of chicken IgY and IgG primary</i>	5	R14203	R14203-1	Bovine serum albumin <i>Used as a negative control</i>	5	R14064	R14064-1
	10	R14204	R14204-1		10	R14065	R14065-1
Rabbit anti-goat IgG (H+L) <i>Used in immunocytochemistry of goat IgG primary</i>	5	R14040	R14040-1		15	R14066	R14066-1
	10	R14041	R14041-1		20	R14067	R14067-1
	15	R14042	R14042-1	Goat F(ab')₂ anti-rabbit IgG (H+L) (human abs) <i>Use of F(ab')₂ fractions can reduce background</i>	5	R14006	R14006-1
	20	R14043	R14043-1		10	R14007	R14007-1
Rabbit anti-goat IgG (H+L) (human abs) <i>For use with goat IgG primary</i>	5	R14213	R14213-1	Goat F(ab')₂ anti-mouse IgG (H) <i>Use of F(ab')₂ fractions can reduce background with mouse IgG primary</i>	5	R14248	R14248-1
	10	R14214	R14214-1		10	R14249	R14249-1
Donkey anti-sheep IgG (H+L) <i>Used in immunocytochemistry of sheep IgG primary</i>	5	R14044	R14044-1	Goat F(ab')₂ anti-mouse IgM (mu chain specific) <i>Use of F(ab')₂ fractions can reduce background with mouse IgM primary</i>	5	R14252	R14252-1
	10	R14045	R14045-1		10	R14253	R14253-1
	15	R14046	R14046-1	Goat F(ab')₂ anti-mouse IgG+IgM (H+L) (human abs) <i>Use of F(ab')₂ fractions can reduce background with mouse IgG and IgM primary</i>	5	R14022	R14022-1
	20	R14047	R14047-1		10	R14023	R14023-1
Goat anti-horseradish peroxidase <i>Used as a secondary corroborative probe</i>	5	R14218	R14218-1	Protein A <i>Mimics a secondary antibody by binding to the FC part of the primary</i>	5	R14048	R14048-1
	10	R14219	R14219-1		10	R14049	R14049-1
Protein G <i>Mimics a secondary antibody by binding to the FC part of the primary</i>	5	R14052	R14052-1		15	R14050	R14050-1
	10	R14053	R14053-1		20	R14051	R14051-1
	15	R14054	R14054-1		(H) heavy chain specific (H+L) heavy and light chain specific (human abs) absorbed against human serum proteins (mouse abs) absorbed against mouse serum proteins (rat abs) absorbed against rat serum proteins		
	20	R14055	R14055-1				

Gold conjugates for LM

Light microscope immunocytochemistry with gold conjugates follows similar principles to the well established indirect labelling methods using peroxidase labelled secondary antibodies, but offers several advantages. No hazardous materials are involved, the procedure is simpler and yields a much more intense stain than conventional peroxidase and PAP techniques. The natural red colour of the gold in transmitted light gives a stain which is easily visualised. Of particular importance is the ability to label sections from the same block for both light and electron microscopy.

Gold conjugate particles are available in 2 nm and 5 nm particle diameters and produce excellent labelling in wax, resin embedded or frozen sections. With silver enhancement a dense black reaction is produced which is clearly localised and can be counterstained with a variety of tissue stains.

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Goat anti-rabbit IgG (H+L) (human abs)	5	R14085	R14085-1
Goat anti-mouse IgG (H) (human abs)	5	R14087	R14087-1
Goat anti-mouse IgG (H+L) (human abs)	5	R14168	R14168-1
Goat anti-mouse IgM (mu chain specific)	5	R14088	R14088-1
Goat anti-mouse IgG+IgM (H+L) (human abs)	5	R14089	R14089-1
Goat anti-mouse IgG (rat abs)	5	R14178	R14178-1
Goat anti-rat IgG (H+L) (human abs)	5	R14091	R14091-1
Goat anti-rat IgG (H+L) (mouse abs)	5	R14184	R14184-1
Goat anti-human IgG gamma chain specific	5	R14092	R14092-1
Goat anti-human IgG (H+L)	5	R14196	R14196-1
Goat anti-human IgM (mu chain specific)	5	R14200	R14200-1
Goat anti-guinea pig IgG (H+L)	5	R14093	R14093-1
Rabbit anti-chicken IgG (H+L)	5	R14205	R14205-1

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Rabbit anti-goat IgG (H+L)	5	R14095	R14095-1
Rabbit anti-goat IgG (H+L) (human abs)	5	R14215	R14215-1
Donkey anti-sheep IgG (H+L)	5	R14096	R14096-1
Goat anti-horseradish peroxidase	5	R14220	R14220-1
Protein A	5	R14097	R14097-1
Protein G	5	R14098	R14098-1
Protein AG	5	R14229	R14229-1
Cationic colloidal gold (poly-L-lysine)	5	R14100	R14100-1
Bovine serum albumin (negative control)	5	R14101	R14101-1
Goat F(ab')₂ anti-rabbit IgG (H+L) (human abs)	5	R14086	R14086-1
Goat F(ab')₂ anti-mouse IgG (H) (human abs)	5	R14250	R14250-1
Goat F(ab')₂ anti-mouse IgM (mu chain specific)	5	R14254	R14254-1
Goat F(ab')₂ anti-mouse IgG+IgM (H+L) (human abs)	5	R14090	R14090-1

Gold conjugates for *in situ* hybridisation

Specific DNA or RNA fragments may be readily identified within cells after *in situ* hybridisation with appropriate DNA/RNA probes. During hybridisation at the cellular level, the labelled probe binds to complementary cellular target DNA/RNA, then is detected with a suitably labelled antibody. Immunogold techniques have become popular for both light and electron microscope detection of DNA/RNA in cells, due to their sensitivity and high resolution. Tissue sections can be counterstained with all normal stains.

Gold conjugates for blotted proteins

Blotting techniques have been established for a number of years, and the use of immunogold conjugates provides a safe, sensitive and economical alternative to traditional labels.

Gold conjugates are an extremely sensitive method for identification of immobilised proteins and nucleic acids on membranes such as nitrocellulose or PVDF. This technique can be used for determination of protein bands transferred from electrophoretic gels, DNA fragments hybridised on membranes with nucleic acid probes, and identification and quantitation of proteins and macromolecules blotted directly in small quantities.

Blotting grade gold conjugates are available with either 2 or 20 nm gold particles.

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 0.25 ml
Streptavidin	2 EM grade	R14280	R14280-1
<i>For the indirect detection of a biotinylated primary antibody</i>	5 EM grade	R14056	R14056-1
	10 EM grade	R14057	R14057-1
	15 EM grade	R14058	R14058-1
	5 LM grade	R14099	R14099-1
Goat anti-biotin <i>For the indirect detection of a biotinylated primary antibody</i>	5 EM grade	R14036	R14036-1
	10 EM grade	R14037	R14037-1
	15 EM grade	R14038	R14038-1
	20 EM grade	R14039	R14039-1
	5 LM grade	R14094	R14094-1
Monoclonal anti-biotin <i>For the indirect detection of a biotinylated primary antibody</i>	2 EM grade	R14281	R14281-1

Conjugated protein	Gold particle size (nm)	Quantity 1 ml	Quantity 2 ml
Goat anti-rabbit IgG (H+L) (human abs)	20	R14102	R14102-2
Goat anti-mouse IgG (H) (human abs)	20	R14103	R14103-2
Goat anti-mouse IgG (H+L) (human abs)	20	R14269	R14269-2
Goat anti-mouse IgM (mu chain specific)	20	R14104	R14104-2
Goat anti-mouse IgG+IgM (H+L) (human abs)	20	R14105	R14105-2
Goat anti-rat IgG (H+L) (human abs)	20	R14106	R14106-2
Goat anti-human IgG (gamma chain specific)	20	R14107	R14107-2
Goat anti-human IgG (H+L)	20	R14268	R14268-2
Goat anti-guinea pig IgG (H+L)	20	R14108	R14108-2
Goat anti-biotin	20	R14109	R14109-2
Rabbit anti-goat IgG (H+L)	20	R14110	R14110-2
Donkey anti-sheep IgG (H+L)	20	R14111	R14111-2
Protein A	20	R14112	R14112-2
Protein G	20	R14113	R14113-2
Protein AG	20	R14267	R14267-2
Streptavidin	20	R14114	R14114-2
Cationic colloidal gold (poly-L-lysine)	20	R14232	R14232-2
Bovine serum albumin	20	R14245	R14245-2

2 nm gold conjugates

These ultra small gold conjugates may be used effectively for EM, LM and immunoblotting applications. A major advantage of the ultra small gold particles is the increased penetration of antibodies into cells and tissues without the need for permeabilisation. Together with silver enhancing techniques they provide an extremely sensitive low background labelling with increased penetration into tissue sections and through cell membranes.

Specialist gold and silver nanoparticles

R14282	Raman specification gold colloid. 150 µl
R14283	Raman specification silver colloid. 150 µl

Customised conjugation service

A customised conjugation service is available, offering bespoke production of high quality conjugates optimised for your sensitivity specifications. Antibodies, proteins and macromolecules can be labelled with gold particles of all sizes.

Please contact us with your requirements.

PROTOGOLD®

PROTOGOLD colloidal gold solution has been specifically developed to give intense staining of proteins blotted onto membranes by dot blots or electrophoretic gels. The gold particles are negatively charged and bind selectively to the blotted proteins with very low background signal, producing a dark red stain. The gold signal does not fade and, when used in conjunction with the BL silver enhancing kit, a further amplification of the signal by x10 - 100 can be achieved. This method is more sensitive than Coomassie Blue or silver staining, and the kit will stain over 20 blots with high intensity.

Benefits of using PROTOGOLD include fast quantitative protein determination, with visible staining of the proteins occurring within minutes; greater sensitivity, allowing detection of picogram quantities of protein; and high resolution, offering crisp images. It is easy to use, with no destaining necessary, allowing staining to be continued indefinitely without over staining.

R14084	PROTOGOLD kit. 500 ml
R14257	Test strips. 50

Silver enhancing kits

These kits offer ease of use and sensitivity for the amplification of immunogold labels in electron and light microscopy and blotting applications. Enhancement occurs during the reduction of silver from one solution (the enhancer) by another (the initiator) in the presence of gold particles. This reduction causes silver to build up on the surface of gold particles, with no diffusion of the signal due to the discrete nature of the silver growth on the gold particles. Enhancement is rapid and easily controllable; the reaction can be monitored in bright light and terminated by washing in tap water, with signal amplifications of x10 - 100 readily achieved.

Two silver enhancing kits are available for different applications, an EM/LM kit and a blotting kit.

R14115	LM and EM kit*. 2 x 15 ml
R14116	Blotting kit**. 2 x 250 ml
R14117	Test strips. 10

* Sufficient for 300 slides

** Sufficient for 20 - 30 blots

Blocking reagents

Non-specific labelling of specimens can occur during immunolabelling procedures. The source of this background labelling must be determined by careful and systematic use of controls, and eliminated to allow proper analysis of the specimen. Background labelling can occur either in the specimen or in the incubating solutions and in either case the background can be substantially reduced by the careful use of blocking reagents.

R14120	Tween® 20. 10 ml
R14121	BSA (powder) (fatty acid free). 10 g
R14122	Gelatin (fish) (45 %). 10 ml

Biomount™

Biomount mounting medium reduces fading of immunogold/silver signals in sections on glass slides, and can be used with both resin and wax embedded tissue sections. It is miscible with xylene, and may be applied following dehydration. Labelling will retain its intensity and contrast.

R1354	Biomount. 100 ml <i>Flammable, harmful</i>
--------------	--

BioBond™

For light and electron microscopy the choice of specimen preparation is critical for the preservation of antigens in the sample. Having prepared the tissue sample for immunolabelling it is imperative to perform the incubation to maximise the specific signal and minimise the background.

Under some incubation conditions tissue sections can be removed from the slide. Some adhesives are not suitable for use with immunogold labelling because of the increased background caused by the attraction of gold particles to the adhesive on the slide.

BioBond is a special adhesive which gives a good bonding between the slide and tissue section for subsequent incubation. It is particularly effective for use with severe incubation conditions such as those used for *in situ* hybridisation. It is suitable for many types of tissue specimen including paraffin wax or resin sections, cell smears, cytospin or cryostat sections.

Sufficient to coat at least 1000 slides.

R1355	BioBond tissue adhesive. 20 ml
--------------	--------------------------------