



Epoxy embeddings sandwiched between ACLAR® sheets



ACLAR® pulled away from the flat specimen embedding.



Cut out the area of interest and glue it on a 00 cylinder for sectioning.



Disc Punches are available in four sizes:

5/16" (7.9mm)


3/8" (9.5mm)

7/16 (11.1mm)


1/2" (12.7mm)



For comparative TEM and SEM procedures, cut an ACLAR® circle in half and compare the two after fixation (Kingsley3); use a blunt needle to mark.



Press between glass plates



Fresh material is cut on a vibratory slicer, treated with HRP and placed on ACLAR® cut into a slide shape - observe under LM - if OK, process for TEM on the slide.

ACLAR® Embedding Film

Physical & Chemical Resistance Data

Physical Data

Density	2.12
Thickness	7.8 mil (0.198mm)
Clarity	Clear
Water Absorption	nil
Water Vapor Transmission Rate @ 100°F (37.7°C) / 100% RH	0.003gm / 100in ² / day (0.047gm / m ² / day)
Dimensional Stability, 10 min @ 300°F (149°C)	≤2%
Dimensional Change, 10 min. @ 300°F (149°C)	< 2%
Thermal Conductivity	4.7 x 10 ⁻⁴ cal-cm/cm ² sec °C
Crystalline Melting Point	395.6 - 399.2°F (202 - 204°C)
Flammability	Nonflammable

Chemical Resistance Data

Acetone	No effect
100% Ethyl alcohol	No effect
Liquid Nitrogen	Remains flexible
Osmium Tetroxide	No effect
Propylene Oxide	No effect
Butyl Alcohol	None
Carbon Tetrachloride	Slightly flexible
1,2-Dichloroethane	None
Ethyl Acetate	Very flexible
Ethyl Ether	Very flexible
Ethylene Oxide	Very flexible
Formic Acid	None
Gasoline	None
All Acids (HCl, H ₂ SO ₄ ..)	None
Methanol	None
Toluene	Slightly flexible
Plastisolve	None