


Certificate of Conformance for Critical Dimension Standard





Product Number: CDMS-1C

Product Description: 2.5x2.5mm Critical Dimension Magnification Standard

Product Serial Number: CD-PG01-0000

The accuracy of these products was determined by reference comparison to working standards traceable to the National Institute of Standards and Technology (NIST).

Line	Certified average pitch	Number of lines	Certified distance between first and last line (1σ)	Total expanded uncertainty (3σ)
2.0mm	2.00mm	2	2.00mm ± 0.05%	± 0.15%
1.0mm	1.00mm	2	1.00mm ± 0.03%	± 0.09%
0.5mm	0.500mm	2	0.500mm ± 0.03%	± 0.09%
0.25mm	0.250mm	2	0.250mm ± 0.03%	± 0.09%
10µm	10.00µm	9	80.00µm ± 0.03%	± 0.09%
5.0µm	5.00µm	10	45.00µm ± 0.03%	± 0.09%
2.0µm	2.00µm	16	30.00μm ± 0.03%	± 0.09%
1.0µm	1.00µm	17	16.00µm ± 0.03%	± 0.09%

The certified average pitch is derived from the stated certified length that was determined using 10 measurements (taken center-to-center) over the stated number of lines (i.e., length divided by the number of lines minus one). The total expanded uncertainty includes both Type A and Type B uncertainties corrected for sample size using an appropriate Student t-factor.

Equipment used:

Instrument	Model number	Serial #	NIST Certified CD/Recalibration	Resolution	Repeatability
FE-SEM	FEI Quanta 3D FEG	D8894	CD-PG01-0211/June 2013	1.2nm	0.03%

Certified by	Signature	Date

This certificate shall not be reproduced without the permission of AISthesis Products, Inc. P.O. Box 1167, Ashland, Oregon 97520 Tel: 541.482.2186 E-mail: CD@aisthesisproducts.com