

Part #	Description	Unit	Price €
AGA0100	Objective Aperture Plate	each	766,80
AGA0100P	Platinum Aperture 1.85mm dia x 0.1mm Hole 5µm	each	312,50
AGA0101P	Platinum Aperture 1.85mm dia x 0.1mm Hole 10µm	each	181,80
AGA0102P	Platinum Aperture 1.85mm dia x 0.1mm Hole 20µm	each	165,60
AGA0103P	Platinum Aperture 1.85mm dia x 0.1mm Hole 25µm	each	162,60
AGA0104P	Platinum Aperture 1.85mm dia x 0.1mm Hole 30µm	each	162,80
AGA0105P	Platinum Aperture 1.85mm dia x 0.1mm Hole 40µm	each	152,80
AGA0106P	Platinum Aperture 1.85mm dia x 0.1mm Hole 50µm	each	152,80
AGA0107P	Platinum Aperture 1.85mm dia x 0.1mm Hole 60µm	each	152,80
AGA0108P	Platinum Aperture 1.85mm dia x 0.1mm Hole 70µm	each	152,80
AGA0109P	Platinum Aperture 1.85mm dia x 0.1mm Hole 100µm	each	143,90
AGA0110P	Platinum Aperture 1.85mm dia x 0.1mm Hole 150µm	each	176,30
AGA0111P	Platinum Aperture 1.85mm dia x 0.1mm Hole 200µm	each	145,80
AGA0111T	Tantalum Aperture for Wehnelt - 200µm Hole	each	187,90
AGA0112P	Platinum Aperture 1.85mm dia x 0.1mm Hole 250µm	each	150,90
AGA0113P	Platinum Aperture 1.85mm dia x 0.1mm Hole 300µm	each	150,90
AGA0113T	Tantalum Aperture for Wehnelt - 300µm Hole	each	190,30
AGA0114P	Platinum Aperture 1.85mm dia x 0.1mm Hole 400µm	each	150,90
AGA0115P	Platinum Aperture 1.85mm dia x 0.1mm Hole 500µm	each	150,90
AGA0115T	Tantalum Aperture for Wehnelt - 500µm Hole	each	187,90
AGA0116P	Platinum Aperture 1.85mm dia x 0.1mm Hole 600µm	each	150,90
AGA0117P	Platinum Aperture 1.85mm dia x 0.1mm Hole 750µm	each	150,90
AGA0117T	Tantalum Aperture for Wehnelt - 750µm Hole	each	238,70
AGA0118P	Platinum Aperture 1.85mm dia x 0.1mm Hole 1000µm	each	150,90
AGA0119P	Platinum Aperture 1.85mm dia x 0.1mm Hole 1250µm	each	150,90
AGA0200P	Platinum Aperture 2mm dia x 0.6mm Hole 5µm	each	319,00
AGA0201M	Molybdenum Aperture 2mm dia x 0.6mm Hole 10µm	each	200,20
AGA0201P	Platinum Aperture 2mm dia x 0.6mm Hole 10µm	each	201,20
AGA0202M	Molybdenum Aperture 2mm dia x 0.6mm Hole 20µm	each	114,80
AGA0202P	Platinum Aperture 2mm dia x 0.6mm Hole 20µm	each	132,00
AGA0203M	Molybdenum Aperture 2mm dia x 0.6mm Hole 25µm	each	108,70
AGA0203P	Platinum Aperture 2mm dia x 0.6mm Hole 25µm	each	108,90
AGA0204M	Molybdenum Aperture 2mm dia x 0.6mm Hole 30µm	each	108,60
AGA0204P	Platinum Aperture 2mm dia x 0.6mm Hole 30µm	each	100,80
AGA0205M	Molybdenum Aperture 2mm dia x 0.6mm Hole 40µm	each	99,60
AGA0205P	Platinum Aperture 2mm dia x 0.6mm Hole 40µm	each	95,70
AGA0206M	Molybdenum Aperture 2mm dia x 0.6mm Hole 50µm	each	98,10
AGA0206P	Platinum Aperture 2mm dia x 0.6mm Hole 50µm	each	100,60
AGA0207M	Molybdenum Aperture 2mm dia x 0.6mm Hole 60µm	each	96,00
AGA0207P	Platinum Aperture 2mm dia x 0.6mm Hole 60µm	each	93,30
AGA0208M	Molybdenum Aperture 2mm dia x 0.6mm Hole 70µm	each	96,00
AGA0208P	Platinum Aperture 2mm dia x 0.6mm Hole 70µm	each	96,00
AGA0209M	Molybdenum Aperture 2mm dia x 0.6mm Hole 100µm	each	96,00
AGA0209P	Platinum Aperture 2mm dia x 0.6mm Hole 100µm	each	106,10
AGA0209T	Tantalum Aperture 2mm dia x 0.6mm Hole 100µm	each	96,00
AGA0210M	Molybdenum Aperture 2mm dia x 0.6mm Hole 150µm	each	107,20
AGA0210P	Platinum Aperture 2mm dia x 0.6mm Hole 150µm	each	107,20
AGA0211M	Molybdenum Aperture 2mm dia x 0.6mm Hole 200µm	each	96,00
AGA0211P	Platinum Aperture 2mm dia x 0.6mm Hole 200µm	each	95,20
AGA0211T	Tantalum Aperture 2mm dia x 0.6mm Hole 200µm	each	96,00
AGA0212M	Molybdenum Aperture 2mm dia x 0.6mm Hole 250µm	each	96,00
AGA0212P	Platinum Aperture 2mm dia x 0.6mm Hole 250µm	each	96,00
AGA0213M	Molybdenum Aperture 2mm dia x 0.6mm Hole 300µm	each	127,30
AGA0213P	Platinum Aperture 2mm dia x 0.6mm Hole 300µm	each	96,00
AGA0214M	Molybdenum Aperture 2mm dia x 0.6mm Hole 400µm	each	96,00
AGA0214P	Platinum Aperture 2mm dia x 0.6mm Hole 400µm	each	96,00
AGA0214T	Tantalum Aperture 2mm dia x 0.6mm Hole 400µm	each	96,70
AGA0215M	Molybdenum Aperture 2mm dia x 0.6mm Hole 500µm	each	96,00

Part #	Description	Unit	Price €
AGA0215P	Platinum Aperture 2mm dia x 0.6mm Hole 500µm	each	96,00
AGA0216M	Molybdenum Aperture 2mm dia x 0.6mm Hole 600µm	each	96,00
AGA0216P	Platinum Aperture 2mm dia x 0.6mm Hole 600µm	each	96,00
AGA0217M	Molybdenum Aperture 2mm dia x 0.6mm Hole 750µm	each	93,30
AGA0217P	Platinum Aperture 2mm dia x 0.6mm Hole 750µm	each	99,10
AGA0218M	Molybdenum Aperture 2mm dia x 0.6mm Hole 1000µm	each	97,60
AGA0218P	Platinum Aperture 2mm dia x 0.6mm Hole 1000µm	each	97,60
AGA0219M	Molybdenum Aperture 2mm dia x 0.6mm Hole 1250µm	each	97,60
AGA0219P	Platinum Aperture 2mm dia x 0.6mm Hole 1200µm	each	97,60
AGA0220P	Platinum Aperture 2mm dia x 0.25mm Hole 5µm	each	296,40
AGA0221P	Platinum Aperture 2mm dia x 0.25mm Hole 200µm	each	94,20
AGA0222P	Platinum Aperture 2mm dia x 0.25mm Hole 20µm	each	118,60
AGA0223P	Platinum Aperture 2mm dia x 0.25mm Hole 300µm	each	94,60
AGA0224P	Platinum Aperture 2mm dia x 0.25mm Hole 400µm	each	94,60
AGA0225P	Platinum Aperture 2mm dia x 0.25mm Hole 40µm	each	94,60
AGA0226P	Platinum Aperture 2mm dia x 0.25mm Hole 50µm	each	94,60
AGA0227P	Platinum Aperture 2mm dia x 0.25mm Hole 60µm	each	94,60
AGA0228P	Platinum Aperture 2mm dia x 0.25mm Hole 70µm	each	94,60
AGA0229P	Platinum Aperture 2mm dia x 0.25mm Hole 100µm	each	96,70
AGA0230P	Platinum Aperture 2mm dia x 0.25mm Hole 150µm	each	94,20
AGA0231P	Platinum Aperture 2mm dia x 0.25mm Hole 10µm	each	204,40
AGA0232P	Platinum Aperture 2mm dia x 0.25mm Hole 250µm	each	96,70
AGA0233P	Platinum Aperture 2mm dia x 0.25mm Hole 25µm	each	108,40
AGA0234P	Platinum Aperture 2mm dia x 0.25mm Hole 30µm	each	111,10
AGA0235P	Platinum Aperture 2mm dia x 0.25mm Hole 500µm	each	94,20
AGA0236P	Platinum Aperture 2mm dia x 0.25mm Hole 600µm	each	96,70
AGA0237P	Platinum Aperture 2mm dia x 0.25mm Hole 750µm	each	96,70
AGA0238P	Platinum Aperture 2mm dia x 0.25mm Hole 1000µm	each	96,70
AGA0239P	Platinum Aperture 2mm dia x 0.25mm Hole 1250µm	each	94,40
AGA0300P	Platinum Aperture 3.04mm dia x 0.25mm Hole 5µm	each	331,10
AGA0301M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 10µm	each	210,40
AGA0301P	Platinum Aperture 3.04mm dia x 0.25mm Hole 10µm	each	255,70
AGA0302M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 20µm	each	119,50
AGA0302P	Platinum Aperture 3.04mm dia x 0.25mm Hole 20µm	each	119,50
AGA0303M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 25µm	each	112,70
AGA0303P	Platinum Aperture 3.04mm dia x 0.25mm Hole 25µm	each	116,30
AGA0304M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 30µm	each	116,30
AGA0304P	Platinum Aperture 3.04mm dia x 0.25mm Hole 30µm	each	112,70
AGA0305M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 40µm	each	106,50
AGA0305P	Platinum Aperture 3.04mm dia x 0.25mm Hole 40µm	each	106,50
AGA0306M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 50µm	each	102,60
AGA0306P	Platinum Aperture 3.04mm dia x 0.25mm Hole 50µm	each	103,20
AGA0307M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 60µm	each	99,60
AGA0307P	Platinum Aperture 3.04mm dia x 0.25mm Hole 60µm	each	103,20
AGA0308M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 70µm	each	102,60
AGA0308P	Platinum Aperture 3.04mm dia x 0.25mm Hole 70µm	each	103,20
AGA0309M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 100µm	each	97,60
AGA0309P	Platinum Aperture 3.04mm dia x 0.25mm Hole 100µm	each	97,60
AGA0310M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 150µm	each	97,60
AGA0310P	Platinum Aperture 3.04mm dia x 0.25mm Hole 150µm	each	122,80
AGA0311M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 200µm	each	121,30
AGA0311P	Platinum Aperture 3.04mm dia x 0.25mm Hole 200µm	each	97,90
AGA0312M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 250µm	each	97,60
AGA0312P	Platinum Aperture 3.04mm dia x 0.25mm Hole 250µm	each	97,60
AGA0313M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 300µm	each	101,10
AGA0313P	Platinum Aperture 3.04mm dia x 0.25mm Hole 300µm	each	97,60
AGA0314M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 400µm	each	118,40
AGA0314P	Platinum Aperture 3.04mm dia x 0.25mm Hole 400µm	each	97,60

Part #	Description	Unit	Price €
AGA0315M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 500um	each	96,10
AGA0315P	Platinum Aperture 3.04mm dia x 0.25mm Hole 500um	each	100,20
AGA0316M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 600µm	each	97,60
AGA0316P	Platinum Aperture 3.04mm dia x 0.25mm Hole 600µm	each	97,60
AGA0317M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 750µm	each	97,60
AGA0317P	Platinum Aperture 3.04mm dia x 0.25mm Hole 750µm	each	118,40
AGA0318M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 1000µm	each	97,60
AGA0318P	Platinum Aperture 3.04mm dia x 0.25mm Hole 1000µm	each	97,60
AGA0319M	Molybdenum Aperture 3.04mm dia x 0.25mm Hole 1250µm	each	103,50
AGA0319P	Platinum Aperture 3.04mm dia x 0.25mm Hole 1250µm	each	90,00
AGA0349	Top Hat Aperture Platinum/Iridium 3mm dia 20µm Hole	each	457,10
AGA0350	Top Hat Aperture Platinum/Iridium 3mm dia 30µm Hole	each	453,00
AGA0351	Top Hat Aperture Platinum/Iridium 3mm dia 50µm Hole	each	423,30
AGA0352	Top Hat Aperture Platinum/Iridium 3mm dia 70µm Hole	each	331,10
AGA0353	Top Hat Aperture Platinum/Iridium 3mm dia 100µm Hole	each	360,90
AGA0354	Top Hat Aperture Platinum/Iridium 3mm dia 150µm Hole	each	297,80
AGA0355	Top Hat Aperture Platinum/Iridium 3mm dia 200µm Hole	each	297,80
AGA0356	Top Hat Aperture Platinum/Iridium 3mm dia 300µm Hole	each	329,00
AGA0357	Top Hat Aperture Platinum/Iridium 3mm dia 500um Hole	each	360,90
AGA0358	Top Hat Aperture Platinum/Iridium 3mm dia 750µm Hole	each	297,80
AGA0359	Top Hat Aperture Platinum/Iridium 3mm dia 1000µm Hole	each	297,80
AGA0400P	Platinum Aperture 4mm dia x 0.2mm Hole 5µm	each	395,50
AGA0401M	Molybdenum Aperture 4mm dia x 0.2mm Hole 10µm	each	180,50
AGA0401P	Platinum Aperture 4mm dia x 0.2mm Hole 10µm	each	180,50
AGA0402M	Molybdenum Aperture 4mm dia x 0.2mm Hole 20µm	each	171,00
AGA0402P	Platinum Aperture 4mm dia x 0.2mm Hole 20µm	each	171,00
AGA0403M	Molybdenum Aperture 4mm dia x 0.2mm Hole 25µm	each	143,10
AGA0403P	Platinum Aperture 4mm dia x 0.2mm Hole 25µm	each	143,10
AGA0404M	Molybdenum Aperture 4mm dia x 0.2mm Hole 30µm	each	143,10
AGA0404P	Platinum Aperture 4mm dia x 0.2mm Hole 30µm	each	137,90
AGA0405M	Molybdenum Aperture 4mm dia x 0.2mm Hole 40µm	each	143,10
AGA0405P	Platinum Aperture 4mm dia x 0.2mm Hole 40µm	each	137,90
AGA0406M	Molybdenum Aperture 4mm dia x 0.2mm Hole 50µm	each	123,60
AGA0406P	Platinum Aperture 4mm dia x 0.2mm Hole 50µm	each	123,60
AGA0407M	Molybdenum Aperture 4mm dia x 0.2mm Hole 60µm	each	123,60
AGA0407P	Platinum Aperture 4mm dia x 0.2mm Hole 60µm	each	118,40
AGA0408M	Molybdenum Aperture 4mm dia x 0.2mm Hole 70µm	each	123,60
AGA0408P	Platinum Aperture 4mm dia x 0.2mm Hole 70µm	each	123,60
AGA0409M	Molybdenum Aperture 4mm dia x 0.2mm Hole 100µm	each	104,80
AGA0409P	Platinum Aperture 4mm dia x 0.2mm Hole 100µm	each	105,30
AGA0410M	Molybdenum Aperture 4mm dia x 0.2mm Hole 150µm	each	104,80
AGA0410P	Platinum Aperture 4mm dia x 0.2mm Hole 150µm	each	105,30
AGA0411M	Molybdenum Aperture 4mm dia x 0.2mm Hole 200µm	each	104,80
AGA0411P	Platinum Aperture 4mm dia x 0.2mm Hole 200µm	each	138,00
AGA0412M	Molybdenum Aperture 4mm dia x 0.2mm Hole 250µm	each	104,80
AGA0412P	Platinum Aperture 4mm dia x 0.2mm Hole 250µm	each	104,80
AGA0413M	Molybdenum Aperture 4mm dia x 0.2mm Hole 300µm	each	119,10
AGA0413P	Platinum Aperture 4mm dia x 0.2mm Hole 300µm	each	104,80
AGA0414M	Molybdenum Aperture 4mm dia x 0.2mm Hole 400µm	each	104,80
AGA0414P	Platinum Aperture 4mm dia x 0.2mm Hole 400µm	each	104,80
AGA0415M	Molybdenum Aperture 4mm dia x 0.2mm Hole 500um	each	106,50
AGA0415P	Platinum Aperture 4mm dia x 0.2mm Hole 500um	each	106,50
AGA0416M	Molybdenum Aperture 4mm dia x 0.2mm Hole 600µm	each	102,10
AGA0416P	Platinum Aperture 4mm dia x 0.2mm Hole 600µm	each	106,50
AGA0417M	Molybdenum Aperture 4mm dia x 0.2mm Hole 750µm	each	106,50
AGA0417P	Platinum Aperture 4mm dia x 0.2mm Hole 750µm	each	106,50
AGA0418M	Molybdenum Aperture 4mm dia x 0.2mm Hole 1000µm	each	133,80
AGA0418P	Platinum Aperture 4mm dia x 0.2mm Hole 1000µm	each	112,70

Part #	Description	Unit	Price €
AGA0419M	Molybdenum Aperture 4mm dia x 0.2mm Hole 1250µm	each	112,70
AGA0419P	Platinum Aperture 4mm dia x 0.2mm Hole 1250µm	each	112,70
AGA042	Filaments for Hitachi H/S/X Series	10	663,20
AGA0420P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 5µm	each	426,30
AGA0421P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 10µm	each	297,10
AGA0422P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 20µm	each	241,30
AGA0423P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 25µm	each	185,40
AGA0424P	Platinum/Iridium Aperture 4mm dia x 0.6mm Hole 30µm	each	162,10
AGA0429P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 100µm	each	144,30
AGA042C	Cartridge Filaments for Hitachi	10	1.353,70
AGA0430P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 150µm	each	135,60
AGA0431P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 200µm	each	144,30
AGA0432P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 250µm	each	135,60
AGA0433P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 300µm	each	135,60
AGA0434P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 400µm	each	144,30
AGA0438P	Platinum/Iridium Aperture 4mm dia x 0.6mm Thick Hole 1000µm	each	152,30
AGA043X	Filaments for Zeiss	10	379,40
AGA045	Filaments for JEOL GC Type (Box of 12)	12	371,30
AGA047	Agar Filaments for TESCAN VEGA series	10	379,40
AGA047T	Agar Filaments for TESCAN VEGA series	10	#N/B
AGA053	Refilamenting service - Siemens Type	each	#N/B
AGA054	Agar Filaments for Cambridge/LEO/Zeiss & AEI	10	#N/B
AGA054A	Agar Tungsten/Rhenium Filaments	10	#N/B
AGA054L	Agar Filaments for LEO Microscope 0.15mm wire	10	175,10
AGA054X	Agar Filaments for Cambridge/LEO/Zeiss & AEI	10	119,80
AGA0601M	Molybdenum Aperture 6.35mm dia x 0.125mm Hole 20µm	each	246,50
AGA0601P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 20µm	each	246,50
AGA0606M	Molybdenum Aperture 6.35mm dia x 0.125mm Hole 50µm	each	174,70
AGA0606P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 50µm	each	178,70
AGA0607M	Molybdenum Aperture 6.35mm dia x 0.125mm Hole 60µm	each	174,70
AGA0607P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 60µm	each	178,70
AGA0608M	Molybdenum Aperture 6.35mm dia x 0.125mm Hole 70µm	each	166,30
AGA0608P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 70µm	each	166,30
AGA0609M	Molybdenum Aperture 6.35mm dia X0.125mm Hole 100µm	each	149,60
AGA0609P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 100µm	each	149,60
AGA0610M	Molybdenum Aperture 6.35mm dia X0.125mm Hole 150µm	each	149,60
AGA0610P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 150µm	each	153,30
AGA0611M	Molybdenum Aperture 6.35mm dia X0.125mm Hole 200µm	each	149,60
AGA0611P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 200µm	each	153,30
AGA0613M	Molybdenum Aperture 6.35mm dia X0.125mm Hole 300µm	each	149,60
AGA0613P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 300µm	each	153,30
AGA0615M	Molybdenum Aperture 6.35mm dia x 0.02mm Hole 500um	each	140,50
AGA0618M	Molybdenum Aperture 6.35mm dia X0.125mm Hole 1000µm	each	135,30
AGA0618P	Platinum/Iridium Aperture 6.35mm dia x 0.125mm Hole 1000µm	each	138,60
AGA062	Refilamenting service - Hitachi (H, S and X Series)	each	76,00
AGA064	Refilamenting service - ISI/ABT/Topcon	each	90,10
AGA0690-6006	LaB6 Filament 60/06 on FEI/Philips Base (not XL30)	each	1.438,10
AGA0690-6006	LaB6 Filament 60/06 on FEI/Philips Base (not XL30)	each	1.438,10
AGA0690-9015	LaB6 Filament 90/15 on FEI/Philips Base (not XL30)	each	1.718,90
AGA0690-9020	LaB6 Filament 90/20 on FEI/Philips Base (not XL30)	each	1.438,10
AGA0691-9020	LaB6 Filament 90/20 on FEI XL30 Base	each	1.438,10
AGA0692-6006	LaB6 Filament 60/06 on AEI Base for Camscan, Cameca, Electro	each	1.276,80
AGA0692-9015	LaB6 Filament 90/15 on AEI Base for Camscan, Cameca, Electro	each	1.701,50
AGA0692-9020	LaB6 Filament 90/20 on AEI Base for Camscan, Cameca, Electro	each	1.660,80
AGA0693-6006	LaB6 Filament 60/06 for JEOL K Base	each	1.324,20
AGA0693-9015	LaB6 Filament 90/15 on JEOL K Base	each	1.763,70
AGA0693-9020	LaB6 Filament 90/20 on JEOL K Base	each	2.091,20
AGA0694-6006	LaB6 Filament 60/06 on Zeiss Base (please specify make & mod	each	1.083,60

Part #	Description	Unit	Price €
AGA0694-9015	LaB6 Filament 90/15 on Zeiss Base (please specify make & mod	each	1.210,20
AGA0694-9020	LaB6 Filament 90/20 on Zeiss Base (please specify make & mod	each	1.691,90
AGA0695-6006	LaB6 Filament 60/06 on ISI/ABT/Topcon two pin Base	each	1.083,60
AGA0695-9015	LaB6 Filament 90/15 on ISI/ABT/Topcon two pin Base	each	1.083,60
AGA0695-9020	LaB6 Filament 90/20 on ISI/ABT/Topcon two pin Base	each	1.083,60
AGA0698-6006	LaB6 Filament 60/06 on Hitachi S Base	each	1.083,60
AGA0698-9015	LaB6 Filament 90/15 on Hitachi S Base	each	1.302,30
AGA0698-9020	LaB6 Filament 90/20 on Hitachi S Base	each	1.083,60
AGA07201	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 100µm	each	147,60
AGA07202	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 200µm	each	147,60
AGA07210	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 10µm	each	173,70
AGA07220	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 20µm	each	173,70
AGA07230	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 30µm	each	173,70
AGA07240	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 40µm	each	147,60
AGA07250	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 50µm	each	147,60
AGA07260	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 60µm	each	147,60
AGA07270	Thin Foil Aperture 2mm dia, 0.25mm thick, Hole 70µm	each	147,60
AGA07301	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 100µm	each	147,60
AGA07301	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 100µm	each	147,60
AGA07302	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 200µm	each	147,60
AGA07302	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 200µm	each	147,60
AGA07310	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 10µm	each	173,70
AGA07320	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 20µm	each	173,70
AGA07330	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 30µm	each	173,70
AGA07340	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 40µm	each	147,60
AGA07350	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 50µm	each	147,60
AGA07360	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 60µm	each	147,60
AGA07370	Thin Foil Aperture 3.04mm dia x 0.25mm thick, Hole 70µm	each	147,60
AGA0801	Gilded Copper Aperture Strip 5 30 40 60µm Holes	each	o.a.
AGA0802	Gilded Copper Aperture Strip 50 50 30 30µm Holes	each	o.a.
AGA0803	Gilded Copper Aperture Strip 5 20 50 100µm Holes	each	o.a.
AGA080G	Gilded Copper Aperture Strip 5 10 20 55µm Holes	each	o.a.
AGA086	Agar Filaments for FEI	10	296,60
AGA089G	Gilded Copper Aperture Strip 5 10 40µm Holes	each	o.a.
AGA0902M	Molybdenum Aperture 10mm dia x 0.1mm Hole 20µm	each	382,20
AGA0902P	Platinum Aperture 10mm dia x 0.1mm Hole 20µm	each	323,50
AGA0904P	Platinum Aperture 10mm dia x 0.1mm Hole 30µm	each	365,60
AGA0906M	Molybdenum Aperture 10mm dia x 0.1mm Hole 50µm	each	244,10
AGA0906P	Platinum Aperture 10mm dia x 0.1mm Hole 50µm	each	206,80
AGA0907P	Platinum Aperture 10mm dia x 0.1mm Hole 80µm	each	194,50
AGA0908M	Molybdenum Aperture 10mm dia x 0.1mm Hole 70µm	each	247,10
AGA0908P	Platinum Aperture 10mm dia x 0.1mm Hole 70µm	each	206,80
AGA0909M	Molybdenum Aperture 10mm dia x 0.1mm Hole 100µm	each	235,70
AGA0909P	Platinum Aperture 10mm dia x 0.1mm Hole 100µm	each	#N/B
AGA0910M	Molybdenum Aperture 10mm dia x 0.1mm Hole 150µm	each	196,00
AGA0910P	Platinum Aperture 10mm dia x 0.1mm Hole 150µm	each	174,70
AGA0911M	Molybdenum Aperture 10mm dia x 0.1mm Hole 200µm	each	196,00
AGA0911P	Platinum Aperture 10mm dia x 0.1mm Hole 200µm	each	180,20
AGA0913M	Molybdenum Aperture 10mm dia x 0.1mm Hole 300µm	each	196,00
AGA0913P	Platinum Aperture 10mm dia x 0.1mm Hole 300µm	each	174,70
AGA0914M	Molybdenum Aperture 10mm dia x 0.1mm Hole 400µm	each	185,40
AGA0914P	Platinum Aperture 10mm dia x 0.1mm Hole 400µm	each	174,70
AGA0916M	Molybdenum Aperture 10mm dia x 0.1mm Hole 600µm	each	193,10
AGA0916P	Platinum Aperture 10mm dia x 0.1mm Hole 600µm	each	162,10
AGA0918M	Molybdenum Aperture 10mm dia x 0.1mm Hole 1000µm	each	138,30
AGA0918P	Platinum Aperture 10mm dia x 0.1mm Hole 1000µm	each	135,30
AGA092	Agar Filaments for JEOL (K Type)	6	341,00
AGA0920M	Molybdenum Aperture 10mm dia x 0.1mm Hole 1500µm	each	135,30

Part #	Description	Unit	Price €
AGA095-50	9.53mm Stainless Steel Disc - 50µm Hole.	10	419,70
AGA095-635	9.53mm Stainless Steel Disc - 635µm Hole.	10	419,70
AGA095-76	9.53mm Stainless Steel Disc - 76µm Hole.	10	419,70
AGA09606P	Platinum Aperture 10mm dia x 0.6mm 50µm	each	321,20
AGA09609P	Platinum Aperture 10mm dia x 0.6mm 100µm	each	253,60
AGA09611P	Platinum Aperture 10mm dia x 0.6mm 200µm	each	256,80
AGA09614P	Platinum Aperture 10mm dia x 0.6mm 400µm	each	o.a.
AGA09616P	Platinum Aperture 10mm dia x 0.6mm 600µm	each	236,00
AGA09618P	Platinum Aperture 10mm dia x 0.6mm 1000µm	each	236,00
AGA097	Hitachi Multihole Aperture Strip	each	o.a.
AGA098	ISI/ABT Multihole Aperture Strip	each	o.a.
AGA099	JEOL Multi-Hole Aperture Strip	each	o.a.
AGA1005M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 50µm	each	250,40
AGA1010M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 100µm	each	199,90
AGA1015M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 150µm	each	200,00
AGA1020M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 200µm	each	199,90
AGA1030M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 300µm	each	203,50
AGA1050M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 500µm	each	179,30
AGA1051M	Molybdenum Spray Aperture 10.5 x 0.25mm Hole 500µm	each	177,60
AGA1060T	Tantalum Aperture 12.69mm X 0.25mm	each	171,70
AGA1075M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 75µm	each	251,60
AGA1080M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 800µm	each	171,60
AGA1100M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 1000µm	each	180,40
AGA1101M	Molybdenum Spray Aperture 10.5 x 0.25mm Hole 1mm	each	158,50
AGA1125M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 125µm	each	175,00
AGA1200M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 2000µm	each	180,40
AGA1200T	Tantalum Spray Aperture 12.69mm		153,50
AGA1201M	Molybdenum Spray Aperture 10.5 x 0.25mm Hole 2mm	each	156,40
AGA1202P	Platinum Aperture 12mm dia x 0.1mm Hole 20µm	each	319,20
AGA1203P	Platinum Aperture 12mm dia x 0.1mm Hole 25µm	each	319,20
AGA1204P	Platinum Aperture 12mm dia x 0.1mm Hole 30µm	each	319,20
AGA1205P	Platinum Aperture 12mm dia x 0.1mm Hole 40µm	each	319,20
AGA1206P	Platinum Aperture 12mm dia x 0.1mm Hole 50µm	each	320,80
AGA1207P	Platinum Aperture 12mm dia x 0.1mm Hole 60µm	each	206,80
AGA1208P	Platinum Aperture 12mm dia x 0.1mm Hole 70µm	each	206,80
AGA1209P	Platinum Aperture 12mm dia x 0.1mm Hole 100µm	each	174,70
AGA1210P	Platinum Aperture 12mm dia x 0.1mm Hole 150µm	each	173,10
AGA1211P	Platinum Aperture 12mm dia x 0.1mm Hole 200µm	each	173,10
AGA1212P	Platinum Aperture 12mm dia x 0.1mm Hole 250µm	each	173,10
AGA1213P	Platinum Aperture 12mm dia x 0.1mm Hole 300µm	each	173,10
AGA1214P	Platinum Aperture 12mm dia x 0.1mm Hole 400µm	each	173,10
AGA1217P	Platinum Aperture 12mm dia x 0.1mm Hole 700µm	each	156,40
AGA1218P	Platinum Aperture 12mm dia x 0.1mm Hole 1000µm	each	150,70
AGA1219P	Platinum Aperture 12mm dia x 0.1mm Hole 1200µm	each	150,70
AGA1220M	Molybdenum Aperture 19.9mm x 0.25mm 1000µm Hole	each	178,90
AGA1220P	Platinum Aperture 19.9mm x 0.25mm 1000µm Hole	each	201,10
AGA1220T	Tantalum Aperture 19.9mm x 0.25mm 1000µm Hole	each	192,90
AGA1221M	Molybdenum Aperture 19.9mm x 0.25mm 1200µm Hole	each	178,90
AGA1221P	Platinum Aperture 19.9mm x 0.25mm 1200µm Hole	each	201,10
AGA1221T	Tantalum Aperture 19.9mm x 0.25mm 1200µm Hole	each	201,10
AGA1222M	Molybdenum Aperture 19.9mm x 0.25mm 1500µm Hole	each	190,10
AGA1222P	Platinum Aperture 19.9mm x 0.25mm 1500µm Hole	each	201,10
AGA1222T	Tantalum Aperture 19.9mm x 0.25mm 1500µm Hole	each	201,10
AGA1223M	Molybdenum Aperture 19.9mm x 0.25mm 2000µm Hole	each	203,00
AGA1223P	Platinum Aperture 19.9mm x 0.25mm 2000µm Hole	each	201,10
AGA1223T	Tantalum Aperture 19.9mm x 0.25mm 2000µm Hole	each	201,10
AGA1500M	Molybdenum Spray Aperture 12.68 x 0.25mm Hole 1500µm	each	160,30
AGABT711	Analog Balometer with 2ft x 2ft hood	each	#N/B

Part #	Description	Unit	Price €
AGABT713	Analog Balometer with 16in x 16in hood	each	#N/B
AGAD-22C-00-X4	Black Gel-Box, 50.8x50.8x6.4mm, Gel X4	50	609,00
AGAD-23C-00-X4	Gel-Box 2" x 3" x 3/8" X4 Gel, black.	50	525,70
AG-ADAM-CAL-F1	Adam Calibration Certificate for F1 Weights (per weight)	each	#N/B
AG-ADAM-CAL-F1WS1	ADAM Calibration Certificate for F1 1g-100g Weight Set	each	#N/B
AG-ADAM-CAL-F1WS2	ADAM Calibration Certificate for F1 1g-200g Weight Set	each	#N/B
AG-ADAM-CAL-F1WS3	ADAM Calibration Certificate for F1 1g-2kg Weight Set	each	#N/B
AG-ADAM-CAL-F1WS5	ADAM Calibration Certificate for F1 1g-500g Weight Set	each	#N/B
AG-ADAM-CAL-M1	Adam Calibration Certificate for M1 Weights (per weight)	each	#N/B
AG-ADAM-CAL-M1WS1	ADAM Calibration Certificate for M1 1g-100g Weight Set	each	#N/B
AG-ADAM-CAL-M1WS2	ADAM Calibration Certificate for M1 1g-200g Weight Set	each	#N/B
AG-ADAM-CAL-M1WS3	ADAM Calibration Certificate for M1 1g-2kg Weight Set	each	#N/B
AG-ADAM-CAL-M1WS5	ADAM Calibration Certificate for M1 1g-500g Weight Set	each	#N/B
AGAFL2888001	Panda Flow Verification Collar	each	#N/B
AGAFL7013501	Pitot 8mm x 0.48m	each	#N/B
AGAFL7013502	Pitot 8mm x 0.8m(31.25in)	each	#N/B
AGAFL7013503	Pitot 8mm x 1.0m(39.5in)	each	#N/B
AGAFL7013506	Pitot 8mm x 0.48m(19in) gland	each	#N/B
AGAFL7048101	Pitot Static Tube (4mm x 305mm long)	each	#N/B
AGAFL7048104	Pitot 4mm x 305mm C/W gland (JJ)	each	#N/B
AGAFL71549801	Smoke cap holder	each	#N/B
AGAFL71805301	Telescopic Pitot Assy	each	#N/B
AGAFL7253801	Pitot 9.5 x 1.22m	each	#N/B
AGAFL7253802	Pitot 9.5 x 1.52m	each	#N/B
AGAFL72852201	Air Cone Kit Airflow Instrument	each	#N/B
AGAFL72852204	Air Cone Round only - Airflow	each	#N/B
AGAFL72900002	LCA501 + Rectangular Air Cone Kit	each	#N/B
AGAFL81571601	PVC Tubing grey 2m Length	each	#N/B
AGAFL81571803	PVC blue Tubing 9m length	each	#N/B
AGAFL81571903	PVC red Tubing 9m length	each	#N/B
AGAFL82859201	Thermocouple Probe	each	#N/B
AGAFL82859401	Barbed Reducer 6 x 4	each	#N/B
AGAFL9004167	Smoke pellets	6	#N/B
AGAG3003M	Molybdenum aperture 30mm dia x 0.6mm thick, 30µm hole	each	o.a.
AGAG3003P	Platinum aperture 30mm dia x 0.6mm thick, 30µm hole	each	o.a.
AGAG3005M	Molybdenum aperture 30mm dia x 0.6mm thick, 50µm hole	each	o.a.
AGAG3005P	Platinum aperture 30mm dia x 0.6mm thick, 50µm hole	each	o.a.
AGAG3010M	Molybdenum aperture 30mm dia x 0.6mm thick, 100µm hole	each	o.a.
AGAG3010P	Platinum aperture 30mm dia x 0.6mm thick, 100µm hole	each	o.a.
AGAG3050M	Molybdenum aperture 30mm dia x 0.6mm thick, 500µm hole	each	o.a.
AGAG3050P	Platinum aperture 30mm dia x 0.6mm thick, 500µm hole	each	o.a.
AGAG3100M	Molybdenum aperture 30mm dia x 0.6mm thick, 1000µm hole	each	o.a.
AGAG3100P	Platinum aperture 30mm dia x 0.6mm thick, 1000µm hole	each	o.a.
AGAG32003M	Molybdenum aperture 30mm dia x 0.25mm thick, 30µm hole	each	o.a.
AGAG32005M	Molybdenum aperture 30mm dia x 0.25mm thick, 50µm hole	each	o.a.
AGAG32010M	Molybdenum aperture 30mm dia x 0.25mm thick, 100µm hole	each	o.a.
AGAG32050M	Molybdenum aperture 30mm dia x 0.25mm thick, 500µm hole	each	o.a.
AGAG32100M	Molybdenum aperture 30mm dia x 0.25mm thick, 1000µm hole	each	o.a.
AGAJ200P	Platinum Aperture 2mm dia x 0.1mm Hole 5µm	each	304,00
AGAJ201M	Molybdenum Aperture 2mm dia x 0.1mm Hole 10µm	each	177,80
AGAJ201P	Platinum Aperture 2mm dia x 0.1mm Hole 10µm	each	192,30
AGAJ202M	Molybdenum Aperture 2mm dia x 0.1mm Hole 20µm	each	175,90
AGAJ202P	Platinum Aperture 2mm dia x 0.1mm Hole 20µm	each	164,30
AGAJ203M	Molybdenum Aperture 2mm dia x 0.1mm Hole 25µm	each	158,50
AGAJ203P	Platinum Aperture 2mm dia x 0.1mm Hole 25µm	each	158,50
AGAJ204M	Molybdenum Aperture 2mm dia x 0.1mm Hole 30µm	each	158,50
AGAJ204P	Platinum Aperture 2mm dia x 0.1mm Hole 30µm	each	179,70
AGAJ205M	Molybdenum Aperture 2mm dia x 0.1mm Hole 40µm	each	130,60

Part #	Description	Unit	Price €
AGAJ205P	Platinum Aperture 2mm dia x 0.1mm Hole 40µm	each	135,90
AGAJ206M	Molybdenum Aperture 2mm dia x 0.1mm Hole 50µm	each	135,90
AGAJ206P	Platinum Aperture 2mm dia x 0.1mm Hole 50µm	each	135,90
AGAJ207M	Molybdenum Aperture 2mm dia x 0.1mm Hole 60µm	each	130,60
AGAJ207P	Platinum Aperture 2mm dia x 0.1mm Hole 60µm	each	135,90
AGAJ208M	Molybdenum Aperture 2mm dia x 0.1mm Hole 70µm	each	135,90
AGAJ208P	Platinum Aperture 2mm dia x 0.1mm Hole 70µm	each	150,10
AGAJ209M	Molybdenum Aperture 2mm dia x 0.1mm Hole 100µm	each	116,50
AGAJ209P	Platinum Aperture 2mm dia x 0.1mm Hole 100µm	each	116,50
AGAJ210M	Molybdenum Aperture 2mm dia x 0.1mm Hole 150µm	each	128,90
AGAJ210P	Platinum Aperture 2mm dia x 0.1mm Hole 150µm	each	116,50
AGAJ211M	Molybdenum Aperture 2mm dia x 0.1mm Hole 200µm	each	116,50
AGAJ211P	Platinum Aperture 2mm dia x 0.1mm Hole 200µm	each	116,50
AGAJ212M	Molybdenum Aperture 2mm dia x 0.1mm Hole 250µm	each	116,50
AGAJ212P	Platinum Aperture 2mm dia x 0.1mm Hole 250µm	each	116,50
AGAJ213M	Molybdenum Aperture 2mm dia x 0.1mm Hole 300µm	each	116,50
AGAJ213P	Platinum Aperture 2mm dia x 0.1mm Hole 300µm	each	116,50
AGAJ214M	Molybdenum Aperture 2mm dia x 0.1mm Hole 400µm	each	116,50
AGAJ214P	Platinum Aperture 2mm dia x 0.1mm Hole 400µm	each	116,50
AGAJ215M	Molybdenum Aperture 2mm dia x 0.1mm Hole 500um	each	139,70
AGAJ215P	Platinum Aperture 2mm dia x 0.1mm Hole 500um	each	139,70
AGAJ216M	Molybdenum Aperture 2mm dia x 0.1mm Hole 600µm	each	139,70
AGAJ216P	Platinum Aperture 2mm dia x 0.1mm Hole 600µm	each	139,70
AGAJ217M	Molybdenum Aperture 2mm dia x 0.1mm Hole 750µm	each	139,70
AGAJ217P	Platinum Aperture 2mm dia x 0.1mm Hole 750µm	each	139,70
AGAJ218M	Molybdenum Aperture 2mm dia x 0.1mm Hole 1000µm	each	139,70
AGAJ218P	Platinum Aperture 2mm dia x 0.1mm Hole 1000µm	each	139,70
AGAJ219M	Molybdenum Aperture 2mm dia x 0.1mm Hole 1250µm	each	139,70
AGAJ219P	Platinum Aperture 2mm dia x 0.1mm Hole 1250µm	each	139,70
AGAP500	AirPro Probe Communications Handle	each	643,90
AGAP800	AirPro Micromanometer	each	748,00
AGAPFS-A-AL	AirPro Advanced Features-1Y License	each	127,90
AGAPFS-A-LTL	AirPro Advanced - LT License	each	394,70
AGAPFS-P-AL	AirPro Professional - 1 Y License	each	266,20
AGAPFS-P-LTL	AirPro Professional - LT License	each	664,70
AGAS2004P	Platinum apertures 2mm dia x 0.6mm thick, 40 um square	each	o.a.
AGAS2005P	Platinum apertures 2mm dia x 0.6mm thick, 50µm square	each	o.a.
AGAS2010P	Platinum apertures 2mm dia x 0.6mm thick, 100µm square	each	o.a.
AGAS2050P	Platinum apertures 2mm dia x 0.6mm thick, 500um square	each	o.a.
AGAS2060P	Platinum apertures 2mm dia x 0.6mm thick, 600µm square	each	o.a.
AGAS3004P	Platinum aperture 3.04mm dia x 0.25mm thick, 40µm square	each	449,10
AGAS3005P	Platinum aperture 3.04mm dia x 0.25mm thick, 50µm square	each	449,10
AGAS3010P	Platinum aperture 3.04mm dia x 0.25mm thick, 100µm square	each	449,10
AGAS3050P	Platinum aperture 3.04mm dia x 0.25mm thick, 500um square	each	449,10
AGAS3060P	Platinum aperture 3.04mm dia x 0.25mm thick, 600µm square	each	449,10
AGAV301M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 10µm	each	255,70
AGAV302M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 20µm	each	119,50
AGAV303M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 25µm	each	116,30
AGAV304M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 30µm	each	116,30
AGAV305M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 40µm	each	106,50
AGAV306M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 50µm	each	103,20
AGAV307M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 60µm	each	103,20
AGAV308M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 70µm	each	103,20
AGAV309M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 100µm	each	97,60
AGAV310M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 150µm	each	97,60
AGAV311M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 200µm	each	97,60
AGAV312M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 250µm	each	97,60
AGAV313M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 300µm	each	97,60

Part #	Description	Unit	Price €
AGAV314M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 400µm	each	97,60
AGAV315M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 500µm	each	97,60
AGAV316M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 600µm	each	97,60
AGAV317M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 750µm	each	97,60
AGAV318M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 1000µm	each	97,60
AGAV319M	Molybdenum Aperture Exact 3mm x 0.25mm Hole 1250µm	each	92,80
AGAZ300P	Platinum Aperture 3mm dia x 0.1mm Hole 5µm	each	304,00
AGAZ301P	Platinum Aperture 3mm dia x 0.1mm Hole 10µm	each	198,50
AGAZ302P	Platinum Aperture 3mm dia x 0.1mm Hole 20µm	each	135,30
AGAZ303P	Platinum Aperture 3mm dia x 0.1mm Hole 25µm	each	150,20
AGAZ304P	Platinum Aperture 3mm dia x 0.1mm Hole 30µm	each	140,10
AGAZ305P	Platinum Aperture 3mm dia x 0.1mm Hole 40µm	each	121,70
AGAZ306P	Platinum Aperture 3mm dia x 0.1mm Hole 50µm	each	143,20
AGAZ307P	Platinum Aperture 3mm dia x 0.1mm Hole 60µm	each	121,70
AGAZ308P	Platinum Aperture 3mm dia x 0.1mm Hole 70µm	each	127,10
AGAZ309P	Platinum Aperture 3mm dia x 0.1mm Hole 100µm	each	116,30
AGAZ310P	Platinum Aperture 3mm dia x 0.1mm Hole 150µm	each	116,90
AGAZ311P	Platinum Aperture 3mm dia x 0.1mm Hole 200µm	each	117,10
AGAZ312P	Platinum Aperture 3mm dia x 0.1mm Hole 250µm	each	116,30
AGAZ313P	Platinum Aperture 3mm dia x 0.1mm Hole 300µm	each	116,90
AGAZ314P	Platinum Aperture 3mm dia x 0.1mm Hole 400µm	each	116,90
AGAZ315P	Platinum Aperture 3mm dia x 0.1mm Hole 500µm	each	133,10
AGAZ316P	Platinum Aperture 3mm dia x 0.1mm Hole 600µm	each	133,10
AGAZ317P	Platinum Aperture 3mm dia x 0.1mm Hole 750µm	each	134,00
AGAZ318P	Platinum Aperture 3mm dia x 0.1mm Hole 1000µm	each	116,50
AGAZ319P	Platinum Aperture 3mm dia x 0.1mm Hole 1250µm	each	116,50
AGB10015	Screw Holder for Irregular Samples	each	#N/B
AGB10020	V-Block Holder Small	each	#N/B
AGB10022	V-Block Holder Large	each	#N/B
AGB10024	Double Clamp Holder	each	#N/B
AGB10028	Ceramic Disc Holder with 3 Discs	each	#N/B
AGB10030	Mounting Plate 1 1/2" x 1" (38mm x 25mm)	each	#N/B
AGB10031	Mounting Plate 1 1/2" (38mm) dia.	each	#N/B
AGB10032	Mounting Plate 1" (25mm) dia.	each	#N/B
AGB10034	Mounting Stubs 3/4" (19mm) dia.	5	#N/B
AGB10035	PCB Holder	each	#N/B
AGB10040	Rotational Holder	each	#N/B
AGB10047	Vacuum Chuck Complete with Vacuum Pump	each	#N/B
AGB10132	Extec Spacer for Double Mounting 25.4mm (1")	each	#N/B
AGB10133	Extec Spacer for Double Mounting 31.75mm (1.25")	each	#N/B
AGB10135	Extec Spacer for Double Mounting 38mm (1.5")	each	#N/B
AGB10146	Extec Air Filter	each	#N/B
AGB10148	Extec Water Filter Kit	each	#N/B
AGB10149	Water Filter Cartridge for AGB10148 Water Filter Kit	3	#N/B
AGB10150	Precision Diamond and Abrasive Wheel Saw with Digital Reading	each	#N/B
AGB10153	Splash Safety Shield	each	#N/B
AGB10164	Extec Labpress 40 - Advanced Automatic Mounting Press	each	#N/B
AGB10167	Extec Mould Assembly 25.4mm (1")	each	#N/B
AGB10168	Extec Mould Assembly 31.75mm (1.25")	each	#N/B
AGB10169	Extec Mould Assembly 38mm (1.5")	each	#N/B
AGB10204	Lid for Labpol Duo AGB10248	each	#N/B
AGB10206	Platen Guard 203mm for Labpol Duo	each	#N/B
AGB10211	Paper Holding Band 250mm	each	#N/B
AGB10248	Labpol Duo 8 - Twin Grinding/Polishing Machine	each	#N/B
AGB10249	Paper Holding Band 203mm	each	#N/B
AGB10263	Extec Labpol 8-12S Advanced Convertible Grinder/Polisher (220v)	each	#N/B
AGB10264	Extec Stainless Steel 6 Sample Holder Tear Shape for Use with Extec Semi-Aut	each	#N/B
AGB10265	Semi-Automatic head with central and variable (0-250N) pressure. Sample hol	each	#N/B

Part #	Description	Unit	Price €
AGB10266	Extec Stainless Steel 6 Sample Holder Tear Shape for 10265 Central Pressure S	each	#N/B
AGB10269	Small Levelling Disc 3mm for AGB10264	each	#N/B
AGB10270	EXTEC Dispensing Unit for Extec Labpol 8-12 Advanced 3 unit with Peristaltic P	each	#N/B
AGB10271	High Viscosity Nozzle for AGB10270 Extec Dispensing Unit	each	#N/B
AGB10272	Stirrer for AGB10270 Extec Dispensing Unit	each	#N/B
AGB10331	Labcut 300E Advanced Bench Top Abrasive Cutting Machine	each	#N/B
AGB10338	Extec Large Capacity Vise Set	each	#N/B
AGB10339	Extec Base Cabinet with drawer and manual pocket	each	#N/B
AGB10342	Extec Recirculating Coolant System 16 gal (~60 l) capacity with 100u filter	each	#N/B
AGB10364	Aluminium Oxide Abrasive Cut-Off Wheels for tool steel or stainless steel gen	10	#N/B
AGB10365	Silicon Carbide Abrasive Cut-Off Wheels for non-metallic or non-ferrous mater	10	#N/B
AGB10390	Dressing Stick for Abrasive Cut-Off Wheels	each	#N/B
AGB10428	Paper Holding Band 200mm	each	#N/B
AGB10429	Platen Guard 8" / 200mm	each	#N/B
AGB10431	Aluminum Platen 203mm	each	#N/B
AGB10432	Aluminum Platen 200mm	each	#N/B
AGB10435	Platen Guard 250mm	each	#N/B
AGB10436	Aluminum Platen 250mm	each	#N/B
AGB10437	Aluminum Platen 300mm	each	#N/B
AGB10438	Paper Holding Band 300mm	each	#N/B
AGB10439	Paper Holding Band 305mm	each	#N/B
AGB10440	Aluminum Platen 305mm	each	#N/B
AGB10449	Lid for Advanced Labpol Grinder / Polisher	each	#N/B
AGB10451	Extec Labpol 12 1.0 Advanced Automatic Grinder/Polisher	each	#N/B
AGB10452	Extec Labpol 12 2.0 Advanced Automatic Grinder/Polisher	each	#N/B
AGB10453	Extec Labpol 12 3.0 Advanced Automatic Grinder/Polisher	each	#N/B
AGB10462	Extec Oscillating Head for Central and Combination Pressure	each	#N/B
AGB10463	Large Levelling Disc 3mm for AGB10266	each	#N/B
AGB10464	Extec Electronic Dispensing Unit 4 units w/ Peristaltic Pump	each	#N/B
AGB10502	Aluminium Oxide Abrasive Cut-Off Wheels for hard ferrous materials, 5in(127r	10	#N/B
AGB10508	Aluminium Oxide Abrasive Cut-Off Wheels Precision Series for med-soft ferrou	10	#N/B
AGB10512	Silicon Carbide Abrasive Cut-Off Wheels Precision Series for non-ferrous mate	10	#N/B
AGB10602	Aluminium Oxide Abrasive Cut-Off Wheels Precision Series for hard ferrous ma	10	#N/B
AGB10608	Aluminium Oxide Abrasive Cut-Off Wheels Precision Series for med-soft ferrou	10	#N/B
AGB10612	Silicon Carbide Abrasive Cut-Off Wheels Precision Series for non-ferrous mate	10	#N/B
AGB12043	AC Dressing Stick 1" x 1/2" x 3"	each	#N/B
AGB12045	DW Dressing Stick 1" x 1/2" x 3"	each	#N/B
AGB12050	Universal Cutting Fluid	0,95l	#N/B
AGB12052	Universal Cutting Fluid	3,8l	#N/B
AGB12065	Water Soluble Cutting Fluid	0,95l	#N/B
AGB12067	Water Soluble Cutting Fluid	3,8l	#N/B
AGB4001	NLS-MT vibrating blade micro-slicer - manual	each	#N/B
AGB4002	NLS-AT Vibrating Blade Tissue Micro-Slicer - Auto	each	#N/B
AGB4531N	Microslicer DTK-1000N, 220volts, 50/60Hz	each	#N/B
AGB4532N	Microslicer DTK-Zero1N, 220volts, 50/60Hz	each	#N/B
AGB5001	Cryo-EM Puck Starter Kit	each	#N/B
AGB5001-M	Cryo-EM Locking Puck Starter Kit	each	#N/B
AGB5002-G2	Cryo-EM Pucks, generation 2.0 (12 hole)		#N/B
AGB5003	Cryo-EM Pucks - set of seven (standard, 12 hole)	each	#N/B
AGB5003-G2	Cryo-EM Pucks Gen 2.0- set of seven (standard, 12 hole)		#N/B
AGB5004	Shelved Storage Cane holds 10 CryoEM pucks (HC34 Compatible)	each	#N/B
AGB5005	Shelved Puck Shipping Cane (ALS Style)	each	#N/B
AGB5006	Angled Cryo-Tongs - for puck manipulation	each	#N/B
AGB5007	Protective Storage Case	each	#N/B
AGB5008	Cryo-EM Storage Cane Starter Kit	each	#N/B
AGB5009	Cryo-EM Full Kit - 6 Full Canes	each	#N/B
AGB5114	TRUEScience SMART pH Cap Kit with Electrode and Android Tablet - White SM	each	#N/B
AGB600	Magnetic Slim-line Stirrer Plate, 15-1500rpm	800ml	#N/B

Part #	Description	Unit	Price €
AGB6202	UVC3 Cryo Chamber, 110-240VAC, 50/60Hz + 3x capsule holders	each	#N/B
AGB6203	Pop-up Isolator glovebag	each	#N/B
AGB6299	Auto Polisher SMARTLAM 3.0 Individual Pressure		#N/B
AGB6300	Manual Polisher SMARTLAM 2.0 220V		#N/B
AGB6301	Micro Cutting Machine CUTLAM MICRO 1.1 220V		#N/B
AGB6302	Micro Cutting Machine CUTLAM MICRO 2.0		#N/B
AGB6303	Micro Cutting Machine CUTLAM MICRO 3.0 Y & Z Axis Automatic		#N/B
AGB6304	Manual Cutting Machine CUTLAM 1.1		#N/B
AGB6305	Cutting Machine CUTLAM 3.1 with Y and Z Axis Motorized		#N/B
AGB6305-C	Cabinet for Cutting Machine CUTLAM 3.1		#N/B
AGB6305-L	Positioning Laser for CUTLAM 3.1		#N/B
AGB6306	Cutting Machine CUTLAM 4.0 Z auto & X axis manual		#N/B
AGB6306-L	Positioning Laser for CUTLAM 4.0		#N/B
AGB6307	Cutting machine CUTLAM 4.0 Z auto & Y axis manual		#N/B
AGB6308	Cutting Machine CUTLAM 4.0 Automatic Descent of the Wheel		#N/B
AGB6309	Cutting Machine CUTLAM 5.0 Y and Z axis automatic		#N/B
AGB6309-L	Positioning Laser CUTLAM 5.0		#N/B
AGB6310	Cutting Machine CUTLAM 5.0 X,Y and Z axis automatic		#N/B
AGB6311	Cutting Machine CUTLAM 5.0, 3Axis auto X Long Stoke		#N/B
AGB6312	Automatic Mounting Press PRESSLAM 1.1 D25,4 to 50mm		#N/B
AGB6312-C	Cooling Unit for PRESSLAM 1.1		#N/B
AGB6313	Automatic Polisher MASTERLAM 1.0 Central Pressure		#N/B
AGB6314	Auto Polisher MASTERLAM 1.0 Central Pressure + Oscillating		#N/B
AGB6315	Auto Polisher MASTERLAM 1.1 D400mm Centrale Press + Oscill.		#N/B
AGB6316	Auto Polisher MASTERLAM 2.0 Individual Pressure		#N/B
AGB6317	Auto Polisher MASTERLAM 3.0 Combo Pressure		#N/B
AGB6320	Set of Flange D35mm for Cutting Wheel D75 to 100mm		#N/B
AGB6321	Screw V-Shape Sample Holder max D24mm max		#N/B
AGB6322	Screw V-Shape Sample Holder max D12mm max		#N/B
AGB6323	Double V-Block Screw Holder for Samples D12mm max		#N/B
AGB6324	Screw Holder for Irregular Samples		#N/B
AGB6325	Cabinet for Metallographic Machine		#N/B
AGB6326	Power Supply Cable-Chinese Socket		#N/B
AGB6327	Small Vice, Screw Clamping for Samples D12mm max		#N/B
AGB6329	Set of Flange D75mm for Cutting Wheel up to D200mm		#N/B
AGB6331	Adaptor Ring for Dint 25,4/12,7 PVC		#N/B
AGB6332	Medium Vice, Screw Clamping for Samples D24mm max		#N/B
AGB6334	Set of Flange D50mm for Cutting Wheel D125 to 150mm		#N/B
AGB6335	Cross-Feed table (X-axis) Manual, Stroke 55mm, acc 0,01mm		#N/B
AGB6336	Double Vice (Left & Right) with Screw Tightening System		#N/B
AGB6337	Small Vice Stroke 53mm, Right Side, Fast Clamping		#N/B
AGB6338	Small Vice Stroke 53mm, Left Side, Fast Clamping		#N/B
AGB6339	Universal Vice 54 (Vertical Clamping)		#N/B
AGB6340	Cross-Feed table (X-axis) Motorised, Stroke 50mm, acc 0,01		#N/B
AGB6341	Positioning Laser for CUTLAM Micro 3.0		#N/B
AGB6342	Adaptor Ring for ID 22/12,7 Steel		#N/B
AGB6343	Adaptor Ring for ID 22/12,7 Plastic		#N/B
AGB6345	Double Vice (Left & Right) with Fast Tightening System		#N/B
AGB6346	Universal Vice 61 (Vertical Clamping)		#N/B
AGB6347	Cross-Feed table (X-axis) Manual, Stroke 100mm, acc 0,01mm		#N/B
AGB6348	Simple Vice for Cross Table 60CT20520 (Fast Clamping)		#N/B
AGB6349	Adaptor Ring for ID32/25,4 Steel		#N/B
AGB6350	Adaptor Ring for ID32/25,4 plastic		#N/B
AGB6354	Simple Vice (Left Side) for Cross-feed Table of CUTLAM 3.1		#N/B
AGB6355	Manual Cross-feed Table for CUTLAM 3.1		#N/B
AGB6356	Automatic Cross-feed table for CUTLAM 3.1 (X Axis)		#N/B
AGB6359	Double Vice with Fast Tightening System High Capacity		#N/B
AGB6360	Filtration/Regulation Kit for Compressed Air - 10µm		#N/B

Part #	Description	Unit	Price €
AGB6364	Option Cutting Wheel D400mm for CUTLAM 5.0		#N/B
AGB6365	UV Light in Cooling Tank, to avoid Moistures		#N/B
AGB6366	Stirring System in Cooling Tank, to avoid Moistures		#N/B
AGB6372	Heating/Cooling Unit for Mould Assembly D25,4 to 40mm		#N/B
AGB6373	Heating/Cooling Unit for Mould Assembly D50mm et 2p		#N/B
AGB6374	Mould Assembly D25,4mm For Single Mounting		#N/B
AGB6375	Mould Assembly D30mm		#N/B
AGB6376	Mould Assembly D1-1/4p		#N/B
AGB6377	Mould Assembly D1-1/2p		#N/B
AGB6378	Mould Assembly D40mm		#N/B
AGB6379	Mould Assembly D50mm		#N/B
AGB6380	Mould Assembly D2p		#N/B
AGB6381	Spacer D25,4mm		#N/B
AGB6382	Spacer D30mm		#N/B
AGB6383	Spacer D40mm		#N/B
AGB6384	Spacer D50mm		#N/B
AGB6385	Spacer D1-1/4p		#N/B
AGB6386	Spacer D1-1/2p		#N/B
AGB6387	Spacer D2p		#N/B
AGB6389	Automatic Dispenser DISTILAM - 4 Pumps + Plate Cleaning		#N/B
AGB6390	Recirculation/Sedimentation Tank with Pump		#N/B
AGB6392	Cover for SMARTLAM 2.0 or MASTERLAM 1.0, 2.0 & 3.0		#N/B
AGB6393	Bowl reducer for Plate D250mm		#N/B
AGB6394	Aluminium Working Plate D250mm		#N/B
AGB6395	Aluminium Working Plate D300mm		#N/B
AGB6396	Aluminium Working Plate D10p		#N/B
AGB6397	Aluminium Working Plate D12p		#N/B
AGB6399	Sample Holder CP D119mm, 6 Tear Shape Cavities D20 to 32mm		#N/B
AGB6400	Sample Holder CP D160mm, 6 Tear Shape Cavities D20-40mm		#N/B
AGB6401	Sample Holder CP D160mm with 9 Cavities D30mm		#N/B
AGB6402	Sample Holder CP D160mm with 6 Cavities D40mm		#N/B
AGB6403	Sample Holder CP D160mm with 5 Cavities D50mm		#N/B
AGB6404	Sample Holder CP D160mm with 3Rectangular Cavities 40x70mm		#N/B
AGB6405	Sample Holder CP D160mm with 6Rectangular Cavities 25x34mm		#N/B
AGB6406	Levelling Disc 3mm for Sample Holder D119mm		#N/B
AGB6407	Levelling Disc 3mm for Sample Holder D160mm		#N/B
AGB6408	Sample Holder CP D160mm with 5 Cavities D2p		#N/B
AGB6409	Sample Holder CP D160mm (Without Cavity, for Large Section		#N/B
AGB6410	Universal Sample Holder D110mm with 6 Cavities D28mm Mold		#N/B
AGB6411	Universal Sample Holder D160mm with 6 Cavities D36mm Mold		#N/B
AGB6412	Universal Sample Holder D160mm with 10 Cavities D28mm Mold		#N/B
AGB6413	Universal Sample Holder D145mm,3 Molds (6 Holes D32,35&40mm		#N/B
AGB6414	Universal Sample Holder D145mm, 3 Molds (6 hol. 20,25&30mm)		#N/B
AGB6415	Universal Sample Holder D145mm with 3 Cavities D50mm Mold		#N/B
AGB6416	Universal Sample Holder D145mm, 2 Molds (3 x30x55		#N/B
AGB6418	Clamping Ring (Metal) for Abrasive Papers D250mm		#N/B
AGB6419	Clamping Ring (Metal) for Abrasive Papers D254mm (10p)		#N/B
AGB6420	Clamping Ring (Metal)for Abrasive Papers D300mm		#N/B
AGB6421	Clamping Ring (Metal) for Abrasive Papers D305mm (12p)		#N/B
AGB6422	Transparent Safety Canopy for MASTERLAM 1.0, 2.0, 3.0		#N/B
AGB6423	FAS Disc D250mm SA		#N/B
AGB6424	FAS Disc D300mm SA		#N/B
AGB6425	FAS Disc D250mm MG		#N/B
AGB6426	FAS Disc D300mm MG		#N/B
AGB6427	FMS Fixing System D250mm		#N/B
AGB6428	FMS Fixing System D300mm		#N/B
AGB6429	Splash Guard Device for MASTERLAM and SMARTLAM		#N/B
AGB6432	Transparent Safety Canopy for MASTERLAM 1.1		#N/B

Part #	Description	Unit	Price €
AGB6433	Aluminium Working Plate D400mm		#N/B
AGB6434	Sample Holder CP D210mm (Without Cavity, for Large Section)		#N/B
AGB6436	Sample Holder IP 6 Cavities D2p		#N/B
AGB6437	Sample Holder IP D160mm 6 Cavities D50mm		#N/B
AGB6438	Distance Rings D30mm (6 Units)		#N/B
AGB6439	Distance Rings D32mm (6 Units)		#N/B
AGB6440	Distance Rings D40mm (6 Units)		#N/B
AGB6441	Distance Rings 1p (6 Units)		#N/B
AGB6442	Distance Rings 1 1/4p (6 Units)		#N/B
AGB6443	Distance Rings 1 1/2p (6 Units)		#N/B
AGB6447	SPRINGLAM Central Pressure Head Semi-Auto Polishing		#N/B
AGB6448	Dosing Unit MM909M Connected to the Machine		#N/B
AGB6449	Clamping Ring (Metal) for Abrasive Paper D200mm		#N/B
AGB6450	Clamping Ring (Metal) for Abrasive Papers D203mm (8p)		#N/B
AGB6451	Clamping Ring (Metal) for Abrasive Papers D230mm		#N/B
AGB6453	Aluminium Working Plate D200mm		#N/B
AGB6454	Aluminum Working Plate D230mm		#N/B
AGB6455	Aluminium Working Plate D8p		#N/B
AGB6457	Bowl reducer for Plate D200mm and 8p		#N/B
AGB6458	Bowl reducer for Plate D230mm		#N/B
AGB6459	FAS Disc D200mm SA		#N/B
AGB6460	FAS Disc D200mm MG		#N/B
AGB6461	FMS Fixing System D200mm		#N/B
AGB6462	Arm to Hold Lapping/Polishing Ring		#N/B
AGB6463	PVC Conditioning Ring OD140mm ID106mm		#N/B
AGB6465	Steel Weight 1Kg, with Stainless Treatment D107mm		#N/B
AGB6467	Distance Rings D1p (3 units)		#N/B
AGB6468	Distance Rings D30mm (3 units)		#N/B
AGB6469	Distance Rings D32mm (3 units)		#N/B
AGB6470	Distance Rings D40mm (3 units)		#N/B
AGB6471	Distance Rings D1-1/4p (3 units)		#N/B
AGB6472	Distance Rings 1-1/2p (3 units)		#N/B
AGB6801	Circulaire 550 Non-Ducted Fume Cabinet, 550x600x1121 mm		#N/B
AGB6801-A	Circulaire 550 Non-Ducted Fume Cabinet with Alarm		#N/B
AGB6801-B1	Base Stand (550) *		#N/B
AGB6801-F1	Solvents Filter for Circulaire 550		#N/B
AGB6801-F10	Educational Filter for Circulaire 550		#N/B
AGB6801-F12	Prefilter Replacements (550), Pk10		#N/B
AGB6801-F2	Mineral Acids Filter for Circulaire 550		#N/B
AGB6801-F3	Alkali Filter for Circulaire 550		#N/B
AGB6801-F4	Ammonia Filter for Circulaire 550		#N/B
AGB6801-F5	Cyanide Filter for Circulaire 550		#N/B
AGB6801-F6	Diethyl Ether Filter for Circulaire 550		#N/B
AGB6801-F7	Aldehydes Filter for Circulaire 550		#N/B
AGB6801-F8	Sulphurs Filter for Circulaire 550		#N/B
AGB6801-F9	Custom Filter for Circulaire 550		#N/B
AGB6801-T1	GRP Spillage Tray (550)		#N/B
AGB6802	Circulaire 650 Non-Ducted Fume Cabinet, 650x600x1121 mm		#N/B
AGB6802-A	Agar Scientific 650 Fume Cabinet with Filter Alarm	each	#N/B
AGB6802-A	Circulaire 650 Non-Ducted Fume Cabinet with Alarm		#N/B
AGB6802-B1	Base Stand (650) *		#N/B
AGB6802-F1	Solvents Filter for Circulaire 650		#N/B
AGB6802-F10	Educational Filter for Circulaire 650		#N/B
AGB6802-F12	Prefilter Replacements (650), Pk10		#N/B
AGB6802-F2	Mineral Acids Filter for Circulaire 650		#N/B
AGB6802-F3	Alkali Filter for Circulaire 650		#N/B
AGB6802-F4	Ammonia Filter for Circulaire 650		#N/B
AGB6802-F5	Cyanide Filter for Circulaire 650		#N/B

Part #	Description	Unit	Price €
AGB6802-F6	Diethyl Ether Filter for Circulaire 650		#N/B
AGB6802-F7	Aldehydes Filter for Circulaire 650		#N/B
AGB6802-F8	Sulphurs Filter for Circulaire 650		#N/B
AGB6802-F9	Custom Filter for Circulaire 650		#N/B
AGB6802-T1	GRP Spillage Tray (650)		#N/B
AGB6803	Circulaire 900 Non-Ducted Fume Cabinet, 900x600x1121 mm		#N/B
AGB6803-A	Circulaire 900 Non-Ducted Fume Cabinet with Alarm		#N/B
AGB6803-B1	Base Stand (900) *		#N/B
AGB6803-F1	Solvents Filter for Circulaire 900		#N/B
AGB6803-F10	Educational Filter for Circulaire 900		#N/B
AGB6803-F12	Prefilter Replacements (900), Pk10		#N/B
AGB6803-F2	Mineral Acids Filter for Circulaire 900		#N/B
AGB6803-F3	Alkali Filter for Circulaire 900		#N/B
AGB6803-F4	Ammonia Filter for Circulaire 900		#N/B
AGB6803-F5	Cyanide Filter for Circulaire 900		#N/B
AGB6803-F6	Diethyl Ether Filter for Circulaire 900		#N/B
AGB6803-F7	Aldehydes Filter for Circulaire 900		#N/B
AGB6803-F8	Sulphurs Filter for Circulaire 900		#N/B
AGB6803-F9	Custom Filter for Circulaire 900		#N/B
AGB6803-T1	GRP Spillage Tray (900)		#N/B
AGB701	Glass microspray	each	#N/B
AGB7019	Nickel Bursting Disc for CDP	each	#N/B
AGB702	Embedding Oven, 230 V	each	#N/B
AGB703	Embedding Oven, 110 V	each	#N/B
AGB705	Vacuum Oven Lid	each	#N/B
AGB7150	Palladium Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7159	Palladium Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7164	Gold/Palladium Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7165	Platinum Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7166	Platinum/Palladium Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7167	Nickel Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7168	Silver Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7169	Platinum/Palladium Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7180	Gold target replating service	each	#N/B
AGB7184	Gold Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7185	Gold/Palladium Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7186	Platinum Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7187	Nickel Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7188	Silver Foil Target, 60mm dia x 0.1mm, Coater Type 5	each	#N/B
AGB7190	Gold Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7191	Gold/Palladium Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7192	Platinum Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7193	Nickel Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7194	Silver Foil Target, 58mm (nominal) dia x 0.1mm, Coater Type 4	each	#N/B
AGB7195	Gold Foil Target, 40mm dia x 0.1mm, Coater Type 6	each	#N/B
AGB7196	Gold/Palladium Foil Target, 40mm dia x 0.1mm, Coater Type 6	each	#N/B
AGB7197	Gold Foil Target, 20mm dia x 0.1mm, Coater Type 7	each	#N/B
AGB7198	Gold/Palladium Foil Target, 20mm dia x 0.1mm, Coater Type 7	each	#N/B
AGB7199	Gold Disc Target, 54mm dia x 0.1mm, Coater Type 3	each	#N/B
AGB7199-2	Gold Disc Target, 54mm dia x 0.2mm, Coater Type 3	each	#N/B
AGB7221	Sample Table for 6 x 12.5mm pin stubs	each	#N/B
AGB7222	Sample Table for 4 x 10mm dia stubs	each	#N/B
AGB7225	Sample Table for 3 x M4 and 3 x 12.5mm pin stubs	each	#N/B
AGB7226	Sample Table for 4 x 25mm mount	each	#N/B
AGB7227	Sample Table for 1 x 25mm mount	each	#N/B
AGB7228	Sample Table for 1 x 32mm mount	each	#N/B
AGB7229	Sample Table for 4 x 12.5mm stubs	each	#N/B
AGB7230	Agar SEM Turbo Carbon Coater	each	#N/B

Part #	Description	Unit	Price €
AGB7230H3	Sample Holder for 1 x 32mm Mount	each	#N/B
AGB7231	Rotary-Planetary-Tilting (R-P-T) Stage for Turbo Coater	each	#N/B
AGB7232	Agar TEM Turbo Carbon Coater	each	#N/B
AGB7234	High Resolution Sputter Coater	each	#N/B
AGB7235	Adaptor for 32mm mount to 25mm mount	each	#N/B
AGB7238	Glass Chamber 150mm x 150mm for Agar Turbo	each	#N/B
AGB7238R	Replacement Chamber O-rings Ø150mm	2	#N/B
AGB7240	Film Thickness Monitor for Agar SEM Turbo Coater	each	#N/B
AGB7241	Apiezon W40 Wax	250g	#N/B
AGB7242	KJLC 704 Silicone Fluid	500ml	#N/B
AGB7246	Apiezon Oil AP201	4l	#N/B
AGB7247	WACKER Silicone Paste.	90ml	#N/B
AGB7250	Film Thickness Monitor for TEM Turbo Carbon Coater	each	#N/B
AGB7251	Rotary-Tilting (R-T) Stage for SEM/TEM	each	#N/B
AGB7259	QuickStick temporary mounting wax	20	#N/B
AGB7267	Ultrasonic Cleaner	2l	#N/B
AGB7268	Ultrasonic Cleaner	2l	#N/B
AGB7269	Ultrasonic Cleaner	6l	#N/B
AGB7270	Apiezon L Grease	50g	#N/B
AGB7270A	Apiezon L Grease	25g	#N/B
AGB7271	Apiezon H Grease	25g	#N/B
AGB7272	Apiezon M Grease	100g	#N/B
AGB7272-1	Apiezon M Grease	25g	#N/B
AGB7272-2	Apiezon M Grease	4kg	#N/B
AGB7273	Apiezon N Grease	25g	#N/B
AGB7274	Apiezon Q Compound	1kg	#N/B
AGB7275	Apiezon W Wax	500g	#N/B
AGB7275-1	Apiezon W Wax block	1kg	#N/B
AGB7277	Apiezon T Grease	25g	#N/B
AGB7278	Silicone High Vacuum Grease	150g	#N/B
AGB7282	Apiezon W100 Wax	250g	#N/B
AGB7288	Silicone 704 Oil SS	500ml	#N/B
AGB7289	Santovac 5 Oil	100ml	#N/B
AGB7290	Santovac 5 Oil	500ml	#N/B
AGB7292	* High Vacuum Leak Sealant Spray	454g	#N/B
AGB7293	* High Vacuum Leak Sealant Bottle	80g	#N/B
AGB7294A	Pfeiffer A119 Rotary pump oil	1l	#N/B
AGB7299	Fomblin GRRt15 Grease	100g	#N/B
AGB7300KIT	Combined Coater System	each	#N/B
AGB7309	Ultrasonic Cleaner	9l	#N/B
AGB7311	Crystalbond 555 Adhesive Stick	68g	#N/B
AGB7312	Crystalbond 509-3 Adhesive Stick	90g	#N/B
AGB7313	Crystalbond 555-HMP Adhesive Stick	68g	#N/B
AGB7314	Crystalbond 590 Adhesive Stick	227g	#N/B
AGB7316	Wafer-Mount 559, 10in x 10in, AdhesiveSheet.	2	#N/B
AGB7317	Wafer-Mount 562, 8in x 10in, AdhesiveSheet.	2	#N/B
AGB7321	PLANOTEC X-ray beam finder	each	#N/B
AGB7322	* Crystalbond 509-S Stripper	950ml	#N/B
AGB7323	* Crystalbond 590-S Stripper	454g	#N/B
AGB7324	Crystalbond 562-S Stripper	450	#N/B
AGB7340	Agar Sputter Coater - Manual	each	#N/B
AGB7341	Agar Sputter Coater - Automatic	each	#N/B
AGB7341SE	Agar Sputter Coater with larger chamber - Automatic	each	#N/B
AGB7348	Manual Film Thickness Monitor	each	#N/B
AGB7349	Terminating Film Thickness Monitor	each	#N/B
AGB7351	Gold Edwards S150B Target, 73mm disc, Scancoat 6	each	#N/B
AGB7355A	Gold Edwards S150A plated Target	each	#N/B
AGB7356	Gold/Palladium Edwards S150B Target, 73mm disc, Scancoat 6	each	#N/B

Part #	Description	Unit	Price €
AGB7357	Platinum Edwards S150B Target, 73mm disc, Scancoat 6	each	#N/B
AGB7361	Easiglow Glow Discharge Unit	each	#N/B
AGB7361-CRYO	Easiglow Glow Discharge System - Cryo-EM Kit		#N/B
AGB7361-DUO	easiGlow Dual Chamber Glow Discharge System	each	#N/B
AGB7363	GD4 Vacuum Pumping System, 3.6m ³ /hr. Inc connection kit. 230	each	#N/B
AGB7365	Rotary Tilting Stage	each	#N/B
AGB7366	Pumping System for Coater, including connection kit	each	#N/B
AGB7367A	Agar Automatic Carbon Coater	each	#N/B
AGB7367SE	Agar Automatic Carbon Coater 150mm x 165-250mm	each	#N/B
AGB7369	Spare Carbon Source for Coater	each	#N/B
AGB7370	Gold Annular Target, OD: 82mm ID: 60mm x 0.1mm, Coater Type 9	each	#N/B
AGB7371	Gold/Palladium Annular Target, OD: 82mm ID: 60mm x 0.1mm, Coater Type 9	each	#N/B
AGB7372	Platinum Annular Target, OD: 82mm ID: 60mm x 0.1mm, Coater Type 9	each	#N/B
AGB7374	Silver Annular Target, OD: 82mm ID: 60mm x 0.1mm, Coater Type 9	each	#N/B
AGB7375	Spare 'O' ring kit for 120mm dia chambers	each	#N/B
AGB7375A	Kit of 'O' rings for 150mm dia chambers	each	#N/B
AGB7375B	Replacement Springs for Carbon Coater	each	#N/B
AGB7375D	Carbon Rod Slider	each	#N/B
AGB7378	Replacement Glass Chamber for Agar Coater	each	#N/B
AGB7379	Dual Pumping System Upgrade Kit	each	#N/B
AGB7386	Film Thickness Monitor Upgrade Kit for Standard Stage	each	o.a.
AGB7386-RTP	Film Thickness Monitor Upgrade Kit for Rotary-Planetary-Tilt	each	839,40
AGB7390	Gold Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7390-2	Gold Disc Target, 57mm dia x 0.2mm, Coater Type 1	each	o.a.
AGB7391	Gold/Palladium Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7391-2	Gold/Palladium Disc Target, 57mm dia x 0.2mm, Coater Type 1	each	o.a.
AGB7392	Platinum Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7392-2	Platinum Disc Target, 57mm dia x 0.2mm, Coater Type 1	each	o.a.
AGB7393	Nickel Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7394	Silver Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7395	Platinum/Palladium Disc Target, 57mm dia x 0.1mm, Coater Type 1	each	o.a.
AGB7395-2	Platinum/Palladium Disc Target, 57mm dia x 0.2mm, Coater Type 1	each	o.a.
AGB7396	Chromium Disc Target, 57mm dia x 3.2mm, Coater Type 2	each	o.a.
AGB7397	Palladium Disc Target, 57mm dia x 0.1mm Coater Type 1	each	o.a.
AGB7398	Palladium Annular Target, OD: 82mm ID: 60mm x 0.1mm Coater	each	o.a.
AGB7414	Sonorex RK52 Ultrasonic Cleaner 1.8L	each	805,30
AGB7415	Sonorex RK100 Ultrasonic Cleaner 3L	each	853,40
AGB7425	Quartz crystals	3	215,30
AGB7436	SC22 Magnetic Field Cancelling System AC/HH/SEM/FrK	each	o.a.
AGB7441	Enamelled Dewar with Handle and Lid - 2 Litre	each	219,00
AGB7442	Enamelled Dewar with Handle and Lid - 1 Litre	each	233,10
AGB7443	Stainless Steel Dewar 2 litre	each	253,80
AGB7446	Stainless Steel Dewar 1 Litre	each	220,80
AGB7448	Replacement glass Insert for AGB7441 or AGB7443 2 Litre Dewa	each	283,90
AGB7449	Replacement glass Insert for AGB7442, AGB7446 or AGB7589 1 L	each	172,10
AGB7452	Stainless Steel Dewar 1 Litre	each	491,60
AGB7452L	Lid for Stainless Steel Dewar 1 Litre	each	97,00
AGB7453	Wide Neck Stainless Dewar 1 Litre	each	465,20
AGB7453L	Lid for Wide Neck Stainless Dewar 1 Litre	each	93,40
AGB7454	Stainless Steel Dewar 2 Litre	each	691,30
AGB7454L	Lid for Stainless Steel Dewar 2 Litre	each	104,90
AGB7455	Stainless Steel Dewar 3 Litre	each	1.042,30
AGB7455L	Lid for Stainless Steel Dewar 3 Litre	each	244,80
AGB7456	Stainless Steel Dewar 6 Litre	each	991,90
AGB7457	Aluminium Dewar 2 Litre	each	o.a.
AGB7458	Aluminium Dewar 5 Litre	each	o.a.
AGB7465	Aluminium Dewar 10 Litre	each	o.a.
AGB7475	Liquid Nitrogen Carrier 1 Litre	each	355,60

Part #	Description	Unit	Price €
AGB7476	Liquid Nitrogen Carrier 2 Litre	each	423,50
AGB7477	Liquid Nitrogen Carrier 4 Litre	each	359,50
AGB7530	Tote Bag, Bio 2	each	111,00
AGB7531	Cryoguard with integral Data logger & temprature display	each	1.919,40
AGB7532	Cryoguard LN2 low level alarm, Mains and Battery	each	1.355,70
AGB7535	Goblet SINGLE 35mm	each	7,80
AGB7535-6	Goblets 35mm	6	40,80
AGB7536	Goblet SINGLE 65mm	each	10,30
AGB7536-10	Goblets 65mm	10	91,50
AGB7537	Bio 3 Cryo Refridgerator	each	891,20
AGB7540	Cryo Sleeves - Cardboard	100	179,50
AGB7541	Bio 22/36 Lid/Neck Plate	each	102,50
AGB7543	Bio 12 Lid/Neck Plate	each	102,50
AGB7544	Neck Plug, Bio 2	each	123,70
AGB7545	Neck Plug, Bio 3	each	159,30
AGB7546	Neck Plug, Bio 8, 10, 20 & 35	each	159,30
AGB7547	Neck Plug, Bio 12	each	175,80
AGB7548	Neck Plug, Bio 21 & 34	each	199,90
AGB7549	Neck Plug, Bio 22 & 36	each	202,50
AGB7550	Liquid Withdrawal Device, 63mm neck Cryolab 35		1.224,80
AGB7554	Filling Manifold for Cryolab 25	each	882,20
AGB7556	Roller Base for Cryolab 20, 35 & 50	each	525,70
AGB7557	Dermatology Dipper 10 ml Ladle	each	151,40
AGB7558	Lifter for Goblets in the Bio 36 - (other sizes available on request)	each	13,30
AGB7559	Aluminium Cryo Canes - 29cm Long	100	251,50
AGB7560	Cryolab 35 Aluminium Dewar	each	1.348,00
AGB7561	Cryolab 50 Aluminum Dewar	each	1.597,80
AGB7562	Cryolab 25 Dispenser Shelf for 25 Litre	each	181,20
AGB7563	Neck Plug, Cryolab 2	each	123,70
AGB7564	Neck Plug, Cryolab 3-20	each	108,60
AGB7565	Neck Plug, Cryolab 25	each	111,40
AGB7566	Neck Plug, Cryolab 35 & 50	each	174,20
AGB7567	Liquid Withdrawal Device for Cryolab 10	each	1.091,70
AGB7568	Liquid Withdrawal Device, 51mm neck Cryolab 20,25	each	1.184,90
AGB7569	Liquid Withdrawal Device, 63mm neck Cryolab 50	each	1.224,80
AGB7570	Bio 3 Cryo Refridgerator	each	988,70
AGB7571	Bio 8 Cryo Refridgerator	each	1.169,60
AGB7572	Bio 10 Cryo Refridgerator	each	1.215,90
AGB7573	Bio 12 Cryo Refridgerator	each	1.383,90
AGB7574	Bio 20 Cryo Refridgerator	each	1.364,50
AGB7575	Bio 21 Cryo Refridgerator	each	1.774,20
AGB7576	Bio 34 Cryo Refridgerator	each	1.861,40
AGB7577	Bio 35 Cryo Refridgerator	each	1.693,90
AGB7578	Bio 36 Cryo Refridgerator	each	2.128,00
AGB7581	Bio22 Cryo Refrigerator	each	1.774,20
AGB7582	Roller Base for 34,35 & 36L refrigerators	each	555,90
AGB7584	Replacement glass Insert for AGB7586 7 Litre Dewar	each	926,90
AGB7585	Stainless Steel Dewar 4.5 Litre	each	324,60
AGB7586	Stainless Steel Dewar 7 Litre	each	1.015,10
AGB7587	Replacement glass Insert for AGB7585 or AGB7598 4.5 Litre De	each	335,50
AGB7588	Replacement glass Insert for AGB7599 10 Litre Dewar	each	1.011,20
AGB7590	Enamelled Dewar - 0.2 Litres	each	177,50
AGB7591	Enamelled Dewar - 0.7 Litres	each	187,40
AGB7592	Enamelled Dewar - 0.75 Litre	each	181,10
AGB7593	Enamelled Dewar - 1.1 Litres	1	168,80
AGB7593R	Replacement Glass Insert for Enamelled Dewar - 1.1 Litres	each	187,30
AGB7594	Enamelled Dewar - 1.5 Litres	each	203,10
AGB7596	Enamelled Dewar - 2 Litres	each	219,00

Part #	Description	Unit	Price €
AGB7597	Enamelled Dewar - 3 Litres	each	228,40
AGB7598	Enamelled Dewar - with Handle and Lid 4.5 Litres	each	279,50
AGB7599	Enamelled Dewar - with Handle and Lid 10 Litres	each	693,60
AGB7606	Compact Embedding Oven - 6 litre	each	1.335,30
AGB7608	General Purpose Oven - 18 litre	each	1.383,00
AGB7610	Fomblin Y VAC 06/6 Pump Oil (532ml)	each	412,60
AGB7611	Fomblin Y VAC 06/6 Pump Oil (2660ml)	each	o.a.
AGB7634	Edwards L9 Diffusion Pump Oil (1 Litre)	1l	1.039,70
AGB7638	Lion S Oil (1 Litre)	1l	647,30
AGB7639	Lion S Oil	100ml	#N/B
AGB7640	Ultragrade 19 Rotary Pump Oil	4l	113,70
AGB7641	Ultragrade 19 Rotary Pump Oil	1l	57,60
AGB7642	Ultragrade 20 Rotary Pump Oil	4l	156,20
AGB7643	Ultragrade 20 Rotary Pump Oil	1l	57,50
AGB7644	Ultragrade 15 Rotary Pump Oil	4l	156,90
AGB7645	Ultragrade 15 Rotary Pump Oil	1l	56,00
AGB7646	Pfeiffer P3 Rotary Pump Oil	500ml	48,00
AGB7650	Plano P47 Scintillator Disc 7.7mm	each	146,30
AGB7651	Plano P47 Scintillator Disc 13.8mm	each	165,00
AGB7652	Plano P47 Scintillator Disc 10mm	each	146,30
AGB7653	Plano P47 Scintillator Disc 12mm	each	146,30
AGB7654	Plano P47 Scintillator Disc 19mm	each	146,30
AGB7655	Plano P47 Scintillator Disc 18mm	each	146,30
AGB7656	Plano P47 Scintillator Disc 19.8mm x 1mm	each	167,50
AGB7657	Plano P47 Scintillator Disc 6.2mm	each	#N/B
AGB7658	Plano P47 Scintillator Disc 19mm x 7.3mm	each	146,30
AGB7659	Plano P47 Scintillator Disc 20mm	each	146,30
AGB7661	Plano P47 Scintillator Disc 9mm x 3mm	each	166,10
AGB7662	Plano P47 Scintillator Disc 8.8mm x 2mm	each	170,90
AGB7663	Plano P47 Scintillator Disc 16.4mm x 0.17mm	each	146,30
AGB7665	Plano P47 Scintillator Disc 12mm x 0.17mm	each	146,30
AGB7666	Plano P47 Scintillator Disc 6mm	each	146,30
AGB7667	Plano P47 Scintillator Disc 12.4mm x 3.2mm	each	146,30
AGB7668	Plano P47 Scintillator Disc 20mm x 2mm	each	173,10
AGB7680	YAG Scintillator Disc 8.8mm x 2mm	each	983,30
AGB7681	YAG Scintillator Disc 10mm x 1mm	each	811,90
AGB7682	YAG Scintillator Disc 13.8mm x 1mm	each	916,40
AGB7683	YAG Scintillator Disc 18mm x 1mm	each	811,90
AGB7684	YAG Scintillator Disc 20mm x 1mm	each	1.255,10
AGB7732	6MHz Quartz Crystals	10	203,50
AGB7736	Dual Pumping System for Sputter/Carbon Coater	each	o.a.
AGB7753	Serial Section Grids Copper	50	64,10
AGB7753N	Serial Section Grids Nickel	50	66,60
AGB776	Liquid Nitrogen Dewar 20 Litre	each	o.a.
AGB777	Roller Base for 20 Litre Liquid Nitrogen Dewar	each	o.a.
AGB778A	Liquid Nitrogen Dewar 25 Litre	each	o.a.
AGB7830	Slicemaster HS-1 vertical slicer	each	o.a.
AGB7831	Slicemaster HK-1 angled slicer	each	o.a.
AGB7832	Slicemaster HW-1 multi-angle slicer	each	o.a.
AGB7833	Standard sample holder for HS-1 and HW-1	each	403,20
AGB7834	Sample holder for thick samples, sample size: up to 8 mm for HW-1	each	546,30
AGB7835	Sample tray sheet (PTFE), for HS-1/HW-1, Sample holder	10	182,10
AGB7836	Spare blade, 0.245 mm, for Slice Master	100	308,20
AGB7837	Blade holder with spare double-edged blade for HS-1/HW-1	10	252,10
AGB7838	Spare double-edged blade, 0.1mm 100 pcs	100	294,20
AGB7839	Jig for spare blade, replacement	each	96,10
AGB7882	Platinum Annular Target, OD: 115mm ID: 92mm x 0.1mm, Coater Type 10	each	o.a.
AGB7883	Nickel Annular Target, OD: 115mm ID: 92mm x 0.1mm, Coater Type 10	each	o.a.

Part #	Description	Unit	Price €
AGB7884	Silver Annular Target, OD: 115mm ID: 92mm x 0.1mm, Coater Type 10	each	o.a.
AGB7885	Gold Foil Target, 40mm dia x 0.1mm, Coater Type 8	4	o.a.
AGB7886	Platinum Foil Target, 40mm dia x 0.1mm, Coater Type 8	4	o.a.
AGB7887	Nickel Foil Target, 40mm dia x 0.1mm, Coater Type 8	4	o.a.
AGB7888	Gold/Palladium Foil Target, 40mm dia x 0.1mm, Coater Type 8	4	o.a.
AGB7889	Silver Foil Target, 40mm dia x 0.1mm, Coater Type 8	4	o.a.
AGB7902	Small Sample Holder with Lid 15mm dia x 10mm High	each	44,70
AGB7903	Large Sample Holder with Lid 23mm dia x 18mm High	each	59,00
AGB793	Glass Vials 10ml with Lids	± 190	94,50
AGB794	Glass Vials 15ml with Lids	± 162	96,10
AGB7962	Heat Pen	each	#N/B
AGB7963	Spare Filament for B7962 Heat Pen	each	69,30
AGB797	Roller Base for 25 Litre Stainless Steel Dewar	each	o.a.
AGB798T	Trolley/Pouring Frame for 25 Litre Dewar	each	o.a.
AGB8001	InForce 50 Nanomechanical Actuator	each	2.671,90
AGB8002	Thin Film Method Pack (Hard on Soft)	each	4.453,20
AGB8003	Polymer Method Pack	each	8.040,50
AGB8004	User Method Development for InView	each	12.864,70
AGB8005	InView Simulation Software	each	556,70
AGB8105	Ph Meter	each	#N/B
AGB8106	Spare Electrode for Ph Meter	each	#N/B
AGB8107	Buffer Solution Ph4 (460ml)	each	#N/B
AGB8108	Buffer Solution Ph7 (460ml)	each	#N/B
AGB8109	Calibration Screwdriver	each	#N/B
AGB8168D	Cutter, 3mm dia for Ultrasonic Disc Cutter	each	o.a.
AGB8180	Bakelite Moulding Powder, pink	2,25kg	141,10
AGB8182	Acrylic Moulding Powder, Transparent	2,25kg	164,00
AGB8187	* EpoThin 2 Resin	950ml	190,20
AGB8188	* EpoThin 2 Hardener	470ml	154,00
AGB8198	* Epo-Kwick Resin Kit	each	696,10
AGB8209	Silicon Carbide Paper Discs 203mm P080	each	32,60
AGB8211	Silicon Carbide Paper Discs 203mm P120	25	25,30
AGB8212	Silicon Carbide Paper Discs 203mm P240	25	22,30
AGB8214	Silicon Carbide Paper Discs 203mm P400	25	23,20
AGB8215	Silicon Carbide Paper Discs 203mm P600	25	23,20
AGB8216	Silicon Carbide Paper Discs 203mm P800	25	18,90
AGB8217	Silicon Carbide Paper Discs 203mm P1200	25	44,80
AGB8218	Silicon Carbide Paper Discs 203mm P2500	25	39,00
AGB8219	Silicon Carbide Paper Discs 203mm P4000	25	55,30
AGB8220A	8" Nylon Cloth PSA	10	188,60
AGB8221	Composite Disc Adhesive Backed 203mm Dia.	10	170,70
AGB8222	Napped Cloth Adhesive Backed 203mm dia	10	119,70
AGB8225	Alumina Suspension Medium (1µm)	175ml	76,50
AGB8226	Alumina Suspension Fine (0.3µm)	175ml	76,50
AGB8227	Alumina Suspension Extra Fine (0.05µm)	180ml	76,00
AGB8233	Specimen levelling press	each	476,60
AGB8234	Unidisc Mk2 disc punch	each	o.a.
AGB8300	YAG Scintillator Disc 6mm x 1mm	each	864,70
AGB8301	YAG Scintillator Disc 6.2mm x 1mm	each	811,90
AGB8302	YAG Scintillator Disc 7.7mm x 1mm	each	811,90
AGB8303	YAG Scintillator Disc 9mm x 3mm	each	o.a.
AGB8304	YAG Scintillator Disc 12mm x 1mm	each	811,90
AGB8305	YAG Scintillator Disc 12mm x 0.17mm	each	922,90
AGB8306	YAG Scintillator Disc 12.4mm x 3.2mm	each	893,70
AGB8307	YAG Scintillator Disc 16.4mm x 0.17mm	each	986,30
AGB8308	YAG Scintillator Disc 19mm x 1mm	each	811,90
AGB8309	YAG Scintillator Disc 19mm x 7.3mm	each	1.177,00
AGB8310	YAG Scintillator Disc 19.8mm x 1mm In Metal Ring	each	o.a.

Part #	Description	Unit	Price €
AGB8311	YAG Scintillator Disc 20mm x 2mm	each	866,60
AGB8320	YAP Scintillator Disc 6mm x 1mm	each	946,20
AGB8321	YAP Scintillator Disc 6.2mm x 1mm	each	#N/B
AGB8322	YAP Scintillator Disc 7.7mm x 1mm	each	888,40
AGB8323	YAP Scintillator Disc 8.8mm x 2mm	each	946,20
AGB8324	YAP Scintillator Disc 9mm x 3mm	each	o.a.
AGB8325	YAP Scintillator Disc 10mm x 1mm	each	946,20
AGB8326	YAP Scintillator Disc 12mm x 1mm	each	888,40
AGB8327	YAP Scintillator Disc 12mm x 0.17mm	each	942,60
AGB8328	YAP Scintillator Disc 12.4mm x 3.2mm	each	969,90
AGB8329	YAP Scintillator Disc 13.7mm x 1mm	each	888,40
AGB8330	YAP Scintillator Disc 16.4mm x 0.17mm	each	942,60
AGB8331	YAP Scintillator Disc 18mm x 1mm	each	888,40
AGB8332	YAP Scintillator Disc 19mm x 1mm	each	888,40
AGB8333	YAP Scintillator Disc 19mm x 7.3mm	each	1.550,40
AGB8334	YAP Scintillator Disc 19.8mm x 1mm In Metal Ring	each	o.a.
AGB8335	YAP Scintillator Disc 20mm x 1mm	each	888,40
AGB8336	YAP Scintillator Disc 20mm x 2mm	each	942,60
AGB8359-10	Storage cabinet with 10 shelves	each	346,10
AGB8359-5	Storage cabinet with 5 shelves	each	321,30
AGB8373	127mm Metal Bonded Diamond Wafer Blade High Concentration	each	657,40
AGB8374	152mm Metal Bonded Diamond Wafer Blade High Concentration	each	858,30
AGB8376	127mm Metal Bonded Diamond Wafer Blade Low Concentration	each	645,90
AGB8377	152mm Metal Bonded Diamond Wafer Blade Low Concentration	each	862,60
AGB8378	203mm Metal Bonded Diamond Wafer Blade Low Concentration	each	1.353,00
AGB8379	127mm Metal Bonded Series 20 High Concentration Wafer Blade	each	657,40
AGB8380	178mm Metal Bonded Series 20 High Concentration Wafer Blade	each	1.083,70
AGB8402	NL50 bench-top nanoparticle deposition system		o.a.
AGB8461	Perspex Storage Box	each	o.a.
AGB8462	Wire Saw Rewinding Device	each	o.a.
AGB8463	3-Axis Goniometer Head	each	o.a.
AGB8464	2-Axis Goniometer Head	each	o.a.
AGB8465	1-Axis Goniometer Head	each	o.a.
AGB8466	Sample Positioning Micrometer Table	each	o.a.
AGB8467	90° divider with 75mm Holder	each	o.a.
AGB8469	Diamond Wire 0.22mm dia x 30 metres	each	o.a.
AGB8470	Diamond Wire 0.13mm dia x 30 metres	each	o.a.
AGB8471	Diamond Wire 0.3mm dia x 30 metres	each	o.a.
AGB8499	Single wheel Polisher/Grinder, 203mm (8") variable speed, 50	each	o.a.
AGB8500	WELL Vertical Diamond Wire Saw	each	o.a.
AGB8500B-KIT	Base Kit- WELL Vertical Diamond Wire Saw, 3500	each	o.a.
AGB8500-KIT	Premium Kit - WELL Vertical Diamond Wire Saw, 3500	each	o.a.
AGB8503	Drum w/ 10 meters of diamond 0.3mm - 20um Pos 3 (Green)	each	o.a.
AGB8503-3	Drum w/ 10 meters of diamond 0.13mm - 20um Pos 3 (Green)	each	420,00
AGB8508-AR	Cutting liquid concentrate 5% (anti-rust)	500ml	92,20
AGB8511	Aluminium Oxide Wheel, 250mm dia.	10	o.a.
AGB8512	Silicon Carbide Wheel, 250mm dia.	10	o.a.
AGB8532-KIT	WELL Horizontal Diamond Wire Saw, 3400 Kit	each	o.a.
AGB8539	Adhesive Silicon Carbide Discs P2500 203mm	25	60,50
AGB8541	Adhesive Silicon Carbide Discs P80 203mm	25	42,10
AGB8542	Adhesive Silicon Carbide Discs P120 203mm	25	43,60
AGB8543	Adhesive Silicon Carbide Discs P180 203mm	25	39,80
AGB8544	Adhesive Silicon Carbide Discs P240 203mm	25	34,90
AGB8545	Adhesive Silicon Carbide Discs P320 203mm	25	36,60
AGB8546	Adhesive Silicon Carbide Discs P400 203mm	25	34,50
AGB8547	Adhesive Silicon Carbide Discs P600 203mm	25	36,20
AGB8548	Adhesive Silicon Carbide Discs P800 203mm	25	34,70
AGB8549	Adhesive Silicon Carbide Discs P1200 203mm	25	38,30

Part #	Description	Unit	Price €
AGB8570	Premium WELL Vertical Diamond Wire Saw 6500 KIT	each	o.a.
AGB8580	Premium WELL Vertical Diamond Wire Saw 4500 KIT	each	o.a.
AGB8640	35mm Aluminium oxide polishing strip, 60µm.	10	8,30
AGB8642	35mm Aluminium oxide polishing strip, 30µm.	10	11,90
AGB8643	35mm Aluminium oxide polishing strip, 12µm.	10	11,70
AGB8644	35mm Aluminium oxide polishing strip, 9µm.	10	8,30
AGB8645	35mm Aluminium oxide polishing strip, 3µm.	10	11,90
AGB8646	35mm Aluminium oxide polishing strip, 1µm.	10	11,90
AGB8647	35mm Aluminium oxide polishing strip, 0.3µm.	10	11,90
AGB8652	35mm Silicon carbide polishing strip, 30µm.	10	8,30
AGB8653	35mm Silicon carbide polishing strip, 15µm.	10	8,10
AGB8654	35mm Silicon carbide polishing strip, 9µm.	10	11,90
AGB8655	35mm Silicon carbide polishing strip, 5µm.	10	11,90
AGB8661	35mm Diamond polishing strip, 15µm	each	24,60
AGB8662	35mm Diamond polishing strip, 9µm	each	23,90
AGB8663	35mm Diamond polishing strip, 6µm	each	24,60
AGB8664	35mm Diamond polishing strip, 3µm	each	24,60
AGB8665	35mm Diamond polishing strip, 1µm	each	24,60
AGB8666	35mm Diamond polishing strip, 0.5µm	each	23,90
AGB8670	70mm Aluminium oxide polishing strip, 60µm.	10	15,10
AGB8672	70mm Aluminium oxide polishing strip, 30µm.	10	15,10
AGB8673	70mm Aluminium oxide polishing strip, 12µm.	10	15,50
AGB8674	70mm Aluminium oxide polishing strip, 9µm.	10	15,10
AGB8675	70mm Aluminium oxide polishing strip, 3µm.	10	15,50
AGB8676	70mm Aluminium oxide polishing strip, 1µm.	10	15,50
AGB8677	70mm Aluminium oxide polishing strip, 0.3µm.	10	20,00
AGB8681	70mm Silicon carbide polishing strip, 30µm.	10	15,10
AGB8682	70mm Silicon carbide polishing strip, 15µm.	10	18,70
AGB8683	70mm Silicon carbide polishing strip, 9 µm.	10	15,10
AGB8684	70mm Silicon carbide polishing strip, 5µm.	10	18,70
AGB8691	70mm Diamond polishing strip, 15µm	each	48,80
AGB8692	70mm Diamond polishing strip, 9µm	each	48,80
AGB8693	70mm Diamond polishing strip, 6µm	each	66,20
AGB8694	70mm Diamond polishing strip, 3µm	each	66,20
AGB8695	70mm Diamond polishing strip, 1µm	each	48,80
AGB8696	70mm Diamond polishing strip, 0.5µm	each	48,80
AGB8700	215 x 279mm Aluminium oxide polishing sheet, 60µm.	each	4,20
AGB8702	215 x 279mm Aluminium oxide polishing sheet, 30µm.	each	4,10
AGB8703	215 x 279mm Aluminium oxide polishing sheet, 12µm.	each	11,60
AGB8704	215 x 279mm Aluminium oxide polishing sheet, 9µm.	each	4,10
AGB8705	215 x 279mm Aluminium oxide polishing sheet, 3µm.	each	5,80
AGB8706	215 x 279mm Aluminium oxide polishing sheet, 1µm.	each	4,10
AGB8707	215 x 279mm Aluminium oxide polishing sheet, 0.3µm.	each	6,10
AGB8712	215 x 279mm Silicon carbide polishing sheet, 30µm.	each	11,00
AGB8713	215 x 279mm Silicon carbide polishing sheet, 15µm.	each	11,00
AGB8715	215 x 279mm Silicon carbide polishing sheet, 5µm.	each	11,00
AGB8721	215 x 279mm Diamond polishing sheet, 15µm	each	134,80
AGB8722	215 x 279mm Diamond polishing sheet, 9µm	each	134,80
AGB8723	215 x 279mm Diamond polishing sheet, 6µm	each	131,50
AGB8724	215 x 279mm Diamond polishing sheet, 3µm	each	134,80
AGB8725	215 x 279mm Diamond polishing sheet, 1µm	each	128,00
AGB8726	215 x 279mm Diamond polishing sheet, 0.5µm	each	134,80
AGB8740	Ultraprep Diamond Lapping Film 203mm 30µm	each	#N/B
AGB8741	Ultraprep Diamond Lapping Film 203mm 15µm	each	416,60
AGB8742	Ultraprep Diamond Lapping Film 203mm 9µm (Pk of 5)	each	#N/B
AGB8743	Ultraprep Diamond Lapping Film 203mm 6µm	each	416,60
AGB8744	Ultraprep Diamond Lapping Film 203mm 3µm	each	409,70
AGB8746	Ultraprep Diamond Lapping Film 203mm 0.5µm	each	479,60

Part #	Description	Unit	Price €
AGB8747	Ultraprep Diamond Lapping Film 203mm 0.1µm	each	404,80
AGB8748	Mixed Pack of Ultraprep Diamond Lapping Film	each	413,10
AGB8754	600 Grit P1200 Silicon Carbide Powder	2,3kg	239,80
AGB8755	1000 Grit P2500 Silicon Carbide Powder	2,3kg	459,00
AGB8762	600 Grit P1200 Alumina Powder	2,3kg	#N/B
AGB8768	12.5µm Alumina Powder	2,3kg	126,00
AGB8769	9.5µm Alumina Powder	2,3kg	126,00
AGB8770	5.0µm Alumina Powder	2,3kg	188,90
AGB8771	3.0µm Alumina Powder	2,3kg	206,80
AGB8780	MiroMet Cerium Oxide Suspension	180ml	127,80
AGB8782H-2	* EpoxiCure Hardener	460ml	102,20
AGB8782R-2	* EpoxiCure Resin	1,9l	204,80
AGB8784	EpoMet G Moulding Powder	1,85kg	177,80
AGB8784F	EpoMet Fine Grain Moulding Powder	1,85kg	172,90
AGB8785	* ProbeMet Conductive Moulding Powder	1,85kg	259,70
AGB8786	* Samplkwick Acrylic Cold mtg Kit Powder 450g + Hardener	450g	189,40
AGB8786H	Samplkwick Acrylic Cold Mounting Hardener	360ml	61,20
AGB8786P	Samplkwick Acrylic Cold Mounting Powder	450g	131,70
AGB8787	TransOptic	2,3kg	156,90
AGB8788	KonductoMet, 450g	450g	97,70
AGB8789	Varidur10 Fast Curing Acrylic Kit (Powder 1kg + Hardener 500ml)	each	177,00
AGB8789H	Varidur10 Hardener	500ml	64,20
AGB8789P	Varidur10 Powder Semi Transparent	1kg	159,60
AGB8790H	* EpoFix Hardener 130ml	130ml	153,40
AGB8790R	* EpoFix Resin	1l	190,90
AGB8800	* Scandiplast Embedding Compound	1kg	181,90
AGB8801	* Hardener for Scandiplast Resin	35ml	23,50
AGB8802	* Scandiplex 82a Resin	1kg	176,30
AGB8804	Scan-dia Embedding Set	each	115,20
AGB8805	Scandia Varioplex	200g	91,60
AGB8808	Aequidir Hardness Equaliser Soft	100g	117,10
AGB8809	Aequidir Hardness Equaliser Medium	100g	117,10
AGB8810	Aequidir Hardness Equaliser Hard	100g	117,10
AGB8814	Scan-Dia Plastic Clips	100	94,00
AGB8815	Scandiquick Powder	1kg	182,50
AGB8816	Scandiquick Liquid	750ml	138,00
AGB8820	Scandiform 25mm Embedding Moulds	each	130,40
AGB8821	Scandiform 32mm Embedding Moulds	each	130,40
AGB8822	Scandiform 38mm Embedding Moulds	each	154,60
AGB8823	Scandiform 50mm Embedding Moulds	each	232,00
AGB8824	Scandiform 25 x 57 x 20mm Moulds	each	206,60
AGB8825	Scandia Paper / Plastic Cups	50	37,90
AGB8828	Demotec 15 Plus Powder	1kg	186,80
AGB8828-10	Demotec 15 Plus Powder	10kg	1.042,20
AGB8828-5	Demotec 15 Plus Powder	5kg	623,40
AGB8829	* Demotec 15 Plus Liquid	1l	207,90
AGB8829-05	* Demotec 15 Plus Liquid	500ml	131,00
AGB8829-5	* Demotec 15 Plus Liquid	5l	664,90
AGB8830	Demotec 70 Powder	1kg	243,10
AGB8831	* Demotec 70 Liquid	1l	174,80
AGB8832	Silicone Rubber Moulds 1inch Diameter	10	129,00
AGB8832-125	Silicone Rubber Moulds 1.25inch Diameter	10	136,00
AGB8832-150	Silicone Rubber Moulds 1.5inch Diameter	10	163,30
AGB8832-200	Silicone Rubber Moulds 2inch Diameter	each	163,30
AGB8833	Silicone Rubber Moulds 2.25inchx1inch	each	305,20
AGB8834	Silicone Rubber Mould 3inch x 2inch	each	323,90
AGB8835	Silicone Rubber Moulds Blue 1.25inch Dia.	10	156,90
AGB8836	2-Part Plastic Moulds 1inch Diameter	12	76,40

Part #	Description	Unit	Price €
AGB8836-125	2-Part Plastic Moulds 1.25inch Diameter	10	85,40
AGB8836-150	2-Part Plastic Moulds 1.5inch Diameter	10	91,70
AGB8837	Metal mounting clips.	100	67,80
AGB8856	PetroBond - Thin Section Bonding Tool	each	o.a.
AGB8958	ELMO Glow Discharge system for TEM grids	each	o.a.
AGB8959	easiSlicer	each	o.a.
AGB8965	Steady Temp Pro Digital	each	8.264,90
AGB9001	Precision Diamond Wire Saw with Digital Control	each	o.a.
AGB9002	Low speed 4" diamond saw with digital readout	each	o.a.
AGB9002B	Low speed 6" diamond saw with digital readout	each	o.a.
AGB9003	Multi-Purpose Grinder/Polisher	each	o.a.
AGB9200	Gyratory Shaker Stuart mini SSM3 5-70 rpm, 3-12? Throw, 230V 50/60Hz	each	1.139,10
AGB9201	Gyratory Shaker Stuart mini SSM4 5-70 rpm, 7? Throw, 230V 50/60Hz	each	897,20
AGB9202	Gyratory Shaker Stuart SSL3 5-70 rpm, 3-12? Throw, 230V 50/60Hz	each	1.527,20
AGB9203	Gyratory Shaker Stuart SSL4 5-70 rpm, 7? Throw, 230V 50/60Hz	each	1.498,00
AGB9204	Tier system (2 platforms plus 8 bar) (SSL3/1)	each	331,90
AGB9205	Tier system (2 platforms plus 8 bar) (SSM3/1)	each	264,70
AGB9230-3	8" Dia. Diamond Lapping Film (PSA) 3 microns grit	each	154,10
AGB9304	8" Backing Plates for PSA Diamond Plate, Sand Paper and Po	each	239,70
AGB9305-1200	SiC Sand Discs (PSA), Waterproof, 8" Diameter, 1200	10	99,10
AGB9305-600	SiC Sand Discs (PSA), Waterproof, 8" Diameter, 600	10	91,80
AGB9305-800	SiC Sand Discs (PSA), Waterproof, 8" Diameter, 800	10	99,10
AGB9320-05	8" Dia. Diamond Lapping Film (PSA) 0.5 microns grit	each	157,80
AGB9320-1	8" Dia. Diamond Lapping Film (PSA) 1 micron grit	each	154,10
AGB9320-6	8" Dia. Diamond Lapping Film (PSA) 6 microns grit	each	150,40
AGB9330	8" Poromeric Polishing Pad (PSA) for Final Polishing	each	39,20
AGB9340	8" Polyamide Foam polishing Pad	each	91,80
AGB9350-05	Diamond Compound Polishing Paste, 0.5 Microns, 10 Grams Syri	each	52,90
AGB9350-1	Diamond Compound Polishing Paste, 1.0 Micron, 10 Grams Syrin	each	50,90
AGB9350-3	Diamond Compound Polishing Paste, 3.0 Microns, 10 Grams Syri	each	50,90
AGB9350-5	Diamond Compound Polishing Paste, 5.0 Microns, 10 Grams Syri	each	47,00
AGB9600	KTR Chiller, low pressure pump, 208-230V 50/60Hz	each	#N/B
AGB9601	KTR Chiller, 2 litre high pressure pump, 208-230V 50/60Hz	each	#N/B
AGB9602	KTR Chiller, 3.5 litre high pressure pump, 208-230V 50/60Hz	each	#N/B
AGB9603	KTC Chiller, LED Display, 2 l high pressure pump, 208-230V 5	each	#N/B
AGB9604	KTC Chiller, LED Display, 3.5 l high pressure pump, 208-230V	each	#N/B
AGB9605	KTC Chiller, LED Display, Low Pressure pump, 208-230V 50/60H	each	#N/B
AGB9606	KTD Digital controlled Chiller, 2 l high pressure pump, 208-	each	#N/B
AGB9607	KTD Digital controlled Chiller, 3.5 l high pressure pump, 20	each	#N/B
AGB9608	KTD Digital controlled Chiller, low pressure pump, 208-230V	each	#N/B
AGB9609	KT1 Chiller with 5l/min pump	each	#N/B
AGB9610	KT1 Chiller with 10l/min pump	each	#N/B
AGB9611	KT1 Chiller with 17l/min pump	each	#N/B
AGB9612	K1 Chiller with 5l/min pump	each	#N/B
AGB9613	K1 Chiller with 10l/min pump	each	#N/B
AGB9614	K1 Chiller with 17l/min pump	each	#N/B
AGB9615	K3 Chiller with 5l/min pump	each	#N/B
AGB9616	K3 Chiller with 10l/min pump	each	#N/B
AGB9617	K3 Chiller with 17l/min pump	each	#N/B
AGB9618	K4 Chiller with 5l/min pump	each	#N/B
AGB9619	K4 Chiller with 10l/min pump	each	#N/B
AGB9620	K4 Chiller with 17l/min pump	each	#N/B
AGB9621	K6 Chiller with 5l/min pump, 208/230V 50/60Hz	each	#N/B
AGB9622	K6 Chiller with 10l/min pump, 208/230V 50/60Hz	each	#N/B
AGB9623	K6 Chiller with 17l/min pump, 208/230V 50/60Hz	each	#N/B
AGB9624	K9 Chiller with 10l/min pump, 380-420V 50/60Hz 3 phase	each	#N/B
AGB9625	K9 Chiller with 17l/min pump, 380-420V 50/60Hz 3 phase	each	#N/B
AGB9626	Chiller with P17 pump 380-420V 50/60Hz Three Phase	each	#N/B

Part #	Description	Unit	Price €
AGB9627	A4 in 5l Container, temp range -15C to +90C	each	#N/B
AGB9628	A6 in 5l Container, temp range -45C to +90C	each	#N/B
AGB9629	K12 Chiller with 17l/min pump, 380-420V 50/60Hz Three Phase	each	#N/B
AGB9630	K12 Chiller with 25l/min pump, 380-420V 50/60Hz Three Phase	each	#N/B
AGBCR-3407	BCR-3407 Styrofoam Cryogen Jug	each	391,70
AGBT1040078	Software, FlowMeter Setup Pack	each	#N/B
AGBT200AG-E	AirGard Fume Hood Monitor, Int	each	#N/B
AGBT4040	4040 Flowmeter with LCD, 300L/min Air, O2, N2 22mm Tube	each	#N/B
AGBT4043	4043 Flowmeter with LCD 200L/min Air, O2, N2 1/2" Tube	each	#N/B
AGBT4045	4045 FlowMeter with LCD, 300L/min Air, O2, N2 3/4" Tube	each	#N/B
AGBT405AG-E	AirGard Fume Hood Monitor, Int	each	#N/B
AGBT4073	4073 Certifier FA and FA Plus Oxygen Sensor Kit	each	#N/B
AGBT4075	4075 Certifier FA, Flow Analyzer Low Flow Kit	each	#N/B
AGBT4080	4080 Certifier FA Plus Flow Analyzer System: Inc 4088 Interface Module & 408	each	#N/B
AGBT4081	4081 Certifier FA Plus, High Flow Module Kit	each	#N/B
AGBT4082	4082 Certifier FA Plus, Low Flow Module Kit	each	#N/B
AGBT4088	4088 Certifier FA Plus, Interface Module	each	#N/B
AGBT4140	4140 Flowmeter with LCD 20L/min Air, O2, N2 1/4" Tube	each	#N/B
AGBT41403	Flowmeter with LCD - 0.01-20L/Min NAO Air O2 NA 1/4 Tube	each	#N/B
AGBT4143	4143 Flowmeter with LCD 20L/min Air, O2, N2 3/8" Tube	each	#N/B
AGBT4199	Flowmeter Battery Pack	each	#N/B
AGBT5200-1	Flowmeter Kit, 30 L/Min, Air, N2, 3/8 Inch 1/4 Inch and 10mm	each	#N/B
AGBT5200-10MM	5200 10mm Tube Ends (2)	each	#N/B
AGBT5200-2	Flowmeter Kit, 30 L/min, Air, O2, N2, 3/8 inch 1/4 inch 10mm	each	#N/B
AGBT5200-3	Flowmeter Kit, 30 L/min, Air, CO2, N2, 3/8 inch 1/4 inch 10m	each	#N/B
AGBT5200-6	Flowmeter Kit, 30 L/min, O2, 3/8 inch 1/4 inch and 10mm tube	each	#N/B
AGBT5200-7	Flowmeter Kit, 30 L/min, CO2, 3/8 inch 1/4 inch and 10mm tub	each	#N/B
AGBT5200-C	5200 Collars (2)	each	#N/B
AGBT5200-EC	5200 End Caps (2)	each	#N/B
AGBT5200-IF-025IN	5200 Inlet Filter 1/4 Inch - Low Pressure	each	#N/B
AGBT5200-IF-0375IN	5200 Inlet Filter 3/8 Inch - Low Pressure	each	#N/B
AGBT5200-LPMK	5200 Low Pressure Measurement Kit	each	#N/B
AGBT5203-2	Flowmeter +WIDE ACC Kit, 30 L/min, Air, O2, N2, 3/8 inch 1/4	each	#N/B
AGBT5210-1	Flowmeter +DP Kit, 30 L/min, Air, N2, 3/8 inch 1/4 inch 10mm	each	#N/B
AGBT5210-2	Flowmeter +DP Kit, 30 L/min, Air, N2, 3/8 inch 1/4 inch 10mm	each	#N/B
AGBT5210-3	Flowmeter +DP Kit, 30 L/min, Air, CO2, N2, 3/8 inch 1/4 inch	each	#N/B
AGBT5210-6	Flowmeter +DP Kit, 30 L/min, O2, 3/8 inch 1/4 inch and 10mm	each	#N/B
AGBT5210-7	Flowmeter +DP Kit, 30 L/min, CO2, 3/8 inch 1/4 inch and 10mm	each	#N/B
AGBT5220-2	Flowmeter +DP +%RH Kit, 30 L/min, Air, O2, N2, 3/8 1/4 inch	each	#N/B
AGBT5230-2	Flowmeter +DP +%RH +HIGH ACC Kit, 30 L/min, Air, O2, N2, 3/8	each	#N/B
AGBT53000	5300 Standard Accessory Kit	each	#N/B
AGBT5300-05IN	5300 1/2 Inch Tube Ends (2)	each	#N/B
AGBT5300-075IN	5300 3/4 Inch Tube Ends (2)	each	#N/B
AGBT5300-1	Flowmeter Kit, 300 L/min, Air, N2, 22mm and 1/2 inch tube en	each	#N/B
AGBT5300-2	Flowmeter Kit, 300 L/min, Air, O2, N2, 22mm and 1/2 inch tub	each	#N/B
AGBT5300-3	Flowmeter Kit, 300 L/min, Air, CO2, N2, 22mm and 1/2 inch tu	each	#N/B
AGBT5300-4	Flowmeter Kit, 300 L/min, O2, 22mm and 1/2 inch tube ends	each	#N/B
AGBT5300-5	Flowmeter Kit, 300 L/min, CO2, 22mm and 1/2 tube ends	each	#N/B
AGBT5303-2	Flowmeter +WIDE ACC Kit, 300 L/min, Air, O2, N2, 22mm and 1/	each	#N/B
AGBT5310-1	Flowmeter + DP Kit, 300 L/min, Air, N2, 22mm and 1/2 inch tu	each	#N/B
AGBT5310-2	Flowmeter +DP Kit, 300 L/min, Air, O2, N2, 22mm and 1/2 inch	each	#N/B
AGBT5310-3	Flowmeter +DP Kit, 300 L/min, Air, CO2, N2, 22mm and 1/2 inc	each	#N/B
AGBT5310-4	Flowmeter + DP Kit, 300 L/Min, O2, 22mm and 1/2 inch tube en	each	#N/B
AGBT5310-5	Flowmeter + DP Kit, 300 L/min, CO2, 22mm and 1/2 inch tube e	each	#N/B
AGBT5320-2	Flowmeter +DP +%RH Kit, 300 L/min, Air, O2, N2, 22mm 1/2 inc	each	#N/B
AGBT5330-2	Flowmeter +DP +%RH +HIGH ACC Kit, 300 L/min, Air, O2, N2, 22	each	#N/B
AGBTCL-5200	Cal of 5200-2, 5200-3, 5200-6, 5200-7, 5203-2	each	#N/B
AGBTCL-5200-AO	Cal of 5200-1 (Air Only)	each	#N/B

Part #	Description	Unit	Price €
AGBTCL-5210	Cal of 5210-2, 5210-3, 5210-4, 5210-6, 5210-7, 5210-8	each	#N/B
AGBTCL-5210-AO	Cal of 5210-1 (Air Only)	each	#N/B
AGBTCL-5220	Cal of 5220-2	each	#N/B
AGBTCL-5230	Cal of 5230-2	each	#N/B
AGBTCL-5300	Cal of 5300-2, 5300-3, 5300-4, 5300-5, 5303-2	each	#N/B
AGBTCL-5300-AO	Cal of 5300-1 (Air Only)	each	#N/B
AGBTCL-5310	Cal of 5310-2, 5310-3, 5310-4, 5310-5	each	#N/B
AGBTCL-5310-AO	Cal of 5310-1 (Air only)	each	#N/B
AGBTCL-5320	Cal of 5320-2	each	#N/B
AGBTCL-5330	Cal of 5330-2	each	#N/B
AGC350	Cryo-Gel, 120g	each	48,70
AGC351	CytoFix Spray (Non-Aerosol, pump bottle)	110ml	37,20
AGC352	CytoFix Spray (Aerosol, can)	110ml	37,20
AGC353	LV CryoOil	1,5ml	148,90
AGC401L-40	Chem Neo Nitrile 300 Gloves Red L (Pk40)		#N/B
AGC401M-40	Chem Neo Nitrile 300 Gloves Red M (Pk40)		#N/B
AGC401S-40	Chem Neo Nitrile 300 Gloves Red S (Pk40)		#N/B
AGC401XL-40	Chem Neo Nitrile 300 Gloves Red X L (Pk40)		#N/B
AGC401XS-40	Chem Neo Nitrile 300 Gloves Red XS (Pk40)		#N/B
AGC401XXL-40	Chem Neo Nitrile 300 Gloves Red XX L (Pk40)		#N/B
AGC402L-150	Ecoshield 250 Nitrile Gloves Green L (Pk150)		#N/B
AGC402M-150	Ecoshield 250 Nitrile Gloves Green M (Pk150)		#N/B
AGC402S-150	Ecoshield 250 Nitrile Gloves Green S (Pk150)		#N/B
AGC402XL-150	Ecoshield 250 Nitrile Gloves Green X L (Pk150)		#N/B
AGC402XS-150	Ecoshield 250 Nitrile Gloves Green XS (Pk150)		#N/B
AGC403L-20	Bio Orange 300 Nitrile Gloves L (Pk20)		#N/B
AGC403L-50	Bio Orange 300 Nitrile Gloves L (Pk50)		#N/B
AGC403L-90	Bio Orange 260 Nitrile Gloves L (Pk90)		#N/B
AGC403M-20	Bio Orange 300 Nitrile Gloves M (Pk20)		#N/B
AGC403M-50	Bio Orange 300 Nitrile Gloves M (Pk50)		#N/B
AGC403M-90	Bio Orange 260 Nitrile Gloves M (Pk90)		#N/B
AGC403S-20	Bio Orange 300 Nitrile Gloves S (Pk20)		#N/B
AGC403S-50	Bio Orange 300 Nitrile Gloves S (Pk50)		#N/B
AGC403S-90	Bio Orange 260 Nitrile Gloves S (Pk90)		#N/B
AGC403XL-20	Bio Orange 300 Nitrile Gloves XL (Pk20)		#N/B
AGC403XL-50	Bio Orange 300 Nitrile Gloves X L (Pk50)		#N/B
AGC403XL-90	Bio Orange 260 Nitrile Gloves X L (Pk90)		#N/B
AGC403XS-50	Bio Orange 300 Nitrile Gloves XS (Pk50)		#N/B
AGC403XS-90	Bio Orange 260 Nitrile Gloves XS (Pk90)		#N/B
AGC404L-100	Microtouch Coated Latex Powder Free Glove - Large Pk100		53,50
AGC404M-100	Microtouch Coated Latex Powder Free Glove - Medium Pk100		53,50
AGC404S-100	Microtouch Coated Latex Powder Free Glove - Small Pk100		53,50
AGC404XS-100	Microtouch Coated Latex Powder Free Glove - ExtraSmall Pk100		32,30
AGC5120	Buffer Solution pH4.00 ± 0.01 @ 20°C 250ml Red		17,30
AGC5121	Buffer Solution pH7.00 ± 0.01 @ 20°C 250ml Yellow		17,30
AGC5122	Buffer Solution pH10.00 ± 0.01 @ 20°C 250ml Blue		17,30
AGC5123	Buffer Solution pH4.00 (Red) ± 0.01 @ 20°C 500ml		26,40
AGC5124	Buffer Solution pH7.00 (Yellow) ± 0.01 @ 20°C 500ml		27,60
AGC5125	Buffer Solution pH10.00 (Blue) ± 0.01 @ 20°C 500ml		27,60
AGC5126	Buffer Solution pH4.00 (Red) ± 0.01 @ 20°C 1L		40,60
AGC5127	Buffer Solution pH7.00 (Yellow) ± 0.01 @ 20°C 1L		40,50
AGC5128	Buffer Solution pH10.00 (Blue) ± 0.01 @ 20°C 1L		40,60
AGC8000	Powder-Free Vinyl Gloves Large	100	#N/B
AGC8000A	Vinyl Gloves Powder-Free Large (Pack of 100)	100	23,10
AGC8001	Powder-Free Vinyl Gloves Medium	100	24,70
AGC8002	Powder-Free Vinyl Gloves Small	100	31,50
AGC8004	Latex Gloves Small	100	36,30
AGC8005	Latex Gloves Medium	100	36,30

Part #	Description	Unit	Price €
AGC8006	Latex Gloves Large	100	36,30
AGC801	Selvyt cloth	each	7,00
AGC8010	* Falcon Dust-Off Plus Special Application Complete	186ml	#N/B
AGC8011	* Falcon Dust-Off Plus Special Application Refill	186ml	#N/B
AGC8012	* Falcon Dust-Off Plus Special App. Refills 12/Box	12 x 186ml	#N/B
AGC8013	* Falcon New Formulation Dust-Off Plus Complete	300ml	61,50
AGC8014	* Falcon New Formulation Dust-Off Plus Refill	300ml	23,20
AGC8015	* Falcon New Formulation Dust-Off+ Refills 12/Case	12 x 300ml	176,50
AGC8016	* Puff duster	each	24,80
AGC8017	Nozzle only for Dust-Off Plus	each	44,80
AGC8018	* Falcon Dust-Off Eco 225ml Disposable Duster	225ml	27,10
AGC8020	Powder Free Latex Gloves Small	100	63,20
AGC8021	Powder Free Latex Gloves Medium	100	63,20
AGC8022	Powder Free Latex Gloves Large	100	63,20
AGC8024	* Dust-Off Jumbo	530ml	24,10
AGC8024-12	* Dust-Off Jumbo 530ml Case of 12	12 x 530ml	216,40
AGC8024L	Silvershield Norfoil Gloves Large (Pair)	each	#N/B
AGC8024M	Silvershield Norfoil Gloves Medium (Pair)	each	27,90
AGC8025M	Ansell Neoprene Gloves Medium	100	31,20
AGC8025S	Ansell Neoprene Gloves Small	100	31,20
AGC8026	Monitor Wipes	200	51,50
AGC8026L	Kevlar Lined Gloves Large	10	136,30
AGC8026S	Kevlar Lined Gloves Small	10	136,30
AGC8026XL	Kevlar Lined Gloves Extra Large	10	136,30
AGC8027L	Sterling Nitrile Gloves Large	150	83,10
AGC8027M	Sterling Nitrile Gloves Medium	150	88,30
AGC8027S	Sterling Nitrile Gloves Small	150	83,10
AGC8027XS	Sterling Nitrile Gloves X-Small	150	83,10
AGC8028L	Cotton Gloves Large	12 pair	16,10
AGC8028M	Cotton Gloves Medium	12 pair	21,80
AGC8028S	Cotton Gloves Small	12 pair	16,10
AGC8030A	Facemask	each	17,20
AGC8031	Epoxy Hand Cleaner	450ml	36,30
AGC8032A	DEB Protect	100ml	64,00
AGC8034A	Deb Restore	100ml	95,90
AGC8038	Microfibre cleaning cloth	each	12,80
AGC8039	Azowipe bactericidal wipes	each	46,10
AGC804	Cocktail Sticks	1000	5,10
AGC8043	* Dust-Off XL 300ml	300ml	13,00
AGC8043-12	* Dust-Off XL 300ml Case of 12	12 x 300ml	99,40
AGC8044	* Dust-Off JR 109ml	109ml	7,90
AGC8044-24	* Dust-Off JR 109ml	24 x 109ml	150,80
AGC8045	* Label Off 50	each	33,10
AGC8046	* Kontakt IPA	each	21,50
AGC8047	* Kontakt 40	each	12,80
AGC8048	* Anti-static cleaning foam	each	9,50
AGC8050	Cryo Apron Small	each	289,70
AGC8051	Cryo Apron Medium	each	315,10
AGC8052	Cryo Apron Large	each	332,30
AGC8053	Cryo Apron Extra Large	each	348,70
AGC806	Cotton Tipped Applicators	100	7,40
AGC8061	Monitor Wipes	80	29,50
AGC8062	Winchester Bottle Carrier	each	238,40
AGC8064	Polythene Ampoule Snappers 11-20ml	100	31,00
AGC8065	*ardous spill clean-up sheets	each	#N/B
AGC8068	* 3M Citrus cleaner, adhesive remover spray	each	47,00
AGC807	Wenol metal polish	each	23,10
AGC8070	Pikal professional polishing paste,	250g	67,50

Part #	Description	Unit	Price €
AGC8072	Speck Grabber Pro	each	65,60
AGC8074	Mopitup spill sheets	each	#N/B
AGC8075	Labtek Liquid	each	85,00
AGC808	Polythene Ampoule Snappers 1ml-4ml	100	24,90
AGC809	Polythene Ampoule Snappers 5ml-10ml	100	29,50
AGC810	Steel Blue Wipes 52cm Wide	each	27,70
AGC811	Steel Blue Wipes 25cm Wide	each	10,60
AGC813A	Kimwipes	196	21,20
AGC814	Service Cloths	76	236,20
AGC815-1	Professional Tissues	100	6,20
AGC816	Medical Wipes	80	2,30
AGC817	White Dispenser for 520mm rolls	each	96,00
AGC818	White Dispenser for 260mm rolls	each	78,20
AGC819	Disposable Gloves Large	100	24,60
AGC819A	Disposable Gloves Large Box 100	100	28,20
AGC819A	Disposable Gloves Large Box 100	100	28,20
AGC820	Disposable Gloves Medium	100	23,90
AGC821	Disposable Gloves Small	100	19,10
AGC821A	Disposable Gloves Small Box 100	100	29,10
AGC833	Lens Cleaning Tissue 20cm x 30cm	100	113,80
AGC835	Gloves Pseudo Leather Size 7 (pair)	each	6,10
AGC837	Gloves Pseudo Leather Size 9 (pair)	each	6,10
AGC838	Gloves Pseudo Leather Size 10 (pair)	each	6,10
AGC839	15cm Bamboo Splints	100	6,60
AGC843	Twill Jean Cloths	10	10,30
AGC844	Twill Jean Cloths	100	94,20
AGC847	Decon 90	each	144,20
AGC848	Decon Neutracon	each	195,40
AGC849	Procure barrier cream	each	14,00
AGC850	Blower	each	18,20
AGC850A	Bakelite blower	each	68,40
AGC860	Fine Nylon Gloves - Small (pair)	each	5,40
AGC861	Fine Nylon Gloves - Medium (pair)	each	6,30
AGC862	Fine Nylon Gloves - Large (pair)	each	6,80
AGC864	Mid Arm Length Cryoglove (Medium)	each	319,40
AGC865	Elbow Length Cryogloves (Medium)	each	357,00
AGC866	1/8inch Hyprez Lapping Sticks	10	19,60
AGC867	3/16inch Hyprez Lapping Sticks	10	19,30
AGC868	1/4inch Hyprez Lapping Sticks	10	20,10
AGC869L	Shoulder Length Cryogloves (Large)	each	364,40
AGC870	Poly-Jean cloth	each	91,50
AGC874	Rubber Gloves - Small (pair)	each	7,70
AGC875	Rubber Gloves - Medium (pair)	each	7,70
AGC876	Rubber Gloves - Large (pair)	each	5,50
AGC879	Nitrile Gloves Large (Pair)	each	31,30
AGC880	Heat Resistant Gloves Mens.(Pair)	each	56,40
AGC881	Heat Resistant Gloves Womens.(Pair)	each	56,40
AGC882	Safety hot hand protector	each	92,00
AGC885	Safety Spectacles Frisco Clear PC Lens Blue Frame	each	42,00
AGC886	Goggles with Ventilation	each	37,70
AGC886V	Visor for AGC886 Goggles	each	26,40
AGC889	Ampoule snappers	each	48,90
AGC892	Spillage Absorption Granules	2,5kg	337,80
AGC899	Polystyrene Cleaning Strips	6	16,00
AGC901	Oxidation Removal from Metals	1l	20,20
AGC902	Oxidation Removal from Metals	5l	66,20
AGC903	Flux remover & PCBA Cleaner	1l	18,40
AGC904	Flux remover & PCBA Cleaner	2l	66,20

Part #	Description	Unit	Price €
AGC905	Flux remover & PCBA Cleaner	25l	217,00
AGC906	Sensitive Metals Cleaner	1l	18,40
AGC907	Sensitive Metals Cleaner	5l	66,20
AGC908	Glass & Optical Lenses Cleaner Inc. Ceramics	1l	20,20
AGC909	Glass & Optical Lenses Cleaner Inc. Ceramics	5l	66,20
AGCL1-PITOT	Check Pitots with Certificate of Conform	each	#N/B
AGCL-4040	Calibration of Flow Meter 4040	each	#N/B
AGCL-4043	Calibration of Flow Meter 4043	each	#N/B
AGCL-4045	Calibration of Flow Meter 4045	each	#N/B
AGCL-4081	Calibration of Certifier FA+ Model 4081	each	#N/B
AGCL-4082	Calibration of Certifier FA+ Model 4082	each	#N/B
AGCL-4140	Calibration of Flow Meter 4140	each	#N/B
AGCL-4143	Calibration of Flow Meter 4143	each	#N/B
AGCL-5725	Calibrate Model QA5725	each	#N/B
AGCL-5815	Calibrate Model 5815	each	#N/B
AGCL-5825	Calibrate Model 5825	each	#N/B
AGCL-6200	Calibration of Balometer LoFlo 6200	each	#N/B
AGCL-7515	Calibration of Model QA7515	each	#N/B
AG-CL7525	Calibration of Model QA7525	each	#N/B
AGCL-7545	Calibration of Model QA7545	each	#N/B
AGCL-7545CO	Calibration of Model QA7545 with CO Repl	each	#N/B
AGCL-800220	Calibration of Humidity / Temp Probe	each	#N/B
AGCL-8380	Calibrate Accubalance QA8380	each	#N/B
AGCL-8455	Calibration of Air Velocity Transducer 8	each	#N/B
AGCL-8465	Calibration of Air Velocity Transducer 8	each	#N/B
AGCL-8475	Calibration of Air Velocity Transducer 8	each	#N/B
AGCL-8715	Calibrate QA8715 Micromanometer	each	#N/B
AGCL-9515	Calibrate Model QA9515	each	#N/B
AGCL-9535	Calibrate Model QA9535	each	#N/B
AGCL-9535-A	Calibrate Model QA9535-A	each	#N/B
AGCL-9545	Calibrate Model QA9545	each	#N/B
AGCL-9545-A	Calibrate Model QA9545-A	each	#N/B
AGCL-9565-P	Calibrate Model QA9565-P (No Probe)	each	#N/B
AGCL-9565-X	Calibrate Model QA9565-X (No Probe)	each	#N/B
AGCL-960	Calibrate Model 960	each	#N/B
AGCL-962	Calibrate Model 962	each	#N/B
AGCL-964	Calibrate Model 964	each	#N/B
AGCL-966	Calibrate Model 966	each	#N/B
AGCL-980	Calibrate Model 980	each	#N/B
AGCL-982	Calibrate Model 982	each	#N/B
AGCL-982CO	Calibrate Model 982CO	each	#N/B
AGCL-984	Calibrate Model 984	each	#N/B
AGCL-985	Calibrate Model 985	each	#N/B
AGCL-986	Calibrate Model 986	each	#N/B
AGCL-ABT	Calibrate ABT Analog Balometer	each	#N/B
AGCL-HM675	Calibration of QAHM675	each	#N/B
AG-CL-HM685	Calibration of QAHM685	each	#N/B
AGCL-LCA501	Calibration of Model QALCA501	each	#N/B
AGCL-PANDA	Calibration of PANDA Leakage Tester	each	#N/B
AGCL-PH730	Calibrate Pro-Hood PH730	each	#N/B
AGCL-PH731	Calibrate QAPH731	each	#N/B
AGCL-PVM610	Calibration of Model QAPVM610	each	#N/B
AGCL-PVM620	Calibration of Model QAPVM620	each	#N/B
AGCQT1501	1500g Core Portable Balance	each	#N/B
AGCQT2000	2000g Core Portable Balance	each	#N/B
AGCQT202	200g Core Portable Balance	each	#N/B
AGCQT251	250g Core Portable Balance	each	#N/B
AGCQT2601	2600g Core Portable Balance	each	#N/B

Part #	Description	Unit	Price €
AGCQT5000	5000g Core Portable Balance	each	#N/B
AGCQT601	600g Core Portable Balance	each	#N/B
AGDCT2000	2000g Dune Compact Balance	each	#N/B
AGDCT201	200g Dune Compact Balance	each	#N/B
AGDCT5000	5000g Dune Compact Balance	each	#N/B
AG-E2-100G	E2 100g	each	#N/B
AG-E2-1G-200G	E2 1g-200g Weight set	each	#N/B
AG-E2-200G	E2 200g	each	#N/B
AGE401-1	Gold Wire 0.2mm dia 1m Length	each	118,00
AGE401-3	Gold Wire 0.2mm dia 3m Length	each	324,30
AGE402-1	Gold/Palladium Wire 0.2mm dia 1m Length	each	324,60
AGE402-3	Gold/Palladium Wire 0.2mm dia 3m Length	each	o.a.
AGE403-1	80% Platinum/Iridium Wire 0.2mm dia 1m Length	each	241,20
AGE403-3	80% Platinum/Iridium Wire 0.2mm dia 3m Length	each	321,20
AGE404	Platinum Wire 0.1mm dia 1m Length	each	92,50
AGE404-1	Platinum Wire 0.1mm dia 3m Length	each	273,10
AGE404-2	Platinum Wire 0.2mm dia 1m Length	each	151,80
AGE404-3	Platinum Wire 0.2mm dia 3m Length	each	406,80
AGE404-5	Platinum Wire 0.2mm dia 25m Length	each	o.a.
AGE405	Aluminium Wire 0.25mm dia 3m Length	each	58,40
AGE406	Silver Wire 0.25mm dia 3m Length	each	41,80
AGE407-1	Tungsten Wire 1mm dia x 10m Length	each	114,90
AGE407-3	Tungsten Wire 0.37mm dia 10m Length	each	42,70
AGE407-4	Tungsten Wire 0.125mm dia 25m Length	each	40,60
AGE407-5	Tungsten Wire 0.2mm dia x 10 Metre Length	each	121,30
AGE407-6	Tungsten Wire 0.25mm dia x 10 Metre Length	each	126,10
AGE408	Molybdenum Sheet	each	112,40
AGE409-2	Carbon Rods 6.3mm x 150mm	10	25,30
AGE412-1	80% Platinum/Palladium Wire 0.2mm dia 1m Length	each	118,30
AGE412-3	80% Platinum/Palladium Wire 0.2mm dia 3m Length	each	311,40
AGE413	Tantalum Sheet	each	104,30
AGE415	Spectrographically Pure Carbon Rods 3.05mm x 300mm	10	158,90
AGE416	Spectrographically Pure Carbon Rods 3.05mm x 300mm	10	137,80
AGE417	Spectrographically Pure Carbon Rods 4.6mm x 300mm	10	169,30
AGE418	Spectrographically Pure Carbon Rods 6.15mm x 300mm	10	201,00
AGE419E	Carbon Fibre 0.27g/m (5 Metre Length)	each	39,10
AGE420	Tungsten Filament 40 Amp	10	74,10
AGE428	Carbon Fibre 1.55g/m 2.5mm dia x 3m Length	3m	91,80
AGE429	Shaped Carbon Rods 6.15mm x 55mm	10	58,20
AGE430	Shaped Carbon Rods 4.6mm x 55mm	10	55,50
AGE431	Shaped Carbon Rods 3.05mm x 55mm	10	44,60
AGE432	Indium Wire	each	179,10
AGE435	Chromium chips	100g	295,50
AGE436-1	80% Platinum/Iridium Wire 0.25mm dia 1m Length	1m	411,60
AGE436-3	80% Platinum/Iridium Wire 0.25mm dia 3m Length	3m	1.054,60
AGE437	75mm x 0.05mm Platinum Boat	each	278,10
AGE438	75mm x 0.05mm Tantalum Boats	each	88,80
AGE439	32mm x 0.05mm Tungsten Boats	each	91,70
AGE441	Shaped Carbon Rods (55mm long) for Agar Coaters	10	57,60
AGE443	Carbon Fibre 0.7g/m 1.0mm dia	1m	38,00
AGE444	Nichrome Wire - 0.25mm x 30m	30m	23,50
AGE445	Carbon Fibre 0.4g/m 1mm dia	1m	33,30
AGE446	Molybdenum boats for aperture cleaning	each	45,30
AGE447-1	Palladium Wire 0.2mm dia	1m	275,70
AGE447-25	Palladium Wire 0.2mm dia	25m	o.a.
AGE447-3	Palladium Wire 0.2mm dia	3m	174,60
AGE448	Carbon Fibre for Balzers/Bal-Tec Coaters	3,5m	146,10
AGE450	29.64 x 14.30mm, High Carbon Crucibles	5	o.a.

Part #	Description	Unit	Price €
AGE451	29.64 x 14.30mm, Vitrified Carbon Crucibles	5	o.a.
AGE452	50.80 x 26.97mm, High Carbon Crucibles	5	o.a.
AGE453	50.80 x 26.97mm, Vitrified Carbon Crucibles	5	o.a.
AGE454	48.01 x 17.78mm, High Carbon Crucibles	5	o.a.
AGE455	48.01 x 17.78mm, Vitrified Carbon Crucibles	5	o.a.
AGE456	22.48 x 15.11µm High Carbon Crucibles	5	o.a.
AGEBL1023E	1020g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL1023I	1020g Eclipse Precision Balance, Pan Size 160mm dia -Interna	each	#N/B
AGEBL104E	100g Eclipse Precision Balance Pan Size 90mm dia	each	#N/B
AGEBL104I	100g Eclipse Precision Balance, Pan Size 90mm dia -Internall	each	#N/B
AGEBL12001E	12000g Eclipse Precision Balance Pan Size 400x300mm	each	#N/B
AGEBL16001E	16000g Eclipse Precision Balance Pan Size 400x300mm	each	#N/B
AGEBL1602E	1600g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL1602I	1600g Eclipse Precision Balance, Pan Size 160mm dia-Internal	each	#N/B
AGEBL1623E	1620g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL1623I	1620g Eclipse Precision Balance, Pan Size 160mm dia -Interna	each	#N/B
AGEBL164E	160g Eclipse Precision Balance Pan Size 90mm dia	each	#N/B
AGEBL164I	160g Eclipse Precision Balance, Pan Size 90mm dia -Internall	each	#N/B
AGEBL214E	210g Eclipse Precision Balance Pan Size 90mm dia	each	#N/B
AGEBL214I	210g Eclipse Precision Balance, Pan Size 90mm dia-Internally	each	#N/B
AGEBL22001E	22000g Eclipse Precision Balance Pan Size 400x300mm	each	#N/B
AGEBL223E	220g Eclipse Precision Balance Pan Size 120mm dia	each	#N/B
AGEBL223I	220g Eclipse Precision Balance, Pan Size 120mm dia-Internall	each	#N/B
AGEBL254E	250g Eclipse Precision Balance Pan Size 90mm dia	each	#N/B
AGEBL254I	250g Eclipse Precision Balance, Pan Size 90mm dia-Internally	each	#N/B
AGEBL2602E	2600g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL2602I	2600g Eclipse Precision Balance, Pan Size 160mm dia-Internal	each	#N/B
AGEBL314E	310g Eclipse Precision Balance Pan Size 90mm dia	each	#N/B
AGEBL314I	310g Eclipse Precision Balance, Pan Size 90mm dia -Internall	each	#N/B
AGEBL32001E	32000g Eclipse Precision Balance Pan Size 400x300mm	each	#N/B
AGEBL3602E	3600g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL3602I	3600g Eclipse Precision Balance, Pan Size 160mm dia-Internal	each	#N/B
AGEBL4201E	4200g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL423E	420g Eclipse Precision Balance Pan Size 120mm dia	each	#N/B
AGEBL423I	420g Eclipse Precision Balance, Pan Size 120mm dia-Internall	each	#N/B
AGEBL4602E	4600g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL4602I	4600g Eclipse Precision Balance, Pan Size 160mm dia-Internal	each	#N/B
AGEBL6201E	6200g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL6202E	6200g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL6202I	6200g Eclipse Precision Balance, Pan Size 160mm dia-Internal	each	#N/B
AGEBL623E	620g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL623I	620g Eclipse Precision Balance, Pan Size 160mm dia-Internall	each	#N/B
AGEBL8201E	8200g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL823E	820g Eclipse Precision Balance Pan Size 160mm dia	each	#N/B
AGEBL823I	820g Eclipse Precision Balance, Pan Size 160mm dia-Internall	each	#N/B
AGEMT2-PBH	Stereo microscope for bulk analysis	each	1.521,20
AGF1100-100	Pyramidal silicon etched AFM probe Al BS coated with 4 lever	100	2.155,40
AGF1100-15	Pyramidal silicon etched AFM probe Al BS coated with 4 lever	15	500,40
AGF1100-200	Pyramidal silicon etched AFM probe Al BS coated with 4 lever	200	3.708,00
AGF1100-400	Pyramidal silicon etched AFM probe Al BS coated with 4 lever	400	5.561,90
AGF1100-50	Pyramidal silicon etched AFM probe Al BS coated with 4 lever	50	1.238,70
AG-F1-100G	F1 100g	each	#N/B
AG-F1-10G	F1 10g	each	#N/B
AGF1110-100	Pyramidal silicon etched AFM probe no Al coated with 4 lever	100	2.155,40
AGF1110-15	Pyramidal silicon etched AFM probe no Al coated with 4 lever	15	500,40
AGF1110-200	Pyramidal silicon etched AFM probe no Al coated with 4 lever	200	3.708,00
AGF1110-400	Pyramidal silicon etched AFM probe no Al coated with 4 lever	400	5.561,90
AGF1110-50	Pyramidal silicon etched AFM probe no Al coated with 4 lever	50	1.238,70

Part #	Description	Unit	Price €
AGF1120-100	Pyramidal silicon etched AFM probe Pt coated with 4 lever	100	2.155,40
AGF1120-15	Pyramidal silicon etched AFM probe Pt coated with 4 lever	15	500,40
AGF1120-50	Pyramidal silicon etched AFM probe Pt coated with 4 lever	50	1.238,70
AGF1140-100	Pyramidal silicon etched AFM probe Al BS coated 125x25x2.1 µm	100	1.937,90
AGF1140-15	Pyramidal silicon etched AFM probe Al BS coated 125x25x2.1 µm	15	414,60
AGF1140-200	Pyramidal silicon etched AFM probe Al BS coated 125x25x2.1 µm	200	3.151,80
AGF1140-400	Pyramidal silicon etched AFM probe Al BS coated 125x25x2.1 µm	400	4.727,60
AGF1140-50	Pyramidal silicon etched AFM probe Al BS coated 125x25x2.1 µm	50	1.052,90
AGF1150-100	Pyramidal silicon etched AFM probe Cr-Au coated 125x25x2.1 µm	100	1.939,80
AGF1150-15	Pyramidal silicon etched AFM probe Cr-Au coated 125x25x2.1 µm	15	462,00
AGF1150-50	Pyramidal silicon etched AFM probe Cr-Au coated 125x25x2.1 µm	50	1.114,90
AGF1160-100	Pyramidal silicon etched AFM probe no Al coated 125x25x2.1 µm	100	1.832,10
AGF1160-15	Pyramidal silicon etched AFM probe no Al coated 125x25x2.1 µm	15	414,60
AGF1160-200	Pyramidal silicon etched AFM probe no Al coated 125x25x2.1 µm	200	3.151,80
AGF1160-400	Pyramidal silicon etched AFM probe no Al coated 125x25x2.1 µm	400	4.727,60
AGF1160-50	Pyramidal silicon etched AFM probe no Al coated 125x25x2.1 µm	50	1.052,90
AGF1170-100	Pyramidal silicon etched AFM probe Pt coated 125x25x2.1 µm	100	1.939,80
AGF1170-15	Pyramidal silicon etched AFM probe Pt coated 125x25x2.1 µm	15	462,00
AGF1170-50	Pyramidal silicon etched AFM probe Pt coated 125x25x2.1 µm	50	1.114,90
AGF1180-100	Pyramidal silicon etched AFM probe AL BS coated 125x30x4.0 µm	100	1.965,30
AGF1180-15	Pyramidal silicon etched AFM probe AL BS coated 125x30x4.0 µm	15	414,60
AGF1180-200	Pyramidal silicon etched AFM probe AL BS coated 125x30x4.0 µm	200	3.151,80
AGF1180-400	Pyramidal silicon etched AFM probe AL BS coated 125x30x4.0 µm	400	4.727,60
AGF1180-50	Pyramidal silicon etched AFM probe AL BS coated 125x30x4.0 µm	50	1.052,90
AGF1190-100	Pyramidal silicon etched AFM probe Cr-Au coated 125x30x4.0 µm	100	1.939,80
AGF1190-15	Pyramidal silicon etched AFM probe Cr-Au coated 125x30x4.0 µm	15	462,00
AGF1190-50	Pyramidal silicon etched AFM probe Cr-Au coated 125x30x4.0 µm	50	1.114,90
AG-F1-1G	F1 1g	each	#N/B
AG-F1-1G-100G	F1 1g-100g Weight set	each	#N/B
AG-F1-1G-200G	F1 1g-200g Weight set	each	#N/B
AG-F1-1G-2KG	F1 1g-2kg Weight set	each	#N/B
AG-F1-1G-500G	F1 1g-500g Weight set	each	#N/B
AG-F1-1KG	F1	1kg	#N/B
AGF1200-15	Pyramidal silicon etched AFM probe Cr-Au BS coated 125x30x4.0 µm	15	500,40
AGF1200-50	Pyramidal silicon etched AFM probe Cr-Au BS coated 125x30x4.0 µm	50	1.238,70
AG-F1-200G	F1 200g	each	#N/B
AG-F1-20G	F1 20g	each	#N/B
AGF1210-100	Pyramidal silicon etched AFM probe no Al coated 125x30x4.0 µm	100	1.832,10
AGF1210-15	Pyramidal silicon etched AFM probe no Al coated 125x30x4.0 µm	15	414,60
AGF1210-200	Pyramidal silicon etched AFM probe no Al coated 125x30x4.0 µm	200	3.151,80
AGF1210-400	Pyramidal silicon etched AFM probe no Al coated 125x30x4.0 µm	400	4.727,60
AGF1210-50	Pyramidal silicon etched AFM probe no Al coated 125x30x4.0 µm	50	1.052,90
AGF1220-100	Pyramidal silicon etched AFM probe Pt coated 125x30x4.0 µm	100	1.939,80
AGF1220-15	Pyramidal silicon etched AFM probe Pt coated 125x30x4.0 µm	15	462,00
AGF1220-50	Pyramidal silicon etched AFM probe Pt coated 125x30x4.0 µm	50	1.114,90
AGF1230-100	Pyramidal silicon etched AFM probe Al BS coated 225x37.5x7.0 µm	100	1.832,10
AGF1230-15	Pyramidal silicon etched AFM probe Al BS coated 225x37.5x7.0 µm	15	418,40
AGF1230-200	Pyramidal silicon etched AFM probe Al BS coated 225x37.5x7.0 µm	200	3.151,80
AGF1230-400	Pyramidal silicon etched AFM probe Al BS coated 225x37.5x7.0 µm	400	4.727,60
AGF1230-50	Pyramidal silicon etched AFM probe Al BS coated 225x37.5x7.0 µm	50	1.052,90
AGF1240-100	Pyramidal silicon etched AFM probe Cr-Au coated 225x37.5x7.0 µm	100	1.939,80
AGF1240-15	Pyramidal silicon etched AFM probe Cr-Au coated 225x37.5x7.0 µm	15	462,00
AGF1240-50	Pyramidal silicon etched AFM probe Cr-Au coated 225x37.5x7.0 µm	50	1.114,90
AGF1250-100	Pyramidal silicon etched AFM probe no Al coated 225x37.5x7.0 µm	100	1.832,10
AGF1250-15	Pyramidal silicon etched AFM probe no Al coated 225x37.5x7.0 µm	15	414,60
AGF1250-200	Pyramidal silicon etched AFM probe no Al coated 225x37.5x7.0 µm	200	3.151,80
AGF1250-400	Pyramidal silicon etched AFM probe no Al coated 225x37.5x7.0 µm	400	4.727,60
AGF1250-50	Pyramidal silicon etched AFM probe no Al coated 225x37.5x7.0 µm	50	1.052,90

Part #	Description	Unit	Price €
AGF1260-100	Pyramidal silicon etched AFM probe Al BS coated 450x50x2.0 µm	100	1.832,10
AGF1260-15	Pyramidal silicon etched AFM probe Al BS coated 450x50x2.0 µm	15	414,60
AGF1260-200	Pyramidal silicon etched AFM probe Al BS coated 450x50x2.0 µm	200	3.151,80
AGF1260-400	Pyramidal silicon etched AFM probe Al BS coated 450x50x2.0 µm	400	4.727,60
AGF1260-50	Pyramidal silicon etched AFM probe Al BS coated 450x50x2.0 µm	50	1.100,70
AGF1270-100	Pyramidal silicon etched AFM probe Cr-Au coated 450x50x2.0 µm	100	1.939,80
AGF1270-15	Pyramidal silicon etched AFM probe Cr-Au coated 450x50x2.0 µm	15	462,00
AGF1270-50	Pyramidal silicon etched AFM probe Cr-Au coated 450x50x2.0 µm	50	1.114,90
AGF1280-100	Pyramidal silicon etched AFM probe no Al coated 450x50x2.0 µm	100	1.832,10
AGF1280-15	Pyramidal silicon etched AFM probe no Al coated 450x50x2.0 µm	15	414,60
AGF1280-200	Pyramidal silicon etched AFM probe no Al coated 450x50x2.0 µm	200	3.151,80
AGF1280-400	Pyramidal silicon etched AFM probe no Al coated 450x50x2.0 µm	400	4.727,60
AGF1280-50	Pyramidal silicon etched AFM probe no Al coated 450x50x2.0 µm	50	1.052,90
AGF1290-100	Pyramidal silicon etched AFM probe Pt coated 450x50x2.0 µm	100	1.939,80
AGF1290-15	Pyramidal silicon etched AFM probe Pt coated 450x50x2.0 µm	15	462,00
AGF1290-50	Pyramidal silicon etched AFM probe Pt coated 450x50x2.0 µm	50	1.266,00
AG-F1-2G	F1 2g	each	#N/B
AG-F1-2KG	F1 2kg	each	#N/B
AGF1300-100	Pyramidal silicon etched AFM probe Al BS coated 225x27.5x3.0 µm	100	1.832,10
AGF1300-15	Pyramidal silicon etched AFM probe Al BS coated 225x27.5x3.0 µm	15	414,60
AGF1300-200	Pyramidal silicon etched AFM probe Al BS coated 225x27.5x3.0 µm	200	3.151,80
AGF1300-400	Pyramidal silicon etched AFM probe Al BS coated 225x27.5x3.0 µm	400	4.727,60
AGF1300-50	Pyramidal silicon etched AFM probe Al BS coated 225x27.5x3.0 µm	50	1.100,70
AGF1310-15	Magnetic Noncontact silicon probe Al BS coated 225x27.5x3.0 µm	15	510,90
AGF1310-50	Magnetic Noncontact silicon probe Al BS coated 225x27.5x3.0 µm	50	1.238,70
AGF1320-100	Pyramidal silicon etched AFM probe Cr-Au coated 225x27.5x3.0 µm	100	1.939,80
AGF1320-15	Pyramidal silicon etched AFM probe Cr-Au coated 225x27.5x3.0 µm	15	462,00
AGF1320-50	Pyramidal silicon etched AFM probe Cr-Au coated 225x27.5x3.0 µm	50	1.227,70
AGF1330-15	Pyramidal silicon etched AFM probe Cr-Au BS coated 225x27.5x3.0 µm	15	500,40
AGF1330-50	Pyramidal silicon etched AFM probe Cr-Au BS coated 225x27.5x3.0 µm	50	1.238,70
AGF1340-100	Pyramidal silicon etched AFM probe no Al coated 225x27.5x3.0 µm	100	1.832,10
AGF1340-15	Pyramidal silicon etched AFM probe no Al coated 225x27.5x3.0 µm	15	414,60
AGF1340-200	Pyramidal silicon etched AFM probe no Al coated 225x27.5x3.0 µm	200	3.151,80
AGF1340-400	Pyramidal silicon etched AFM probe no Al coated 225x27.5x3.0 µm	400	4.727,60
AGF1340-50	Pyramidal silicon etched AFM probe no Al coated 225x27.5x3.0 µm	50	1.052,90
AGF1350-100	Pyramidal silicon etched AFM probe Pt coated 225x27.5x3.0 µm	100	1.939,80
AGF1350-15	Pyramidal silicon etched AFM probe Pt coated 225x27.5x3.0 µm	15	462,00
AGF1350-50	Pyramidal silicon etched AFM probe Pt coated 225x27.5x3.0 µm	50	1.266,00
AGF1360-100	Pyramidal silicon etched AFM probe Al BS coated 125x22.5x1.0 µm	100	1.832,10
AGF1360-15	Pyramidal silicon etched AFM probe Al BS coated 125x22.5x1.0 µm	15	414,60
AGF1360-200	Pyramidal silicon etched AFM probe Al BS coated 125x22.5x1.0 µm	200	3.151,80
AGF1360-400	Pyramidal silicon etched AFM probe Al BS coated 125x22.5x1.0 µm	400	4.727,60
AGF1360-50	Pyramidal silicon etched AFM probe Al BS coated 125x22.5x1.0 µm	50	1.052,90
AGF1370-100	Pyramidal silicon etched AFM probe Cr-Au coated 125x22.5x1.0 µm	100	1.939,80
AGF1370-15	Pyramidal silicon etched AFM probe Cr-Au coated 125x22.5x1.0 µm	15	462,00
AGF1370-50	Pyramidal silicon etched AFM probe Cr-Au coated 125x22.5x1.0 µm	50	1.114,90
AGF1380-100	Pyramidal silicon etched AFM probe no Al coated 125x22.5x1.0 µm	100	1.832,10
AGF1380-15	Pyramidal silicon etched AFM probe no Al coated 125x22.5x1.0 µm	15	414,60
AGF1380-200	Pyramidal silicon etched AFM probe no Al coated 125x22.5x1.0 µm	200	3.151,80
AGF1380-400	Pyramidal silicon etched AFM probe no Al coated 125x22.5x1.0 µm	400	4.727,60
AGF1380-50	Pyramidal silicon etched AFM probe no Al coated 125x22.5x1.0 µm	50	1.052,90
AGF1390-100	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1390-15	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	15	414,60
AGF1390-200	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1390-400	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1390-50	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1400-100	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	100	1.939,80
AGF1400-15	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	15	462,00

Part #	Description	Unit	Price €
AGF1400-50	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	50	1.114,90
AGF1410-15	Pyramidal silicon etched AFM probe Cr-Au BS coated with 3 lever	15	500,40
AGF1410-50	Pyramidal silicon etched AFM probe Cr-Au BS coated with 3 lever	50	1.238,70
AGF1420-100	Pyramidal silicon etched AFM probe no Al coated with 3 lever	100	1.832,10
AGF1420-15	Pyramidal silicon etched AFM probe no Al coated with 3 lever	15	414,60
AGF1420-200	Pyramidal silicon etched AFM probe no Al coated with 3 lever	200	3.151,80
AGF1420-400	Pyramidal silicon etched AFM probe no Al coated with 3 lever	400	4.727,60
AGF1420-50	Pyramidal silicon etched AFM probe no Al coated with 3 lever	50	1.052,90
AGF1430-100	Pyramidal silicon etched AFM probe Pt coated with 3 lever	100	1.939,80
AGF1430-15	Pyramidal silicon etched AFM probe Pt coated with 3 lever	15	462,00
AGF1430-50	Pyramidal silicon etched AFM probe Pt coated with 3 lever	50	1.114,90
AGF1440-100	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1440-15	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	15	414,60
AGF1440-200	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1440-400	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1440-50	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1450-15	Magnetic Noncontact silicon probe Al BS coated with 3 lever	15	500,40
AGF1450-50	Magnetic Noncontact silicon probe Al BS coated with 3 lever	50	1.238,70
AGF1460-100	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	100	1.939,80
AGF1460-15	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	15	462,00
AGF1460-50	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	50	1.114,90
AGF1470-15	Pyramidal silicon etched AFM probe Cr-Au BS coated with 3 lever	15	505,10
AGF1470-50	Pyramidal silicon etched AFM probe Cr-Au BS coated with 3 lever	50	1.238,70
AGF1480-100	Pyramidal silicon etched AFM probe no Al coated with 3 lever	100	1.832,10
AGF1480-15	Pyramidal silicon etched AFM probe no Al coated with 3 lever	15	414,60
AGF1480-200	Pyramidal silicon etched AFM probe no Al coated with 3 lever	200	3.151,80
AGF1480-400	Pyramidal silicon etched AFM probe no Al coated with 3 lever	400	4.727,60
AGF1480-50	Pyramidal silicon etched AFM probe no Al coated with 3 lever	50	1.052,90
AGF1490-100	Pyramidal silicon etched AFM probe Pt coated with 3 lever	100	1.939,80
AGF1490-15	Pyramidal silicon etched AFM probe Pt coated with 3 lever	15	462,00
AGF1490-50	Pyramidal silicon etched AFM probe Pt coated with 3 lever	50	1.114,90
AGF1500-100	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1500-15	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	15	414,60
AGF1500-200	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1500-400	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1500-50	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	50	1.052,90
AG-F1-500G	F1 500g	500g	#N/B
AG-F1-50G	F1 50g	50g	#N/B
AGF1510-100	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	100	1.939,80
AGF1510-15	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	15	462,00
AGF1510-50	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	50	1.114,90
AGF1520-100	Pyramidal silicon etched AFM probe no Al coated with 3 lever	100	1.832,10
AGF1520-15	Pyramidal silicon etched AFM probe no Al coated with 3 lever	15	414,60
AGF1520-200	Pyramidal silicon etched AFM probe no Al coated with 3 lever	200	3.151,80
AGF1520-400	Pyramidal silicon etched AFM probe no Al coated with 3 lever	400	4.727,60
AGF1520-50	Pyramidal silicon etched AFM probe no Al coated with 3 lever	50	1.052,90
AGF1530-100	Pyramidal silicon etched AFM probe Pt coated with 3 lever	100	1.939,80
AGF1530-15	Pyramidal silicon etched AFM probe Pt coated with 3 lever	15	462,00
AGF1530-50	Pyramidal silicon etched AFM probe Pt coated with 3 lever	50	1.114,90
AGF1540-100	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1540-15	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	15	414,60
AGF1540-200	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1540-400	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1540-50	Pyramidal silicon etched AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1550-100	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	100	1.939,80
AGF1550-15	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	15	462,00
AGF1550-50	Pyramidal silicon etched AFM probe Cr-Au coated with 3 lever	50	1.114,90
AGF1560-100	Pyramidal silicon etched AFM probe no Al coated with 3 lever	100	1.832,10

Part #	Description	Unit	Price €
AGF1560-15	Pyramidal silicon etched AFM probe no Al coated with 3 lever	15	414,60
AGF1560-200	Pyramidal silicon etched AFM probe no Al coated with 3 lever	200	3.151,80
AGF1560-400	Pyramidal silicon etched AFM probe no Al coated with 3 lever	400	4.727,60
AGF1560-50	Pyramidal silicon etched AFM probe no Al coated with 3 lever	50	1.052,90
AGF1570-100	Tipless Noncontact AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1570-15	Tipless Noncontact AFM probe Al BS coated with 3 lever	15	414,60
AGF1570-200	Tipless Noncontact AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1570-400	Tipless Noncontact AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1570-50	Tipless Noncontact AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1580-100	Tipless Noncontact AFM probe Cr Au coated with 3 lever	100	1.939,80
AGF1580-15	Tipless Noncontact AFM probe Cr Au coated with 3 lever	15	462,00
AGF1580-50	Tipless Noncontact AFM probe Cr Au coated with 3 lever	50	1.114,90
AGF1590-100	Tipless Noncontact AFM probe no Al coated with 3 lever	100	1.832,10
AGF1590-15	Tipless Noncontact AFM probe no Al coated with 3 lever	15	414,60
AGF1590-200	Tipless Noncontact AFM probe no Al coated with 3 lever	200	3.151,80
AGF1590-400	Tipless Noncontact AFM probe no Al coated with 3 lever	400	4.727,60
AGF1590-50	Tipless Noncontact AFM probe no Al coated with 3 lever	50	1.052,90
AG-F1-5G	F1 5g	5g	#N/B
AG-F1-5KG	F1 5kg	each	#N/B
AGF1600-100	Tipless Noncontact AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1600-15	Tipless Noncontact AFM probe Al BS coated with 3 lever	15	414,60
AGF1600-200	Tipless Noncontact AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1600-400	Tipless Noncontact AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1600-50	Tipless Noncontact AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1610-100	Tipless Noncontact AFM probe Cr Au coated with 3 lever	100	1.939,80
AGF1610-15	Tipless Noncontact AFM probe Cr Au coated with 3 lever	15	462,00
AGF1610-50	Tipless Noncontact AFM probe Cr Au coated with 3 lever	50	1.114,90
AGF1620-100	Tipless Noncontact AFM probe no Al coated with 3 lever	100	1.832,10
AGF1620-15	Tipless Noncontact AFM probe no Al coated with 3 lever	15	414,60
AGF1620-200	Tipless Noncontact AFM probe no Al coated with 3 lever	200	3.151,80
AGF1620-400	Tipless Noncontact AFM probe no Al coated with 3 lever	400	4.727,60
AGF1620-50	Tipless Noncontact AFM probe no Al coated with 3 lever	50	1.052,90
AGF1630-100	Tipless Noncontact AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1630-15	Tipless Noncontact AFM probe Al BS coated with 3 lever	15	414,60
AGF1630-200	Tipless Noncontact AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1630-400	Tipless Noncontact AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1630-50	Tipless Noncontact AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1640-100	Tipless Noncontact AFM probe Cr Au coated with 3 lever	100	1.939,80
AGF1640-15	Tipless Noncontact AFM probe Cr Au coated with 3 lever	15	462,00
AGF1640-50	Tipless Noncontact AFM probe Cr Au coated with 3 lever	50	1.114,90
AGF1650-100	Tipless Noncontact AFM probe no Al coated with 3 lever	100	1.832,10
AGF1650-15	Tipless Noncontact AFM probe no Al coated with 3 lever	15	414,60
AGF1650-200	Tipless Noncontact AFM probe no Al coated with 3 lever	200	3.151,80
AGF1650-400	Tipless Noncontact AFM probe no Al coated with 3 lever	400	4.727,60
AGF1650-50	Tipless Noncontact AFM probe no Al coated with 3 lever	50	1.052,90
AGF1660-100	Tipless Noncontact AFM probe Al BS coated with 3 lever	100	1.832,10
AGF1660-15	Tipless Noncontact AFM probe Al BS coated with 3 lever	15	414,60
AGF1660-200	Tipless Noncontact AFM probe Al BS coated with 3 lever	200	3.151,80
AGF1660-400	Tipless Noncontact AFM probe Al BS coated with 3 lever	400	4.727,60
AGF1660-50	Tipless Noncontact AFM probe Al BS coated with 3 lever	50	1.052,90
AGF1670-100	Tipless Noncontact AFM probe Cr Au coated with 3 lever	100	1.939,80
AGF1670-15	Tipless Noncontact AFM probe Cr Au coated with 3 lever	15	462,00
AGF1670-50	Tipless Noncontact AFM probe Cr Au coated with 3 lever	50	1.114,90
AGF1680-100	Tipless Noncontact AFM probe no Al coated with 3 lever	100	1.832,10
AGF1680-15	Tipless Noncontact AFM probe no Al coated with 3 lever	15	414,60
AGF1680-200	Tipless Noncontact AFM probe no Al coated with 3 lever	200	3.151,80
AGF1680-400	Tipless Noncontact AFM probe no Al coated with 3 lever	400	4.727,60
AGF1680-50	Tipless Noncontact AFM probe no Al coated with 3 lever	50	1.052,90

Part #	Description	Unit	Price €
AGF1690-5	High Res silicon probe Cr Au coated 125x25x2.1 µm	5	843,70
AGF1700-5	High Resolution silicon probe Cr Au coated 125x30x4.0 µm	5	675,30
AGF1710-5	High Resolution silicon probe Cr Au coated 225x37.5x7.0 µm	5	843,70
AGF1720-15	Hardened DLC silicon probe Al BS coated with 4 lever	15	500,40
AGF1720-50	Hardened DLC silicon probe Al BS coated with 4 lever	50	1.238,70
AGF1730-15	Hardened DLC silicon probe Al BS coated 125x25x2.1 µm	15	500,40
AGF1730-50	Hardened DLC silicon probe Al BS coated 125x25x2.1 µm	50	1.364,10
AGF1740-15	Hardened DLC silicon probe Al BS coated 125x30x4.0 µm	15	500,40
AGF1740-50	Hardened DLC silicon probe Al BS coated 125x30x4.0 µm	50	1.238,70
AGF1750-15	Hardened DLC silicon probe Al BS coated 225x37.5x7.0 µm	15	500,40
AGF1750-50	Hardened DLC silicon probe Al BS coated 225x37.5x7.0 µm	50	1.238,70
AGF1760-15	Hardened DLC silicon probe Al BS coated 450x50x2.0 µm	15	500,40
AGF1760-50	Hardened DLC silicon probe Al BS coated 450x50x2.0 µm	50	500,40
AGF1770-15	Hardened DLC silicon probe Al BS coated 225x27.5x27.5 µm	15	500,40
AGF1770-50	Hardened DLC silicon probe Al BS coated 225x27.5x27.5 µm	50	1.238,70
AGF1780-15	Hardened DLC silicon probe Al BS coated with 3 lever	15	500,40
AGF1780-50	Hardened DLC silicon probe Al BS coated with 3 lever	50	1.238,70
AGF1790-15	Hardened DLC silicon probe Al BS coated with 3 lever	15	500,40
AGF1790-50	Hardened DLC silicon probe Al BS coated with 3 lever	50	1.238,70
AGF1800-1-1	HOPG grades	each	565,90
AGF1800-2-1	HOPG grades	each	846,40
AGF1810-1-1	HOPG grades	each	176,60
AGF1820-1-1	HOPG grades	each	103,10
AGF1820-1-5	HOPG grades	each	323,80
AGF1820-2-1	HOPG grades	each	125,20
AGF1820-2-5	HOPG grades	each	551,10
AGF1830-1	Sample for characterisation of tip shape mounted on a round metal plate 12mm	each	290,50
AGF1830-2	Sample for characterisation of tip shape, without metal plate	each	281,50
AGF1840-1	Calibration gratings mounted on a round metal plate 12mm	each	281,90
AGF1840-2	Calibration gartings, without metal plate	each	281,50
AGF1850	Piezoresistive 300x100µm self sensing silicon cantilever, 6x5mm PCB, 30-65kHz	each	#N/B
AGF1851	Piezoresistive 300x100µm self sensing silicon cantilever, 6x5mm PCB 65-95kHz;	each	#N/B
AGF1852	Piezoresistive 300x100µm self sensing silicon cantilever, 30-65kHz (chip only, \	each	#N/B
AGF1853	Piezoresistive 300x100µm self sensing silicon cantilever, 65-95kHz (chip only, \	each	#N/B
AGF1855	Piezoresistive 70x30µm self sensing silicon cantilever, 6x5mm PCB 600-1300kt	each	#N/B
AGF1857	Piezoresistive 70x30µm self sensing silicon cantilever, 600-1300kHz (chip only,	each	#N/B
AGF1858	Self-sensing cantilever starter kit plus	each	#N/B
AGF1859	Self-sensing cantilever starter kit	each	#N/B
AGF1860	AFM Cantilever Pre-Amplifier	each	#N/B
AGF1861	Counter part PCB (with 10 pin connector)	each	#N/B
AGF1862	Counter part PCB soldered to a flat flex cable (with 10 pin connector)	each	#N/B
AGF1863	Output cable for pre-amplifier	each	#N/B
AGF1864	Piezoresistive 100x48um self sensing silicon cantilever, with PCB 250-550kHz	each	#N/B
AGF1865	Piezoresistive 100x48um self sensing silicon cantilever, 250-550kHz (chip only,	each	#N/B
AGF1866	Piezoresistive 70x30µm self-sensing silicon cantilever, with 6x5mm PCB with 1	each	#N/B
AGF1867	Piezoresistive 70x30µm self-sensing silicon cantilever, (chip only, without PCB)	each	#N/B
AGF7001	SPM Specimen Discs, 12mm	50	#N/B
AGF7003	SPM Specimen Discs, 15mm	50	#N/B
AGF7005	Disc Gripper 12mm	each	#N/B
AGF7006	Disc Gripper 15mm	each	#N/B
AGF7007	Disc Carrier for ten 12mm discs	each	#N/B
AGF7008	Disc Carrier for eight 15mm discs	each	#N/B
AGF7009	Dispenser for discs	each	#N/B
AGF7010	Magnetic Pickup Tool	each	#N/B
AGF7011	12mm dia Round Glass Coverslips, 0.13 - 0.17mm thick	100	#N/B
AGF7012	15mm dia Round Glass Coverslips, 0.13 - 0.17mm thick	100	#N/B
AGF7013	Mica Disk 10mm V1	10	#N/B
AGF7018	Extended gold calibration kit	each	#N/B

Part #	Description	Unit	Price €
AGF7019	14mm dia Mica Discs 0.1mm thick	10	#N/B
AGF7020	Gold Standard Kit	each	#N/B
AGF7021	Grating TGXYZ01 Step Height 20nm	each	#N/B
AGF7022	Grating TGXYZ02 Step Height 100nm	each	#N/B
AGF7023	Grating TGXYZ03 Step Height 500nm	each	#N/B
AGF7026	Grating Type TGX01 3 Micron Pitch	each	#N/B
AGF7033	Grating Type TGX11 10 Micron Pitch	each	#N/B
AGF7034	Set for vertical calibration	each	#N/B
AGF7035	Set of gratings and test structures	each	#N/B
AGF7036	Set for tip characterisation	each	#N/B
AGF7037	11mm x 11mm Gold Arrandee Substrates	50	#N/B
AGF7038	12mm x 12mm Gold Arrandee Substrates	50	#N/B
AGF7040	1inch x 1.5inch Mica 200nm Gold	each	#N/B
AGF7041	1inch x 3inch Mica 200nm Gold	each	#N/B
AGF7042	4inch x 525µm Si Wafers 100nm Gold	3	#N/B
AGF7043	4inch x 525µm Si Wafers 100nm Gold	12	#N/B
AGF7044	Aluminosilicate Slides 10nm Gold	each	#N/B
AGF7045	Aluminosilicate Slides 50nm Gold	each	#N/B
AGF7046	Aluminosilicate Slides 100nm Gold	each	#N/B
AGF7047	22mm Coverslips Coated 10nm Gold	12	#N/B
AGF7048	22mm Coverslips Coated 50nm Gold	12	#N/B
AGF7049	15mm diameter Coverslips Coated 10nm Gold	24	#N/B
AGF7050	AFM 144nm reference standards, certified, non-traceable, mounted on disc	each	#N/B
AGF7051	AFM 144nm reference standards, certified (certificate provided), traceable, m	each	#N/B
AGF7052	AFM 300nm reference standards, certified, non-traceable, mounted on disc	each	#N/B
AGF7054	145nm AFM Reference Standard on 12mm Disc	each	#N/B
AGF7054U	145nm AFM Reference Standard, Certified, Unmounted	each	#N/B
AGF7055	292nm AFM Reference Standard on 12mm Disc	each	#N/B
AGF705-5	High aspect ratio AFM cantilever, Au coating 300 kHz 26 N/m	each	#N/B
AGF7055T	292nm Traceable AFM Reference Standard 12mm Disc	each	#N/B
AGF706-10	AFM cantilever, Au coating 300 kHz 26 N/m	10	#N/B
AGF706-100	AFM cantilever, Au coating 300 kHz 26 N/m	100	#N/B
AGF706-24	AFM cantilever, Au coating 300 kHz 26 N/m	24	#N/B
AGF706-50	AFM cantilever, Au coating 300 kHz 26 N/m	50	#N/B
AGF707-10	AFM cantilever, Al reflective coating 300 kHz 26 N/m	10	#N/B
AGF707-100	AFM cantilever, Al reflective coating 300 kHz 26 N/m	100	#N/B
AGF707-24	AFM cantilever, Al reflective coating 300 kHz 26 N/m	24	#N/B
AGF707-50	AFM cantilever, Al reflective coating 300 kHz 26 N/m	50	#N/B
AGF708-10	AFM cantilever, Au reflective coating 300 kHz 26 N/m	10	#N/B
AGF708-100	AFM cantilever, Au reflective coating 300 kHz 26 N/m	100	#N/B
AGF708-24	AFM cantilever, Au reflective coating 300 kHz 26 N/m	24	#N/B
AGF708-50	AFM cantilever, Au reflective coating 300 kHz 26 N/m	50	#N/B
AGF709-10	AFM cantilever, no coating 300 kHz 26 N/m	10	#N/B
AGF709-100	AFM cantilever, no coating 300 kHz 26 N/m	100	#N/B
AGF709-24	AFM cantilever, no coating 300 kHz 26 N/m	24	#N/B
AGF709-50	AFM cantilever, no coating 300 kHz 26 N/m	50	#N/B
AGF710-5	High res AFM cantilever, Au coating 300 kGz 26 N/m	5	#N/B
AGF711-10	AFM cantilever, Au reflective coating 135 kHz 9 N/m	10	#N/B
AGF711-100	AFM cantilever, Au reflective coating 135 kHz 9 N/m	100	#N/B
AGF711-24	AFM cantilever, Au reflective coating 135 kHz 9 N/m	24	#N/B
AGF711-50	AFM cantilever, Au reflective coating 135 kHz 9 N/m	50	#N/B
AGF712-10	AFM cantilever, Au coating 70 kHz 2 N/m	10	#N/B
AGF712-100	AFM cantilever, Au coating 70 kHz 2 N/m	100	#N/B
AGF712-24	AFM cantilever, Au coating 70 kHz 2 N/m	24	#N/B
AGF712-5	High aspect ratio AFM cantilever, Au coating 70 kHz 2 N/m	5	#N/B
AGF712-50	AFM cantilever, Au coating 70 kHz 2 N/m	50	#N/B
AGF713-10	AFM cantilever, magnetic tip coating 70 kHz 2 N/m	10	#N/B
AGF713-100	AFM cantilever, magnetic tip coating 70 kHz 2 N/m	100	#N/B

Part #	Description	Unit	Price €
AGF713-24	AFM cantilever, magnetic tip coating 70 kHz 2 N/m	24	#N/B
AGF713-50	AFM cantilever, magnetic tip coating 70 kHz 2 N/m	50	#N/B
AGF714-10	AFM cantilever, Al reflective coating 70 kHz 2 N/m	10	#N/B
AGF714-100	AFM cantilever, Al reflective coating 70 kHz 2 N/m	100	#N/B
AGF714-24	AFM cantilever, Al reflective coating 70 kHz 2 N/m	24	#N/B
AGF714-50	AFM cantilever, Al reflective coating 70 kHz 2 N/m	50	#N/B
AGF715-10	AFM cantilever, Au reflective coating 70 kHz 2 N/m	10	#N/B
AGF715-100	AFM cantilever, Au reflective coating 70 kHz 2 N/m	100	#N/B
AGF715-24	AFM cantilever, Au reflective coating 70 kHz 2 N/m	24	#N/B
AGF715-50	AFM cantilever, Au reflective coating 70 kHz 2 N/m	50	#N/B
AGF716-10	AFM cantilever, no coating 70 kHz 2 N/m	10	#N/B
AGF716-100	AFM cantilever, no coating 70 kHz 2 N/m	100	#N/B
AGF716-24	AFM cantilever, no coating 70 kHz 2 N/m	24	#N/B
AGF716-50	AFM cantilever, no coating 70 kHz 2 N/m	50	#N/B
AGF717-10	AFM cantilever, Pt coating 70 kHz 2 N/m	10	#N/B
AGF717-100	AFM cantilever, Pt coating 70 kHz 2 N/m	100	#N/B
AGF717-24	AFM cantilever, Pt coating 70 kHz 2 N/m	24	#N/B
AGF717-50	AFM cantilever, Pt coating 70 kHz 2 N/m	50	#N/B
AGF718-5	High Res AFM cantilever, Au coating 70 kHz 2 N/m	5	#N/B
AGF719-10	Mixed AFM cantilevers Au coating 17/150/75 kHz 0.3/9/2.5 N/m	10	#N/B
AGF719-100	Mixed AFM cantilevers Au coating 17/150/75 kHz 0.3/9/2.5 N/m	100	#N/B
AGF719-24	Mixed AFM cantilevers Au coating 17/150/75 kHz 0.3/9/2.5 N/m	24	#N/B
AGF719-50	Mixed AFM cantilevers Au coating 17/150/75 kHz 0.3/9/2.5 N/m	50	#N/B
AGF720-10	Mixed AFM cantilevers, Al reflective coating, 17/150/75 kHz, 0.3/9/2.5 N/m	10	#N/B
AGF720-100	Mixed AFM cantilevers, Al reflective coating, 17/150/75 kHz, 0.3/9/2.5 N/m	100	#N/B
AGF720-24	Mixed AFM cantilevers, Al reflective coating, 17/150/75 kHz, 0.3/9/2.5 N/m	24	#N/B
AGF720-50	Mixed AFM cantilevers, Al reflective coating, 17/150/75 kHz, 0.3/9/2.5 N/m	50	#N/B
AGF721-10	Mixed AFM cantilevers, no coating 17/150/75 kHz, 0.3/9/2.5 N/m	10	#N/B
AGF721-100	Mixed AFM cantilevers, no coating 17/150/75 kHz, 0.3/9/2.5 N/m	100	#N/B
AGF721-24	Mixed AFM cantilevers, no coating 17/150/75 kHz, 0.3/9/2.5 N/m	24	#N/B
AGF721-50	Mixed AFM cantilevers, no coating 17/150/75 kHz, 0.3/9/2.5 N/m	50	#N/B
AGF722-10	Mixed AFM cantilevers, Au coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	10	#N/B
AGF722-100	Mixed AFM cantilevers, Au coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	100	#N/B
AGF722-24	Mixed AFM cantilevers, Au coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	24	#N/B
AGF722-50	Mixed AFM cantilevers, Au coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	50	#N/B
AGF723-10	Mixed AFM cantilevers, no coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	10	#N/B
AGF723-100	Mixed AFM cantilevers, no coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	100	#N/B
AGF723-24	Mixed AFM cantilevers, no coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	24	#N/B
AGF723-50	Mixed AFM cantilevers, no coating 17/75/150/1200 kHz, 0.3/2.5/9/100 N/m	50	#N/B
AGF724-10	High freq. cantilever, Al reflective coating 1200kHz 100 N/m	10	#N/B
AGF724-100	High freq. cantilever, Al reflective coating 1200kHz 100 N/m	100	#N/B
AGF724-24	High freq. cantilever, Al reflective coating 1200kHz 100 N/m	24	#N/B
AGF724-5	High freq. cantilever, Al reflective coating 1200kHz 100 N/m	5	#N/B
AGF724-50	High freq. cantilever, Al reflective coating 1200kHz 100 N/m	50	#N/B
AGF725-10	High freq. cantilever, Au reflective coating 1200kHz 100 N/m	10	#N/B
AGF725-100	High freq. cantilever, Au reflective coating 1200kHz 100 N/m	100	#N/B
AGF725-24	High freq. cantilever, Au reflective coating 1200kHz 100 N/m	24	#N/B
AGF725-5	High freq. cantilever, Au reflective coating 1200kHz 100 N/m	5	#N/B
AGF725-50	High freq. cantilever, Au reflective coating 1200kHz 100 N/m	50	#N/B
AGF726-10	High freq. cantilever, no coating 1200kHz 100N/m	10	#N/B
AGF726-100	High freq. cantilever, no coating 1200kHz 100N/m	100	#N/B
AGF726-24	High freq. cantilever, no coating 1200kHz 100N/m	24	#N/B
AGF726-5	High freq. cantilever, no coating 1200kHz 100N/m	5	#N/B
AGF726-50	High freq. cantilever, no coating 1200kHz 100N/m	50	#N/B
AGFIBMS1	Micro-Sections prepared on your sample, with FIB SE images	each	#N/B
AGFIBMS2	Custom options - fast turnaround	each	#N/B
AGFIBMS3	Custom options - sacrificial metal top coat	each	#N/B
AGFIBMS4	Custom options - wider or deeper sections	each	#N/B

Part #	Description	Unit	Price €
AGFIBMS5	Custom options - additional slices	each	#N/B
AGFIBMS6	Custom options - very specific location	each	#N/B
AGFIBMS7	Custom options - in-situ Chemical decoration	each	#N/B
AGFIBMS8	Custom options - Report and description	each	#N/B
AGFIBMS9	Custom options - additional images	each	#N/B
AGFIBXTEM1	FIB prepared ultra-thin TEM sectioning - standard section	each	#N/B
AGFIBXTEM2	Custom options - fast turnaround	each	458,00
AGFIBXTEM3	Custom options - very specific	each	95,60
AGFIBXTEM4	Custom options - thinner sub r	each	95,60
AGFIBXTEM6	Custom options - amorphous skin reduction step - low kV clean	each	95,60
AGFIBXTEM7	Custom options - wedge shaped	each	95,60
AGFIBXTEM8	Custom options - Report with p	each	170,20
AGG1001	Metal Cryo Grid Box	each	54,80
AGG1092450	Microtome stub	100	156,90
AGG1092450-10	Microtome stub	10	19,40
AGG1092450-20	Microtome stub	20	37,30
AGG1092450-50	Microtome stub	50	85,40
AGG1201	Auxiliary Display PAD7 for Pioneer balances	each	91,80
AGG1202	Density Kit for Solids (Pioneer balances)	each	544,40
AGG1203	Sinker Glass for Density Determination (Pioneer Balances)	each	354,00
AGG1204	Printer SF40A for Pioneer balances	each	513,70
AGG1205	SF40A Ink Ribbon Cassette for Pioneer balances	each	9,90
AGG1206	SF40A Paper Roll, 57.5 mm (2 Rolls) for Pioneer balances	each	13,10
AGG1207	Security Device for Pioneer balances	each	180,60
AGG1208	In-Use Cover for Pioneer balances	each	25,50
AGG1209	Dust Cover (for 0,1mg and 1mg Pioneer models)	each	39,40
AGG1210	Stand-alone Ionizer ION-100A for Pioneer balances	each	1.159,20
AGG1211	RS232 Cable (9-pin) for Pioneer balances	each	98,30
AGG1212	Interface Cable, USB (Type A to B) for Pioneer balances	each	21,40
AGG2002	Athene 200 Mesh Thin Bar Grids Copper 3.05mm	100	42,30
AGG2002G	Athene 200 Mesh Thin Bar Grids Au+Cu 3.05mm	100	67,10
AGG2002N	Athene 200 Mesh Thin Bar Grids Nickel 3.05mm	100	44,60
AGG2002P	Athene 200 Mesh Thin Bar Grids Pt+Cu 3.05mm	100	67,10
AGG2003	Athene 300 Mesh Thin Bar Grids Copper 3.05mm	100	46,20
AGG2003G	Athene 300 Mesh Thin Bar Grids Au+Cu 3.05mm	100	71,10
AGG2003N	Athene 300 Mesh Thin Bar Grids Nickel 3.05mm	100	30,00
AGG2003P	Athene 300 Mesh Thin Bar Grids Pt+Cu 3.05mm	100	71,10
AGG2004	Athene 400 Mesh Thin Bar Grids Copper 3.05mm	100	49,30
AGG2004G	Athene 400 Mesh Thin Bar Grids Au+Cu 3.05mm	100	74,60
AGG2004N	Athene 400 Mesh Thin Bar Grids Nickel 3.05mm	100	52,00
AGG2004P	Athene 400 Mesh Thin Bar Grids Pt+Cu 3.05mm	100	74,60
AGG201	Athene Old 150 Grids Copper 3.05mm	100	30,90
AGG2010A	Agar Grids 50 Parallel Bar Gold 3.05mm	50	63,00
AGG2010C	Agar Grids 50 Parallel Bar Copper 3.05mm	100	25,40
AGG2010N	Agar Grids 50 Parallel Bar Nickel 3.05mm	100	35,80
AGG2010PD	Agar Grids 50 Parallel Bar Cu/Pd 3.05mm	100	35,70
AGG2011A	Agar Grids 75 Parallel Bar Gold 3.05mm	50	63,00
AGG2011C	Agar Grids 75 Parallel Bar Copper 3.05mm	100	25,40
AGG2011N	Agar Grids 75 Parallel Bar Nickel 3.05mm	100	35,80
AGG2011PD	Agar Grids 75 Parallel Bar Cu/Pd 3.05mm	100	35,70
AGG2012A	Agar Grids 100 Parallel Bar Gold 3.05mm	50	63,00
AGG2012C	Agar Grids 100 Parallel Bar Copper 3.05mm	100	25,40
AGG2012N	Agar Grids 100 Parallel Bar Nickel 3.05mm	100	35,80
AGG2012PD	Agar Grids 100 Parallel Bar Cu/Pd 3.05mm	100	35,70
AGG2013A	Agar Grids 150 Parallel Bar Gold 3.05mm	50	65,60
AGG2013C	Agar Grids 150 Parallel Bar Copper 3.05mm	100	27,90
AGG2013N	Agar Grids 150 Parallel Bar Nickel 3.05mm	100	38,80
AGG2013PD	Agar Grids 150 Parallel Bar Cu/Pd 3.05mm	100	38,70

Part #	Description	Unit	Price €
AGG2014A	Agar Grids 200 Parallel Bar Gold 3.05mm	50	65,60
AGG2014C	Agar Grids 200 Parallel Bar Copper 3.05mm	100	27,90
AGG2014N	Agar Grids 200 Parallel Bar Nickel 3.05mm (Pk 100)	100	38,80
AGG2014PD	Agar Grids 200 Parallel Bar Cu/Pd 3.05mm	100	38,70
AGG2015A	Agar Grids 300 Parallel Bar Gold 3.05mm	50	65,60
AGG2015C	Agar Grids 300 Parallel Bar Copper 3.05mm	100	28,20
AGG2015N	Agar Grids 300 Parallel Bar Nickel 3.05mm	100	38,80
AGG2015PD	Agar Grids 300 Parallel Bar Cu/Pd 3.05mm	100	38,70
AGG2016A	Agar Grids 400 Parallel Bar Gold 3.05mm	50	69,50
AGG2016C	Agar Grids 400 Parallel Bar Copper 3.05mm	100	30,00
AGG2016N	Agar Grids 400 Parallel Bar Nickel 3.05mm	100	39,10
AGG2016PD	Agar Grids 400 Parallel Bar Cu/Pd 3.05mm	100	40,10
AGG2017A	Agar Grids 50 Bar + Single Gold 3.05mm	50	63,00
AGG2017C	Agar Grids 50 Bar + Single Copper 3.05mm	100	25,40
AGG2017N	Agar Grids 50 Bar + Single Nickel 3.05mm	100	35,80
AGG2017PD	Agar Grids 50 Bar + Single Cu/Pd 3.05mm	100	35,70
AGG2018A	Agar Grids 75 Bar + Single Gold 3.05mm	50	63,00
AGG2018C	Agar Grids 75 Bar + Single Copper 3.05mm	100	25,40
AGG2018N	Agar Grids 75 Bar + Single Nickel 3.05mm	100	35,80
AGG2018PD	Agar Grids 75 Bar + Single Cu/Pd 3.05mm	100	49,00
AGG2019A	Agar Grids 100 Bar + Single Gold 3.05mm	50	86,40
AGG2019C	Agar Grids 100 Bar + Single Copper 3.05mm	100	25,40
AGG2019N	Agar Grids 100 Bar + Single Nickel 3.05mm	100	35,80
AGG2019PD	Agar Grids 100 Bar + Single Cu/Pd 3.05mm	100	49,00
AGG201G	Athene Old 150 Grids Au on Cu 3.05mm	100	56,10
AGG201N	Athene Old 150 Grids Nickel 3.05mm	100	22,60
AGG201P	Athene Old 150 Grids Pt on Cu 3.05mm	100	56,10
AGG202	Athene Old 200 Grids Copper 3.05mm	100	30,90
AGG2020A	Agar Grids 150 Bar + Single Gold 3.05mm	50	92,40
AGG2020C	Agar Grids 150 Bar + Single Copper 3.05mm	100	27,90
AGG2020N	Agar Grids 150 Bar + Single Nickel 3.05mm	100	38,80
AGG2020PD	Agar Grids 150 Bar + Single Cu/Pd 3.05mm	100	53,10
AGG2021A	Agar Grids 200 Bar + Single Gold 3.05mm	50	65,60
AGG2021C	Agar Grids 200 Bar + Single Copper 3.05mm	100	27,90
AGG2021N	Agar Grids 200 Bar + Single Nickel 3.05mm	100	38,80
AGG2021PD	Agar Grids 200 Bar + Single Cu/Pd 3.05mm	100	53,10
AGG2022A	Agar Grids 300 Bar + Single Gold 3.05mm	50	65,60
AGG2022C	Agar Grids 300 Bar + Single Copper 3.05mm	100	27,90
AGG2022N	Agar Grids 300 Bar + Single Nickel 3.05mm	100	38,80
AGG2022PD	Agar Grids 300 Bar + Single Cu/Pd 3.05mm	100	53,10
AGG2023A	Agar Grids 400 Bar + Single Gold 3.05mm	50	91,50
AGG2023C	Agar Grids 400 Bar + Single Copper 3.05mm	100	30,00
AGG2023N	Agar Grids 400 Bar + Single Nickel 3.05mm	100	39,10
AGG2023PD	Agar Grids 400 Bar + Single Cu/Pd 3.05mm	100	54,90
AGG202G	Athene Old 200 Grids Au on Cu 3.05mm	100	56,10
AGG202N	Athene Old 200 Grids Nickel 3.05mm	100	33,60
AGG202P	Athene Old 200 Grids Pt on Cu 3.05mm	100	56,10
AGG203	Athene Old 300 Grids Copper 3.05mm	100	30,90
AGG203G	Athene Old 300 Grids Au on Cu 3.05mm	100	56,10
AGG203N	Athene Old 300 Grids Nickel 3.05mm	100	23,30
AGG203P	Athene Old 300 Grids Pt on Cu 3.05mm	100	56,10
AGG204	Athene Old 400 Grids Copper 3.05mm	100	31,60
AGG204G	Athene Old 400 Grids Au on Cu 3.05mm	100	56,50
AGG204N	Athene Old 400 Grids Nickel 3.05mm	100	34,00
AGG204P	Athene Old 400 Grids Pt on Cu 3.05mm	100	56,50
AGG205	Athene New 100 Grids Copper 3.05mm	100	20,80
AGG2050A	Agar Grids 50 Mesh Gold 3.05mm	50	63,00
AGG2050C	Agar Grids 50 Mesh Copper 3.05mm	100	27,20

Part #	Description	Unit	Price €
AGG2050N	Agar Grids 50 Mesh Nickel 3.05mm	100	35,80
AGG2050PD	Agar Grids 50 Mesh Copper/Palladium 3.05mm	100	35,70
AGG205G	Athene New 100 Grids Au on Cu 3.05mm	100	56,10
AGG205N	Athene New 100 Grids Nickel 3.05mm	100	22,60
AGG205P	Athene New 100 Grids Pt on Cu 3.05mm	100	104,60
AGG206	Athene New 200 Grids Copper 3.05mm	100	30,90
AGG206G	Athene New 200 Grids Au on Cu 3.05mm	100	56,10
AGG206N	Athene New 200 Grids Nickel 3.05mm	100	22,60
AGG206P	Athene New 200 Grids Pt on Cu 3.05mm	100	56,10
AGG207	Athene New 300 Grids Copper 3.05mm	100	20,80
AGG2075A	Agar Grids 75 Mesh Gold 3.05mm	50	63,00
AGG2075C	Agar Grids 75 Mesh Copper 3.05mm	100	25,80
AGG2075N	Agar Grids 75 Mesh Nickel 3.05mm	100	35,80
AGG2075PD	Agar Grids 75 Mesh Copper/Palladium 3.05mm	100	35,70
AGG207G	Athene New 300 Grids Au on Cu 3.05mm	100	56,10
AGG207N	Athene New 300 Grids Nickel 3.05mm	100	33,60
AGG207P	Athene New 300 Grids Pt on Cu 3.05mm	100	58,70
AGG209	Athene M50 Grids Copper 3.05mm	100	34,20
AGG209G	Athene M50 Grids Au on Cu 3.05mm	100	59,20
AGG209N	Athene M50 Grids Nickel 3.05mm	100	36,60
AGG209P	Athene M50 Grids Pt on Cu 3.05mm	100	59,20
AGG210	Athene M75 Grids Copper 3.05mm	100	34,20
AGG2100A	Agar Grids 100 Mesh Gold 3.05mm	50	63,00
AGG2100C	Agar Grids 100 Mesh Copper 3.05mm	100	25,80
AGG2100N	Agar Grids 100 Mesh Nickel 3.05mm	100	35,80
AGG2100PD	Agar Grids 100 Mesh Copper/Palladium 3.05mm	100	35,70
AGG210G	Athene M75 Grids Au on Cu 3.05mm	100	59,20
AGG210N	Athene M75 Grids Nickel 3.05mm	100	36,60
AGG210P	Athene M75 Grids Pt on Cu 3.05mm	100	59,20
AGG211	Athene D/50 Grids Copper 3.05mm	100	31,10
AGG211G	Athene D/50 Grids Au on Cu 3.05mm	100	71,10
AGG211N	Athene D/50 Grids Nickel 3.05mm	100	48,60
AGG211P	Athene D/50 Grids Pt on Cu 3.05mm	100	71,10
AGG212	Athene D/100 Grids Copper 3.05mm	100	31,10
AGG212G	Athene D/100 Grids Au on Cu 3.05mm	100	71,10
AGG212N	Athene D/100 Grids Nickel 3.05mm	100	48,60
AGG212P	Athene D/100 Grids Pt on Cu 3.05mm	100	71,10
AGG213	Athene D/200 Grids Copper 3.05mm	100	46,20
AGG2130	Athene D1000 Folding 1mm Hole Grids + Latch Cu	100	76,80
AGG2130G	Athene D1000 Folding 1mm Hole Grids + Latch Cu/Au	100	95,40
AGG2130N	Athene D1000 Folding 1mm Hole Grids + Latch Ni	100	79,10
AGG2130P	Athene D1000 Folding 1mm Hole Grids + Latch Cu/Pt	100	102,50
AGG213G	Athene D/200 Grids Au on Cu 3.05mm	100	71,10
AGG213N	Athene D/200 Grids Nickel 3.05mm	100	48,60
AGG213P	Athene D/200 Grids Pt on Cu 3.05mm	100	71,10
AGG214	Athene Hexagon 100 Mesh Copper 3.05mm	100	30,10
AGG2140A	400 x 100 Mesh Grids Gold 3.05mm	50	65,60
AGG2140C	400 x 100 Mesh Grids Copper 3.05mm	100	27,90
AGG2140N	400 x 100 Mesh Grids Nickel 3.05mm	100	61,60
AGG2140PD	400 x 100 Mesh Grids Cu/Pd 3.05mm	100	53,10
AGG2145	Athene Hexagon 150 Mesh Copper 3.05mm	100	30,10
AGG2145G	Athene Hexagon 150 Mesh Au on Cu 3.05mm	100	57,60
AGG2145N	Athene Hexagon 150 Mesh Nickel 3.05mm	100	32,60
AGG2145P	Athene Hexagon 150 Mesh Pt on Cu 3.05mm	100	55,20
AGG214G	Athene Hexagon 100 Mesh Au on Cu 3.05mm	100	57,60
AGG214N	Athene Hexagon 100 Mesh Nickel 3.05mm	100	32,60
AGG214P	Athene Hexagon 100 Mesh Pt on Cu 3.05mm	100	55,20
AGG215	Athene Aei Grids Copper 3.05mm	100	34,80

Part #	Description	Unit	Price €
AGG2150A	Agar Grids 150 Mesh Gold 3.05mm	50	65,60
AGG2150C	Agar Grids 150 Mesh Copper 3.05mm	100	27,90
AGG2150N	Agar Grids 150 Mesh Nickel 3.05mm	100	38,80
AGG2150PD	Agar Grids 150 Mesh Copper/Palladium 3.05mm	100	38,70
AGG215G	Athene AEI Grids Au on Cu 3.05mm	100	59,90
AGG215N	Athene AEI Grids Nickel 3.05mm	100	37,50
AGG215P	Athene AEI Grids Platinised 3.05mm	100	59,90
AGG2175C	Agar Grids 175 Mesh Copper 3.05mm	100	27,90
AGG2175N	Agar Grids 175 Mesh Nickel 3.05mm	100	38,80
AGG2175PD	Agar Grids 175 Mesh Copper/Palladium 3.05mm	100	53,10
AGG218	Athene Robertson Grids Copper 3.05mm	100	34,80
AGG218G	Athene Robertson Grids Au on Cu 3.05mm	100	59,90
AGG218N	Athene Robertson Grids Nickel 3.05mm	100	19,90
AGG218P	Athene Robertson Grids Pt on Cu 3.05mm	100	59,90
AGG219	Athene 483 Grids Copper 3.05mm	100	37,70
AGG219G	Athene 483 Grids Au on Cu 3.05mm	100	62,80
AGG219N	Athene 483 Grids Nickel 3.05mm	100	40,30
AGG219P	Athene 483 Grids Pt on Cu 3.05mm	100	62,80
AGG2200A	Agar Grids 200 Mesh Gold 3.05mm	50	65,60
AGG2200C	Agar Grids 200 Mesh Copper 3.05mm	100	28,10
AGG2200N	Agar Grids 200 Mesh Nickel 3.05mm	100	39,00
AGG2200PD	Agar Grids 200 Mesh Copper/Palladium 3.05mm	100	41,90
AGG220-1	Athene 2mm x 12.5µm Slots Copper 3.05mm	100	30,70
AGG220-3	Athene 2mm x 25µm Slots Copper 3.05mm	100	30,70
AGG220-4	Athene 2mm x 60µm Slots Copper 3.05mm (Pk 100)	100	30,70
AGG220-5	Athene 2mm x 125µm Slots Copper 3.05mm	100	30,70
AGG220-6	Athene 2mm x 500um Slots Copper 3.05mm	100	30,70
AGG220-7	Athene 2mm x 1000µm Slots Copper 3.05mm	100	22,50
AGG220-8	Athene 2mm x 1000µm Slot + Bar Copper 3.05mm	100	43,80
AGG220G1	Athene 2mm x 12.5µm Slots Au on Cu 3.05mm	100	56,50
AGG220G3	Athene 2mm x 25µm Slots Au on Cu 3.05mm	100	56,50
AGG220G4	Athene 2mm x 60µm Slots Au on Cu 3.05mm	100	56,50
AGG220G5	Athene 2mm x 125µm Slots Au on Cu 3.05mm	100	56,50
AGG220G6	Athene 2mm x 500um Slots Au on Cu 3.05mm	100	56,50
AGG220G7	Athene 2mm x 1000µm Slots Au on Cu 3.05mm	100	56,50
AGG220G8	Athene 2mm x 1000µm Slot + Bar Cu/Au 3.05mm	100	67,80
AGG220N1	Athene 2mm x 12.5µm Slots Nickel 3.05mm	100	22,50
AGG220N3	Athene 2mm x 25µm Slots Nickel 3.05mm	100	33,40
AGG220N4	Athene 2mm x 60µm Slots Nickel 3.05mm	100	33,40
AGG220N5	Athene 2mm x 125µm Slots Nickel 3.05mm	100	33,40
AGG220N6	Athene 2mm x 500um Slots Nickel 3.05mm	100	22,50
AGG220N7	Athene 2mm x 1000µm Slots Nickel 3.05mm	100	34,00
AGG220N8	Athene 2mm x 1000µm Slot + Bar Nickel 3.05mm	100	31,30
AGG220P1	Athene 2mm x 12.5µm Slots Pt on Cu 3.05mm	100	55,80
AGG220P3	Athene 2mm x 25µm Slots Pt on Cu 3.05mm	100	55,80
AGG220P4	Athene 2mm x 60µm Slots Pt on Cu 3.05mm	100	55,80
AGG220P5	Athene 2mm x 125µm Slots Pt on Cu 3.05mm	100	55,80
AGG220P6	Athene 2mm x 500um Slots Pt on Cu 3.05mm	100	55,80
AGG220P7	Athene 2mm x 1000µm Slots Pt on Cu3.05	100	55,80
AGG220P8	Athene 2mm x 1000µm Slot + Bar Cu/Pt 3.05mm	100	68,50
AGG220T6	Athene Thick 2mm x 500um Slots Copper 3.05mm	100	46,20
AGG220T7	Athene Thick 2mm x 1000µm Slot Copper 3.05mm	100	46,20
AGG220TG6	Athene Thick 2mmx500um Slots Au on Cu 3.05mm	100	71,10
AGG220TG7	Athene Thick 2mmx1mm Slots Au on Cu 3.05mm	100	71,10
AGG220TN6	Athene Thick 2mm x 500um Slots Nickel 3.05mm	100	32,70
AGG220TN7	Athene Thick 2mm x 1000µm Slots Nickel 3.05mm	100	48,60
AGG220TP6	Athene Thick 2mmx500um Slots Pt on Cu 3.05mm	100	73,60
AGG220TP7	Athene Thick 2mmx1mm Slots Pt on Cu 3.05mm	100	73,60

Part #	Description	Unit	Price €
AGG221	Athene Sjostrand Grids Copper 3.05mm	100	34,80
AGG221G	Athene Sjostrand Grids Au on Cu 3.05mm	100	59,90
AGG221N	Athene Sjostrand Grids Nickel 3.05mm	100	37,50
AGG221P	Athene Sjostrand Grids Pt on Cu 3.05mm	100	59,90
AGG2220A	Agar Grids 200 Mesh Gold Centre mark 3.05mm	50	69,20
AGG2220C	Agar Grids 200 Mesh Copper Centre mark 3.05mm	100	26,40
AGG2220N	Agar Grids 200 Mesh Nickel Centre mark 3.05mm	100	41,50
AGG2250A	Agar Grids 250 Mesh Gold 3.05mm (Tube of 50)		65,60
AGG2250A	Agar Grids 250 Mesh Gold 3.05mm	50	65,60
AGG2250C	Agar Grids 250 Mesh Copper 3.05mm (Tube of 100)		27,90
AGG2250C	Agar Grids 250 Mesh Copper 3.05mm	100	27,90
AGG2250N	Agar Grids 250 Mesh Nickel 3.05mm (Tube of 100)		38,80
AGG2250N	Agar Grids 250 Mesh Nickel 3.05mm	100	38,80
AGG2250PD	Agar Grids 250 Mesh Copper/Palladium 3.05mm	100	38,70
AGG225-1	Athene 100µm Hole Grid Copper 3.05mm	100	22,50
AGG225-2	Athene 200µm Hole Grid Copper 3.05mm	100	22,50
AGG225-3	Athene 375µm Hole Grid Copper 3.05mm	100	30,70
AGG225-4	Athene 600µm Hole Grid Copper 3.05mm	100	30,70
AGG225-5	Athene 800µm Hole Grid Copper 3.05mm	100	30,70
AGG225-6	Athene 1000µm Hole Grid Copper 3.05mm	100	30,70
AGG225-7	Athene 2000µm Hole Grid Copper 3.05mm	100	30,70
AGG225G1	Athene 100µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G2	Athene 200µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G3	Athene 375µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G4	Athene 600µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G5	Athene 800µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G6	Athene 1000µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225G7	Athene 2000µm Hole Grids Au on Cu 3.05mm	100	55,80
AGG225N1	Athene 100µm Hole Grids Nickel 3.05mm	100	33,40
AGG225N2	Athene 200µm Hole Grids Nickel 3.05mm	100	33,40
AGG225N3	Athene 375µm Hole Grids Nickel 3.05mm	100	33,40
AGG225N4	Athene 600µm Hole Grids Nickel 3.05mm	100	33,40
AGG225N5	Athene 800µm Hole Grids Nickel 3.05mm	100	#N/B
AGG225N6	Athene 1000µm Hole Grids Nickel 3.05mm	100	33,40
AGG225N7	Athene 2000µm Hole Grids Nickel 3.05mm	100	33,40
AGG225P1	Athene 100µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P2	Athene 200µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P3	Athene 375µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P4	Athene 600µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P5	Athene 800µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P6	Athene 1000µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG225P7	Athene 2000µm Hole Grids Pt on Cu 3.05mm	100	55,80
AGG226	Athene R20m Grids Copper 3.05mm	100	37,70
AGG226G	Athene R20m Grids Au on Cu 3.05mm	100	62,80
AGG226N	Athene R20m Grids Nickel 3.05mm	100	40,30
AGG226P	Athene R20m Grids Platinised 3.05mm	100	62,80
AGG227	Athene Polyslot Grids Copper 3.05mm	100	37,70
AGG227G	Athene Polyslot Grids Au on Cu 3.05mm	100	62,80
AGG227N	Athene Polyslot Grids Nickel 3.05mm	100	40,30
AGG227P	Athene Polyslot Grids Pt on Cu 3.05mm	100	62,80
AGG230	Folding Grids 50/50 Mesh Copper 3.05mm	100	60,20
AGG2300A	Agar Grids 300 Mesh Gold 3.05mm	50	68,70
AGG2300C	Agar Grids 300 Mesh Copper 3.05mm	100	28,10
AGG2300N	Agar Grids 300 Mesh Nickel 3.05mm	100	39,00
AGG2300PD	Agar Grids 300 Mesh Copper/Palladium 3.05mm	100	38,70
AGG230A	Folding Grids 50/50 Mesh Gold 3.05mm	50	87,80
AGG230N	Folding Grids 50/50 Mesh Nickel 3.05mm	100	66,50
AGG230PD	Folding Grids 50/50 Mesh Cu/Pd 3.05mm	100	95,50

Part #	Description	Unit	Price €
AGG231	Folding Grids 100/100 Mesh Copper 3.05mm	100	60,20
AGG231A	Folding Grids 100/100 Mesh Gold 3.05mm	50	87,80
AGG231PD	Folding Grids 100/100 Mesh Cu/Pd 3.05mm	100	95,50
AGG233	Finder Grids H7 Copper 3.05mm	100	50,80
AGG233A	Finder Grids H7 Gold 3.05mm	100	294,00
AGG233N	Finder Grids H7 Nickel 3.05mm	100	59,80
AGG234	Folding Grids 100/200 Mesh Copper 3.05mm	100	60,20
AGG234A	Folding Grids 100/200 Mesh Gold 3.05mm	50	87,80
AGG234N	Folding Grids 100/200 Mesh Nickel 3.05mm	100	66,50
AGG234PD	Folding Grids 100/200 Mesh Cu/Pd 3.05mm	100	95,50
AGG236	Folding Grids 50/100 Mesh Copper 3.05mm	100	60,20
AGG236A	Folding Grids 50/100 Mesh Gold 3.05mm	50	113,00
AGG236N	Folding Grids 50/100 Mesh Nickel 3.05mm	100	66,50
AGG236PD	Folding Grids 50/100 Mesh Cu/Pd 3.05mm	100	95,50
AGG2375A	300 x 75 Mesh Grids Gold 3.05mm	50	65,60
AGG2375C	300 x 75 Mesh Grids Copper 3.05mm	100	28,20
AGG2375N	300 x 75 Mesh Grids Nickel 3.05mm	100	38,80
AGG2375PD	300 x 75 Mesh Grids Cu/Pd 3.05mm	100	53,10
AGG2394C	12mm Numbered Finder Grids Copper	20	189,90
AGG2394N	12mm Numbered Finder Grids Nickel	20	198,90
AGG2395C	2x1 Slot Grids for Compustage Cu	100	37,70
AGG2396	Tomography grid	each	62,80
AGG2400A	Agar Grids 400 Mesh Gold 3.05mm	50	69,50
AGG2400C	Agar Grids 400 Mesh Copper 3.05mm	100	30,40
AGG2400N	Agar Grids 400 Mesh Nickel 3.05mm	100	39,30
AGG2400PD	Agar Grids 400 Mesh Copper/Palladium 3.05mm	100	40,10
AGG2405A	Agar Grids Hexagonal 50 Mesh Gold 3.05mm	50	63,00
AGG2405C	Agar Grids Hexagonal 50 Mesh Copper 3.05mm	100	25,40
AGG2405N	Agar Grids Hexagonal 50 Mesh Nickel 3.05mm	100	36,30
AGG2405PD	Agar Grids Hexagonal 50 Mesh Cu/Pd 3.05mm	100	35,70
AGG2410A	Agar Grids Hexagonal 100 Mesh Gold 3.05mm	50	63,00
AGG2410C	Agar Grids Hexagonal 100 Mesh Copper 3.05mm	100	25,40
AGG2410N	Agar Grids Hexagonal 100 Mesh Nickel 3.05mm	100	35,80
AGG2410PD	Agar Grids Hexagonal 100 Mesh Cu/Pd 3.05mm	100	35,70
AGG2415A	Agar Grids Hexagonal 150 Mesh Gold 3.05mm	50	92,40
AGG2415C	Agar Grids Hexagonal 150 Mesh Copper 3.05mm	100	27,90
AGG2415N	Agar Grids Hexagonal 150 Mesh Nickel 3.05mm	100	38,80
AGG2415PD	Agar Grids Hexagonal 150 Mesh Cu/Pd 3.05mm	100	38,70
AGG2430A	Agar Grids Hexagonal 300 Mesh Gold 3.05mm	50	65,60
AGG2430C	Agar Grids Hexagonal 300 Mesh Copper 3.05mm	100	28,20
AGG2430N	Agar Grids Hexagonal 300 Mesh Nickel 3.05mm	100	38,80
AGG2430PD	Agar Grids Hexagonal 300 Mesh Cu/Pd 3.05mm	100	38,70
AGG243A	Fine Wire Mesh Gold 1500 Mesh	each	337,70
AGG243C	Fine Wire Mesh Copper 1500 Mesh	each	195,60
AGG243N	Fine Wire Mesh Nickel 1500 Mesh	each	107,30
AGG2440A	Agar Grids Hexagonal 400 Mesh Gold 3.05mm	50	69,50
AGG2440C	Agar Grids Hexagonal 400 Mesh Copper 3.05mm	100	30,00
AGG2440PD	Agar Grids Hexagonal 400 Mesh Cu/Pd 3.05mm	100	40,10
AGG245	Finder Grids HF15 Copper / Rhodium 3.05mm	100	50,10
AGG2450A	Agar Grids Hexagonal 200 Mesh Gold 3.05mm	50	65,60
AGG2450C	Agar Grids Hexagonal 200 Mesh Copper 3.05mm	100	28,20
AGG2450N	Agar Grids Hexagonal 200 Mesh Nickel 3.05mm	100	38,80
AGG2450PD	Agar Grids Hexagonal 200 Mesh Cu/Pd 3.05mm	100	38,70
AGG245A	Finder Grids HF15 Gold 3.05mm	100	294,00
AGG245N	Finder Grids HF15 Nickel 3.05mm	100	57,50
AGG246	Finder Grids H2 Copper 3.05mm	100	50,10
AGG2460AL	Aluminium Grids 50 Square Mesh 3.05mm	25	63,80
AGG2460MO	Molybdenum Grids 50 Square Mesh 3.05mm	25	76,20

Part #	Description	Unit	Price €
AGG2460SS	Stainless Steel Grids 50 Square Mesh 3.05mm	25	74,40
AGG2461AL	Aluminium Grids 75 Square Mesh 3.05mm	25	63,80
AGG2461MO	Molybdenum Grids 75 Square Mesh 3.05mm	25	76,60
AGG2461SS	Stainless Steel Grids 75 Square Mesh 3.05mm	25	53,70
AGG2461TI	Titanium Grids 75 Square Mesh 3.05mm	25	127,90
AGG2462AL	Aluminium Grids 100 Square Mesh 3.05mm	25	63,80
AGG2462MO	Molybdenum Grids 100 Square Mesh 3.05mm	25	76,60
AGG2462SS	Stainless Steel Grids 100 Square Mesh 3.05mm	25	53,70
AGG2462TI	Titanium Grids 100 Square Mesh 3.05mm	25	43,10
AGG2463AL	Aluminium Grids 150 Square Mesh 3.05mm	25	63,80
AGG2463MO	Molybdenum Grids 150 Square Mesh 3.05mm	25	79,60
AGG2463SS	Stainless Steel Grids 150 Square Mesh 3.05mm	25	53,70
AGG2463TI	Titanium Grids 150 Square Mesh 3.05mm	25	43,10
AGG2464AL	Aluminium Grids 200 Square Mesh 3.05mm	25	66,00
AGG2464MO	Molybdenum Grids 200 Square Mesh 3.05mm	25	76,60
AGG2464SS	Stainless Steel Grids 200 Square Mesh 3.05mm	25	53,70
AGG2464TI	Titanium Grids 200 Square Mesh 3.05mm	25	141,40
AGG2465AL	Aluminium Grids 300 Square Mesh 3.05mm	25	63,80
AGG2465MO	Molybdenum Grids 300 Square Mesh 3.05mm	25	76,60
AGG2465SS	Stainless Steel Grids 300 Square Mesh 3.05mm	25	53,60
AGG2466MO	Folding Grids 100/100 Mesh Molybdenum 3.05mm	25	141,30
AGG2466SS	Folding Grids 100/100 Mesh S.Steel 3.05mm	25	85,20
AGG2466TI	Folding Grids 100/100 Mesh Titanium 3.05mm	25	85,20
AGG2467MO	Folding Grids 100/200 Mesh Molybdenum 3.05mm	25	141,30
AGG2467SS	Folding Grids 100/200 Mesh S.Steel 3.05mm	25	85,20
AGG2467TI	Folding Grids 100/200 Mesh Titanium 3.05mm	25	85,20
AGG2468SS	Folding Grids 100/300 Mesh S Steel 3.05mm	25	85,20
AGG2468TI	Folding Grids 100/300 Mesh Titanium 3.05mm	25	79,10
AGG2469MO	Molybdenum Grids 400µm Hole 3.05mm	25	90,50
AGG2469SS	Stainless Steel Grids 400µm Hole 3.05mm	25	76,60
AGG246A	Finder Grids H2 Gold 3.05mm	100	294,00
AGG246N	Finder Grids H2 Nickel 3.05mm	100	57,50
AGG247	Finder Grids H6 Copper 3.05mm	100	50,10
AGG2470AL	Aluminium Grids 1000µm Hole 3.05mm	25	89,80
AGG2470MO	Molybdenum Grids 1000µm Hole 3.05mm	25	90,50
AGG2470SS	Stainless Steel Grids 1000µm Hole 3.05mm	25	76,60
AGG2470TI	Titanium Grids 1000µm Hole 3.05mm	25	62,60
AGG2471AL	Aluminium Grids 1500µm Hole 3.05mm	25	92,40
AGG2471MO	Molybdenum Grids 1500µm Hole 3.05mm	25	90,50
AGG2471SS	Stainless Steel Grids 1500µm Hole 3.05mm	25	76,60
AGG2471TI	Titanium Grids 1500µm Hole 3.05mm	25	62,60
AGG2472AL	Aluminium Grids 2mm x 1.2mm Slot 3.05mm	25	92,40
AGG2472MO	Molybdenum Grids 2mmx1.2mm Slot 3.05mm	25	96,70
AGG2472SS	Stainless Grids 2 x 1.2mm Slot 3.05mm	25	76,60
AGG2472TI	Titanium Grids 2 x 1.2mm Slot 3.05mm	25	205,60
AGG2473AL	Aluminium Grids 2mm x 0.42mm Slot 3.05mm	25	92,40
AGG2473MO	Molybdenum Grids 2mm x 0.42mm Slot 3.05mm	25	90,50
AGG2473SS	Stainless Grids 2mm x 0.42mm Slot 3.05mm	25	76,60
AGG2473TI	Titanium Grids 2mm x 0.42mm Slot 3.05mm	25	62,60
AGG2474AL	Aluminium Grids 400 Square Mesh 3.05mm	25	#N/B
AGG2475A	Agar Grids Hexagonal 75 Mesh Gold 3.05mm	50	63,00
AGG2475C	Agar Grids Hexagonal 75 Mesh Copper 3.05mm	100	25,40
AGG2475N	Agar Grids Hexagonal 75 Mesh Nickel 3.05mm	100	35,80
AGG2475PD	Agar Grids Hexagonal 75 Mesh Cu/Pd 3.05mm	100	35,70
AGG247A	Finder Grids H6 Gold 3.05mm	100	294,00
AGG247N	Finder Grids H6 Nickel 3.05mm	100	56,10
AGG2480A	Finder Grids Gold 3.05mm	50	97,80
AGG2480C	Finder Grids Copper 3.05mm	100	47,90

Part #	Description	Unit	Price €
AGG2480N	Finder Grids Nickel 3.05mm	100	108,70
AGG2481A	SEM Finder Grids Gold 10mm	5	301,90
AGG2481C	SEM Finder Grids Copper 10mm (Tube of 10)	10	93,70
AGG2481N	SEM Finder Grids Nickel 10mm (Tube of 10)	10	178,90
AGG2482A	SEM F2 Finder Grids Gold 10mm	5	301,90
AGG2482C	SEM F2 Finder Grids Copper 10mm (Tube of 10)	10	151,50
AGG2482N	SEM F2 Finder Grids Nickel 10mm (Tube of 10)	10	110,60
AGG2483A	SEM F3 Finder Grids Gold 10mm	5	301,90
AGG2483C	SEM F3 Finder Grids Copper 10mm (Tube of 10)	10	151,50
AGG2483N	SEM F3 Finder Grids Nickel 10mm (Tube of 10)	10	178,90
AGG248A	Fine Wire Mesh Gold 1000 Mesh	each	323,50
AGG248C	Fine Wire Mesh Copper 1000 Mesh	each	175,60
AGG248N	Fine Wire Mesh Nickel 1000 Mesh	each	87,40
AGG2495A	Agar Grids 2mm x 1.5mm Slot Gold 3.05mm	50	75,40
AGG2495C	Agar Grids 2mm x 1.5mm Slot Copper 3.05mm	100	36,50
AGG2495N	Agar Grids 2mm x 1.5mm Slot Nickel 3.05mm	100	51,50
AGG2495PD	Agar Grids 2mm x 1.5mm Slot Cu/Pd 3.05mm	100	66,60
AGG249N	Fine Wire Mesh Nickel 2000 Mesh	each	107,30
AGG2500A	Agar Grids 2x1mm Slot Gold 3.05mm	50	75,40
AGG2500C	Agar Grids 2x1mm Slot Copper 3.05mm	100	36,50
AGG2500N	Agar Grids 2x1mm Slot Nickel 3.05mm	100	51,50
AGG2500PD	Agar Grids 2x1mm Slot Cu/Pd 3.05mm	100	48,60
AGG250-1	Mica Sheets 75mm x 25mm x 0.15mm thick	20	93,60
AGG250-2	Mica Sheets 25mm x 25mm x 0.15mm thick	20	57,90
AGG250-3	Mica Sheets 11mm x 11mm x 0.15mm thick	20	85,50
AGG250-4	Mica Sheets 100mm x 20mm, 0.15mm thick	20	119,50
AGG250-5	Mica Sheets 150mm x 150mm x 0.15mm thick	3	172,60
AGG250-7	Mica Sheets Ultraclean 20mm x 20mm x 0.25mm thick	10	90,20
AGG251P	Plain Slides, 76 x 26mm, 1.2-1.5mm	50	15,90
AGG252	Stainless Steel Mesh 150mm x 150mm	3	45,00
AGG2525A	Agar Grids 2mm x 0.75mm Slot Gold 3.05mm	50	75,40
AGG2525C	Agar Grids 2mm x 0.75mm Slot Copper 3.05mm	100	36,50
AGG2525N	Agar Grids 2mm x 0.75mm Slot Nickel 3.05mm	100	51,50
AGG2525PD	Agar Grids 2mm x 0.75mm Slot Cu/Pd 3.05mm	100	48,60
AGG253	Copper Tape - 12.5mm	16m	96,50
AGG253A	Copper Tape - 6.4mm	16m	62,40
AGG254B	Cellulose Acetate 180um, 150mm x 100mm	20	40,00
AGG255	Cellulose Acetate 35µm, 150mm x 100mm	20	24,50
AGG2550A	Agar Grids 2x0.5mm Slot Gold 3.05mm	50	75,40
AGG2550C	Agar Grids 2x0.5mm Slot Copper 3.05mm	100	38,90
AGG2550N	Agar Grids 2x0.5mm Slot Nickel 3.05mm	100	51,50
AGG2550PD	Agar Grids 2 x 0.5mm Slot Cu/Pd 3.05mm	100	48,60
AGG2553A	Agar Grids 1.5mm x 0.3mm Slot Gold 3.05mm	50	75,40
AGG2553C	Agar Grids 1.5mm x 0.3mm Slot Copper 3.05mm	100	38,90
AGG2553N	Agar Grids 1.5mm x 0.3mm Slot Nickel 3.05mm	100	81,80
AGG2553PD	Agar Grids 1.5mm x 0.3mm Slot Cu/Pd 3.05mm	100	40,20
AGG255A	Cellulose Acetate 50µm, 150mm x 100mm	20	27,20
AGG255B	Cellulose Acetate 75µm, 150mm x 100mm	20	24,50
AGG256	White Tile, 150 x 150mm	each	6,10
AGG2560A	Agar Grids 1x0.2mm Slot Gold 3.05mm	50	75,40
AGG2560C	Agar Grids 1x0.2mm Slot Copper 3.05mm	100	38,90
AGG2560N	Agar Grids 1x0.2mm Slot Nickel 3.05mm	100	51,50
AGG2560PD	Agar Grids 1x0.2mm Slot Cu/Pd 3.05mm	100	48,60
AGG2564A	Triple slot grid, Gold 3.05mm	50	75,40
AGG2564C	Triple slot grid, Copper 3.05mm	100	37,80
AGG2564N	Triple slot grid, Nickel 3.05mm	100	52,80
AGG2564PD	Triple slot grid, Cu/Pd 3.05mm	100	50,50
AGG257	Pasteur Pipettes Short 5.75inch (Pack of 250)	250	26,80

Part #	Description	Unit	Price €
AGG258	Rubber Bulbs	10	32,40
AGG259	Pasteur Pipettes Long 9inch	250	35,20
AGG260	PVC tape, 9mm	33m	13,30
AGG2600A	Agar Grids 1000µm Hole Gold 3.05mm	50	75,40
AGG2600C	Agar Grids 1000µm Hole Copper 3.05mm	100	36,50
AGG2600N	Agar Grids 1000µm Hole Nickel 3.05mm	100	51,50
AGG2600PD	Agar Grids 1000µm Hole Cu/Pd 3.05mm	100	66,60
AGG2605A	Agar Grids 1500µm Hole Gold 3.05mm	50	75,40
AGG2605C	Agar Grids 1500µm Hole Copper 3.05mm	100	36,50
AGG2605N	Agar Grids 1500µm Hole Nickel 3.05mm	100	51,50
AGG2605PD	Agar Grids 1500µm Hole Cu/Pd 3.05mm	100	66,60
AGG261	BEEM Dial-a-Grid Storage box	each	9,40
AGG2610A	Agar Grids 100µm Hole Gold 3.05mm	50	75,40
AGG2610C	Agar Grids 100µm Hole Copper 3.05mm	100	41,00
AGG2610N	Agar Grids 100µm Hole Nickel 3.05mm	100	55,80
AGG2610PD	Agar Grids 100µm Hole Cu/Pd 3.05mm	100	76,00
AGG2615A	Agar Grids 150µm Hole Gold 3.05mm	50	75,40
AGG2615C	Agar Grids 150µm Hole Copper 3.05mm	100	41,00
AGG2615N	Agar Grids 150µm Hole Nickel 3.05mm	100	55,80
AGG2615PD	Agar Grids 150µm Hole Cu/Pd 3.05mm	100	76,00
AGG2618A	Agar Grids 200µm Hole Gold 3.05mm	50	97,00
AGG2618C	Agar Grids 200µm Hole Copper 3.05mm	100	36,50
AGG2618N	Agar Grids 200µm Hole Nickel 3.05mm	100	81,80
AGG2618PD	Agar Grids 200µm Hole Cu/Pd 3.05mm	100	66,60
AGG262	Silver tape, 9mm	66m	69,70
AGG2620A	Agar Grids 2000µm Hole Gold 3.05mm	50	75,40
AGG2620C	Agar Grids 2000µm Hole Copper 3.05mm	100	36,50
AGG2620PD	Agar Grids 2000µm Hole Cu/Pd 3.05mm	100	48,60
AGG263	Double Sided Tape	55m	27,80
AGG2630A	Agar Grids 300µm Hole Gold 3.05mm	50	97,00
AGG2630C	Agar Grids 300µm Hole Copper 3.05mm	100	36,50
AGG2630N	Agar Grids 300µm Hole Nickel 3.05mm	100	51,50
AGG2630PD	Agar Grids 300µm Hole Cu/Pd 3.05mm	100	66,60
AGG264	Double Sided Tape - 38mm	55m	55,50
AGG2640A	Agar Grids 400µm Hole Gold 3.05mm	50	75,40
AGG2640C	Agar Grids 400µm Hole Copper 3.05mm	100	60,20
AGG2640N	Agar Grids 400µm Hole Nickel 3.05mm	100	51,50
AGG2640PD	Agar Grids 400µm Hole Cu/Pd 3.05mm	100	66,60
AGG2645A	Agar Grids 500um Hole Gold 3.05mm	50	97,00
AGG2645C	Agar Grids 500um Hole Copper 3.05mm	100	36,50
AGG2645N	Agar Grids 500um Hole Nickel 3.05mm	100	81,80
AGG2645PD	Agar Grids 500um Hole Cu/Pd 3.05mm	100	66,60
AGG265	Aluminium Tape with Adhesive Backing - 12mm	55m	35,60
AGG2650A	Agar Grids 600 Mesh Gold 3.05mm	50	72,00
AGG2650C	Agar Grids 600 Mesh Copper 3.05mm	100	47,70
AGG2650N	Agar Grids 600 Mesh Nickel 3.05mm	100	59,30
AGG2650PD	Agar Grids 600 Mesh Copper/Palladium 3.05mm	100	92,70
AGG2655A	Agar Grids 600 Mesh Thin Bar Gold 3.05mm	50	91,90
AGG2655C	Agar Grids 600 Mesh Thin Bar Copper 3.05mm	100	81,70
AGG2655N	Agar Grids 600 Mesh Thin Bar Nickel 3.05mm	100	106,60
AGG265A	Aluminium Tape with Adhesive Backing - 25mm	55m	59,30
AGG265C	Aluminium Tape with Adhesive Backing - 50mm	55m	93,30
AGG266	Plastic Pipette 1ml	25	17,30
AGG2660A	Agar Grids 600µm Hole Gold 3.05mm	50	97,00
AGG2660C	Agar Grids 600µm Hole Copper 3.05mm	100	36,50
AGG2660N	Agar Grids 600µm Hole Nickel 3.05mm	100	51,50
AGG2660PD	Agar Grids 600µm Hole Cu/Pd 3.05mm	100	48,60
AGG267	Plastic Pipette 2ml	25	6,30

Part #	Description	Unit	Price €
AGG2670A	Agar Grids 600 Mesh Hexagon Gold 3.05mm	50	91,90
AGG2670C	Agar Grids 600 Mesh Hexagon Copper 3.05mm	100	81,70
AGG2670N	Agar Grids 600 Mesh Hexagon Nickel 3.05mm	100	106,60
AGG2672A	Agar Grids 50µm Hole Gold 3.05mm	50	97,00
AGG2672C	Agar Grids 50µm Hole Copper 3.05mm	100	41,00
AGG2672N	Agar Grids 50µm Hole Nickel 3.05mm	100	55,80
AGG2672PD	Agar Grids 50µm Hole Cu/Pd 3.05mm	100	76,00
AGG2675A	Agar Grids 75µm Hole Gold 3.05mm	50	75,40
AGG2675C	Agar Grids 75µm Hole Copper 3.05mm	100	41,00
AGG2675N	Agar Grids 75µm Hole Nickel 3.05mm	100	55,80
AGG2675PD	Agar Grids 75µm Hole Cu/Pd 3.05mm	100	76,00
AGG268	Plastic Pipette 5ml	25	11,00
AGG2680A	Agar Grids 800um Hole Gold 3.05mm	50	97,00
AGG2680C	Agar Grids 800um Hole Copper 3.05mm	100	36,50
AGG2680N	Agar Grids 800um Hole Nickel 3.05mm	100	51,50
AGG2680PD	Agar Grids 800um Hole Cu/Pd 3.05mm	100	66,60
AGG269	Plastic Pipette 10ml	25	14,20
AGG270	Plastic Pipette 20ml	25	20,60
AGG2700A	High Transmission 200 Mesh Grids Gold 3.05mm	50	64,00
AGG2700C	High Transmission 200 Mesh Grids Cu 3.05mm	100	49,30
AGG271	Plastic Pipette 60ml	10	33,40
AGG2710A	Hex 200 Mesh Thin Bar Grids Gold 3.05mm	50	64,00
AGG2710C	Hex 200 Mesh Thin Bar Grids Copper 3.05mm	100	49,30
AGG2710N	Hex 200 Mesh Thin Bar Grids Nickel 3.05mm	100	62,30
AGG272	Coarse Needles 21g x 1.5"	10	6,10
AGG2720A	High Transmission 300 Mesh Grids Gold 3.05mm	50	66,10
AGG2720C	High Transmission 300 Mesh Grids Cu 3.05mm	100	58,70
AGG2720N	High Transmission 300 Mesh Grids Ni 3.05mm	100	73,40
AGG273	Fine Needles 23g x 1.25"	10	5,80
AGG2730A	High Transmission 400 Mesh Grids Au 3.05mm	50	81,10
AGG2730C	High Transmission 400 Mesh Grids Cu 3.05mm	100	74,40
AGG2730N	High Transmission 400 Mesh Grids Ni 3.05mm	100	93,90
AGG2740A	Hex 300 Mesh Thin Bar Grids Gold 3.05mm	50	98,90
AGG2740C	Hex 300 Mesh Thin Bar Grids Copper 3.05mm	100	57,00
AGG2740N	Hex 300 Mesh Thin Bar Grids Nickel 3.05mm	100	73,40
AGG2750A	Hex 460 Mesh Thin Bar Grids Gold 3.05mm	50	81,10
AGG2750C	Hex 460 Mesh Thin Bar Grids Copper 3.05mm	100	74,40
AGG2750N	Hex 460 Mesh Thin Bar Grids Nickel 3.05mm	100	93,90
AGG276	Grid storage box for 100 grids	each	20,10
AGG2760A	Hex 700 Mesh Thin Bar Grids Gold 3.05mm	50	91,90
AGG2760C	Hex 700 Mesh Thin Bar Grids Copper 3.05mm	100	81,70
AGG2761A	F1 Finder Grids Gold 3.05mm	50	76,00
AGG2761C	F1 Finder Grids Copper 3.05mm	100	49,70
AGG2761N	F1 Finder Grids Nickel 3.05mm	100	73,00
AGG2761PD	F1 Finder Grids Copper/Palladium 3.05mm	100	97,00
AGG2762A	F2 Finder Grids Gold 3.05mm	50	76,00
AGG2762C	F2 Finder Grids Copper 3.05mm	100	47,90
AGG2762N	F2 Finder Grids Nickel 3.05mm	100	73,00
AGG2762PD	F2 Finder Grids Copper/Palladium 3.05mm	100	97,00
AGG2763A	F3 Finder Grids Gold 3.05mm	25	81,50
AGG2763C	F3 Finder Grids Copper 3.05mm	25	25,70
AGG2763N	F3 Finder Grids Nickel 3.05mm	25	40,70
AGG2764A	F4 Finder Grids Gold 3.05mm	50	97,80
AGG2764C	F4 Finder Grids Copper 3.05mm	100	47,90
AGG2764N	F4 Finder Grids Nickel 3.05mm	100	73,00
AGG2764PD	F4 Finder Grids Copper/Palladium 3.05mm	100	97,00
AGG276A	Grid storage box for 50 grids	each	13,00
AGG276N	Grid storage box with identification number	each	22,00

Part #	Description	Unit	Price €
AGG276P	Grid storage box	each	47,60
AGG277-0	Storage Box Hinged Lid 1 Compartment	each	22,80
AGG2770A	High Transmission Mixed Mesh Grids Au 3.05mm	50	66,10
AGG2770C	High Transmission Mixed Mesh Grids Cu 3.05mm	100	57,00
AGG2770N	High Transmission Mixed Mesh Grids Ni 3.05mm	100	67,90
AGG277-12	Storage Box Hinged Lid 12 Compartments	each	23,60
AGG277-18	Storage Box Hinged Lid 18 Compartments	each	23,50
AGG277-3	Storage Box Hinged Lid 3 Compartments	each	22,80
AGG277-36	Storage Box Hinged Lid 36 Compartments	each	23,50
AGG277-6	Storage Box Hinged Lid 6 Compartments	each	23,50
AGG278	Box with Sliding Lid, 46 x 34 x 13mm.	50	53,60
AGG2780A	Agar Grids 1000 Mesh Gold 3.05mm	25	129,40
AGG2780C	Agar Grids 1000 Mesh Copper 3.05mm	25	52,20
AGG2780N	Agar Grids 1000 Mesh Nickel 3.05mm	25	77,80
AGG2785A	Agar Grids 1500 Mesh Gold 3.05mm	15	129,40
AGG2785C	Agar Grids 1500 Mesh Copper 3.05mm	15	52,20
AGG2785N	Agar Grids 1500 Mesh Nickel 3.05mm	15	77,80
AGG2786C	Agar Grids 2000 Mesh Copper 3.05mm	10	52,20
AGG2786N	Agar Grids 2000 Mesh Nickel 3.05mm	10	77,80
AGG281A	Storage rack, plastic	each	70,80
AGG282	Slide dispenser	each	122,30
AGG283	Square Plastic Tubes, 8.5 x 8.5 x 23.7mm	100	#N/B
AGG284	Glass Vials 7ml with Lids	± 190	91,60
AGG2840	Set of 12 Tubes Identigrids A-M	50	394,80
AGG2841	Set of 12 Tubes Identigrids N-Z	50	361,90
AGG284B	Spare Caps for Vials	100	8,20
AGG285	Glass Vials 3.5ml (Pack of 666 approx. with Lids)	± 666	233,40
AGG2851	Identigrids (A) Copper 3.05mm	50	39,70
AGG2852	Identigrids (B) Copper 3.05mm	50	37,70
AGG2853	Identigrids (C) Copper 3.05mm	50	37,70
AGG2854	Identigrids (D) Copper 3.05mm	50	37,70
AGG2855	Identigrids (E) Copper 3.05mm	50	37,70
AGG2856	Identigrids (F) Copper 3.05mm	50	37,70
AGG2857	Identigrids (G) Copper 3.05mm	50	37,70
AGG2858	Identigrids (H) Copper 3.05mm	100	37,70
AGG286	Polystyrene Tubes 9mm x 30mm	50	7,30
AGG2860	Identigrids (J) Copper 3.05mm	50	37,70
AGG2861	Identigrids (K) Copper 3.05mm	50	37,70
AGG2862	Identigrids (L) Copper 3.05mm	50	37,70
AGG2863	Identigrids (M) Copper 3.05mm	50	37,70
AGG2864	Identigrids (N) Copper 3.05mm	100	37,70
AGG2865	Identigrids (P) Copper 3.05mm	50	37,70
AGG2866	Identigrids (Q) Copper 3.05mm	50	37,70
AGG2867	Identigrids (R) Copper 3.05mm	50	37,70
AGG2868	Identigrids (S) Copper 3.05mm	50	37,70
AGG2869	Identigrids (T) Copper 3.05mm	50	37,70
AGG2870	Identigrids (U) Copper 3.05mm	50	37,70
AGG2871	Identigrids (V) Copper 3.05mm	50	37,70
AGG2873	Identigrids (X) Copper 3.05mm	50	37,70
AGG2874	Identigrids (Y) Copper 3.05mm	50	37,70
AGG2875	Identigrids (Z) Copper 3.05mm	50	37,70
AGG2876	Identigrids (O) Copper 3.05mm	50	37,70
AGG288	Polythene tubes.	50	35,70
AGG2880C	Particle Analysis Grids Copper	25	129,30
AGG2880N	Particle Analysis Grids Nickel	25	134,00
AGG289	Polystyrene Tubes 13mm x 60mm High	50	11,50
AGG2905C	Square 50 Mesh Tabbed Grids Copper	100	34,50
AGG2905N	Square 50 Mesh Tabbed Grids Nickel	100	38,40

Part #	Description	Unit	Price €
AGG2910C	Square 100 Mesh Tabbed Grids Copper	100	34,50
AGG2910N	Square 100 Mesh Tabbed Grids Nickel	100	38,40
AGG29200	Gelatine Capsules Size 0	100	4,30
AGG29201	Gelatine Capsules Size 1	100	4,20
AGG29202	Gelatine Capsules Size 2	100	3,60
AGG29203	Gelatine Capsules Size 3	100	3,80
AGG29204	Gelatine Capsules Size 4	100	4,20
AGG29208	Gelatine Capsules Size 00	100	7,80
AGG29209	Gelatine Capsules Size 000	100	9,70
AGG2920A	Square 200 Mesh Tabbed Grids Gold.	25	49,70
AGG2920C	Square 200 Mesh Tabbed Grids Copper.	100	34,50
AGG2920N	Square 200 Mesh Tabbed Grids Nickel.	100	38,40
AGG29210	Gelatine Capsules Size 0	1000	44,70
AGG29211	Gelatine Capsules Size 1	1000	43,70
AGG29212	Gelatine Capsules Size 2	1000	44,70
AGG29213	Gelatine Capsules Size 3	1000	35,60
AGG29214	Gelatine Capsules Size 4	1000	33,30
AGG29218	Gelatine Capsules Size 00	1000	79,80
AGG29219	Gelatine Capsules Size 000	1000	95,80
AGG2930A	Square 300 Mesh Tabbed Grids Gold	25	55,20
AGG2930C	Square 300 Mesh Tabbed Grids Copper.	100	34,50
AGG2930N	Square 300 Mesh Tabbed Grid Nickel.	100	38,40
AGG2940A	Square 400 Mesh Tabbed Grids Gold	25	63,80
AGG2940C	Square 400 Mesh Tabbed Grids Copper.	100	38,30
AGG2940N	Square 400 Mesh Tabbed Grids Nickel.	100	41,60
AGG2950C	Square 500 Mesh Tabbed Grids Copper	100	41,60
AGG2950N	Square 500 Mesh Tabbed Grids Nickel	100	51,10
AGG2970C	Slotted 300x75 Mesh Tabbed Grids Copper	100	34,50
AGG2970N	Slotted 300x75 Mesh Tabbed Grids Nickel	100	38,40
AGG2980A	Slot 2mm x 1mm Tabbed Grids Gold	100	65,60
AGG2980C	Slot 2mm x 1mm Tabbed Grids Copper	100	38,40
AGG2980N	Slot 2mm x 1mm Tabbed Grids Nickel	100	47,90
AGG2982	CVD Diamond Hole Grid 1000µm	each	155,10
AGG2983	CVD Diamond Slot Grid 2x1mm	each	156,60
AGG2985	Finder Grids for SEM specimens	each	62,80
AGG299	Transparent Film Bags 2.75inch x 3.75inch	1000	90,10
AGG2998C	LM-SEM Locator Grid Copper	each	161,30
AGG2998N	LM-SEM Locator Grid Nickel	each	176,20
AGG301	Aluminium pin stubs, 12mm	100	36,70
AGG301A	SEM Specimen Stubs 12.5mm dia, 3.2x8mm pin no groove	100	41,90
AGG301B	Brass pin stubs.	10	15,20
AGG301C	Copper pin stubs.	10	102,60
AGG301D	Slotted specimen stub	each	10,70
AGG301E	SEM Specimen Stubs, 12.5mm dia, 45° chamfer	10	21,20
AGG301F	SEM Specimen Stubs Aluminium 12.5mm dia 3.2 x 6mm pin	100	38,70
AGG301F-B	SEM Specimen Stubs, Brass, 12.5mm dia, 3.2 x 6mm pin	100	107,80
AGG301FN	12.5mm Aluminium Stubs Short Pin Numbered 1-10	each	o.a.
AGG301-LT-10	SEM Stub with Leit Carbon tab in tube	10	23,90
AGG301-LT-100	SEM Stub with Leit Carbon tab in tube	100	170,40
AGG301-LT-12	SEM Stub w/ Leit Carbon tab in plastic box	12	29,10
AGG301-LT-50	SEM Stub with Leit Carbon tab in tube	50	96,80
AGG301-LTS-10	SEM Stub w/ Leit Carbon tab in tube - Serialised	10	48,80
AGG301-LTS-100	SEM Stub w/ Leit Carbon tab in tube - Serialised	100	324,70
AGG301-LTS-12	SEM Stubwith Leit tab in plastic box - Serialised	12	54,90
AGG301-LTS-50	SEM Stub w/ Leit Carbon tab in tube - Serialised	50	184,30
AGG301N	12.5mm Aluminium Pin Stubs Numbered 1-10	each	o.a.
AGG301P	Aluminium pin stubs - high purity	100	72,70
AGG301-S-50	SEM Specimen Stubs- Serialised	50	97,90

Part #	Description	Unit	Price €
AGG301-T-10	SEM Specimen Stub with Al core carbon Tab & Tube	each	29,60
AGG301-T-100	SEM Specimen Stub with Al core carbon Tab & Tube	100	197,00
AGG301-T-12	SEM Stub w/ Al tab in plastic box	12	35,00
AGG301-T-50	SEM Specimen Stub with Al core carbon Tab & Tube	50	119,80
AGG301-TS-10	SEM Stub with Al core Tab & Tube -Serialised	10	54,40
AGG301-TS-100	SEM Stub with Al core Tab Tab & Tube -Serialised	100	362,60
AGG301-TS-12	SEM Stub w/ Al tab in plastic box AGG3103 -Serialised	12	60,80
AGG301-TS-50	SEM Stub with Al core Tab & Tube -Serialised	50	205,80
AGG302	* Quick Drying Silver Paint	50g	159,70
AGG3020	SEM Specimen Stubs, 12.5mm dia, 45° chamfer	each	25,80
AGG3020A	SEM Specimen Stubs, 12.5mm dia, 20° chamfer	each	192,20
AGG3021	JEOL 15 x 5 mm stubs.	50	#N/B
AGG3022	JEOL 50 x 5 mm stubs.	50	200,90
AGG3023	JEOL 25 x 5 mm cylinder stubs.	50	97,70
AGG3024	JEOL 25 x 10 mm cylinder stubs.	50	150,00
AGG3025	SEM Specimen Stubs, 25mm dia, threaded pin	each	315,30
AGG3026	Flat pin stub, low profile.	10	43,60
AGG3027	90° pin stub, low profile.	10	48,80
AGG3028	38° pin stub, low profile.	10	64,20
AGG3029	36° pin stub, low profile.	10	57,10
AGG3030	Flat pin stub, low profile.	10	34,60
AGG3031	90° pin stub, low profile.	10	44,00
AGG303-10	CorrStub SEM Specimen Stubs with V-Notch	each	35,20
AGG3032-151H	SEMClip, 15 x 6mm, M4 cylinder mount, 1 clip	each	43,20
AGG3032-181	Semclip 18mm Pin Mount 1 Clip	each	44,70
AGG3032-251	Semclip 25mm Pin Mount 1 Clip	each	45,30
AGG3032-251H	SEMClip, 25 x 6mm, M4 cylinder mount, 1 clip	each	46,80
AGG3032-251J	Semclip 25mm x 8mm Cylinder Mount 1 Clip	each	39,50
AGG3032-252	Semclip 25mm Pin Mount 2 Clips	each	67,60
AGG3032-252H	SEMClip, 25 x 6mm, M4 cylinder mount, 2 clips	each	65,90
AGG3032-252J	Semclip 25mm x 8mm Cylinder Mount 2 Clips	each	57,10
AGG3032-253	Semclip 25mm Pin Mount 3 Clips	each	86,50
AGG3032-253H	SEMClip, 25 x 6mm, M4 cylinder mount, 3 clips	each	86,00
AGG3032-253J	Semclip 25mm x 8mm Cylinder Mount 3 Clips	each	74,30
AGG3032-25A	Semclip 25 x 20mm Cylinder 1 Clip 45°. Angled	each	56,20
AGG3032-321	Semclip 32mm Pin Mount 1 Clip	each	41,60
AGG3032-321H	SEMClip, 32 x 10mm, M4 cylinder mount, 1 clip	each	45,30
AGG3032-321J	Semclip 32mm x 6mm Cylinder Mount 1 Clip	each	41,60
AGG3032-322	Semclip 32mm Pin Mount 2 Clips	each	66,50
AGG3032-322H	SEMClip, 32 x 10mm, M4 cylinder mount, 2 clips	each	65,90
AGG3032-322J	Semclip 32mm x 10mm Cylinder Mount 2 Clips	each	59,00
AGG3032-323	Semclip 32mm Pin Mount 3 Clips	each	87,10
AGG3032-323H	SEMClip, 32 x 10mm, M4 cylinder mount, 3 clips	each	86,50
AGG3032-323J	Semclip 32mm x 10mm Cylinder Mount 3 Clips	each	76,20
AGG3032-381	Semclip 38mm Pin Mount 1 Clip	each	55,60
AGG3032-382	Semclip 38mm Pin Mount 2 Clips	each	73,40
AGG3032-383	Semclip 38mm Pin Mount 3 Clips	each	80,20
AGG3032-501	Semclip 50mm Pin Mount 1 Clip	each	51,80
AGG3032-502	Semclip 50mm Pin Mount 2 Clips	each	79,10
AGG3032-503	Semclip 50mm Pin Mount 3 Clips	each	92,40
AGG3032C	Replacement Clips for Semclip	10	142,60
AGG3032S	Replacement Screws for Semclip Mounts	10	43,00
AGG3033-12	Semclip 12.5mm Pin Mount 90° 1 Clip	each	74,50
AGG3034-25	Semclip 25mm Pin Mount Double 90° 2 Clips	each	80,10
AGG3034-25H	SEMClip, 25 x 16mm, M4 cylinder mount, double 90°, 2 clips	each	77,80
AGG3034-25J	Semclip 25 X16mm Cylinder Mount Double 90° 2 Clips	each	77,80
AGG3035-25	Semclip 25mm Pin Mount 45/90° 2 Clips	each	80,10
AGG3035-25H	SEMClip, 25 x 16mm, M4 cylinder mount, 45/90°, 2 clips	each	80,10

Part #	Description	Unit	Price €
AGG3035-25J	Semclip 25 X 16mm Cylinder Mount 45/90° 2 Clips	each	80,10
AGG303E	Colloidal graphite (Aquadag)	500g	45,70
AGG303E	Aquadag Graphite - ECO	500g	45,70
AGG303-LT-10	CorrStub with Leit Carbon tab in tube	10	48,00
AGG303-LT-12	CorrStub with Leit Carbon tab in plastic box	12	66,60
AGG303-LTS-10	CorrStub with Leit Carbon tab in tube - Serialised	10	67,90
AGG303-LTS-12	CorrStub w/ Leit tab in plastic box -Serialised	12	92,40
AGG303-S-10	CorrStub SEM Specimen Stubs with V-Notch - Serialised (Pk10)	10	56,20
AGG303-T-10	CorrStubs with V-Notch,with Tab & Tube	10	52,60
AGG303-T-12	CorrStub with Al tab in box	12	72,50
AGG303-TS-10	CorrStub with Al core carbon & Tube - Serialised	10	72,40
AGG303-TS-12	CorrStub w/Al tab - Serialised in box	12	92,10
AGG3040	* EPO-TEK H20E epoxy	28g	452,90
AGG3042	* EPO-TEK 353ND, epoxy	4g	21,10
AGG3045-1010	HOPG 0.4um 0.1 Mosaic Spread 10mm x 10mm x 1mm	each	445,40
AGG3045-55	HOPG 0.4um 0.1 Mosaic Spread 5mm x 5mm x 1mm	each	419,50
AGG3046-55	HOPG 0.8um 0.2 Mosaic Spread 5mm x 5mm x 1mm	each	180,00
AGG3047	AFM/SPM Disc Storage Box - 2 x 20mm Strips	each	27,80
AGG3048	AFM/SPM Disc Storage Box - 3 x 12mm Strips	each	16,60
AGG305	SEM Specimen Stubs, 32mm dia, re-entrant base	each	144,90
AGG306	JEOL stubs, aluminium.	50	20,90
AGG306A	SEM Specimen Stubs, 10mm dia, 5mm high, for JEOL heating stage	each	20,70
AGG306B	JEOL stubs, brass.	50	67,30
AGG306D	10mm x 5mm Stub with 1x1mm Slot	each	195,40
AGG306E	SEM Specimen Stubs, 10mm dia, 5mm high	each	#N/B
AGG306N	Plain Stubs for JEOL 10mm Dia. x 10mm Numbered 1-10	each	o.a.
AGG307	SEM Specimen Stubs, 15mm dia, 10mm high	each	#N/B
AGG308	SEM Specimen Stubs, 15mm dia, angled stub, 45° chamfer	each	#N/B
AGG3100	Stub Storage Box for 14 Pin-Type Stubs	each	14,70
AGG3101	Stub Storage Box for 12 x 10mm or 15mm dia cylinder stubs	each	13,70
AGG3102	Stub Storage Box for 14 x 15mm Hitachi Stubs with M4 thread	each	24,30
AGG3103	Storage Box for 12 x 12.5mm Pin Stubs	each	7,10
AGG3104	Storage Box for 12 x 10mm Cylinder Stubs	each	7,10
AGG3105	Storage Box for 3 Re-Entrant Stubs	each	7,10
AGG3106	Storage Box with 10 Tubes	each	17,00
AGG3107	Storage Kit with Carbon Tabs	each	41,70
AGG3108	Storage Kit with Adhesive Tabs	each	41,00
AGG3109	Sticky tabs.	500	72,30
AGG310C	Storage Box for 8 Stubs 1inch Dia	each	37,00
AGG310D	Storage Box for 8 Stereoscan Stubs	each	9,60
AGG311	Storage Box for 4 Stubs	each	10,60
AGG3110	Storage Kit with 10 Spectrotabs	each	46,30
AGG3111	Conductive carbon discs, roll of 250	each	79,20
AGG3112	Glue pen for mounting SEM samples.	each	9,70
AGG3113	SEM-STOR paper box storage for pin mounts.	100	105,00
AGG3113-10	Storage box for EM pin mounts (10 Packs of 100)	100	1.199,30
AGG3114	White Box for 10 Storage Tubes	each	12,00
AGG3115	Leit Tab on Pin Stub In Storage Tube	each	406,50
AGG3120	10ml Polypropylene Beakers	10	34,30
AGG3121	25ml Polypropylene Beakers	10	41,30
AGG3122	50ml Polypropylene Beakers	10	41,30
AGG3123	100ml Polypropylene Beakers	10	41,30
AGG3124	250ml Polypropylene Beakers	10	63,30
AGG3125	500ml Polypropylene Beakers	500ml	80,40
AGG3126	25ml Stacking Beakers	10	27,30
AGG3127	50ml Stacking Beakers	10	29,30
AGG3128	100ml Stacking Beakers	10	31,10
AGG3132	Gel-Pak Boxes Grade 4 89x114x12.5mm	10	#N/B

Part #	Description	Unit	Price €
AGG3133	Membrane Box 55x55x25mm	10	75,70
AGG3135-15	Cleanroom Storage Box - 15 Compartments	each	122,90
AGG3135-32	Cleanroom Storage Box - 32 Compartments	each	126,80
AGG3140-1	Polyethylene Cups 35mm Diameter x 22mm	50	15,00
AGG3140-2	Polyethylene Cups 28mm Diameter x 18mm	50	12,20
AGG3140-3	Polyethylene Cups 22mm Diameter x 16mm	50	16,10
AGG3140-4	Polyethylene Cups 18mm Diameter x 14mm	50	8,90
AGG3152-1	A1 Cutting Mat (600mm x 900mm)	each	89,90
AGG3152-2	A2 Cutting Mat (600mm x 450mm)	each	44,90
AGG3152-3	A3 Cutting Mat (450mm x 300mm)	each	21,20
AGG3152-4	A4 Cutting Mat (300mm x 200mm)	each	14,20
AGG3152-5	A5 Cutting Mat (230mm x 160mm)	each	10,80
AGG3152-6	A6 Cutting Mat (144mm x 105mm)	each	9,90
AGG3160	SEM Specimen Stubs, 12.7mm dia, 45/90° chamfer, 7.8mm pin	each	6,40
AGG3161	SEM Specimen Stubs, 12.7mm dia, low profile 45° chamfer, 9.5mm pin	each	5,10
AGG3162	SEM Specimen Stubs, 12.7mm dia, 45/90° chamfer, 9.5mm pin	each	6,40
AGG3163	SEM Specimen Stubs, 12.7mm dia, low profile 70° chamfer, 9.5mm pin. EBSD.	each	10,00
AGG3164	SEM Specimen Stubs, 25mm dia, 45/90° chamfer, 9.5mm pin	each	17,00
AGG3165	SEM Specimen Stubs, 25mm dia, double 90° chamfer, 9.5mm pin	each	15,00
AGG3166	SEM Specimen Stubs, 10mm dia, 10mm high, 45° chamfer	each	9,50
AGG3167	SEM Specimen Stubs, 15mm dia, 10mm high, 45° chamfer	each	51,00
AGG3168	SEM Specimen Stubs, 15mm dia, 10mm high, angled 70° chamfer for EBSD	each	80,20
AGG3169	SEM Specimen Stubs, 25mm dia, 20mm high, angled 30° chamfer	each	#N/B
AGG3170	SEM Specimen Stubs, 25mm dia, 20mm high, angled 45° chamfer	each	10,20
AGG3171	SEM Specimen Stubs, 25mm dia, 16mm high, angled 45/90° chamfer	each	15,20
AGG3172	SEM Specimen Stubs, 25mm dia, 16mm high, double angled 90° chamfer	each	14,30
AGG318	SEM Specimen Stubs, 32mm dia, 10mm high	each	132,40
AGG3180	Toxic GHS Labels 26x26mm	500	70,40
AGG3181	Corrosive GHS Labels 26x26mm	500	70,40
AGG3182	Dangerous for the Environment GHS Labels 26x26mm	500	70,40
AGG3183	Explosive GHS Labels 26x26mm	500	70,40
AGG3184	Flammable GHS Labels 26x26mm	500	70,40
AGG3185	Gas under Pressure GHS Labels 26x26mm	500	70,40
AGG3186	Irritating GHS Labels 26x26mm	500	70,40
AGG3187	Oxidising GHS Labels 26x26mm	500	70,40
AGG3188	Health Hazard GHS Labels 26x26mm, Roll	500	70,40
AGG319	Pin Stub - 100mm dia x 9.5mm Pin Height	each	167,10
AGG3203	* M-Bond 610 - Single Set	each	104,30
AGG3204-50	Loctite 460 sample bonding adhesive,.	50g	179,30
AGG3205	10mm aluminium SEM Environmental Microscopy stubs	100	117,70
AGG3205-3/8	3/8" Aluminium SEM Environmental Microscopy Stubs	100	117,70
AGG3205C	10mm copper SEM Environmental Microscopy stubs	50	133,00
AGG3205C-3/8	3/8" copper SEM Environmental Microscopy Stubs	100	185,70
AGG3206	10mm aluminium SEM Environmental Microscopy stubs cupped	100	134,40
AGG3206-3/8	3/8" SEM Environmental Microscopy Stubs cupped	100	134,40
AGG3207	* M-Bond 610 - Complete Kit	100	280,60
AGG3208	Flat Topped Coolstage Specimen Stubs	100	347,80
AGG3209	Dished Coolstage Specimen Stubs	100	470,40
AGG321	SEM Specimen Stubs, 12.5mm dia, carbon	each	7,90
AGG3210	11mm x 5.5mm x 4mm Type 1 Fesem Mount for Hitachi	each	80,20
AGG3211	11mm x 5.5mm x 3.5mm Type 2 Fesem Mount for Hitachi	each	95,50
AGG3212	11mm x 5.5mm x 5.5mm Type 3 Fesem Mount for Hitachi	each	108,70
AGG3213	11mm x 5.5mm x 2mm Type 4 Fesem Mount for Hitachi	each	108,70
AGG323A	SEM Specimen Stubs, 10mm dia, 10mm high, Pure Carbon Single		5,00
AGG323A-10	SEM Specimen Stubs, 10mm dia, 10mm high, Pure Carbon Pk10		40,90
AGG323A-100	SEM Specimen Stubs, 10mm dia, 10mm high, Pure Carbon Pk100		297,80
AGG327	SEM Specimen Stubs, 32mm dia, re-entrant base, carbon	each	32,00
AGG329	Ultramicrotome Glass Strips, 406 x 25.4 x 6mm	30	418,90

Part #	Description	Unit	Price €
AGG330	Knife Glass 406mm x 25.4mm x 10mm	15	333,70
AGG3300	* Leit C Conducting Carbon Cement 30ml	each	66,90
AGG3300A	* Leit C Thinners 30ml	30ml	15,90
AGG3301	* Anti-static spray	400ml	50,10
AGG3302	Leit C Plast	15g	80,10
AGG3302A	Leit C Plast Refills	2x15g	94,40
AGG3304	high performance silver paste	each	584,00
AGG3305	Tempfix	15g	73,00
AGG3309	SEM Specimen Stubs, 10mm dia, stub angled 45°	each	17,30
AGG3313	SEM Specimen Stubs, 15mm dia, 6mm high, M4 thread	50	94,40
AGG3313A	SEM Specimen Stubs, 15mm dia, 10mm high, M4 thread, angled 4	each	11,40
AGG3313B	SEM Specimen Stubs, 15mm dia, 12mm high, M4 thread, angled 7	each	11,40
AGG3313D	SEM Specimen Stubs, 15mm dia, 10mm high, M4 thread, angled 4, 45-90° cha	each	15,10
AGG3313N	15mm x 6mm Stubs M4 Thread, Hitachi, Numbered 1-10	each	o.a.
AGG3318	SEM Specimen Stubs, 32mm dia, 10mm high, M4 thread	each	151,30
AGG3319	Membrane Boxes 35mm x 35mm x 18mm	12	29,30
AGG332	Graduated Containers 100ml	100	49,60
AGG3321A	Cylinder stubs (JEOL) 2 slots for 10mm and 12.7mm dia stubs.	each	86,10
AGG3321B	M4 thread stubs (HITACHI) for 3 stubs, up to 32mm dia.	each	86,10
AGG3321C	Pin stubs 3.2mm pin dia. Slot for 12.7mm dia, 2 pin holes	each	86,10
AGG3325	SEM Specimen Stubs, 12.5mm dia, 3.2 x 15mm pin	each	84,50
AGG3327	LKB Trufs	50	403,10
AGG3327A	LKB Trufs	100	99,30
AGG3328	Leica LKB Multiplate, 220V, 25W	each	2.910,30
AGG3329	Magic tape, 12mm	66m	12,70
AGG332A	Graduated Containers 125ml	100	39,00
AGG3330	Syringes for Millipore Filters 20ml	12	150,10
AGG3331	Stainless Steel Holder for 13mm Filters	each	183,90
AGG3332	* Millipore Filters 0.45µm pore size, 13mm dia	100	183,90
AGG3333	* Millipore Filters 0.22µm pore size, 13mm dia	100	196,20
AGG3334	Millipore Cement	20ml	97,70
AGG3336-2	Silicon Wafer 2" (50.8mm)	each	50,50
AGG3336-3	Silicon Wafer 3" (76.2mm)	each	53,30
AGG3336-4P	Silicon Wafer 4" (102mm) P Type	each	82,60
AGG3337-05	Silicon Aperture Frame 0.5mm x 0.5mm	each	250,60
AGG3337-10	Silicon Aperture Frame 1.0mm x 1.0mm	each	250,60
AGG3337-15	Silicon Aperture Frame 0.5mm x 1.5mm	each	250,60
AGG3338	Swinnex Holder 25mm Diameter	each	31,30
AGG3339	Replacement Element for AGG3328 Multiplate	each	354,10
AGG333A	Lids for AGG332A Containers	100	17,10
AGG3342	Conductive Silver Pen	each	72,50
AGG3343	3mm Silicon Nitride Coated Discs	10	204,20
AGG3343H	3mm Silicon Nitride Hydrophilic Coated Discs	10	328,10
AGG3344	Ink pen for SEM	each	6,30
AGG3347A	Carbon Tabs with Aluminium 12mm	100	69,10
AGG3347N	Leit Adhesive Carbon Tabs 12mm	100	21,10
AGG3348A	Carbon Tabs with Aluminium 24mm	50	58,40
AGG3348N	Leit Adhesive Carbon Tabs 25mm	54	39,20
AGG3349	Silver epoxy kit 14g	each	97,30
AGG3350	50ml Tri-Pour Graduated Beakers	25	18,60
AGG3351	100ml Tri-Pour Graduated Beakers	25	21,40
AGG3352	250ml Tri-Pour Graduated Beakers	25	21,20
AGG3353	400ml Tri-Pour Graduated Beakers	25	36,00
AGG3354	1000ml Tri-Pour Graduated Beakers	25	58,30
AGG3355	800ml Tri-Pour Graduated Beakers	25	43,00
AGG3356	Adhesive Carbon Strips 65mm x 20mm	5	22,80
AGG3357A	Carbon Tabs with Aluminium 8mm	100	58,40
AGG3357N	Leit Adhesive Carbon Tabs 9mm	100	21,10

Part #	Description	Unit	Price €
AGG3357S	Leit Tabs 9mm Diameter	98	33,30
AGG3358	Spectrotabs 12mm dia.	120	67,50
AGG3358-25	Spectrotabs 25mm dia.	63	83,30
AGG3359A	Ecosolve NPB	1l	51,50
AGG336	Ultramicrotome Glass Strips, 203 x 50.8 x 6mm	50	248,80
AGG3360	Caps for 50ml Tri-Pour Beakers	25	10,00
AGG3361	Caps for 100ml Tri-Pour Beakers	25	11,80
AGG3362	Caps for 250ml Tri-Pour Beakers	25	#N/B
AGG3363	Caps for 400ml Tri-Pour Beakers	25	21,50
AGG3364	Caps for 1000ml Tri-Pour Beakers	25	32,30
AGG3365	Caps for 800ml Tri-Pour Beakers	25	30,90
AGG3374	Plastic Pasteur Pipettes 3ml. 1ml grads with 0.25ml incs	100	16,80
AGG3374-500	Plastic Pasteur Pipettes 3ml. 1ml grads with 0.25ml incs	500	72,20
AGG3374A	Plastic Pasteur Pipettes 3ml. 3ml grads with 0.5ml incs	100	13,40
AGG3375	Sterile Plastic Pasteur Pipettes 3ml. 1ml grads with 0.25ml incs	120	45,50
AGG3375A	Sterile Plastic Pasteur Pipettes 3ml. 3ml grads with 0.5ml incs	100	34,10
AGG3376	SEM Specimen Stubs, 32mm dia, 5mm high	each	104,20
AGG3377	SEM Specimen Stubs, 25mm dia, 6mm high, M4 thread	50	88,80
AGG3377B	Hitachi Cylinder Specimen Mount, 25 dia x 6mm x M4	10	57,50
AGG3384	SEM Specimen Stubs, 12.5mm dia, 10mm high	each	22,30
AGG3384N	10mm x 12.5mm Stubs for JEOL Numbered 1-10	each	o.a.
AGG3385	SEM Specimen Stubs, 12.5mm dia, 5mm high	each	18,60
AGG3386	SEM Specimen Stubs, 32mm dia, 20mm high	50	204,20
AGG3386-10	SEM Specimen Stubs 32mm dia 20mm high	10	51,30
AGG3387-20	SEM Specimen Stubs, 50mm dia, 10mm high	20	103,50
AGG3388	Silicon wafer, chips 5 x 5 x 0.5mm	36	88,30
AGG3389	HOPG 3.5µm 1.5 Mosaic Spread 10mm x 10mm x 2mm	each	182,60
AGG3389-1212	HOPG 3.5µm 1.5 Mosaic Spread 12mm x 12mm x 2mm	each	173,90
AGG3389-1717	HOPG 3.5µm 1.5 Mosaic Spread 17mm x 17mm x 2mm	each	199,20
AGG3390	Silicon Wafer, chips 5 x 5 mm. (Approx 270 chips per wafer)	each	172,10
AGG3390-10	Silicon wafer, chips 10 x 10 x 0.5mm.	± 55	163,20
AGG3391	Silicon Wafer, chips 5 x 7 mm. (Approx 180 chips	±180	144,30
AGG339-10	DiATOME Histo - Diamond Knife 6mm	each	o.a.
AGG339-11	DiATOME Histo - Diamond Knife 4mm	each	o.a.
AGG339-12	DiATOME Histo - Diamond Knife 8mm	each	o.a.
AGG339-13	DiATOME Ultra 45° Diamond Knife 3.5mm	each	o.a.
AGG339-1335	DiATOME Ultra 35° Diamond Knife 3.5mm	each	o.a.
AGG3392	Gripping Stub	each	91,50
AGG339-2	DiATOME Ultra 45° Diamond Knife 1.5mm	each	o.a.
AGG339-235	DiATOME Ultra 35° Diamond Knife 1.5mm	each	o.a.
AGG3393A	Membrane Box - 78mm x 39mm	each	15,90
AGG339-4	DiATOME Ultra 45° Diamond Knife 2.0mm	each	o.a.
AGG3394-15	Disposable Plastic Syringe 15ml	20	24,60
AGG3394-2	DiATOME Cryo AFM Diamond Knife 2.00mm	each	o.a.
AGG3394-3	DiATOME Cryo AFM Diamond Knife 3.00mm	each	o.a.
AGG3394-30	Disposable Plastic Syringe 30ml	20	30,40
AGG339-435	DiATOME Ultra 35° Diamond Knife 2.0mm	each	o.a.
AGG3394-60	Disposable Plastic Syringe 60ml	20	#N/B
AGG339-5	DiATOME Ultra 45° Diamond Knife 2.5mm	each	o.a.
AGG339-535	DiATOME Ultra 35° Diamond Knife 2.5mm	each	o.a.
AGG3396	Polypropylene cups	each	9,20
AGG3397	Double Sided Copper Tape 12.7mm x 16.4m	each	84,80
AGG339-7	DiATOME Ultra 45° Diamond Knife 3.0mm	each	o.a.
AGG339-735	DiATOME Ultra 35° Diamond Knife 3.0mm	each	o.a.
AGG3398	Parafilm - 100mm x 38m	each	92,10
AGG339-8	DiATOME Ultra 45° Diamond Knife 4.0mm	each	o.a.
AGG339-835	DiATOME Ultra 35° Diamond Knife 4.0mm	each	o.a.
AGG3399-6	DiATOME Histo Jumbo Diamond Knife 6.00mm, 45°	each	o.a.

Part #	Description	Unit	Price €
AGG3399-8	DiATOME Histo Jumbo Diamond Knife 8.00mm	each	o.a.
AGG3400	Kemet Diamond Compound 0.25µm	5g	42,00
AGG3401	Kemet Diamond Compound 1µm	5g	41,30
AGG3402	Kemet Diamond Compound 3µm	5g	43,70
AGG3403	Kemet Diamond Compound 6µm	5g	58,20
AGG3404	Kemet Diamond Compound 14µm	5g	63,60
AGG3407	Kemet CO-42 Cleaning Solution, Trigger Spray	450ml	17,20
AGG3408	* Kemet CO-42 Cleaning Solution	5l	93,30
AGG341	Settlement dish	each	#N/B
AGG3411	Kemet Fluid Diluent, Type OS, Trigger Spray	400ml	22,60
AGG3412	Kemet Fluid Diluent, Type W, Trigger Spray	450ml	23,00
AGG3418	* Conductive Graphite Isopropanol Base	30g	23,90
AGG3418A	* Extender for AGG3418 Conductive Graphite 30ml	30ml	12,30
AGG3420	Carbon disc on 12.5mm stub. Box of 10	10	92,50
AGG3420A	12.5mm Carbon Discs	10	70,00
AGG3421	15mm Carbon Disc	each	60,10
AGG3421A	15mm Carbon Discs	10	165,50
AGG3422	25mm Carbon Disc	each	#N/B
AGG3422A	25mm Carbon Discs	10	#N/B
AGG3423	SEM Specimen Stubs, 15mm dia, 6mm high, M4 thread, carbon	each	9,60
AGG3424	SEM Specimen Stubs, 10mm dia, angled, carbon	each	17,80
AGG3425	SEM Specimen Stubs, 15mm dia, angled, carbon	each	20,20
AGG3427	SEM Specimen Stubs, 12.5mm dia, 10mm high, carbon	each	46,50
AGG3428	Vitreous Carbon Planchette 12.7mm Dia.	each	179,90
AGG3429	Vitreous Carbon Planchette 25mm Dia.	each	219,70
AGG3430	Benchkote Pad 460mm x 570mm	50	194,10
AGG3431	Benchkote Plus 500mm x 600mm	50	189,40
AGG3432	Benchkote Plus 600mm x 50m Roll	each	354,40
AGG3435-1	DiATOME Cryo 35° Diamond Knife 1.5mm	each	o.a.
AGG3435-2	DiATOME Cryo 35° Diamond Knife 2.0mm	each	o.a.
AGG3435-3	DiATOME Cryo 35° Diamond Knife 3.0mm	each	o.a.
AGG3435-4	DiATOME Cryo 35° Diamond Knife 3.5mm	each	o.a.
AGG3435-5	DiATOME Cryo 35° Diamond Knife 4.0mm	each	o.a.
AGG3437-2	DiATOME Ultra AFM Diamond Knife 2.00mm	each	o.a.
AGG3437-3	DiATOME Ultra AFM Diamond Knife 3.00mm	each	o.a.
AGG3438	DiATOME Cryo 25 - Diamond Knife	each	o.a.
AGG3439-2	DiATOME Cryo Immuno 2.0 - Diamond Knife	each	o.a.
AGG3439-3	DiATOME Cryo Immuno 3.0 - Diamond Knife	each	o.a.
AGG3440	Sable Hair Brushes Size 3 (4mm)	3	59,60
AGG3441	Sable Hair Brushes Size 2 (3mm)	3	48,20
AGG3442	Sable Hair Brushes Size 1 (2mm)	3	45,00
AGG3443	Sable Hair Brushes Size 0 (1.5mm)	3	37,80
AGG3444	Sable Hair Brushes Size 00 (1.2mm)	3	33,70
AGG3445	Sable Hair Brushes Size 000 (0.8mm)	3	60,10
AGG3446	Set of Six Sable Hair Brushes, one of each size	each	84,10
AGG3447	Fine drawing pen	each	10,80
AGG3448-1	DiATOME Cryo 45° Diamond Knife 1.5mm	each	o.a.
AGG3448-2	DiATOME Cryo 45° Diamond Knife 2.0mm	each	o.a.
AGG3448-3	DiATOME Cryo 45° Diamond Knife 3.0mm	each	o.a.
AGG3448-4	DiATOME Cryo 45° Diamond Knife 3.5mm	each	o.a.
AGG3448-5	DiATOME Cryo 45° Diamond Knife 4.0mm	each	o.a.
AGG3449-4	DiATOME Histo Cryo Diamond Knife 4mm	each	o.a.
AGG3449-6	DiATOME Histo Cryo Diamond Knife 6mm	each	o.a.
AGG3451-16	Membrane Box 100mm x 100mm x 16mm	each	18,90
AGG3451-25	Membrane Box 100mm x 100mm x 25mm	each	20,30
AGG3451-50	Membrane Box 100mm x 100mm x 50mm	each	19,90
AGG3452A	Membrane Box - 17mm High x 82mm dia with 4 Clips	each	16,00
AGG3453	Storage boxes with membrane - Round. 165 x 42mm	each	25,40

Part #	Description	Unit	Price €
AGG3454	Membrane Box - 80mm High x 165mm Diameter	each	65,30
AGG3455	Membrane Box 125mm x 75mm x 25mm	each	28,00
AGG3456	Membrane Box 175mm x 50mm x 50mm	each	21,00
AGG3457	Membrane Box 175mm x 125mm x 50mm	each	24,10
AGG3458	Membrane Box 150mm x 150mm x 50mm	each	28,60
AGG3459	Membrane Box 200mm x 200mm x 50mm	each	30,80
AGG346	5ml Polythene Bottles	10	#N/B
AGG3462	* Bell Shine Spray 115ml	each	40,90
AGG3465	Specimen pins for cryo-ultramicrotomy -	50	184,20
AGG3466	Freezer Hats Type A	50	297,50
AGG3467	Freezer Hats Type B	50	344,30
AGG3468	DiATOME Static Line II - 230V	each	o.a.
AGG3469	DiATOME Ultra Sonic	each	o.a.
AGG347	10ml Polythene Bottles	10	46,90
AGG348	30ml Polythene Bottles	10	53,90
AGG349	Boxes with lift off lids, 42 x 19 x 17mm.	100	#N/B
AGG350	Glass Vials 1.75ml, with Lids	±880	271,40
AGG350-10	Sample pin for cryo ultramicrotomes	each	26,50
AGG350-100	Sample pin for cryo ultramicrotomes	100	226,00
AGG3504	FIB SEM mould for biological specimens	each	27,50
AGG3505	Micromoulds	10	20,00
AGG350-50	Sample pin for cryo ultramicrotomes	50	121,80
AGG3506	Polypropylene Petri Dishes	100	45,70
AGG3508A	Microwave Capsule Holder with Lid	each	368,60
AGG3509	Microwave Holder for Microcentrifuge Tubes	each	93,20
AGG351	440ml Wash Bottle	each	4,60
AGG3515	Micromoulds	100	174,20
AGG3516	Ultramicrotome Glass Strips, 406 x 25.4 x 8mm	24	345,00
AGG3517	Aspirator bottle	each	42,70
AGG3518	Desiccator cabinet	each	260,60
AGG3519	Hot Vac mini vacuum desiccator	each	1.203,50
AGG352	250ml Wash Bottle	each	4,30
AGG3524	Giammara-Hanker Cast-A-Slide translucent mould	each	23,90
AGG3525	BEEM capsule holder UV transparent	each	26,00
AGG3527	Pestle and mortar 70mm dia	each	28,80
AGG353	Resealable Plastic Bags Plain 3inchx3.25inch	100	3,30
AGG3530	Flat Embedding Mould - 15 Blocks	each	39,30
AGG3531-1	20 Block Flat Embedding Mould, Blue	each	58,80
AGG3532	DWS 24 + 3 Block Embedding Mould	each	44,70
AGG3533	Flat Embedding Mould (20 Cavities & 8 Rectangles)	each	41,10
AGG3535	* Conductive Gold Past	2g	309,90
AGG3537	MICRO liquid cleaner	each	26,90
AGG3538	Grid coating pen	each	52,20
AGG354	Resealable Plastic Bags Plain 4inchx5.5inch	100	5,30
AGG3547	Flat Bottomed Capsules	each	30,20
AGG3549	21 Cavity Blue Silicone Rubber Mould	each	27,40
AGG355	Resealable Plastic Bags Plain 5inchx7.5inch	100	8,10
AGG3552	Plastic embedding stubs	each	91,30
AGG356	Resealable Plastic Bags Plain 8inchx11inch	100	9,90
AGG3560	500ml Wash Bottle - Printed Water	each	8,50
AGG3560A	250ml Wash Bottle - Printed Water	each	7,90
AGG3561	500ml Wash Bottle - Printed Distilled Water	each	9,40
AGG3561A	250ml Wash Bottle - Printed Distilled Water	each	7,90
AGG3562	500ml Wash Bottle - Printed Acetone	each	9,40
AGG3562A	250ml Wash Bottle - Printed Acetone	each	7,90
AGG3563	500ml Wash Bottle - Printed Methanol	each	8,50
AGG3563A	250ml Wash Bottle - Printed Methanol	each	7,90
AGG3564	500ml Wash Bottle - Printed Isopropanol	each	9,50

Part #	Description	Unit	Price €
AGG3564A	250ml Wash Bottle - Printed Isopropanol	each	7,90
AGG3565	500ml Wash Bottle - Printed I.M.S.	each	8,50
AGG3565A	250ml Wash Bottle - Printed I.M.S.	each	7,90
AGG3566	500ml Wash Bottle - Printed Ethanol	each	8,50
AGG3566A	250ml Wash Bottle - Printed Ethanol	each	7,90
AGG3567	Gelatin capsule holder	each	45,80
AGG357	Resealable Plastic Bags Panelled 3inchx3.25inch	100	3,10
AGG3570	Spotting Brush Size 3	each	31,10
AGG3571	Spotting Brush Size 2	each	24,00
AGG3572	Spotting Brush Size 1	each	25,40
AGG3573	Spotting Brush Size 0	each	22,90
AGG3574	Spotting Brush Size 00	each	24,00
AGG3575	Spotting Brush Size 000	each	20,50
AGG3577	DiATOME Cryotrim 20 - Diamond Knife	each	o.a.
AGG3578	DiATOME Cryotrim 45 - Diamond Knife	each	o.a.
AGG358	Resealable Plastic Bags Panelled 4inchx5.5inch	100	5,30
AGG3580	Driplok Wash Bottle 500ml Custom/Diy	each	15,80
AGG3580A	Driplok Wash Bottle 250ml Custom/DIY	each	10,90
AGG3581	Driplok Wash Bottle 500ml Acetone	each	13,50
AGG3581A	Driplok Wash Bottle 250ml Acetone	each	11,50
AGG3582	Driplok Wash Bottle 500ml Methanol	each	12,20
AGG3582A	Driplok Wash Bottle 250ml Methanol	each	11,50
AGG3583	Driplok Wash Bottle 500ml Isopropanol	each	13,00
AGG3583A	Driplok Wash Bottle 250ml Isopropanol	each	11,50
AGG3584	Driplok Wash Bottle 500ml I.M.S.	each	12,20
AGG3584A	Driplok Wash Bottle 250ml I.M.S.	each	11,50
AGG3585	Driplok Wash Bottle 500ml Ethanol	each	12,20
AGG3585A	Driplok Wash Bottle 250ml Ethanol	each	11,50
AGG3586	Driplok Wash Bottle 500ml M.E.K.	each	18,60
AGG359	Resealable Plastic Bags Panelled 5inchx7.5inch	100	8,80
AGG3590	SEM Metallurgical Mount Sample Holder (1.25 inch Dia)	each	112,10
AGG3590A	SEM Metallurgical Mount Sample Holder (1.5 inch Dia)	each	108,20
AGG3590B	SEM Metallurgical Mount Sample Holder (1inch Dia)	each	129,80
AGG3591	Hitachi three stub holder	each	#N/B
AGG3592	Hitachi holder for six pin stubs	each	114,10
AGG3595	SEM Specimen Stubs, 25.4mm square. Gun residue holders	each	98,50
AGG360-1	BEEM Capsules 00	100	14,20
AGG360-2	BEEM Capsules 00	1000	124,00
AGG3608	Ultipette 1-10ul	each	257,00
AGG3608-10	Ultipette - 10 Microlitre	each	168,00
AGG3608-20	Ultipette - 20 Microlitre	each	168,00
AGG3608-5	Ultipette - 5 Microlitre	each	168,00
AGG3608-50	Ultipette - 50 Microlitre	each	168,00
AGG3608-V5	Ultipette 0.5-5ul	each	257,00
AGG3609	Ultipette capillary tips 0.5-10ul	500	84,00
AGG3611	1-10ul Ultipette with Tip Ejector	each	347,80
AGG361-1	BEEM Capsules 00 Conical	100	14,20
AGG3612	Micro-classic pipette controller	each	114,80
AGG361-2	BEEM Capsules 00 Conical	1000	124,00
AGG3613	65mm dia, Disposable Desiccant Cartridges	6	201,90
AGG3614	100mm dia, Disposable Desiccant Cartridges	6	319,40
AGG3615	Oil Mist Filter for AGB7366 Pumping System	each	69,00
AGG3617	Cartridges for 60mm desiccator, 114mm dia	230g	111,00
AGG3618	Cartridges for 200mm desiccator, 197mm dia	680g	243,90
AGG3619	Cartridges for 250mm desiccator, 238mm dia	1360g	295,20
AGG362-1	BEEM Capsules 3	100	15,30
AGG362-2	Beem Capsules 3	1000	120,60
AGG3623	Oil Mist Filter 5m3/hr (83 l/min)	each	103,50

Part #	Description	Unit	Price €
AGG3623A	Adaptor for AGG3623 Oil Mist Filter for Edwards pumps	each	237,10
AGG3624	Oil Mist Filter 15m ³ /hr (250 l/min)	each	351,40
AGG3624A	Adaptor for AGG3624 Oil Mist Filter for Edwards pumps	each	414,60
AGG3626	Storage Tubes for 12.5mm Pin Stubs	10	6,40
AGG3626S	Storage Tubes for 12.5mm pin stubs	1000	294,70
AGG3630	13mm Swin-Lok Filter Holders	each	20,30
AGG3631	25mm Swin-Lok Holder	each	38,10
AGG363-1	BEEM Capsules 00 Hemi-Hyperbola	100	16,00
AGG363-2	BEEM Capsules 00 Hemi-Hyperbola	1000	135,70
AGG3633	47mm Swin-Lok Holder	each	575,20
AGG3634	Specimen grid box for 50 grids	each	17,90
AGG3635	Specimen grid box for 50 grids with unique number	each	19,60
AGG3636WN	Storage Tubes for Pin Type Stubs up to 38mm dia.	10	24,30
AGG364	BEEM capsule holder 00	each	23,40
AGG3643	Storage box for 50 tabbed grids	each	11,30
AGG3644	Tungsten Wire Circlips	25	52,60
AGG3647	Histology Glass 406mm x 38mm x 6mm	20	215,60
AGG3648	* Electrodag 1415	50g	157,80
AGG3649	* Electrodag 1415	150g	403,10
AGG365	BEEM capsule holder 3	each	21,40
AGG3650	Plastic Petri Dishes (sterile).0	500	570,70
AGG3650A	Plastic Petri Dishes (sterile).	10	15,70
AGG3651	Aluminium Dishes 53mm dia x 17mm height	144	33,90
AGG3653	Grid coating plates	each	25,50
AGG3654	BEEM flat embedding mould	each	17,40
AGG3656A	* DAG 580 dry film lubricant	850g	115,40
AGG3657	Pin type adaptor for Cambridge analytical stub	each	183,20
AGG3658	12.5mm dia cylinder stub converter for JEOL	each	18,60
AGG3659	15mm dia cylinder stub converter with M4 thread for Hitachi	each	28,90
AGG3659A	15mm dia cylinder stub converter with M4 thread and 45° an	each	62,40
AGG3660	10mm dia cylinder stub converter for JEOL	each	17,60
AGG3661	15mm cylinder stub converter for ISI/ABT/Topcon	each	72,20
AGG3662	SEM Stub (pin format) for TEM Grids	each	93,80
AGG3662A	Cryo Grid Holder for SEM (Cu)	each	176,50
AGG3679	Polypropylene capsule chambers	each	40,90
AGG3680	Electronic timer - count down/count up	each	39,50
AGG3681	Dual Size SEM Sample Holder 12.5mm	each	96,00
AGG3682	ISI DS 130 Tilting Sample Holder	each	o.a.
AGG3683	Cambridge Carousel for 8 pin stubs	each	91,70
AGG3686	25mm Freeze Drying Pot	3	93,30
AGG3688	Mini desiccator	each	474,00
AGG3689	* Graphit '33' spray.	200ml	35,90
AGG3691	* Quick Drying Silver Paint with Brush	25g	91,30
AGG3692	* Electrodag 1415 with Brush	25g	91,30
AGG3694	BEEM Dial-a-Grid and block holder module	each	10,10
AGG3695	SEM Stub Freeze Dry Holder	each	65,10
AGG3696	Preparation Stand for Pin Type and 10mm dia Stubs	each	183,20
AGG3697	Preparation Stand for Pin Type and 15mm dia Stubs	each	122,70
AGG3698	Specimen Preparation Stand for Pin and 12.5mm Stub	each	167,90
AGG3699	BEEM capsule Chuck, Stainless Steel	each	103,20
AGG3700	Polypropylene Containers 30ml	72	88,20
AGG3701	Polypropylene Containers 60ml	48	75,10
AGG3702	Polypropylene Containers 125ml	36	75,30
AGG3703	Polypropylene Containers 250ml	36	139,80
AGG3704	Polyethylene Vials 1.6ml	100	106,10
AGG3705	Polyethylene Vials 2ml	100	110,30
AGG3706	Polyethylene Vials 8.4ml	100	153,30
AGG3707	Circlips 3mm dia	3	106,50

Part #	Description	Unit	Price €
AGG3707-1	35mm O.D. x 10mm Mortar & Pestle	each	134,40
AGG3707-2	50mm O.D. x 14mm Mortar & Pestle	each	177,50
AGG3707-3	65mm O.D. x 18mm Mortar & Pestle	each	305,60
AGG3708	Circlip Extraction Tool	each	108,20
AGG3708-1	Quartz Disc 1inch x 1/16inch (25mm x 1.5mm)	each	22,40
AGG3708-10	Quartz Disc 4inch x 1/8inch (100mm x 3.2mm)	each	171,90
AGG3708-2	Quartz Disc 1inch x 1/8inch (25mm x 3.2mm)	each	22,40
AGG3708-3	Quartz Disc 2inch x 1/16inch (50mm x 1.5mm)	each	58,50
AGG3708-4	Quartz Disc 2inch x 1/8inch (50mm x 3.2mm)	each	61,90
AGG3708-7	Quartz Disc 3inch x 1/16inch (75mm x 1.5mm)	each	103,30
AGG3708-8	Quartz Disc 3inch x 1/8inch (75mm x 3.2mm)	each	109,10
AGG3708-9	Quartz Disc 4inch x 1/16inch (100mm x 1.5mm)	each	152,90
AGG371	* Lacomit Varnish	500ml	32,90
AGG3710	Jeol 2 Position Single Tilt Holder	each	o.a.
AGG3711	Microcentrifuge Tubes Plain	100	11,60
AGG3711C	Microcentrifuge Tubes Coloured	50	62,30
AGG3712	Glass Scoring Jigs	each	109,80
AGG3713-12	Kapton Plastic Film Tape - 12mm	33m	115,70
AGG3713-19	Kapton Plastic Film Tape - 19mm	33m	168,20
AGG3713-25	Kapton Plastic Film Tape - 25mm	33m	172,30
AGG3714	Embedding mould, polypropylene, 23 x 23 x 8mm	50	33,10
AGG3714-12	Type 92 Polyimide Film Tape - 12mm	33m	61,10
AGG3714-19	Type 92 Polyimide Film Tape - 19mm	33m	96,70
AGG3714-25	Type 92 Polyimide Film Tape - 25mm	33m	143,10
AGG3715	Humidity Indicator Cards 5/10/60% RH	125	53,90
AGG3715A	Humidity Indicator Cards 5/10/15% RH	125	45,70
AGG3716	4 Position Multi Sample Holder for FEI Microscopes	each	o.a.
AGG3717	Single Tilt Single Sample Holder for FEI Microscopes	each	o.a.
AGG3718	Hitachi 3 Postions Multi Sample Holder	each	o.a.
AGG3719	FIB grid storage box	each	24,90
AGG372	* Lacomit Remover	500ml	19,00
AGG3720	3" Diameter SEM stub	each	111,50
AGG3721	* Copper spray	each	63,50
AGG3722	Cardboard Boxes 50mm x 32mm x 11mm	144	160,50
AGG3723	Cardboard Boxes 62mm x 41mm x 18mm	144	150,00
AGG3724	Cardboard Boxes 76mm x 50mm x 32mm	72	92,70
AGG3725	Cardboard Boxes 90mm x 65mm x 35mm	36	100,80
AGG3726	Grid storage box for 96 grids	each	30,70
AGG3727	Round Cryo Grid Box with Lid	each	22,90
AGG3728	Round Cryo Grid Box Base only	each	17,10
AGG3729	Square Cryo Grid Box with Lid	each	41,20
AGG3731	Nalgene cryoware labels	200	99,90
AGG3733	Cryo Grid Box Handling Rod	each	42,60
AGG3734	* Electrodag 502.	30g	24,70
AGG3735	Round Cryo Grid Box with Pin Type Lid	each	42,40
AGG3736	4-hole Cryo-EM storage box, Round shape to fit FEI Vitrobot	each	16,80
AGG373A	* Quadratch instrument cleaner	5l	61,70
AGG3740	Snap-Fit Gelatin Capsules Size 00	1000	35,80
AGG3741	Snap-Fit Gelatin Capsules Size 1	1000	32,00
AGG3742	Easy-Moulds 5.6mm	5	86,50
AGG3743	Easy-Moulds 8mm	5	93,70
AGG3744	Flat Embedding Capsules (Polyethylene)	500	48,30
AGG3745	Pyramid Tip Mould	each	28,40
AGG3746	Lucite Mounting Cylinders	100	96,60
AGG3747	PTFE Mould - 16 Cavities	each	153,20
AGG3748	PTFE Mould - 3 Cavities	each	153,20
AGG3749	Glass Knife Holder for 6-7mm Thickness	each	90,80
AGG3749A	Glass Knife Holder for 8mm Thickness	each	92,90

Part #	Description	Unit	Price €
AGG375	Anti-static gun	each	177,20
AGG3754	Ependoff Microtome Chuck RMC	each	116,50
AGG3756	Grid Staining Matrix System	each	181,70
AGG3757	Staining Vessel	each	150,20
AGG3758	Microwell staining mould	each	28,10
AGG3759	Flat Bottomed Embedding Capsules (Polypropylene)	500	51,60
AGG3760	Peel-A-Way Moulds 22mm x 22mm x 20mm	288	401,80
AGG3761	Peel-A-Way Moulds 22mm x 30mm x 20mm	288	408,20
AGG3762	Peel-A-Way Moulds 22mm x 40mm x 20mm	264	413,00
AGG3763	Peel-A-Way Moulds 22mm Tapered to 8mm	288	#N/B
AGG3764	Peel-A-Way Moulds 22mm Tapered to 12mm	288	362,30
AGG3765	PTFE immunostaining pad	each	106,30
AGG3766	Chien staining pad	each	20,50
AGG3768	Slit Attachment	each	119,80
AGG3769	3 Position Grid Holder, fits standard JEOL 2pos holders	each	121,50
AGG3771	Compact Magnetic Stirrer, 240V	each	477,30
AGG3772	Accessory Retort Rod, 12.5 x 300mm	each	49,00
AGG3773	SEM Stub pin Adaptors 10mm	each	124,80
AGG3774	SEM Stub pin Adaptors 15mm	each	126,40
AGG3775	SEM Stub pin Adaptors 12.5mm	each	128,00
AGG3776	TEM Grid Freeze Drying Sample Holder	each	129,70
AGG3777	Wood Sample Holder boxes	each	131,40
AGG3778	12mm Magnetic Follower	3	16,20
AGG3779	20mm Magnetic Followers	3	18,80
AGG3780	25mm Magnetic Followers	3	23,60
AGG3781	40mm Magnetic Followers	3	27,50
AGG3783	Magnetic Retriever, PTFE, 8mm dia x 350mm	each	22,80
AGG3784	Magnetic Stirrer with Stainless Steel Top Plate, 160 X 160mm	each	441,30
AGG3785	Magnetic Stirrer with Ceramic Top Plate, 160 x 160mm, 230V	each	449,70
AGG3787	Magnetic Stirrer/Hotplate with Ceramic Top Plate, 160 x 160mm	each	642,70
AGG3788	Magnetic Stirrer/Hotplate (Coated Aluminium Silicon) US152	each	628,30
AGG3788C	Magnetic stirrer / hotplate (Glass ceramic) UC152	each	628,30
AGG3790	* Conductive Silver Paint	30g	139,30
AGG3791	* Diluent for AGG3790 + AGG3535	25ml	12,80
AGG3792	Silver adhesive sheet double coated	each	857,80
AGG3793	* Mikrostik non-conductive adhesive	14ml	17,20
AGG3794	Carbon adhesive sheet double coated	10	92,40
AGG3795	PTFE Mould - 11 Cavities	each	161,10
AGG3796	MENCO grid staining pad	each	137,00
AGG3798	Digital magnetic stirrer/hotplate with ceramic top plate, 160 x 160mm, 230V, !	each	1.473,50
AGG3799	Portable Magnetic Stirrer	each	509,20
AGG3799A-EU	9V Mains Adaptor for AGG3799 Portable Stirrer, 230 V AC, European two pin plug	each	70,00
AGG3799A-UK	9V Mains Adaptor for AGG3799 Portable Stirrer, 230 V AC, UK plug	each	70,00
AGG3800	Microtome Chuck, cylindrical, Al, 22mm dia	each	136,40
AGG3801	LKB Knifemaker 7800 Repair Kit	each	172,50
AGG3802	13mm dia 10µm Nuclepore Filters	100	150,00
AGG3803	13mm dia 8µm Nuclepore Filters	100	124,00
AGG3807	13mm dia 1.0µm Nuclepore Filters	100	100,90
AGG3808	13mm dia 0.8µm Nuclepore Filters	100	152,40
AGG3809	LKB Knifemaker 2178 MK2 cutter	each	141,30
AGG381	Dental wax	each	16,60
AGG3811	13mm dia 0.2µm Nuclepore Filters	100	113,70
AGG3812	13mm dia 0.1µm Nuclepore Filters	100	165,00
AGG3816	13mm dia 0.015µm Nuclepore Filters	100	312,00
AGG3821	25mm dia 12µm Nuclepore Filters	100	123,50
AGG3822	25mm dia 10µm Nuclepore Filters	100	118,90
AGG3823	25mm dia 8µm Nuclepore Filters	100	126,60
AGG3824	25mm dia 5µm Nuclepore Filters	100	134,60

Part #	Description	Unit	Price €
AGG3825	25mm dia 3µm Nuclepore Filters	100	118,90
AGG3826	25mm dia 2µm Nuclepore Filters	100	121,30
AGG3827	25mm dia 1µm Nuclepore Filters	100	118,90
AGG3828	25mm dia 0.8µm Nuclepore Filters	100	184,40
AGG3829	25mm dia 0.6µm Nuclepore Filters	100	124,70
AGG3830	25mm dia 0.4µm Nuclepore Filters	100	#N/B
AGG3831	25mm dia 0.2µm Nuclepore Filters	100	145,40
AGG3832	25mm dia 0.1µm Nuclepore Filters	100	245,40
AGG3834	25mm dia 0.05µm Nuclepore Filters	100	182,10
AGG3836	25mm dia 0.015µm Nuclepore Filters	100	334,00
AGG384	Tackiwax	each	10,30
AGG3841	47mm dia 12µm Nuclepore Filters	100	173,50
AGG3842	47mm dia 10µm Nuclepore Filters	100	182,10
AGG3843	47mm dia 8µm Nuclepore Filters	100	285,50
AGG3844	47mm dia 5µm Nuclepore Filters	100	202,70
AGG3845	47mm dia 3µm Nuclepore Filters	100	233,10
AGG3846	47mm dia 2µm Nuclepore Filters	100	182,10
AGG3847	47mm dia 1µm Nuclepore Filters	100	234,80
AGG3848	47mm dia 0.8µm Nuclepore Filters	100	236,00
AGG3849	47mm dia 0.6µm Nuclepore Filters	100	182,10
AGG3850	47mm dia 0.4µm Nuclepore Filters	100	260,20
AGG3851	47mm dia 0.2µm Nuclepore Filters	100	253,80
AGG3852	47mm dia 0.1µm Nuclepore Filters	100	432,40
AGG3853	47mm dia 0.08µm Nuclepore Filters	100	405,20
AGG3854	47mm dia 0.05µm Nuclepore Filters	100	432,40
AGG3856	47mm dia 0.015 Nuclepore Filters	100	363,50
AGG3860	Refillable Dispenser for Polishing Compound	each	12,80
AGG3861	Refillable Dispenser for Lubricating Fluid	each	12,80
AGG3862	Diamond Compound 1/4 Micron	400g	61,50
AGG3863	Diamond Compound 1 Micron	400g	61,50
AGG3864	Diamond Compound 3 Micron	400g	98,30
AGG3865	Diamond Compound 6 Micron	400g	84,40
AGG3866	Diamond Compound 8 Micron	400g	93,70
AGG3867	Diamond Compound 14 Micron	400g	98,30
AGG3868	Diamond Compound 25 Micron	400g	82,30
AGG3869	Diamond Compound 45 Micron	400g	84,40
AGG3870	Lubricating Fluid 2.5 Litre	2,5l	61,90
AGG3871	Plastic Storage Box - 7 Compartments	each	13,00
AGG3872	Plastic Storage Box - 12 Compartment (New Design)	each	25,40
AGG3873	Plastic Storage Box - 9 Compartments	each	20,30
AGG3874	Plastic Storage Box - 18 Compartments	each	18,60
AGG3877	Plastic storage boxes	each	49,50
AGG3880	Thermoplastic Wax (80°C Melting Point)	179g	180,90
AGG3881	Thermoplastic Wax (52°C Melting Point)	100g	94,30
AGG3882	Replicating Tape - 125µm, 19mm x 4.5m	each	78,10
AGG3883	Replicating Tape - 125µm, 38mm x 4.5m	each	111,30
AGG3884A	Replicating Tape - 50µm, 38mm x 1.1m	each	67,40
AGG3885	Storage Box - 14 Compartments	each	19,90
AGG3886	Storage Box - 21 Compartments	each	23,10
AGG3887	Lakeside Cement, 12 bars, approx	250g	206,30
AGG3888	President SEM replication kit, 2 x 53ml	each	#N/B
AGG390	Vacuum Tweezers, complete, 240V	each	664,50
AGG3900	Potassium Bromide Crystal 10mm x 10mm x 50mm	each	140,20
AGG3903	Borosilicate Petri Dish - 55mm Diameter	each	47,50
AGG3904	Borosilicate Petri Dish - 74mm Diameter	each	53,80
AGG3905	Borosilicate Petri Dish - 94mm Diameter	each	62,00
AGG3906	Pyrex Petri Dish - 55mm Diameter	each	31,60
AGG3907	Pyrex Petri Dish - 74mm Diameter	each	34,60

Part #	Description	Unit	Price €
AGG3908	Pyrex Petri Dish - 94mm Diameter	each	41,20
AGG3909	5ml Balance Boats	50	23,20
AGG3909BL	5ml Balance Boats Black	50	22,60
AGG390A	Set of spare nozzles and tubing	each	175,60
AGG3910	30ml Balance Boats	100	9,90
AGG3910B	30ml Balance Boats Black	100	9,90
AGG3911	100ml Balance Boats	50	11,30
AGG3911B	100ml Balance Boats Black	50	11,30
AGG3912	Aluminium Dishes 43mm dia x 12.5mm height	100	59,30
AGG3918	Polyethylene mincing dishes	100	87,90
AGG3919	Swinnex Holder 13mm Diameter	each	17,90
AGG3920	* Millipore Filters 0.45µm pore size, 25mm dia	100	#N/B
AGG3921	* Millipore Filters 0.22µm pore size, 25mm dia	100	195,90
AGG3926	Vacuum Storage Vessel (Small)	each	292,00
AGG3927	Vacuum Storage Vessel (Medium)	each	481,10
AGG3928	Vacuum Storage Vessel (Large)	each	551,20
AGG3934	Anti-static weighing boats.	250	34,10
AGG3935	Carbon Tape with Aluminium 8mm x 20m	each	66,90
AGG3935A	Carbon Tape with Aluminium 12mm x 20m	each	92,80
AGG3935B	Carbon Tape with Aluminium 20mm x 20m	each	154,40
AGG3939	Adhesive Carbon Tape 8mm x 20m	each	57,10
AGG3939A	Adhesive Carbon Tape 12mm x 20m	each	96,30
AGG3939B	Adhesive Carbon Tape 20mm x 20m	each	145,90
AGG3939C	Adhesive Carbon Tape 50mm x 20m	each	261,40
AGG394	Resealable Plastic Bags Plain 57x57mm	100	2,90
AGG3940	Copper Tape with Conductive Adhesive 12mm x 50m	each	86,20
AGG3941	Sodium Chloride Crystal 10mm x 10mm x 50mm	each	138,80
AGG3949	2ml Sample Vials 12mm x 38mm with Caps	288	290,10
AGG395	Protective goggles	each	45,90
AGG3950	Plastic Boxes 25mm x 25mm x 8mm	100	35,80
AGG3951	Plastic Boxes 25mm x 25mm x 16mm	50	19,30
AGG3952	Plastic Boxes 32mm x 32mm x 32mm	16	14,10
AGG3953	Plastic Boxes 73mm x 50mm x 11mm	6	9,20
AGG3953S	Plastic Boxes 73mm x 50mm x 11mm	44	52,70
AGG3954	Plastic Boxes 95mm x 29mm x 8mm	10	9,80
AGG3955	Plastic Boxes 115mm x 89mm x 25mm	6	23,40
AGG3956	Plastic Boxes 190mm x 37mm x 12mm	12	21,70
AGG3957	Plastic Boxes 50mm x 50mm x 19mm	50	34,30
AGG3958	Plastic Boxes 73mm x 50mm x 25mm	10	12,80
AGG3959	Clear round polystyrene containers	20	13,30
AGG396	Bench Protection Sheets 570mm x 460mm	100	199,30
AGG3961	Glass Vials 30ml	± 210	155,10
AGG3962	Lids for 30ml Glass Vials	500	58,80
AGG3964	3.5ml Glass Vials	±544	193,40
AGG3964-1	Lids for AGG3964 Glass Vials	50	45,90
AGG397	Benchkote Roll 46cm	50m	298,30
AGG3970	Storage box	each	#N/B
AGG3971	Mini vials	125	127,80
AGG3971-1000	Mini vials	1000	720,50
AGG3972	Grid vials	each	70,90
AGG3972B	Anti-static plastic vials	100	93,10
AGG3973	Storage Rack - 30 Drawers	each	136,80
AGG3974	Storage Rack - 45 Drawers	each	158,00
AGG3977	Dust cover and base	each	#N/B
AGG3979	Potassium Chloride Crystals 10x10x10mm	5	86,00
AGG398	Parafilm - 50mm x 76m	each	75,60
AGG399	Aluminium pin stubs, 25 mm dia.	50	99,60
AGG3994	Wafer Trays for 1in dia wafers	24	176,90

Part #	Description	Unit	Price €
AGG3995	Wafer Trays for 2in dia wafers	24	202,20
AGG3996	Wafer Trays for 4in dia wafers	24	232,60
AGG3998	Wafer Trays for 3in dia wafers	24	336,90
AGG399B	Brass pin stubs, 25 mm dia.	10	41,40
AGG399F	SEM Specimen Stubs, 25mm dia, 3.2 x 6mm pin	50	99,10
AGG400	Aluminium pin stubs, 32 mm dia.	50	124,40
AGG4001-100	100ml Tri-Cornered Beaker	25	#N/B
AGG4001-1000	1000ml Tri-Cornered Beaker	25	#N/B
AGG4001-250	250ml Tri-Cornered Beaker	25	#N/B
AGG4001-400	400ml Tri-Cornered Beaker	25	#N/B
AGG4001-50	50ml Tri-Cornered Beaker	25	#N/B
AGG4001-800	800ml Tri-Cornered Beaker	25	#N/B
AGG4002	Double Tamper Evident Bottle 200ml white cap, shrink-wrap neck ring	100	397,80
AGG4003	Double Tamper Evident Bottle 200ml white cap, shrink-wrap neck ring	10	58,80
AGG4004	Double Tamper Evident Bottle 200ml white cap	100	397,80
AGG4005	Double Tamper Evident Bottle 200ml white cap	10	58,80
AGG400B	Brass pin stubs, 32 mm dia.	10	188,50
AGG4023	Small sample box	each	12,30
AGG4051	Large Storage Box	each	104,70
AGG410	Micro Pipettes, intraMARK, 1/2/3/4/5 125mm	250	139,70
AGG411	Micro Pipettes, intraMARK, 10 125mm	250	87,00
AGG412	Micro Pipettes, intraMARK, 20 125mm	250	87,00
AGG413	Micro Pipettes, intraMARK, 25 125mm	250	87,00
AGG414	Micro Pipettes, intraMARK, 20/40 125mm	250	117,70
AGG415	Micro Pipettes, intraMARK, 40 125mm	250	87,00
AGG416	Micro Pipettes, intraMARK, 50 125mm	250	87,00
AGG417	Micro Pipettes, intraMARK, 50/100 125mm	250	90,10
AGG418	Micro Pipettes, intraMARK, 100 125mm	250	90,10
AGG419	Micro Pipettes, intraMARK, 200 125mm	200	167,40
AGG420	Pipette aid for Micro Pipettes	each	23,50
AGG4501	25 Place Polycarbonate Box Clear for 0.5-2ml Tubes	8	76,10
AGG4502	25 Place Polycarbonate Box with white base for 0.5-2ml Tubes	8	79,30
AGG4503	81 Place Polycarbonate Storage Box with fixed rack White for 0.5-2ml Tubes	5	113,20
AGG4504	81 Place Polycarbonate Storage Box with fixed rack Blue for 0.5-2ml Tubes	5	127,00
AGG4505	81 Place Polycarbonate Storage Box with fixed rack Green for 0.5-2ml Tubes	5	127,00
AGG4506	81 Place Polycarbonate Storage Box with fixed rack Orange for 0.5-2ml Tubes	5	127,00
AGG4507	81 Place Polycarbonate Storage Box with fixed rack Yellow for 0.5-2ml Tubes	5	127,00
AGG4508	81 Place Polycarbonate Storage Box with fixed rack Mixed Colours for 0.5-2ml	5	93,10
AGG4509	81 Place Polycarbonate Storage Box with removable rack, for tubes up to 12m	5	114,30
AGG4510	100 Place Cryogenic Storage Box, white insert & yellow base for small tubes	5	113,10
AGG4511	100 Place Polycarbonate Storage Box, Yellow base & Insert & Printed Lid for 0.!	5	83,00
AGG4512	100 Place Polycarbonate Storage Box, White base & Insert & Printed Lid for 0.!	5	124,10
AGG4513	100 Place Polycarbonate Completely Clear Box	5	129,10
AGG4514	Cryogenic Storage Box without grids and removable lid	5	40,80
AGG4515	81 Place Polycarbonate Tall Tube Box Blue for 3-5ml tubes	each	17,20
AGG4516	81 Place Polycarbonate Tall Tube Box Green for 3-5ml tubes	each	17,00
AGG4517	81 Place Polycarbonate Tall Tube Box Orange for 3-5ml tubes	each	15,90
AGG4518	81 Place Polycarbonate Tall Tube Box White for 3-5ml tubes	each	15,50
AGG4519	81 Place Polycarbonate Tall Tube Box Yellow for 3-5ml tubes	each	15,90
AGG4520	100 Place Polypropylene Storage Box black	each	13,90
AGG4521	200 Place Polypropylene Storage Box black	each	59,10
AGG4522	0.5ml Blackline Conical Microtube	50	104,70
AGG4523	1.5ml Blackline Conical Tube black	50	31,30
AGG4524	15ml Blackline Conical Tube with screw cap	50	94,90
AGG4525	50ml Blackline Conical Tube with screw cap	50	135,40
AGG4526	1.5ml Conical Microtube amber	50	26,70
AGG4527	1.5ml Conical Microtube blue	50	26,70
AGG4528	1.5ml Conical Microtube green	50	27,50

Part #	Description	Unit	Price €
AGG4529	1.5ml Conical Microtube natural	50	26,70
AGG4530	1.5ml Conical Microtube orange	50	26,10
AGG4531	1.5ml Conical Microtube red	50	28,60
AGG4532	1.5ml Conical Microtube violet	50	26,70
AGG4533	1.5ml Conical Microtube yellow	50	26,70
AGG4536	5ml Prep Tube Conical blue with cap	250	76,60
AGG4537	5ml Prep Tube Conical black with cap	250	83,80
AGG4538	5ml Prep Tube Conical green with cap	250	90,80
AGG4539	5ml Prep Tube Conical red with cap	250	90,80
AGG4540	5ml Prep Tube Conical yellow with cap	250	76,90
AGG4541	5ml Prep Tube Conical natural with cap	250	83,30
AGG4542	5ml Prep Tube Conical natural without cap	250	86,90
AGG4543	5ml Prep Tube Free Standing natural with cap	250	85,10
AGG4544	5ml Prep Tube Free Standing natural without cap	25	98,30
AGG4545	1.5ml Tamper Evident Microtube - Free Standing blue	500	70,80
AGG4546	1.5ml Tamper Evident Microtube - Free Standing green	500	70,80
AGG4547	1.5ml Tamper Evident Microtube - Free Standing natural	500	70,80
AGG4548	1.5ml Tamper Evident Microtube - Free Standing red	500	70,80
AGG4549	1.5ml Tamper Evident Microtube - Free Standing yellow	500	70,80
AGG4550	1.5ml Tamper Evident Microtube - green	500	31,10
AGG4551	1.5ml Tamper Evident Microtube natural	500	69,10
AGG4552	1.5ml Tamper Evident Microtube yellow	500	31,10
AGG4553	1.5ml Tamper Evident Microtube red	500	31,10
AGG4555	5ml Polypropylene Tube with blue cap	250	128,30
AGG4556	5ml Polypropylene Tube with pink cap	250	81,70
AGG4557	5ml Polypropylene Tube with white cap	250	81,70
AGG4558	10ml Polypropylene Tube with blue cap	250	128,30
AGG4559	10ml Polypropylene Tube with pink cap	250	143,90
AGG4560	10ml Polypropylene Tube with white cap	250	143,90
AGG4580	Cryomold Standard 25 x 20 x 5mm	100	62,60
AGG4581	Cryomold Biopsy 10x10x5mm	100	62,60
AGG4582	Cryomold Intermediate 15x15x5mm	100	62,60
AGG4583	Round Cryomold Standard	100	113,80
AGG4584	Round Cryomold Biopsy	100	113,80
AGG5009	Dual-hook system for Cryo shipping cane	each	321,70
AGG5010	Intro to Cryo Consumable Kit	each	414,80
AGG5101-0-10	Gel Box Clear, 25.4 x 25.4 x 6.4mm , AD, Gel X0 (Pk10)		98,40
AGG5101-0-25	Gel-Pak Clear, 25.4 x 25.4 x 6.4mm, AD, X0 Pk25		143,50
AGG5101-4-10	Gel-Pak Box Clear, 25.4 x 25.4 x 6.4mm, AD, Gel X4 (Pk10)		98,40
AGG5101-4-25	Gel-Pak Clear, 25.4 x 25.4 x 6.4mm, AD, X4 Pk25		143,50
AGG5102-0-10	Gel-Pak Box Clear, 50.8 x 50.8 x 6.4mm, AD, Gel X0 (Pk10)		124,80
AGG5102-0-25	Gel-Pak Clear, 50.8 x 50.8 x 6.4mm, AD, X0 Pk25		182,70
AGG5102-4-10	Gel-Pak Box Clear, 50.8 x 50.8 x 6.4mm, AD, Gel X4 (Pk10)		124,80
AGG5102-4-25	Gel-Pak Clear, 50.8 x 50.8 x 6.4mm, AD, X4 Pk25		281,40
AGG5102-8-25	Gel-Pak Clear, 50.8 x 50.8 x 6.4mm, AD, X8 Pk25		328,30
AGG5103-0-20	Gel-Pak Clear, 50.8 x 76.2 x 9.5mm, AD, X0 Pk20		245,30
AGG5103-4-20	Gel-Pak Clear, 50.8 x 76.2 x 9.5mm, AD, X4 Pk20		245,30
AGG5103-8-20	Gel-Pak Clear, 50.8 x 76.2 x 9.5mm, AD, X8 Pk20		245,20
AGG5104-4-10	Gel-Pak Clear, 89 x 114 x 12.7mm, AD, X4 Pk10		224,90
AGG5105-4-25	Gel-Pak Black Conductive, 25.4 x 25.4 x 12.7mm, AD X4 Pk25		174,40
AGG5105-8-25	Gel-Pak Black Conductive, 25.4 x 25.4 x 12.7mm, AD X8 Pk25		191,70
AGG5106-4-20	Gel-Pak Black Conductive, 50.8 x 76.2 x 9.5mm, AD X4 Pk20		291,20
AGG5107-0-25	Gel-Pak Black Conductive, 50.8 x 50.8 x 6.4mm, AD X0 Pk25		356,70
AGG5107-4-25	Gel-Pak Black Conductive, 50.8 x 50.8 x 6.4mm, AD X4 Pk25		356,70
AGG5107-8-25	Gel-Pak Black Conductive, 50.8 x 50.8 x 6.4mm, AD X8 Pk25		393,50
AGG5108-0-10	Gel-Pak BlackCond/CAnti-Static 50.8x50.8x6.4mm X0 (Pk10)		178,50
AGG5108-0-25	Gel-Pak Black Anti-Static, 50.8 x 50.8 x 6.4mm, AD X0 Pk25		396,50
AGG5108-4-20	Gel-Pak Black Anti-Static, 50.8 x 76.2 x 9.5mm, AD X4 Pk20		359,40

Part #	Description	Unit	Price €
AGG5108-4-25	Gel-Pak Black Anti-Static, 50.8 x 50.8 x 6.4mm, AD X4 Pk25		396,50
AGG5108-8-25	Gel-Pak Black Anti-Static, 50.8 x 50.8 x 6.4mm, AD X8 Pk25		436,30
AGG5109-0-10	Gel-Pak Black Anti-static, 63.5x89x12.7mm AD Gel X0 (Pk10)		195,10
AGG5109-0-15	Gel-Pak Black Anti-Static, 63.5 x 89 x 12.7mm, AD X0 Pk15		237,50
AGG5109-4-10	Gel-Pak Black AntiStatic 63.5x89x12mm AD Gel X4 (Pk10)		195,10
AGG5109-4-15	Gel-Pak Black Anti-Static, 63.5 x 89 x 12.7mm, AD X4 Pk15		282,80
AGG5109-8-15	Gel-Pak Black Anti-Static, 63.5 x 89 x 12.7mm, AD X8 Pk15		329,80
AGG5110-0-25	Gel-Pak Clear Anti-Static, 50.8 x 50.8 x 6.4mm, AD X0 Pk25		396,50
AGG5110-4-25	Gel-Pak Clear Anti-Static, 50.8 x 50.8 x 6.4mm, AD X4 Pk25		396,50
AGG5112-0-25	Gel-Tray Black Conductive, 50.8 x 50.8 x 6.4mm, BD X0 Pk25		480,80
AGG5114-4-25	VR Conductive Gel-Pak Tray, 50.8 x 50.8mm X4 Pk25		560,50
AGG5115-4-25	VR Conductive Gel-Pak Tray, 101.6 x 101.6mm X4 Pk25		1.285,50
AGG5116	Clip for 4 x 4" VR Conductive Gel-Pak Tray Pk50		188,20
AGG5117	VP24 Vacuum Base, for 2" and 4" VR Conductive Gel-Pak VR		399,50
AGG5118	Manual vacuum hand pump with VP-24 base for 2" and 4"		521,70
AGG528-244	Metal Probe, 18 ga., bent 25.4mm long	each	7,10
AGG541	Gatan 3View system SEM pin stubs with standard , Al, Pk50		95,90
AGG542	Gatan 3View system SEM pin stubs with large,Al (Pk50)		93,60
AGG6010	Bio Diamond Paste series 100 1.103	10g	145,30
AGG6011	Bio Diamond Paste series 1001.105	10g	180,80
AGG6012	Bio Diamond Stick 1.122	10g	74,30
AGG6013	Bio Diamond Stick 1.123P	10g	81,10
AGG6014	Bio Diamond Paste series 100 2.103	10g	220,40
AGG6015	Bio Diamond Paste series 100 2.105	10g	263,20
AGG6016	Bio Diamond Paste series 100 3.102	10g	169,20
AGG6017	Bio Diamond Paste series 100 3.103	10g	204,30
AGG6018	Bio Diamond Paste series 100 3.105	10g	280,70
AGG6019	Bio Diamond Stick 3.122	10g	90,80
AGG6020	Bio Diamond Stick 3.123P	10g	93,80
AGG6021	Bio Diamond Paste series 100 6.101	10g	180,80
AGG6022	Bio Diamond Paste series 100 6.102	10g	204,30
AGG6023	Bio Diamond Paste series 100 6.103	10g	263,20
AGG6024	Bio Diamond Paste series 100 6.105	10g	333,60
AGG6025	Bio Diamond Stick 6.122	10g	91,40
AGG6026	Bio Diamond Stick 6.123P	10g	98,80
AGG6027	Bio Diamond Paste series 100 8.102	10g	280,70
AGG6028	Bio Diamond Paste series 100 8.103	10g	333,60
AGG6029	Bio Diamond Paste series 100 8.105	10g	389,20
AGG6030	Bio Diamond Stick 9.122	10g	92,50
AGG6031	Bio Diamond Stick 9.123P	10g	106,40
AGG6032	Bio Diamond Paste series 100 14.101	10g	204,30
AGG6033	Bio Diamond Paste series 100 14.102	10g	294,60
AGG6034	Bio Diamond Paste series 100 14.103	10g	389,20
AGG6035	Bio Diamond Paste series 100 14.105	10g	460,70
AGG6036	Bio Diamond Stick 15.122	10g	93,80
AGG6037	Bio Diamond Stick 15.123P	10g	107,90
AGG6038	Bio Diamond Paste series 100 025.102	10g	118,20
AGG6039	Bio Diamond Paste series 100 025.103	10g	177,10
AGG6040	Bio Diamond Stick 1/4.122	10g	62,70
AGG6041	Bio Diamond Stick 1/4.123P	10g	70,10
AGG6042	Bio Diamond Paste series 100 010.102	10g	181,50
AGG6043	Bio Diamond Paste series 100 010.103	10g	209,60
AGG6044	Liquid Bio Diamant NEODIA 1M Refill 1 L		174,20
AGG6045	Liquid Bio Diamant NEODIA 1M Refill 2,5L		391,90
AGG6046	Liquid Bio Diamant NEODIA 1M Refill 400ml		99,30
AGG6047	Liquid Bio Diamant NEODIA 1M 200ml		57,90
AGG6048	Liquid Bio Diamant NEODIA 1P Refill 1 L		260,70
AGG6049	Liquid Bio Diamant NEODIA 1P Refill 2,5L		544,50

Part #	Description	Unit	Price €
AGG6050	Liquid Bio Diamant NEODIA 1P Refill 400ml	400ml	145,70
AGG6051	Liquid Bio Diamant NEODIA 1P 200ml	200ml	84,90
AGG6052	Liquid Bio Diamant NEODIA 1F Refill 1L		159,70
AGG6053	Liquid Bio Diamant NEODIA 1F Refill 2,5L		335,10
AGG6054	Liquid Bio Diamant NEODIA 1F Refill 400ml		82,50
AGG6055	Liquid Bio Diamant NEODIA 2M Refill 1L		174,20
AGG6056	Liquid Bio Diamant NEODIA 2M Refill 2,5L		391,90
AGG6057	Liquid Bio Diamant NEODIA 2M Refill 400ml		99,30
AGG6058	Liquid Bio Diamant NEODIA 2M 200ml		57,90
AGG6059	Liquid Bio Diamant NEODIA 2P Refill 1L		260,70
AGG6060	Liquid Bio Diamant NEODIA 2P Refill 2,5L		544,50
AGG6061	Liquid Bio Diamant NEODIA 2P Refill 400ml		145,70
AGG6062	Liquid Bio Diamant NEODIA 2P 200ml		84,90
AGG6063	Liquid Bio Diamant NEODIA 3M Refill 1L		196,40
AGG6064	Liquid Bio Diamant NEODIA 3M Refill 2,5L		409,90
AGG6065	Liquid Bio Diamant NEODIA 3M Refill 400ml		112,10
AGG6066	Liquid Bio Diamant NEODIA 3M 200ml		65,10
AGG6067	Liquid Bio Diamant NEODIA 3P Refill 1L		295,00
AGG6068	Liquid Bio Diamant NEODIA 3P Refill 2.5L		616,10
AGG6069	Liquid Bio Diamant NEODIA 3P Refill 400ml	400ml	149,80
AGG6070	Liquid Bio Diamant NEODIA 3P 200ml	200ml	96,00
AGG6071	Liquid Bio Diamant NEODIA 3F Refill 1L		159,70
AGG6072	Liquid Bio Diamant NEODIA 3F Refill 2,5L		335,10
AGG6073	Liquid Bio Diamant NEODIA 3F Refill 400ml		82,50
AGG6074	Liquid Bio Diamant NEODIA 6M Refill 1L		213,50
AGG6075	Liquid Bio Diamant NEODIA 6M Refill 2,5L		446,00
AGG6076	Liquid Bio Diamant NEODIA 6M Refill 400ml		121,90
AGG6077	Liquid Bio Diamant NEODIA 6M 200ml		71,00
AGG6078	Liquid Bio Diamant NEODIA 6P Refill 1L		308,70
AGG6079	Liquid Bio Diamant NEODIA 6P Refill 2,5L		644,80
AGG6080	Liquid Bio Diamant NEODIA 6P Refill 400ml		156,80
AGG6081	Liquid Bio Diamant NEODIA 6P 200ml		94,10
AGG6082	Liquid Bio Diamant NEODIA 6F Refill 1L		159,70
AGG6083	Liquid Bio Diamant NEODIA 6F Refill 2,5L		335,10
AGG6084	Liquid Bio Diamant NEODIA 6F Refill 400ml		82,50
AGG6085	Liquid Bio Diamant NEODIA 9M Refill 1L		213,50
AGG6086	Liquid Bio Diamant NEODIA 9M Refill 2,5L		446,00
AGG6087	Liquid Bio Diamant NEODIA 9M Refill 400ml		121,90
AGG6088	Liquid Bio Diamant NEODIA 9M 200ml		71,00
AGG6089	Liquid Bio Diamant NEODIA 9P Refill 1L		308,70
AGG6090	Liquid Bio Diamant NEODIA 9P Refill 2,5L		644,80
AGG6091	Liquid Bio Diamant NEODIA 9P Refill 400ml	400ml	156,80
AGG6092	Liquid Bio Diamant NEODIA 9P 200ml	200ml	94,10
AGG6093	Liquid Bio Diamant NEODIA 9F Refill 1L		172,90
AGG6094	Liquid Bio Diamant NEODIA 9F Refill 2,5L		363,10
AGG6095	Liquid Bio Diamant NEODIA 9F Refill 400ml		89,40
AGG6096	Liquid Bio Diamant NEODIA 14M Refill1L		213,50
AGG6097	Liquid Bio Diamant NEODIA 14M Refill 2,5L		446,00
AGG6098	Liquid Bio Diamant NEODIA 14M Refill 400ml		121,90
AGG6099	Liquid Bio Diamant NEODIA 14M 200ml		71,00
AGG6100	Liquid Bio Diamant NEODIA 14P Refill 1L		308,70
AGG6101	Liquid Bio Diamant NEODIA 14P Refill 2,5L		644,80
AGG6102	Liquid Bio Diamant NEODIA 14P Refill 400ml		156,80
AGG6103	Liquid Bio Diamant NEODIA 14P 200ml		94,10
AGG6104	Liquid Bio Diamant NEODIA 14F Refill 1L		172,90
AGG6105	Liquid Bio Diamant NEODIA 14F Refill 2,5L		363,10
AGG6106	Liquid Bio Diamant NEODIA 14F Refill 400ml		89,40
AGG6107	Liquid Bio Diamant NEODIA 1/4M Refill 1L		192,10

Part #	Description	Unit	Price €
AGG6108	Liquid Bio Diamant NEODIA 1/4M Refill 2,5L		401,30
AGG6109	Liquid Bio Diamant NEODIA 1/4M Refill 400ml		109,70
AGG6110	Liquid Bio Diamant NEODIA 1/4M 200ml		63,80
AGG6111	Liquid Bio Diamant NEODIA 1/4P Refill 1L		288,10
AGG6112	Liquid Bio Diamant NEODIA 1/4P Refill 2,5L		601,70
AGG6113	Liquid Bio Diamant NEODIA 1/4P Refill 400ml		160,80
AGG6114	Liquid Bio Diamant NEODIA 1/4P 200ml		93,50
AGG6115	HAZ Pulmatic Bio Diamant 1.313	75g	67,20
AGG6116	HAZ Pulmatic Bio Diamant 1.323P	75g	76,80
AGG6117	HAZ Pulmatic Bio Diamant 2.313	75g	67,20
AGG6118	HAZ Pulmatic Bio Diamant 2.323P	75g	76,80
AGG6119	HAZ Pulmatic Bio Diamant 3.313	75g	67,20
AGG6120	HAZ Pulmatic Bio Diamant 3.323P	75g	76,80
AGG6121	HAZ Pulmatic Bio Diamant 6.313	75g	87,90
AGG6122	HAZ Pulmatic Bio Diamant 6.323P	75g	94,80
AGG6123	HAZ Pulmatic Bio Diamant 8.313	75g	87,90
AGG6124	HAZ Pulmatic Bio Diamant 8.323P	75g	94,80
AGG6125	HAZ Pulmatic Bio Diamant 10.313	75g	87,90
AGG6126	HAZ Pulmatic Bio Diamant 10.323P	75g	94,80
AGG6127	HAZ Pulmatic Bio Diamant 14.313	75g	94,00
AGG6128	HAZ Pulmatic Bio Diamant 14.323P	75g	107,90
AGG6129	HAZ Pulmatic Bio Diamant 25.313	75g	94,00
AGG6130	HAZ Pulmatic Bio Diamant 25.323P	75g	107,90
AGG6131	HAZ Pulmatic Bio Diamant 40.313	75g	94,00
AGG6132	HAZ Pulmatic Bio Diamant 40.323P	75g	107,90
AGG6133	HAZ Pulmatic Bio Diamant 025.313	75g	67,20
AGG6134	HAZ Pulmatic Bio Diamant 025.323P	75g	76,80
AGG6135	Abrasive Paper Excellence P80 200mm dia NA Pk50	50	114,80
AGG6136	Abrasive Paper Excellence P80 230mm dia NA	50	126,10
AGG6137	Abrasive Paper Excellence P80 250mm dia NA	50	142,20
AGG6138	Abrasive Paper Excellence P80 300mm dia NA	50	162,90
AGG6139	Abrasive Paper Excellence P120 200mm dia NA	50	114,80
AGG6140	Abrasive Paper Excellence P120 230mm dia NA	50	126,10
AGG6141	Abrasive Paper Excellence P120 250mm dia NA	50	142,20
AGG6142	Abrasive Paper Excellence P120 300mm dia NA	50	162,90
AGG6143	Abrasive Paper Excellence P180 200mm dia NA	50	114,80
AGG6144	Abrasive Paper Excellence P180 230mm dia NA	50	126,10
AGG6145	Abrasive Paper Excellence P180 250mm dia NA	50	142,20
AGG6146	Abrasive Paper Excellence P180 300mm dia NA	50	162,90
AGG6147	Abrasive Paper Excellence P240 200mm dia NA	50	114,80
AGG6148	Abrasive Paper Excellence P240 230mm dia NA	50	126,10
AGG6149	Abrasive Paper Excellence P240 250mm dia NA	50	142,20
AGG6150	Abrasive Paper Excellence P240 300mm dia NA	50	162,90
AGG6151	Abrasive Paper Excellence P320 200mm dia NA	50	114,80
AGG6152	Abrasive Paper Excellence P320 230mm dia NA	50	126,10
AGG6153	Abrasive Paper Excellence P320 250mm dia NA	50	142,20
AGG6154	Abrasive Paper Excellence P320 300mm dia NA	50	162,90
AGG6155	Abrasive Paper Excellence P400 200mm dia NA	50	114,80
AGG6156	Abrasive Paper Excellence P400 230mm dia NA	50	126,10
AGG6157	Abrasive Paper Excellence P400 250mm dia NA	50	142,20
AGG6158	Abrasive Paper Excellence P400 300mm dia NA	50	162,90
AGG6159	Abrasive Paper Excellence P600 200mm dia NA	50	114,80
AGG6160	Abrasive Paper Excellence P600 230mm dia NA	50	126,10
AGG6161	Abrasive Paper Excellence P600 250mm dia NA	50	142,20
AGG6162	Abrasive Paper Excellence P600 300mm dia NA	50	162,90
AGG6163	Abrasive Paper Excellence P1200 200mm dia NA	50	114,80
AGG6164	Abrasive Paper Excellence P1200 230mm dia NA	50	126,10
AGG6165	Abrasive Paper Excellence P1200 250mm dia NA	50	142,20

Part #	Description	Unit	Price €
AGG6166	Abrasivo Paper Excellence P1200 300mm dia NA	50	162,90
AGG6167	Abrasivo Paper Excellence P80 200mm dia SA	50	259,30
AGG6168	Abrasivo Paper Excellence P80 230mm dia SA	50	303,60
AGG6169	Abrasivo Paper Excellence P80 250mm dia SA	50	329,40
AGG6170	Abrasivo Paper Excellence P80 300mm dia SA	50	386,30
AGG6171	Abrasivo Paper Excellence P120 200mm dia SA	50	259,30
AGG6172	Abrasivo Paper Excellence P120 230mm dia SA	50	303,60
AGG6173	Abrasivo Paper Excellence P120 250mm dia SA	50	329,40
AGG6174	Abrasivo Paper Excellence P120 300mm dia SA	50	386,30
AGG6175	Abrasivo Paper Excellence P180 200mm dia SA	50	259,30
AGG6176	Abrasivo Paper Excellence P180 230mm dia SA	50	303,60
AGG6177	Abrasivo Paper Excellence P180 250mm dia SA	50	329,40
AGG6178	Abrasivo Paper Excellence P180 300mm dia SA	50	386,30
AGG6179	Abrasivo Paper Excellence P240 200mm dia SA	50	259,30
AGG6180	Abrasivo Paper Excellence P240 230mm dia SA	50	303,60
AGG6181	Abrasivo Paper Excellence P240 250mm dia SA	50	329,40
AGG6182	Abrasivo Paper Excellence P240 300mm dia SA	50	386,30
AGG6183	Abrasivo Paper Excellence P320 200mm dia SA	50	259,30
AGG6184	Abrasivo Paper Excellence P320 230mm dia SA	50	303,60
AGG6185	Abrasivo Paper Excellence P320 250mm dia SA	50	329,40
AGG6186	Abrasivo Paper Excellence P320 300mm dia SA	50	386,30
AGG6187	Abrasivo Paper Excellence P400 200mm dia SA	50	259,30
AGG6188	Abrasivo Paper Excellence P400 230mm dia SA	50	303,60
AGG6189	Abrasivo Paper Excellence P400 250mm dia SA	50	329,40
AGG6190	Abrasivo Paper Excellence P400 300mm dia SA	50	386,30
AGG6191	Abrasivo Paper Excellence P600 200mm dia SA	50	259,30
AGG6192	Abrasivo Paper Excellence P600 230mm dia SA	50	303,60
AGG6193	Abrasivo Paper Excellence P600 250mm dia SA	50	329,40
AGG6194	Abrasivo Paper Excellence P600 300mm dia SA	50	386,30
AGG6195	Abrasivo Paper Excellence P1200 200mm dia SA	50	259,30
AGG6196	Abrasivo Paper Excellence P1200 230mm dia SA	50	303,60
AGG6197	Abrasivo Paper Excellence P1200 250mm dia SA	50	329,40
AGG6198	Abrasivo Paper Excellence P1200 300mm dia SA	50	386,30
AGG6199	FIXLAM Disc 200mm dia SA	5	55,10
AGG6200	FIXLAM Disc 230mm dia SA	5	70,80
AGG6201	FIXLAM Disc 250mm dia SA	5	81,30
AGG6202	FIXLAM Disc 300mm dia SA.	5	103,20
AGG6203	FIXLAM Disc 200mm dia MG	5	169,40
AGG6204	FIXLAM Disc 230mm dia MG	5	178,90
AGG6205	FIXLAM Disc 250mm dia MG	5	197,50
AGG6206	FIXLAM Disc 300mm dia MG	5	213,20
AGG6207	Abrasivo Paper P80 200mm dia NA	100	125,20
AGG6208	Abrasivo Paper P80 230mm dia NA	100	147,30
AGG6209	Abrasivo Paper P80 250mm dia NA	100	176,20
AGG6210	Abrasivo Paper P80 300mm dia NA	100	226,90
AGG6211	Abrasivo Paper P80 305mm dia NA	100	238,20
AGG6212	Abrasivo Paper P80 400mm dia dia NA	100	449,20
AGG6213	Abrasivo Paper P120 200mm dia NA	100	125,20
AGG6214	Abrasivo Paper P120 230mm dia NA	100	147,30
AGG6215	Abrasivo Paper P120 250mm dia NA	100	176,20
AGG6216	Abrasivo Paper P120 300mm dia NA	100	226,90
AGG6217	Abrasivo Paper P120 305mm dia NA	100	238,20
AGG6218	Abrasivo Paper P120 400mm dia dia NA	100	449,20
AGG6219	Abrasivo Paper P180 200mm dia NA	100	125,20
AGG6220	Abrasivo Paper P180 230mm dia NA	100	147,30
AGG6221	Abrasivo Paper P180 250mm dia NA	100	176,20
AGG6222	Abrasivo Paper P180 300mm dia NA	100	226,90
AGG6223	Abrasivo Paper P180 305mm dia NA	100	238,20

Part #	Description	Unit	Price €
AGG6224	Abrasivo Paper P180 400mm dia dia NA	100	449,20
AGG6225	Abrasivo Paper P240 200mm dia NA	100	101,90
AGG6226	Abrasivo Paper P240 230mm dia NA	100	147,30
AGG6227	Abrasivo Paper P240 250mm dia NA	100	176,20
AGG6228	Abrasivo Paper P240 300mm dia NA	100	226,90
AGG6229	Abrasivo Paper P240 305mm dia NA	100	238,20
AGG6230	Abrasivo Paper P240 400mm dia dia NA	100	449,20
AGG6231	Abrasivo Paper P320 200mm dia NA	100	125,20
AGG6232	Abrasivo Paper P320 230mm dia NA	100	147,30
AGG6233	Abrasivo Paper P320 250mm dia NA	100	176,20
AGG6234	Abrasivo Paper P320 300mm dia NA	100	226,90
AGG6235	Abrasivo Paper P320 305mm dia NA	100	238,20
AGG6236	Abrasivo Paper P320 400mm dia dia NA	100	449,20
AGG6237	Abrasivo Paper P400 200mm dia NA	100	101,90
AGG6238	Abrasivo Paper P400 230mm dia NA	100	147,30
AGG6239	Abrasivo Paper P400 250mm dia NA	100	176,20
AGG6240	Abrasivo Paper P400 300mm dia NA	100	226,90
AGG6241	Abrasivo Paper P400 305mm dia NA	100	238,20
AGG6242	Abrasivo Paper P400 400mm dia dia NA Pk100	100	449,20
AGG6243	Abrasivo Paper P600 200mm dia NA	100	101,90
AGG6244	Abrasivo Paper P600 230mm dia NA	100	147,30
AGG6245	Abrasivo Paper P600 250mm dia NA	100	176,20
AGG6246	Abrasivo Paper P600 300mm dia NA	100	226,90
AGG6247	Abrasivo Paper P600 305mm dia NA	100	238,20
AGG6248	Abrasivo Paper P600 400mm dia dia NA	100	449,20
AGG6249	Abrasivo Paper P800 200mm dia NA	100	101,90
AGG6250	Abrasivo Paper P800 230mm dia NA	100	147,30
AGG6251	Abrasivo Paper P800 250mm dia NA	100	176,20
AGG6252	Abrasivo Paper P800 300mm dia NA	100	226,90
AGG6253	Abrasivo Paper P800 305mm dia NA	100	238,20
AGG6254	Abrasivo Paper P800 400mm dia dia NA	100	449,20
AGG6255	Abrasivo Paper P1000 200mm dia NA	100	125,20
AGG6256	Abrasivo Paper P1000 230mm dia NA	100	147,30
AGG6257	Abrasivo Paper P1000 250mm dia NA	100	176,20
AGG6258	Abrasivo Paper P1000 300mm dia NA	100	226,90
AGG6259	Abrasivo Paper P1000 305mm dia NA	100	238,20
AGG6260	Abrasivo Paper P1000 400mm dia dia NA	100	449,20
AGG6261	Abrasivo Paper P1200 200mm dia NA	100	101,90
AGG6262	Abrasivo Paper P1200 230mm dia NA	100	147,30
AGG6263	Abrasivo Paper P1200 250mm dia NA	100	176,20
AGG6264	Abrasivo Paper P1200 300mm dia NA	100	226,90
AGG6265	Abrasivo Paper P1200 305mm dia NA	100	238,20
AGG6266	Abrasivo Paper P1200 400mm dia dia NA	100	449,20
AGG6267	Abrasivo Paper P2400 200mm dia NA	100	313,10
AGG6268	Abrasivo Paper P2400 230mm dia NA	100	422,50
AGG6269	Abrasivo Paper P2400 250mm dia NA	100	556,70
AGG6270	Abrasivo Paper P2400 300mm dia NA	100	794,60
AGG6271	Abrasivo Paper P2400 305mm dia NA	100	794,60
AGG6272	Abrasivo Paper P4000 200mm dia NA	100	375,90
AGG6273	Abrasivo Paper P4000 230mm dia NA	100	415,40
AGG6274	Abrasivo Paper P4000 250mm dia NA	100	556,70
AGG6275	Abrasivo Paper P4000 300mm dia NA	100	794,60
AGG6276	Abrasivo Paper P4000 305mm dia NA	100	794,60
AGG6277	Abrasivo Paper P80 200mm dia SA	100	190,20
AGG6278	Abrasivo Paper P80 230mm dia SA	100	226,50
AGG6279	Abrasivo Paper P80 250mm dia SA	100	271,30
AGG6280	Abrasivo Paper P80 300mm dia SA	100	315,20
AGG6281	Abrasivo Paper P80 305mm dia SA	100	375,00

Part #	Description	Unit	Price €
AGG6282	Abrasivo Paper P80 400mm dia dia SA	100	593,00
AGG6283	Abrasivo Paper P120 200mm dia SA	100	176,20
AGG6284	Abrasivo Paper P120 230mm dia SA	100	226,50
AGG6285	Abrasivo Paper P120 250mm dia SA	100	248,00
AGG6286	Abrasivo Paper P120 300mm dia SA	100	315,20
AGG6287	Abrasivo Paper P120 305mm dia SA	100	328,30
AGG6288	Abrasivo Paper P120 400mm dia dia SA	100	593,00
AGG6289	Abrasivo Paper P180 200mm dia SA	100	176,20
AGG6290	Abrasivo Paper P180 230mm dia SA	100	226,50
AGG6291	Abrasivo Paper P180 250mm dia SA	100	228,60
AGG6292	Abrasivo Paper P180 300mm dia SA	100	315,20
AGG6293	Abrasivo Paper P180 305mm dia SA	100	328,30
AGG6294	Abrasivo Paper P180 400mm dia dia SA	100	593,00
AGG6295	Abrasivo Paper P240 200mm dia SA	100	176,20
AGG6296	Abrasivo Paper P240 230mm dia SA	100	226,50
AGG6297	Abrasivo Paper P240 250mm dia SA	100	228,60
AGG6298	Abrasivo Paper P240 300mm dia SA	100	315,20
AGG6299	Abrasivo Paper P240 305mm dia SA	100	328,30
AGG6300	Abrasivo Paper P240 400mm dia dia SA	100	593,00
AGG6301	Abrasivo Paper P320 200mm dia SA	100	176,20
AGG6302	Abrasivo Paper P320 230mm dia SA	100	226,50
AGG6303	Abrasivo Paper P320 250mm dia SA	100	228,60
AGG6304	Abrasivo Paper P320 300mm dia SA	100	315,20
AGG6305	Abrasivo Paper P320 305mm dia SA	100	328,30
AGG6306	Abrasivo Paper P320 400mm dia dia SA	100	593,00
AGG6307	Abrasivo Paper P400 200mm dia SA	100	176,20
AGG6308	Abrasivo Paper P400 230mm dia SA	100	226,50
AGG6309	Abrasivo Paper P400 250mm dia SA	100	228,60
AGG6310	Abrasivo Paper P400 300mm dia SA	100	315,20
AGG6311	Abrasivo Paper P400 305mm dia SA	100	328,30
AGG6312	Abrasivo Paper P400 400mm dia dia SA	100	593,00
AGG6313	Abrasivo Paper P600 200mm dia SA	100	176,20
AGG6314	Abrasivo Paper P600 230mm dia SA	100	226,50
AGG6315	Abrasivo Paper P600 250mm dia SA	100	228,60
AGG6316	Abrasivo Paper P600 300mm dia SA	100	315,20
AGG6317	Abrasivo Paper P600 305mm dia SA	100	328,30
AGG6318	Abrasivo Paper P600 400mm dia dia SA	100	593,00
AGG6319	Abrasivo Paper P800 200mm dia SA	100	176,20
AGG6320	Abrasivo Paper P800 230mm dia SA	100	226,50
AGG6321	Abrasivo Paper P800 250mm dia SA	100	228,60
AGG6322	Abrasivo Paper P800 300mm dia SA	100	315,20
AGG6323	Abrasivo Paper P800 305mm dia SA	100	328,30
AGG6324	Abrasivo Paper P800 400mm dia dia SA	100	593,00
AGG6325	Abrasivo Paper P1000 200mm dia SA	100	176,20
AGG6326	Abrasivo Paper P1000 230mm dia SA	100	226,50
AGG6327	Abrasivo Paper P1000 250mm dia SA	100	228,60
AGG6328	Abrasivo Paper P1000 300mm dia SA	100	315,20
AGG6329	Abrasivo Paper P1000 305mm dia SA	100	328,30
AGG6330	Abrasivo Paper P1000 400mm dia dia SA	100	593,00
AGG6331	Abrasivo Paper P1200 200mm dia SA	100	176,20
AGG6332	Abrasivo Paper P1200 230mm dia SA	100	226,50
AGG6333	Abrasivo Paper P1200 250mm dia SA	100	228,60
AGG6334	Abrasivo Paper P1200 300mm dia SA	100	315,20
AGG6335	Abrasivo Paper P1200 305mm dia SA	100	328,30
AGG6336	Abrasivo Paper P1200 400mm dia dia SA	100	593,00
AGG6337	Abrasivo Paper P2400 200mm dia SA	100	397,30
AGG6338	Abrasivo Paper P2400 230mm dia SA	100	591,60
AGG6339	Abrasivo Paper P2400 250mm dia SA	100	610,20

Part #	Description	Unit	Price €
AGG6340	Abrasive Paper P2400 300mm dia SA	100	918,10
AGG6341	Abrasive Paper P4000 200mm dia SA	100	466,40
AGG6342	Abrasive Paper P4000 230mm dia SA	100	591,60
AGG6343	Abrasive Paper P4000 250mm dia SA	100	610,20
AGG6344	Abrasive Paper P4000 300mm dia SA	100	918,10
AGG6345	Cut-Off Wheel Blue 150x12,7x0,45mm	5	213,10
AGG6346	Cut-Off Wheel Blue 203x25,4x1mm	10	290,20
AGG6347	Cut-Off Wheel Blue 230x22x1,2mm	10	327,30
AGG6348	Cut-Off Wheel Blue 250x32x1,5mm	10	360,30
AGG6349	Cut-Off Wheel Blue 305x32x1,5mm	10	403,90
AGG6350	Cut-Off Wheel Blue 125x12,7x0,45mm	5	161,10
AGG6351	Cut-Off Wheel Red 150x12,7x0,45mm	5	213,10
AGG6352	Cut-Off Wheel Red 203x25,4x1mm	10	290,20
AGG6353	Cut-Off Wheel Red 230x22x1,2mm	10	327,30
AGG6354	Cut-Off Wheel Red 250x32x1,5mm	10	360,30
AGG6355	Cut-Off Wheel Red 305x32x1,5mm	10	403,90
AGG6356	Cut-Off Wheel Red 125x12,7x0,45mm	5	161,10
AGG6357	Cut-Off Wheel White 150x12,7x0,45mm	5	236,40
AGG6358	Cut-Off Wheel White 203x25,4x1mm	10	290,20
AGG6359	Cut-Off Wheel White 230x22x1,2mm	10	327,30
AGG6360	Cut-Off Wheel White 250x32x1,5mm	10	360,30
AGG6361	Cut-Off Wheel White 305x32x1,5mm	10	403,90
AGG6362	Cut-Off Wheel White 125x12,7x0,45mm	each	178,40
AGG6363	Precision Cut-off Wheel Black 102X12,7x0,3mm	each	505,90
AGG6364	Cutting Wheel CBN Black 152x12,7x0,38mm	each	654,60
AGG6365	Precision Cut-off Wheel Black 127X12,7x0,4mm	each	682,80
AGG6366	Precision Cut-off Wheel Black 152X12,7x0,5mm	each	680,30
AGG6367	Precision Cut-off Wheel Black 203X12,7x0,8mm	each	732,10
AGG6368	Precision Cut-Off Wheel Black 254X32x1mm	each	863,30
AGG6369	Precision Cut-Off Wheel Black 305X32x1,3mm	each	1.065,60
AGG6370	Cut-Off Wheel Green H1 250x32x1,5mm Pk10	10	271,80
AGG6371	Cut-Off Wheel Green H1 305x32x2mm Pk10	10	305,10
AGG6372	Cut-Off Wheel Green H1 350x32x2,5mm Pk10	10	411,60
AGG6373	Cut-Off Wheel Green H1 400x3,0x32mm Pk10	10	515,10
AGG6374	Cut-Off Wheel Excellence Green H1 250x32x1,5mm Pk10	10	304,10
AGG6375	Cut-Off Wheel Excellence Green H1 305x32x2mm Pk10	10	394,50
AGG6376	Cut-Off Wheel Excellence Green H1 350x32x2,5mm Pk10	10	603,40
AGG6377	Cut-Off Wheel Excellence Green H1 400x3,0x32mm Pk10	10	709,80
AGG6378	Cut-Off Wheel Yellow H2 250x32x1,5mm Pk10	10	271,80
AGG6379	Cut-Off Wheel Yellow H2 305x2x32mm Pk10	10	305,10
AGG6380	Cut-Off Wheel Yellow H2 350x32x2,5mm Pk10	10	411,60
AGG6381	Cut-Off Wheel Yellow H2 400x3,0x32mm Pk10	10	515,10
AGG6382	Cut-Off Wheel Excellence Yellow H2 250x32x1,5mm Pk10	10	304,10
AGG6383	Cut-Off Wheel Excellence Yellow H2 305x2x32mm Pk10	10	394,50
AGG6384	Cut-Off Wheel Excellence Yellow H2 350x32x2,5mm Pk10	10	603,40
AGG6385	Cut-Off Wheel Excellence Yellow H2 400x32x3mm Pk10	10	709,80
AGG6386	Cut-Off Wheel White H4 250x32x1,5mm Pk10	10	295,50
AGG6387	Cut-Off Wheel White H4 305x32x2mm Pk10	10	334,60
AGG6388	Cut-Off Wheel White H4 350x32x2,5mm Pk10	10	417,10
AGG6389	Cut-Off Wheel White H4 400x32x3mm Pk10	10	598,20
AGG6390	Cut-Off Wheel Excellence White H4 250x32x1,5mm Pk10	10	312,80
AGG6391	Cut-Off Wheel Excellence White H4 305x32x2mm Pk10	10	413,60
AGG6392	Cut-Off Wheel Excellence White H4 350x32x2,5mm Pk10	10	629,30
AGG6393	Cut-Off Wheel Excellence White H4 400x32x3mm Pk10	10	790,50
AGG6394	Liquid ALPLAN 0,05µ 1L	1l	o.a.
AGG6395	Liquid ALPLAN 0,25V 1L	1l	o.a.
AGG6396	Liquid ALPLAN 1µ 1L	1l	o.a.
AGG6397	Liquid ALPLAN 3µ 1L	1l	o.a.

Part #	Description	Unit	Price €
AGG6399	FAS Disc 230mm dia SA	each	148,50
AGG6403	FAS Disc 230mm dia MG	each	178,90
AGG6406	Grinding Wheel Gr60 355mm dia	each	745,80
AGG6407	Grinding Wheel Gr150 355mm dia	each	756,90
AGG6408	Final Liquid 500g	500g	77,90
AGG6409	Final Liquid 5L	5l	251,60
AGG6410	Liquid L1 1L	1l	53,60
AGG6411	Liquid L2 1L	1l	53,60
AGG6412	Resin 601/2 1L + Varnish 100ml	each	349,20
AGG6413	Green Phenol Powder 602 10Kg	10kg	215,90
AGG6414	Green Phenol Powder 602 2,5Kg	2,5kg	76,50
AGG6415	Cold Resin Epoxy 603 Kit 1,5Kg	1,5kg	194,70
AGG6416	Hot Mounting Resin Copper 2Kg	2kg	164,50
AGG6417	Resin 605 1,5Kg	1,5kg	325,50
AGG6418	Resin 605 15Kg	1,5kg	2.330,20
AGG6419	Resin 607 Kit 1,5Kg	1,5kg	218,90
AGG6420	Catalysor 607 500ml	500ml	125,00
AGG6421	Hardener 607 250ml	250ml	72,20
AGG6422	powder 607 750g	750g	127,10
AGG6423	Resin 608 1,5Kg	1,5kg	295,50
AGG6424	Resin 608 15Kg	1,5kg	1.534,20
AGG6425	Resin 609 1,5Kg	1,5kg	316,00
AGG6426	Resin 609 15Kg	1,5kg	2.226,60
AGG6427	Resin Liquid 605 500ml	500ml	150,00
AGG6428	Resin Liquid 605 5L	5l	870,90
AGG6429	Hot Mounting Resin Transparent 616 10Kg	10kg	571,60
AGG6430	Hot Mounting Resin Transparent 616 2,5Kg	2,5kg	186,90
AGG6431	Black Conductive Phenol Powder 617 1,5Kg	1,5kg	246,20
AGG6432	Resin Liquid 609 500ml	500ml	150,00
AGG6433	Resin Liquid 609 5L	5l	812,70
AGG6434	Red Phenol Powder 622 10Kg	10kg	215,90
AGG6435	Red Phenol Powder 622 2,5Kg	2,5kg	76,50
AGG6436	Flexible Mould 625 25mm dia Pk5	5	100,90
AGG6437	Flexible Mould 630 30mm dia Pk5	5	100,90
AGG6438	Flexible Mould 632 32mm dia Pk5	5	100,90
AGG6439	Black Phenol Powder 632 Loaded with Glass Fibres 10Kg	10kg	428,50
AGG6440	Black Phenol Powder 632 loaded with Glass Fibres 2,5Kg	2,5kg	147,40
AGG6441	Hot Mounting Resin Black Epoxy 10Kg	10kg	612,50
AGG6442	Hot Mounting Resin Black Epoxy 2,5Kg	2,5kg	192,40
AGG6443	Hot Mounting Resin Black Epoxy 634 2,5Kg	2,5kg	248,90
AGG6444	Flexible Mould 640 40mm dia Pk5	5	100,90
AGG6445	Black Phenol Powder 642 10Kg		228,40
AGG6446	Black Phenol Powder 642 2,5Kg		85,30
AGG6447	Flexible Mould 650 50mm dia Pk5		152,80
AGG6448	Flexible Mould 660 60mm dia Pk5		173,20
AGG6449	Kit of 4 Flexible Mould		401,40
AGG6450	Flexible Mould 670 70x30x40mm dia Pk5		396,10
AGG6451	Flexible Mould 670 100x42x30mm dia Pk5		453,30
AGG6452	Flexible Mould 670 120x50x30mm dia Pk5		511,60
AGG6453	Flexible Mould 670 100x70x30mm dia Pk5		521,00
AGG6454	Kit of 3 Flexible Mould		396,10
AGG6455	Mould with Removable Bottom 681 25mm dia Pk5		166,30
AGG6456	Mould with Removable Bottom 681 30mm dia Pk5		182,00
AGG6457	Mould with Removable Bottom 681 32mm dia Pk5		189,30
AGG6458	Mould with Removable Bottom 681 D40mm dia Pk5		228,90
AGG6459	Powder Mould Removal Agent SOLILUB 50g		68,60
AGG6460	Blade Holder 1mm dia Unit 100		183,80
AGG6461	Blade Holder 2mm dia Unit 100		183,80

Part #	Description	Unit	Price €
AGG6462	Blade Holder 3mm dia Unit 100		183,80
AGG6463	Clips 5 Unit 100		163,10
AGG6464	Clips 3 Unit 100		174,40
AGG6465	Cardboard Cups 200ml 200pces Pk200		157,60
AGG6466	Kit 400 Cups + 400 Spatulas		245,50
AGG6467	Resin Liquid 609 1,5Kg Blue		347,50
AGG6468	Resin Liquid 609 1,5Kg Yellow		347,50
AGG6469	Resin Liquid 609 1,5Kg Orange		347,50
AGG6470	Resin Liquid 609 1,5Kg Red		347,50
AGG6471	Resin Liquid 609 1,5Kg Green		347,50
AGG6472	SPHERE VOLUMA Kit 3x250cm3		176,80
AGG6473	SPHERE VOLUMA Red 250cm3		86,80
AGG6474	SPHERE VOLUMA Blue 250cm3		86,80
AGG6475	SPHERE VOLUMA Ivory 250cm3		86,80
AGG6476	Resin 601/2 1L		301,50
AGG6477	Resin 601/2 100ml		92,70
AGG6478	Cold Resin 603 Liquid 1L		152,60
AGG6479	Resin Liquid 608 500ml		151,00
AGG6480	Resin Liquid 608 5L		742,20
AGG6481	Cold Resin 603 Hardener 500 mL		81,10
AGG6482	Powder for Resin 608 1kg		165,20
AGG6483	Powder for Resin 608 10kg		1.008,20
AGG6484	Hot Mounting Resin Transparent 616-2 10Kg		715,10
AGG6485	Hot Mounting Resin Transparent 616-2 2,5Kg		236,30
AGG6486	Grey Composite Powder Phenofree 10Kg		215,90
AGG6487	Grey Composite Powder Phenofree 2,5Kg		76,50
AGG6488	Grey Composite Powder Phenofree 25Kg		474,50
AGG6489	Red Composite Powder Phenofree 10Kg		215,90
AGG6490	Red Composite Powder Phenofree 2,5Kg		76,50
AGG6491	Red Composite Powder Phenofree 25Kg		474,50
AGG6492	White Composite Powder Phenofree 10Kg		215,90
AGG6493	White Composite Powder Phenofree 2,5Kg		76,50
AGG6494	White Composite Powder Phenofree 25Kg		474,50
AGG6495	Lubricant Fluid 704 1L		30,70
AGG6496	Lubricant Fluid 704 5x1L Pk5		125,00
AGG6497	Cutting Fluid 722 5L		201,20
AGG6498	Lubricant Fluid Booster Premium Type 3 Can 5L		97,10
AGG6499	Lubricant Fluid mm dia702 1L		15,00
AGG6500	Lubricant Fluid mm dia702 5L		70,10
AGG6501	Bottle 500ml + Red Sprayer		30,90
AGG6502	Green Sprayer		8,70
AGG6503	Bottle 400ml + Green Sprayer		14,00
AGG6504	Bottle 200ml + White Vaporizer		14,00
AGG6505	Pressurized Mounting Device mm dia 806		1.402,70
AGG6506	Polymerization Lamp mm dia 866		1.318,60
AGG6507	Magnetic Plate 200mm dia		367,80
AGG6509	Magnetic Plate 230mm dia		372,30
AGG6510	FMS Fixing System 230mm dia		147,30
AGG6511	Magnetic Plate 250mm dia		512,70
AGG6513	Magnetic Plate 300mm dia		539,70
AGG6515	FMS Magnetic Plate 400mm dia dia		823,40
AGG6518	Foam Template for Storage of Samples Pk2		55,30
AGG6519	Self Adhesive PlexiGlass Shelf Pk2		89,10
AGG6520	PVC Protection Sheet Pk10		69,80
AGG6521	Cover LAM 50mm dia (Unit 100)		139,70
AGG6522	LAM 15 320x400mm dia 2x24 tissues		100,30
AGG6523	LAM 15 80mm dia 400		100,30
AGG6524	Nital 5% Etchant 1L		129,30

Part #	Description	Unit	Price €
AGG6525	Picral 2-1 Etchant 1L		179,50
AGG6526	Murakami Etchant 1L		209,00
AGG6527	Kalling Etchant 1L		153,90
AGG6528	Oberhoffer Etchant 1L		139,30
AGG6529	CAMEO DISK Silver 200mm dia SA Pk2		137,40
AGG6530	CAMEO DISK Silver 230mm dia SA Pk2		156,00
AGG6531	CAMEO DISK Silver 250mm dia SA Pk2		153,90
AGG6532	CAMEO DISK Silver 300mm dia SA Pk2		175,00
AGG6533	CAMEO DISK Gold 200mm dia SA Pk2		137,40
AGG6534	CAMEO DISK Gold 230mm dia SA Pk2		156,00
AGG6535	CAMEO DISK Gold 250mm dia SA Pk2		153,90
AGG6536	CAMEO DISK Gold 300mm dia SA Pk2		175,00
AGG6537	CAMEO DISK Platinum 1 200mm dia SA		96,10
AGG6538	CAMEO DISK Platinum 1 230mm dia SA		111,10
AGG6539	CAMEO DISK Platinum 1 250mm dia SA		122,90
AGG6540	CAMEO DISK Platinum 1 300mm dia SA		160,30
AGG6541	CAMEO DISK Platinum 2 200mm dia SA		96,10
AGG6542	CAMEO DISK Platinum 2 230mm dia SA		111,10
AGG6543	CAMEO DISK Platinum 2 250mm dia SA		122,90
AGG6544	CAMEO DISK Platinum 2 300mm dia SA		160,30
AGG6545	CAMEO DISK Platinum 3 200mm dia SA		96,10
AGG6546	CAMEO DISK Platinum 3 230mm dia SA		111,10
AGG6547	CAMEO DISK Platinum 3 250mm dia SA		122,90
AGG6548	CAMEO DISK Platinum 3 300mm dia SA		160,30
AGG6549	CAMEO DISK Platinum 0 200mm dia SA		96,10
AGG6550	CAMEO DISK Platinum 0 230mm dia SA		111,10
AGG6551	CAMEO DISK Platinum 0 250mm dia SA		122,90
AGG6552	CAMEO DISK Platinum 0 300mm dia SA		160,30
AGG6553	CAMEO DISK Platinum 4 200mm dia SA		96,10
AGG6554	CAMEO DISK Platinum 4 230mm dia SA		111,10
AGG6555	CAMEO DISK Platinum 4 250mm dia SA		122,90
AGG6556	CAMEO DISK Platinum 4 300mm dia SA		160,30
AGG6557	CAMEO DISK Silver 200mm dia MG Pk2		142,90
AGG6558	CAMEO DISK Silver 230mm dia MG Pk2		149,00
AGG6559	CAMEO DISK Silver 250mm dia MG Pk2		166,00
AGG6560	CAMEO DISK Silver 300mm dia MG Pk2		183,80
AGG6561	CAMEO DISK Gold 200mm dia MG Pk2		142,90
AGG6562	CAMEO DISK Gold 230mm dia MG Pk2		149,00
AGG6563	CAMEO DISK Gold 250mm dia MG Pk2		166,00
AGG6564	CAMEO DISK Gold 300mm dia MG Pk2		183,80
AGG6565	CAMEO DISK Platinum 1 200mm dia MG		98,20
AGG6566	CAMEO DISK Platinum 1 230mm dia MG		119,30
AGG6567	CAMEO DISK Platinum 1 250mm dia MG		130,80
AGG6568	CAMEO DISK Platinum 1 300mm dia MG		168,10
AGG6569	CAMEO DISK Platinum 2 200mm dia MG		98,20
AGG6570	CAMEO DISK Platinum 2 230mm dia MG		119,30
AGG6571	CAMEO DISK Platinum 2 250mm dia MG		130,80
AGG6572	CAMEO DISK Platinum 2 300mm dia MG		168,10
AGG6573	CAMEO DISK Platinum 3 200mm dia MG		98,20
AGG6574	CAMEO DISK Platinum 3 230mm dia MG		119,30
AGG6575	CAMEO DISK Platinum 3 250mm dia MG		130,80
AGG6576	CAMEO DISK Platinum 3 300mm dia MG		168,10
AGG6577	CAMEO DISK Platinum 0 200mm dia MG		98,20
AGG6578	CAMEO DISK Platinum 0 230mm dia MG		119,30
AGG6579	CAMEO DISK Platinum 0 250mm dia MG		130,80
AGG6580	CAMEO DISK Platinum 0 300mm dia MG		168,10
AGG6581	CAMEO DISK Platinum 4 200mm dia MG		98,20
AGG6582	CAMEO DISK Platinum 4 230mm dia MG		119,30

Part #	Description	Unit	Price €
AGG6583	CAMEO DISK Platinum 4 250mm dia MG		130,80
AGG6584	CAMEO DISK Platinum 4 300mm dia MG		168,10
AGG6585	Polishing Cloth Touchlam 2FC1 200mm dia SA Pk5		96,10
AGG6586	Polishing Cloth Touchlam 2FC1 230mm dia SA Pk5		117,70
AGG6587	Polishing Cloth Touchlam 2FC1 250mm dia SA Pk5		125,00
AGG6588	Polishing Cloth Touchlam 2FC1 300mm dia SA Pk5		160,30
AGG6589	Polishing Cloth Touchlam 2FC1 200mm dia MG Pk5		117,10
AGG6590	Polishing Cloth Touchlam 2FC1 230mm dia MG Pk5		152,80
AGG6591	Polishing Cloth Touchlam 2FC1 250mm dia MG Pk5		147,60
AGG6592	Polishing Cloth Touchlam 2FC1 300mm dia MG Pk5		192,90
AGG6593	Polishing Cloth Touchlam 2TS3 200mm dia SA Pk5		84,60
AGG6594	Polishing Cloth Touchlam 2TS3 230mm dia SA Pk5		94,30
AGG6595	Polishing Cloth Touchlam 2TS3 250mm dia SA Pk5		99,30
AGG6596	Polishing Cloth Touchlam 2TS3 300mm dia SA Pk5		133,70
AGG6597	Polishing Cloth Touchlam 2TS3 200mm dia MG Pk5		103,00
AGG6598	Polishing Cloth Touchlam 2TS3 230mm dia MG Pk5		122,60
AGG6599	Polishing Cloth Touchlam 2TS3 250mm dia MG Pk5		128,90
AGG6600	Polishing Cloth Touchlam 2TS3 300mm dia MG Pk5		157,90
AGG6601	Polishing Cloth Touchlam 2TS4 200mm dia SA Pk5		114,00
AGG6602	Polishing Cloth Touchlam 2TS4 230mm dia SA Pk5		133,10
AGG6603	Polishing Cloth Touchlam 2TS4 250mm dia SA Pk5		146,00
AGG6604	Polishing Cloth Touchlam 2TS4 300mm dia SA Pk5		181,10
AGG6605	Polishing Cloth Touchlam 2TS4 200mm dia MG Pk5		148,20
AGG6606	Polishing Cloth Touchlam 2TS4 230mm dia MG Pk5		157,60
AGG6607	Polishing Cloth Touchlam 2TS4 250mm dia MG Pk5		172,40
AGG6608	Polishing Cloth Touchlam 2TS4 300mm dia MG Pk5		217,10
AGG6609	Polishing Cloth Touchlam 2TT1 200mm dia SA Pk5		113,20
AGG6610	Polishing Cloth Touchlam 2TT1 230mm dia SA Pk5		132,70
AGG6611	Polishing Cloth Touchlam 2TT1 250mm dia SA Pk5		146,00
AGG6612	Polishing Cloth Touchlam 2TT1 300mm dia SA Pk5		179,00
AGG6613	Polishing Cloth Touchlam 2TT1 200mm dia MG Pk5		146,90
AGG6614	Polishing Cloth Touchlam 2TT1 230mm dia MG Pk5		157,30
AGG6615	Polishing Cloth Touchlam 2TT1 250mm dia MG Pk5		172,40
AGG6616	Polishing Cloth Touchlam 2TT1 300mm dia MG Pk5		217,10
AGG6617	Polishing Cloth Touchlam 2TT2 200mm dia SA Pk5		124,00
AGG6618	Polishing Cloth Touchlam 2TT2 230mm dia SA Pk5		146,00
AGG6619	Polishing Cloth Touchlam 2TT2 250mm dia SA Pk5		161,00
AGG6620	Polishing Cloth Touchlam 2TT2 300mm dia SA Pk5		183,40
AGG6621	Polishing Cloth Touchlam 2TT2 200mm dia MG Pk5		146,90
AGG6622	Polishing Cloth Touchlam 2TT2 230mm dia MG Pk5		172,80
AGG6623	Polishing Cloth Touchlam 2TT2 250mm dia MG Pk5		177,10
AGG6624	Polishing Cloth Touchlam 2TT2 300mm dia MG Pk5		238,80
AGG6625	Polishing Cloth Touchlam 3FV1 200mm dia SA Pk5		58,10
AGG6626	Polishing Cloth Touchlam 3FV1 230mm dia SA Pk5		69,60
AGG6627	Polishing Cloth Touchlam 3FV1 250mm dia SA Pk5		77,50
AGG6628	Polishing Cloth Touchlam 3FV1 300mm dia SA Pk5		95,80
AGG6629	Polishing Cloth Touchlam 3FV1 200mm dia MG Pk5		69,10
AGG6630	Polishing Cloth Touchlam 3FV1 230mm dia MG Pk5		83,60
AGG6631	Polishing Cloth Touchlam 3FV1 250mm dia MG Pk5		93,00
AGG6632	Polishing Cloth Touchlam 3FV1 300mm dia MG Pk5		114,30
AGG6633	Polishing Cloth Touchlam 3FV2 400mm dia dia SA Pk5		151,90
AGG6634	Polishing Cloth Touchlam 3SA4 200mm dia SA Pk5		66,30
AGG6635	Polishing Cloth Touchlam 3SA4 230mm dia SA Pk5		95,10
AGG6636	Polishing Cloth Touchlam 3SA4 250mm dia SA Pk5		101,20
AGG6637	Polishing Cloth Touchlam 3SA4 300mm dia SA Pk5		124,30
AGG6638	Polishing Cloth Touchlam 3SA4 200mm dia MG Pk5		86,10
AGG6639	Polishing Cloth Touchlam 3SA4 230mm dia MG Pk5		115,80
AGG6640	Polishing Cloth Touchlam 3SA4 250mm dia MG Pk5		131,30

Part #	Description	Unit	Price €
AGG6641	Polishing Cloth Touchlam 3SA4 300mm dia MG Pk5		147,10
AGG6642	Polishing Cloth Touchlam 3SE2 200mm dia SA Pk5		71,70
AGG6643	Polishing Cloth Touchlam 3SE2 230mm dia SA Pk5		95,90
AGG6644	Polishing Cloth Touchlam 3SE2 250mm dia SA Pk5		109,20
AGG6645	Polishing Cloth Touchlam 3SE2 300mm dia SA Pk5		134,00
AGG6646	Polishing Cloth Touchlam 3SE2 200mm dia MG Pk5		92,90
AGG6647	Polishing Cloth Touchlam 3SE2 230mm dia MG Pk5		124,70
AGG6648	Polishing Cloth Touchlam 3SE2 250mm dia MG Pk5		141,80
AGG6649	Polishing Cloth Touchlam 3SE2 300mm dia MG Pk5		158,50
AGG6650	Polishing Cloth Touchlam 3TL1 200mm dia SA Pk5		103,70
AGG6651	Polishing Cloth Touchlam 3TL1 230mm dia SA Pk5		128,50
AGG6652	Polishing Cloth Touchlam 3TL1 250mm dia SA Pk5		145,70
AGG6653	Polishing Cloth Touchlam 3TL1 300mm dia SA Pk5		157,90
AGG6654	Polishing Cloth Touchlam 3TL1 200mm dia MG Pk5		135,20
AGG6655	Polishing Cloth Touchlam 3TL1 230mm dia MG Pk5		152,20
AGG6656	Polishing Cloth Touchlam 3TL1 250mm dia MG Pk5		172,10
AGG6657	Polishing Cloth Touchlam 3TL1 300mm dia MG Pk5		191,40
AGG6658	Polishing Cloth Touchlam 4FV3 200mm dia SA Pk5		50,20
AGG6659	Polishing Cloth Touchlam 4FV3 230mm dia SA Pk5		60,70
AGG6660	Polishing Cloth Touchlam 4FV3 250mm dia SA Pk5		67,90
AGG6661	Polishing Cloth Touchlam 4FV3 300mm dia SA Pk5		92,00
AGG6662	Polishing Cloth Touchlam 4FV3 200mm dia MG Pk5		60,10
AGG6663	Polishing Cloth Touchlam 4FV3 230mm dia MG Pk5		72,90
AGG6664	Polishing Cloth Touchlam 4FV3 250mm dia MG Pk5		81,10
AGG6665	Polishing Cloth Touchlam 4FV3 300mm dia MG Pk5		104,00
AGG6666	Polishing Cloth Touchlam 4MP1 200mm dia SA Pk5		172,80
AGG6667	Polishing Cloth Touchlam 4MP1 230mm dia SA Pk5		209,00
AGG6668	Polishing Cloth Touchlam 4MP1 250mm dia SA Pk5		218,60
AGG6669	Polishing Cloth Touchlam 4MP1 300mm dia SA Pk5		267,90
AGG6670	Polishing Cloth Touchlam 4MP1 200mm dia MG Pk5		193,30
AGG6671	Polishing Cloth Touchlam 4MP1 230mm dia MG Pk5		250,70
AGG6672	Polishing Cloth Touchlam 4MP1 250mm dia MG Pk5		262,20
AGG6673	Polishing Cloth Touchlam 4MP1 300mm dia MG Pk5		321,40
AGG6674	Polishing Cloth Touchlam 4MP2 200mm dia SA Pk5		161,10
AGG6675	Polishing Cloth Touchlam 4MP2 230mm dia SA Pk5		194,80
AGG6676	Polishing Cloth Touchlam 4MP2 250mm dia SA Pk5		203,40
AGG6677	Polishing Cloth Touchlam 4MP2 300mm dia SA Pk5		258,60
AGG6678	Polishing Cloth Touchlam 4MP2 200mm dia MG Pk5		195,40
AGG6679	Polishing Cloth Touchlam 4MP2 230mm dia MG Pk5		253,40
AGG6680	Polishing Cloth Touchlam 4MP2 250mm dia MG Pk5		264,30
AGG6681	Polishing Cloth Touchlam 4MP2 300mm dia MG Pk5		336,10
AGG6682	Filter 100? for CUTLAM 2.02 300mm MG		o.a.
AGG6685	Adaptor Ring for I25,4/22 Steel Pk5		77,20
AGG6686	Adaptor Ring for I25,4/22 PVC Pk5		77,20
AGG6689	Filter 100? for CUTLAM Micro2.0		o.a.
AGG6690	Filter 100? for CUTLAM 4.0		o.a.
AGG6691	Adaptor ring 32-16 steel. Pk5		77,20
AGG6692	Refractometer.		349,50
AGG6693	Low Speed Option for CUTLAM MICRO 2.0		o.a.
AGG6695	Abrasive Stone 6x13x100mm Gr120 for CAMEO P0&P1 Pk2		27,80
AGG6696	Abrasive Stone 6x13x100mm Gr400 for CAMEO P2 Pk2		27,80
AGG6697	Abrasive Stone 6x13x100mm Gr600 for CAMEO P3&P4 Pk2		27,80
AGG6698	Dressing Tool for CAMEO GOLD and SILVER		33,60
AGG700	Rotary Grid Storage Box for 25 grids	each	6,10
AGG700-10	Rotary Grid Storage Boxes for 25 grids	10	50,30
AGG700S	Serialised Rotary Grid Storage Box for 25 Grids	each	11,50
AGG700S-10	Serialised Rotary Grid Storage Box for 25 Grids	10	93,20
AGG701	Rotary Grid Storage Box Kit - Dumont Tweezer N5AC	each	109,30

Part #	Description	Unit	Price €
AGG702	Rotary Grid Storage Box Kit - Dumont HP Tweezers 5 SS	each	84,80
AGG706	Rotary Grid Storage Box for 6 Grids	each	6,10
AGG706-10	Rotary Grid Storage Box for 6 Grids	10	56,60
AGGL29961	Rotating SEM Specimen Holder - Pin Mount	each	282,80
AGGL29961-H	Rotating SEM Specimen Holder - M4 Mount	each	256,00
AGGP29936	Planotec Marker, Copper	20	33,40
AGH100	Cryo Compatible Slide Box, polycarbonate, holds 100 slides,	each	46,40
AGH101	Tall Slide Box, 100 slide capacity	each	43,40
AGH102	Wooden Slide Storage Box for 100 slides 38x76mm	each	114,40
AGH103	Wooden Slide Storage Box for 100 slides 51x76mm	each	143,00
AGH104	Slide-Sette Saver Storage Box	each	57,30
AGH105	Slide-Sette Saver, natural	25	34,40
AGH106	Slide-Sette Saver, black	25	34,40
AGH107	Paraffin Block Mailer	50	137,90
AGH108	LockMailer Microscope Mailer and Staining Jar, without Capin	each	79,90
AGH108-1	Capinserts for LockMailer, blue	100	6,10
AGH108-2	Capinserts for LockMailer, grey	100	6,10
AGH108-3	Capinserts for LockMailer, green	100	6,10
AGH108-4	Capinserts for LockMailer, lila	100	6,10
AGH108-5	Capinserts for LockMailer, orange	100	6,10
AGH109	Mini Slide Box for 25 slides, black	each	6,00
AGH109-1	Mini Slide Box for 25 slides, white	each	6,40
AGH109-2	Mini Slide Box for 25 slides, red	each	6,00
AGH109-3	Mini Slide Box for 25 slides, green	each	6,00
AGH109-4	Mini Slide Box for 25 slides, blue	each	6,00
AGH109-5	Mini Slide Box for 25 slides, yellow	each	6,00
AGH111	Slide Storage System, blue	each	22,90
AGH111-1	Slide Storage System, green	each	22,90
AGH111-2	Slide Storage System, pink	each	22,90
AGH111-3	Slide Storage System, white	each	22,90
AGH111-4	Slide Storage System, yellow	each	22,90
AGH112	Slide Rack accommodates 25 slides	2	101,60
AGH113	Staining Dish holds one Slide Rack	each	130,00
AGH114	Storage Box for 4 slide racks	each	184,00
AGH115	Slide Staining and Storage kit	each	374,00
AGH116	Slide Tray Cabinet	each	547,50
AGH117	Slide Holder Strip	each	26,30
AGH118	* Histo/Cyto-freeze	283g	#N/B
AGH119	* Cryogen Spray	283g	#N/B
AGH120	Handi-Vac Coverslip Pick-Up Tool	each	33,80
AGH121	Pop-Up 20-Slide Holder, white/black	each	16,40
AGH121-1	Pop-Up 20-Slide Holder, white/white	each	16,70
AGH121-2	Pop-Up 20-Slide Holder, white/red	each	16,40
AGH121-3	Pop-Up 20-Slide Holder, white/green	each	16,40
AGH121-4	Pop-Up 20-Slide Holder, white/blue	each	16,40
AGH121-5	Pop-Up 20-Slide Holder, white/yellow	each	16,40
AGHCB1002	1000g Highland Precision Balance	each	#N/B
AGHCB123	120g Highland Precision Balance	each	#N/B
AGHCB1502	1500g Highland Precision Balance	each	#N/B
AGHCB153	150g Highland Precision Balance	each	#N/B
AGHCB3001	3000g Highland Precision Balance	each	#N/B
AGHCB302	300g Highland Precision Balance	each	#N/B
AGHCB602	600g Highland Precision Balance, 0.02g Readability	each	#N/B
AGHCB602H	600g Highland Precision Balance, 0.01g Readability	each	#N/B
AGHM675M	Hydronic Manometer Metric	each	#N/B
AGHM685M	Hydronic Manometer Metric with Flow	each	#N/B
AGJ401	TEM Grid Holder Non-Magnetic	each	o.a.
AGJ403	TEM Grid/Sample Holder - 2 sample pucks, non-magnetic	each	o.a.

Part #	Description	Unit	Price €
AGJ410	Tungsten probe tips with nickel Shank, tip radius 0.5µm, 13° taper angle.	10	450,10
AGJ411	All tungsten probe tips, tip radius <0.5µm, 10-13° taper	10	212,40
AGJ412	AP300 tungsten probe tips, Tip radius 0.5µm, 8-1	20	548,70
AGJ413	AP250 probe tips. Box of 10. Tip radius <0.5µm, 10° taper	10	279,10
AGJ414	Xtreme Access Short-Cut probe tips, 6° taper angle.	10	296,80
AGJ415	AP300 Helios 1200 probe tips, Tip radius 0.5µm,	20	743,30
AGJ420	TEM Lift-Out Grids Copper 3 Post	100	339,80
AGJ421	TEM Lift-Out Grids Molybdenum 3 Post	25	453,10
AGJ422	5 post lift-out grids, copper.	100	380,20
AGJ423	TEM Lift-Out Grids Copper 4 Post	100	348,10
AGJ424	TEM Lift-Out Grids Molybdenum 4 Post	25	551,20
AGJ425	Beryllium half-ring lift-out grids.	10	823,40
AGJ430	Lift-out grid storage box	each	71,00
AGJ440	45° Short-Cut Coupons Top-Side Thinning Copper	20	255,30
AGJ441	45° Short-Cut Coupons Bottom Thinning Copper	20	211,20
AGJ442	Atom probe preparation coupons.	20	236,40
AGJ443	26.5° Coupons Top-Side Thinning Copper	20	240,40
AGJ444	26.5° Coupons Bottom Thinning Copper	20	211,20
AGJ456	AP250 Probe Point Holder	each	939,40
AGJ460	FIB lift-out grids.	25	67,00
AGL2302-1	PELCO Microslides Storage Box, white, 100 slides	each	16,00
AGL2302-1-CS	PELCO Microslides Storage Box, white, 100 slides	12	167,60
AGL2302-2	PELCO Microslides Storage Box, red, 100 slides	each	16,00
AGL2302-2-CS	PELCO Microslides Storage Box, red, 100 slides	12	167,60
AGL2302-3	PELCO Microslides Storage Box, green, 100 slides	each	16,00
AGL2302-3-CS	PELCO Microslides Storage Box, green, 100 slides	12	167,60
AGL2302-4	PELCO Microslides Storage Box, blue, 100 slides	each	16,00
AGL2302-4-CS	PELCO Microslides Storage Box, blue, 100 slides	12	167,60
AGL2302-5	PELCO Microslides Storage Box, yellow, 100 slides	each	16,00
AGL2302-5-CS	PELCO Microslides Storage Box, yellow, 100 slides	12	167,60
AGL2302-6	PELCO Microslides Storage Box, black, 100 slides	each	16,00
AGL2302-6-CS	PELCO Microslides Storage Box, black, 100 slides	12	167,60
AGL2304-1	PELCO Microslides Storage Box, white, 50 slides	each	12,70
AGL2304-1-CS	PELCO Microslides Storage Box, white, 50 slides	12	128,90
AGL2304-2	PELCO Microslides Storage Box, red, 50 slides	each	12,70
AGL2304-2-CS	PELCO Microslides Storage Box, red, 50 slides	12	128,90
AGL2304-3	PELCO Microslides Storage Box, green, 50 slides	each	12,70
AGL2304-3-CS	PELCO Microslides Storage Box, green, 50 slides	12	128,90
AGL2304-4	PELCO Microslides Storage Box, blue, 50 slides	each	12,70
AGL2304-4-CS	PELCO Microslides Storage Box, blue, 50 slides	12	128,90
AGL2304-5	PELCO Microslides Storage Box, yellow, 50 slides	each	12,70
AGL2304-5-CS	PELCO Microslides Storage Box, yellow, 50 slides	12	128,90
AGL2304-6	PELCO Microslides Storage Box, black, 50 slides	each	12,70
AGL2304-6-CS	PELCO Microslides Storage Box, black, 50 slides	12	128,90
AGL2306-1	PELCO Microslides Storage Box, white, 25 slides	each	8,70
AGL2306-1-CS	PELCO Microslides Storage Box, white, 25 slides	12	96,30
AGL2306-2	PELCO Microslides Storage Box, red, 25 slides	each	8,70
AGL2306-2-CS	PELCO Microslides Storage Box, red, 25 slides	12	96,30
AGL2306-3	PELCO Microslides Storage Box, green, 25 slides	each	8,70
AGL2306-3-CS	PELCO Microslides Storage Box, green, 25 slides	12	96,30
AGL2306-4	PELCO Microslides Storage Box, blue, 25 slides	each	8,70
AGL2306-4-CS	PELCO Microslides Storage Box, blue, 25 slides	12	96,30
AGL2306-5	PELCO Microslides Storage Box, yellow, 25 slides	each	8,70
AGL2306-5-CS	PELCO Microslides Storage Box, yellow, 25 slides	12	96,30
AGL2306-6	PELCO Microslides Storage Box, black, 25 slides	each	8,70
AGL2306-6-CS	PELCO Microslides Storage Box, black, 25 slides	12	97,60
AGL4004	THMS600 Freezing-Heating-Stage with Manipulators -196 to +60	each	o.a.
AGL4028	Neutral Density Filter 0.30 ND 1 Stop	each	31,30

Part #	Description	Unit	Price €
AGL4030	Neutral Density Filter 0.60 ND 2 Stops	each	31,30
AGL4032	Neutral Density Filter 0.90 ND 3 Stops	each	31,30
AGL4033	Graticule A 16mm Diameter	each	98,10
AGL4033-23	Graticule A 23mm Diameter	each	92,50
AGL4033-26	Graticule A 26mm Diameter	each	141,30
AGL4033-27	Graticule A 27mm Diameter	each	141,30
AGL4034	Graticule A 21mm Diameter	each	100,50
AGL4035	Graticule B 16mm Diameter	each	113,20
AGL4035-23	Graticule B 23mm Diameter	each	92,90
AGL4036	Graticule B 21mm Diameter	each	113,20
AGL4037	Graticule C 16mm Diameter	each	124,60
AGL4037-23	Graticule C 23mm Diameter	each	124,60
AGL4037-26	Graticule C 26mm Diameter	each	158,90
AGL4038	Graticule C 21mm Diameter	each	124,60
AGL4039	Graticule D 16mm Diameter	each	113,20
AGL4040	Graticule D 21mm Diameter	each	113,20
AGL4041	Graticule E 16mm Diameter	each	98,10
AGL4041-23	Graticule E 23mm Diameter	each	101,10
AGL4041-26	Graticule E 26mm Diameter	each	119,00
AGL4041-27	Graticule E 27mm Diameter	each	126,20
AGL4042	Graticule E 21mm Diameter	each	98,10
AGL4043	Graticule F 16mm Diameter	each	124,60
AGL4044	Graticule F 21mm Diameter	each	124,60
AGL4045	Graticule G 16mm Diameter	each	124,60
AGL4046	Graticule G 21mm Diameter	each	124,60
AGL4047	Graticule H 16mm Diameter	each	124,60
AGL4048	Graticule H 21mm Diameter	each	124,60
AGL4049	Graticule I 16mm Diameter	each	113,20
AGL4050	Graticule I 21mm Diameter	each	113,20
AGL4051	Graticule J 16mm Diameter	each	124,40
AGL4052	Graticule J 21mm Diameter	each	124,40
AGL4053	Graticule K 16mm Diameter	each	98,10
AGL4054	Graticule K 21mm Diameter	each	98,10
AGL4054-26	Graticule K 26mm Diameter	each	141,30
AGL4055	Graticule L 16mm Diameter	each	124,60
AGL4056	Graticule L 21mm Diameter	each	124,60
AGL4057	Graticule M 16mm Diameter	each	218,40
AGL4058	Graticule M 21mm Diameter	each	218,40
AGL4059	Graticule N 16mm Diameter	each	218,40
AGL4060	Graticule N 21mm Diameter	each	218,40
AGL4061	Graticule O 16mm Diameter	each	232,00
AGL4062	Graticule O 21mm Diameter	each	232,00
AGL4063	Graticule P 16mm Diameter	each	232,00
AGL4064	Graticule P 21mm Diameter	each	207,40
AGL4065	Graticule Q 16mm Diameter	each	207,40
AGL4066	Graticule Q 21mm Diameter	each	232,00
AGL4076	The England Finder	each	257,60
AGL4077AP	Refill Ink for AGL4077P Purple Slide Marker	20ml	21,20
AGL4077AR	Refill Ink for AGL4077R Red Slide Marker	20ml	28,00
AGL4077P	Slide Marker (Purple)	each	275,70
AGL4077R	Slide Marker (Red)	each	281,50
AGL4078	Stage Micrometer 1mm Scale x 0.01mm Divisions	each	180,90
AGL4078NPL	NPL Certified Stage Micrometer, 1mm Scale x 0.01mm Sub-Divisions	each	o.a.
AGL4078S	UKAS Certified Stage Micrometer, 1mm Scale x 0.01mm Sub-Divisions	each	o.a.
AGL4079	Stage Micrometer 1mm Scale x 0.01mm for Reflected Light	each	211,10
AGL4079NPL	NPL Certified Stage Micrometer, 1mm Scale x 0.01mm Sub-Divisions, for refl	each	o.a.
AGL4079S	UKAS Certified Stage Micrometer, 1mm Scale x 0.01mm Sub-Divisions, for refl	each	o.a.
AGL4081	Recalibration of UKAS Certified Stage Micrometer 1 x 0.01mm	each	297,90

Part #	Description	Unit	Price €
AGL4081NPL	Stage Micrometer 10mm Scale x 0.1mm for Reflected Light	each	o.a.
AGL4081S	NPL Certified Stage Micrometer, 10mm Scale x 0.1mm Sub-Divisions, for reflex	each	o.a.
AGL4086	UKAS Certified Stage Micrometer, 10mm Scale x 0.1mm Sub-Divisions, for refl	each	19,30
AGL4090	Cavity Slide Single Depression	50	151,70
AGL4091A	Cavity Slide Two Depressions	10	43,70
AGL4102	Plastic Plate with 12 Cavities	each	8,30
AGL4103	Melinex film sheets	5	38,90
AGL4108	Coplin jar	each	19,60
AGL4109	Plastic Coplin jar	each	28,80
AGL4109-R	Plastic Coplin jar (round)	each	28,80
AGL4111	Plastic Staining Rack	each	#N/B
AGL4113	Blotting paper, 100 sheets	each	24,50
AGL4114	Two slide mailer	each	1,60
AGL4115	9 Slide Tray	10	28,90
AGL4116	Slide Tray Box	each	55,00
AGL4121	Dissecting Board, Cork	each	31,50
AGL4124	Dissecting Pins	10	114,30
AGL4125	Processing Cassettes White	1000	#N/B
AGL4126	Processing Cassettes Green	1000	122,90
AGL4127	Processing Cassettes Pink	1000	122,90
AGL4128	Processing Cassettes Yellow	1000	122,90
AGL4129	Processing Cassettes Blue	1000	122,90
AGL4130	Processing Cassettes Orange	1000	122,90
AGL4131	Processing Cassettes Grey	1000	122,90
AGL4132	Pelleted Paraffin Wax	1kg	35,80
AGL4133	Pelleted Paraffin Wax	10kg	275,30
AGL4136	Paraffin section mounting bath	each	o.a.
AGL4146	Steel Knife C-Profile 160mm x 34mm x 10mm	each	620,60
AGL4147	Tungsten Carbide Knife D-Profile 160mm x 34mm x 10mm	each	o.a.
AGL4157	Ralph Knifemaker (Histology) 25mm	each	3.471,70
AGL4158	Ralph Knifemaker for 25mm wide glass	each	o.a.
AGL4160	Slide labels, 19 X 24mm.	1000	19,50
AGL4163	Staining Trough Black Plastic	each	52,80
AGL4164	Filter paper	each	29,30
AGL4166	TK drop bottle, 50ml	each	#N/B
AGL4180	Filing cabinet for cassettes	each	402,10
AGL4181	Fibrowax	1kg	#N/B
AGL4184	Polyester Wax	500g	276,00
AGL4185	Glycerin Albumen 100ml	100ml	163,00
AGL4186	Gelatin 100g	each	71,40
AGL4193	Plastic Cover Slips 22mm x 22mm	100	9,40
AGL4197M	Super Mini Pap Pen (Clear).	7,5ml	86,20
AGL4197S	Super Pap Pen.	16,5ml	116,40
AGL4198	Fro-Tissuer Pen	each	163,10
AGL4199	Para-Tissuer Pen	each	163,10
AGL4201	Stage Micrometer 10mm Scale x 0.1mm Divisions	each	205,40
AGL4201NPL	NPL Certified Stage Micrometer, 10mm Scale x 0.1mm Sub-Divisions	each	1.192,00
AGL4201S	UKAS Certified Stage Micrometer, 10mm Scale x 0.1mm Sub-Divisions	each	o.a.
AGL4202	Stage Micrometer 0.1mm Scale x 0.002mm Divisions	each	324,60
AGL4202NPL	NPL Certified Stage Micrometer, 0.1mm Scale x 0.002mm Sub-Divisions	each	1.527,10
AGL4202S	UKAS Certified Stage Micrometer, 0.1mm Scale x 0.002mm Sub-Divisions	each	o.a.
AGL4203	Stage Micrometer 0.1inch Scale x 0.001inch Divisions	each	204,00
AGL4203NPL	NPL Certified Stage Micrometer, 0.1" Scale x 0.001" Sub-Divisions	each	1.204,40
AGL4203S	UKAS Certified Stage Micrometer, 0.1" Scale x 0.001" Sub-Divisions	each	o.a.
AGL4204	Stage Micrometer 2mm Scale 0.01mm Divisions	each	205,80
AGL4207	UKAS Certified Stage Micrometer 2mm Scale 0.01mm OBSOLETE WHEN SOLD C	each	480,80
AGL4207NPL	Stage Micrometer 20mm Scale x 10 Micron Divisions	each	1.285,50
AGL4207S	NPL Certified Stage Micrometer 20mm Scale x 10µm Divs	each	o.a.

Part #	Description	Unit	Price €
AGL4208	UKAS Certified Stage Micrometer 20mm Scale x 10µm Divs	each	480,90
AGL4208S	Calibration Standard Slide	each	o.a.
AGL4222A	UKAS Certified Calibration Standard Slide- Pattern A (CONCENTRIC CIRCLES) 1,	50	#N/B
AGL4242	Plain Slides, 76x26mm, 0.8mm - 1.0mm Thick	50	12,70
AGL4243	Twin Frost Slides 76 x 26mm 0.8mm - 1.0mm Thick	50	36,70
AGL4244	Plain Slides, 76x26mm, 1.0mm - 1.2mm Thick	50	11,30
AGL4245	Twin Frost Slides 76 x 26mm 1.0mm - 1.2mm Thick	50	37,30
AGL4247A	Coverglasses for Haemocytometer	10	25,80
AGL4247-B	Burker - Double cell	each	169,10
AGL4247-F	Fuchs Rosenthal - Double cell	each	148,20
AGL4247-I	Improved Neubauer - Double cell	each	178,70
AGL4247-N	Neubauer - Double cell	each	148,20
AGL4247-TH	Thoma - Double cell	each	148,20
AGL4248-B	Burker - Single cell	each	151,90
AGL4248-F	Fuchs Rosenthal - Single cell	each	151,90
AGL4248-I	Improved Neubauer - Single cell	each	151,90
AGL4248-N	Neubauer - Single cell	each	148,20
AGL4250	Five slide mailer	each	2,90
AGL4257	Coverslip staining racks to hold 25 coverslips	each	86,70
AGL4261	Filing Cabinet for 6000 Slides	each	439,60
AGL4265	Filing Cabinet for 1500 Slides	each	#N/B
AGL4266	Filing Cabinet for 3000 Slides	each	423,40
AGL4267	Spring Spacer for AGL4265 Slide Cabinet (240mm)	each	25,70
AGL4268	Drawer Divider for AGL4265 Slide Cabinet	each	25,50
AGL4274	Filing Cabinet for 3000 Cassettes, 480 x 480 x 275	each	483,10
AGL4283	Set of 16 High Dispersion Liquids (HD 16). 7.4ml bottles	each	o.a.
AGL4284	Set of 31 High Dispersion Liquids (HD 1/2). 7.4ml bottles	each	o.a.
AGL4285	Set of 6 High Dispersion Liquids for the identification of a	each	827,20
AGL4295	Quick-Stick R.I.=1.662 30ml	30ml	281,00
AGL4295-1	Quick-Stick R.I.=1.539 30ml	30ml	395,50
AGL4295-2	Quick-Stick R.I.=1.582 30ml	30ml	281,00
AGL4295-3	Quick-Stick R.I.=1.605 30ml	30ml	395,50
AGL4295-4	Quick-Stick R.I.=1.686 30ml	30ml	371,30
AGL4295-5	Quick-Stick R.I.=1.704 30ml	30ml	371,30
AGL4299	Top cover for 480mm Unit	each	132,00
AGL4300	Base stand for 480mm Unit	each	138,60
AGL4320	Graticule A 19mm Diameter	each	98,10
AGL4321	Graticule B 19mm Diameter	each	113,20
AGL4322	Graticule C 19mm Diameter	each	124,60
AGL4323	Graticule D 19mm Diameter	each	113,20
AGL4324	Graticule E 19mm Diameter	each	98,10
AGL4325	Graticule F 19mm Diameter	each	111,20
AGL4326	Graticule G 19mm Diameter	each	124,60
AGL4327	Graticule H 19mm Diameter	each	111,20
AGL4328	Graticule I 19mm Diameter	each	113,20
AGL4329	Graticule J 19mm Diameter	each	111,20
AGL4330	Graticule K 19mm Diameter	each	98,10
AGL4331	Graticule L 19mm Diameter	each	124,60
AGL4332	Graticule M 19mm Diameter	each	195,20
AGL4333	Graticule N 19mm Diameter	each	195,20
AGL4334	Graticule O 19mm Diameter	each	232,00
AGL4335	Graticule P 19mm Diameter	each	207,40
AGL4336	Graticule Q 19mm Diameter	each	207,40
AGL4338-1	Rectangular coverglasses, 24 x 60mm, 1 Thickness	10	23,40
AGL4340	Super Frost Slides - Blue	50	22,30
AGL4341	Super Frost Slides - Green	50	22,30
AGL4342	Super Frost Slides - Orange	50	22,30
AGL4343	Super Frost Slides - Pink	50	12,40

Part #	Description	Unit	Price €
AGL4344	Super Frost Slides - Yellow	50	22,30
AGL4345	Polysine microscope slides.	72	74,80
AGL4351	Circular Thermanox Coverslips, 22mm dia	100	199,80
AGL4352	Rectangular Thermanox Coverslips, 22mm x 60mm	100	246,40
AGL4353	Circular Thermanox Coverslips, 15mm dia	100	302,50
AGL4354	Rectangular Thermanox Coverslips, 10.5mm x 22mm	100	310,30
AGL4361	Circular Thermanox Coverslips, 25mm dia	100	199,80
AGL4362	Rectangular Thermanox Coverslips, 24mm x 30mm	100	199,80
AGL4363	Metal Slides, 75 x 25mm.	5	147,00
AGL4364	Plain Glass Slides 28mm x 48mm x 1mm	100	102,40
AGL4365	Glass Slides 46mm x 27mm x 1.2-1.4mm	100	204,10
AGL4365-72	Glass Slides 46mm x 27mm x 1.2-1.4mm (Pack of 72)		116,30
AGL4367	Plain Glass Slides 24mm x 24mm x 1mm	100	#N/B
AGL4368-1	Standard Slide Sheet Labels 24mm x 24mm	1000	44,80
AGL4368-2	Pathology Slide Sheet Labels 24mm x 24mm	1000	41,90
AGL4368-3	Standard Slide Sheet Labels 12mm x 24mm	2000	45,30
AGL4368-4	Pathology Slide Sheet Labels 12mm x 24mm	2000	45,30
AGL4369-1	Standard Slide Roll Labels 24mm x 24mm	1000	41,90
AGL4369-3	Standard Slide Roll Labels 12mm x 24mm	1000	22,60
AGL4370	Tissue capture pen	each	171,60
AGL4371	Securline Lab Marker, extra fine tip, black	each	21,10
AGL4374	Mcllwain Tissue Chopper	each	o.a.
AGL4375	BioWave Pro+ Laboratory Microwave System	each	o.a.
AGL4375-10	Electron Microscopy Tissue Processing Application Kit	each	o.a.
AGL4375-20	Immunolabelling Kit	each	o.a.
AGL4375-30	Formaldehyde Fixation and Decalcification Kit	each	o.a.
AGL4375-40	Paraffin Tissue Processing Kit	each	o.a.
AGL4376	BioWave Pro+ Microwave tissue processing system	each	o.a.
AGL4380-1	102mm x 76mm Large Glass Slides	36	249,20
AGL4380-2	102mm x 83mm Large Glass Slides	36	213,00
AGL4380-3	127mm x 102mm Large Glass Slides	36	295,20
AGL4380-4	152mm x 114mm Large Glass Slides	36	295,70
AGL4380-5	178mm x 127mm Large Glass Slides	36	299,40
AGL4381-1	76mm x 83mm Large Cover Glasses	36	253,70
AGL4381-2	76mm x 89mm Large Cover Glasses	36	263,20
AGL4381-3	95mm x 114mm Large Cover Glasses	36	339,60
AGL4381-4	105mm x 140mm Large Cover Glasses	36	421,40
AGL4381-5	114mm x 159mm Large Cover Glasses	36	444,50
AGL4382	Paraplast Plus	1kg	44,90
AGL4384	Slide drying hotplate, 240V	each	#N/B
AGL4384S	Slimline Slide drying hotplate, 240V	each	#N/B
AGL4385	Trimming Board 12in (305mm) Square	each	74,50
AGL4390	Sedgewick Rafter Counting Chamber - Plastic	each	59,20
AGL4391	Sedgewick Rafter Counting Chamber - Glass	each	#N/B
AGL4392	Cover Glass for AGL4390 & AGL4391 Counting Chambers	each	17,40
AGL4393	Howard Cell for Fruit Juices	each	331,10
AGL4394	Howard Cell coverglass	each	109,00
AGL4395	Standard Set RF-1/2 for general use, 76 Liquids ranging from	each	3.094,50
AGL4396	Standard Set RF-1/5 for general use, 31 Liquids ranging from	each	1.865,10
AGL4397A-1	Immersion Oil Type A	30ml	115,40
AGL4397A-2	Immersion Oil Type A	120ml	#N/B
AGL4397A-3	Immersion Oil Type A	480ml	260,00
AGL4397B-1	Immersion Oil Type B	30ml	123,90
AGL4397B-2	Immersion Oil Type B	120ml	170,70
AGL4397B-3	Immersion Oil Type B	480ml	260,00
AGL4397DF-1	Immersion Oil Type LDF	30ml	154,80
AGL4397DF-2	Immersion Oil Type LDF	120ml	228,50
AGL4397DF-3	Immersion Oil Type DF	480ml	330,90

Part #	Description	Unit	Price €
AGL4397FF-1	Immersion Oil Type FF	30ml	143,70
AGL4397FF-2	Immersion Oil Type FF	120ml	228,50
AGL4397FF-3	Immersion Oil Type FF	480ml	350,60
AGL4398-37	Water Matching Liquid at 37°C 1 oz		268,10
AGL4453	Square Quartz Coverslip 19mm x 19mm x 0.5mm Thick	each	26,30
AGL4454	Square Quartz Coverslip 19mm x 19mm x 0.25mm Thick	each	50,70
AGL4455	Circular Quartz Coverslip 25.4mm Diameter x 0.25mm Thick	each	67,80
AGL4456	Circular Quartz Coverslip 18mm dia. x 0.25mm Thick	each	43,50
AGL4458	200mm x 315mm Aclar Sheets	10	170,80
AGL4459	200mm x 315mm Aclar Sheets	25	340,60
AGL4460	Circular Quartz Coverslips 10mm dia.x 1mm	5	95,00
AGL4461	Circular Quartz Coverslips 25mm dia x 1.25mm	5	318,40
AGL4462	Circular Quartz Coverslips 22mm dia x 0.25mm	5	204,90
AGL4463	Circular Quartz Coverslips 10mm dia x 0.25mm	5	548,60
AGL4464	Circular Quartz Coverslips 20mm dia x 1mm	5	108,50
AGL4465	Square Quartz Coverslips 22mm dia x 22mm x 0.25mm	5	204,90
AGL4466	Circular Quartz Coverslips 10mm dia.x 0.1-0.15mm	5	194,50
AGL4470	Quartz Slides 75 x 25 x 1mm USE AG26011	5	280,30
AGL4470-1	Quartz Slide 25mm x 25mm x 1mm	each	26,80
AGL4470-2	Quartz Slide 50mm x 25mm x1mm	each	39,10
AGL4471	Quartz Slide 75 x 25 x 1mm	each	94,50
AGL4472	Staining jar for coverglasses	each	53,00
AGL4490	Dissecting Board, Autoclavable	each	81,90
AGL4492	Insect Pins 0.25 x 38mm	100	15,90
AGL4493	Insect Pins 0.1 x 12mm	50	46,30
AGL4494	Insect Pins 0.2 x 12mm	50	41,40
AGL4498	Sapphire Blade	each	o.a.
AGL4504	Microwave Polymerisation Tray	each	56,40
AGL4505	Spring Spacer for AGL4266 Slide Cabinet (480mm)	each	#N/B
AGL4506	Drawer Divider for AGL4266 Slide Cabinet	each	28,60
AGL4507	Cardboard Boxes for 50 x 25 x 80mm Blocks	100	232,70
AGL4508	Cardboard Boxes for 63 x 50 x 80mm Blocks	100	180,10
AGL4509	Ralph Knifemaker (Histology) 38mm	each	#N/B
AGL4512	Cardboard Wallets for 3 Slides	100	58,40
AGL4514	Porcelain Plate with 12 Cavities	each	84,00
AGL4516	Disposable Dissecting Boards 297 x 420 x 5mm	10	301,80
AGL4520	Micro-Ruler MR-1	each	o.a.
AGL4526	Economy Slide Box for 25 Large Slides	each	23,70
AGL4545	Plastic Slide Storage Box for 100 slides, White	each	14,10
AGL4546	Plastic Slide Storage Box for 100 slides, Red	each	14,10
AGL4547	Plastic Slide Storage Box for 100 slides, Green	each	14,10
AGL4548	Plastic Slide Storage Box for 100 slides, Blue	each	14,10
AGL4549	Plastic Slide Storage Box for 100 slides, Yellow	each	14,10
AGL4550	Triangular tungsten carbide knives	3	o.a.
AGL4552	Plastic Slide Storage Box for 100 slides, Black	each	14,10
AGL4553	Slide Saver Boxes - 25 Slide Capacity, removable lid	2	27,10
AGL4554	Slide Saver Box - 100 Slide Capacity, hinged lid	each	23,20
AGL4555	Immunostain moisture chamber	each	121,10
AGL4556	Microscope slide tray	each	177,50
AGL4558	Polypropylene staining jar	each	47,80
AGL462226-0	Rectangular coverglasses, 22 x 26mm, 0 Thickness	100	31,60
AGL462226-1	Rectangular coverglasses, 22 x 26mm, 1 Thickness	100	22,70
AGL462226-3	Rectangular coverglasses, 22 x 26mm, 3 Thickness	100	#N/B
AGL462226-5	Rectangular coverglasses, 22 x 26mm, 5 Thickness	100	#N/B
AGL462232-1	Rectangular coverglasses, 22 x 32mm, 1 Thickness	100	8,00
AGL462232-15	Rectangular coverglasses, 22 x 32mm, 1.5 Thickness	100	23,70
AGL462232-2	Rectangular coverglasses, 22 x 32mm, 2 Thickness	100	#N/B
AGL462232-3	Rectangular coverglasses, 22 x 32mm, 3 Thickness	100	40,10

Part #	Description	Unit	Price €
AGL462240-1	Rectangular coverglasses, 22 x 40mm, 1 Thickness	100	11,40
AGL462240-2	Rectangular coverglasses, 22 x 40mm, 2 Thickness	100	#N/B
AGL462240-3	Rectangular coverglasses, 22 x 40mm, 3 Thickness	100	47,00
AGL462240-4	Rectangular coverglasses, 22 x 40mm, 4 Thickness	100	#N/B
AGL462240-5	Rectangular coverglasses, 22 x 40mm, 5 Thickness	100	137,60
AGL462250-00	Rectangular coverglasses, 22 x 50mm, 00 Thickness	100	#N/B
AGL462250-1	Rectangular coverglasses, 22 x 50mm, 1 Thickness	100	13,60
AGL462250-15	Rectangular coverglasses, 22 x 50mm, 1.5 Thickness	100	15,20
AGL462250-2	Rectangular coverglasses, 22 x 50mm, 2 Thickness	100	#N/B
AGL462250-3	Rectangular coverglasses, 22 x 50mm, 3 Thickness	100	#N/B
AGL462264-15	Rectangular coverglasses, 22 x 64mm, 1.5 Thickness	100	42,40
AGL462264-2	Rectangular coverglasses, 22 x 64mm, 2 Thickness	100	40,70
AGL462264-3	Rectangular coverglasses, 22 x 64mm, 3 Thickness	100	#N/B
AGL462432-1	Rectangular coverglasses, 24 x 32mm, 1 Thickness	100	#N/B
AGL462432-3	Rectangular coverglasses, 24 x 32mm, 3 Thickness	100	#N/B
AGL462432-5	Rectangular coverglasses, 24 x 32mm, 5 Thickness	100	#N/B
AGL462440-1	Rectangular coverglasses, 24 x 40mm, 1 Thickness	100	13,80
AGL462440-15	Rectangular coverglasses, 24 x 40mm, 1.5 Thickness	100	14,30
AGL462440-3	Rectangular coverglasses, 24 x 40mm, 3 Thickness	100	50,10
AGL462440-5	Rectangular coverglasses, 24 x 40mm, 5 Thickness	100	55,60
AGL462450-0	Rectangular coverglasses, 24 x 50mm, 0 Thickness	100	14,40
AGL462450-1	Rectangular coverglasses, 24 x 50mm, 1 Thickness	100	14,10
AGL462450-15	Rectangular coverglasses, 24 x 50mm, 1.5 Thickness	100	15,20
AGL462450-2	Rectangular coverglasses, 24 x 50mm, 2 Thickness	100	42,90
AGL462450-3	Rectangular coverglasses, 24 x 50mm, 3 Thickness	100	#N/B
AGL462450-5	Rectangular coverglasses, 24 x 50mm, 5 Thickness	100	#N/B
AGL463564-1	Rectangular coverglasses, 35 x 64mm, 1 Thickness	100	61,40
AGL463564-15	Rectangular coverglasses, 35 x 64mm, 1.5 Thickness	100	64,90
AGL463564-3	Rectangular coverglasses, 35 x 64mm, 3 Thickness	100	#N/B
AGL464864-15	Rectangular coverglasses, 48 x 64mm, 1.5 Thickness	100	83,80
AGL465764-0	Rectangular coverglasses, 57 x 64mm, 0 Thickness	100	#N/B
AGL46R10-15	Circular coverglasses, 10mm dia, 1.5 Thickness	100	22,30
AGL46R10-2	Circular coverglasses, 10mm dia, 2 Thickness	100	#N/B
AGL46R13-0	Circular coverglasses, 13mm dia, 0 Thickness	100	45,50
AGL46R13-1	Circular coverglasses, 13mm dia, 1 Thickness	100	15,40
AGL46R13-15	Circular coverglasses, 13mm dia, 1.5 Thickness	100	17,50
AGL46R13-2	Circular coverglasses, 13mm dia, 2 Thickness	100	19,10
AGL46R19-0	Circular coverglasses, 19mm dia, 0 Thickness	100	18,60
AGL46R19-1	Circular coverglasses, 19mm dia, 1 Thickness	100	18,90
AGL46R19-2	Circular coverglasses, 19mm dia, 2 Thickness	100	30,00
AGL46R22-0	Circular coverglasses, 22mm dia, 0 Thickness	100	23,40
AGL46R22-1	Circular coverglasses, 22mm dia, 1 Thickness	100	23,20
AGL46R22-15	Circular coverglasses, 22mm dia, 1.5 Thickness	100	27,40
AGL46R22-2	Circular coverglasses, 22mm dia, 2 Thickness	100	43,50
AGL46R24-1	Circular coverglasses, 24mm dia, 1 Thickness	100	23,00
AGL46R5-1	Circular coverglasses, 5mm dia, 1 Thickness	100	16,80
AGL46S10-0	Square coverglasses, 10 x 10mm, 0 Thickness	100	32,10
AGL46S10-15	Square coverglasses, 10 x 10mm, 1.5 Thickness	100	27,70
AGL46S10-2	Square coverglasses, 10 x 10mm, 2 Thickness	100	39,50
AGL46S13-0	Square coverglasses, 13 x 13mm, 0 Thickness	100	#N/B
AGL46S13-00	Square coverglasses, 13 x 13mm, 00 Thickness	100	#N/B
AGL46S13-2	Square coverglasses, 13 x 13mm, 2 Thickness	100	#N/B
AGL46S13-4	Square coverglasses, 13 x 13mm, 4 Thickness	100	#N/B
AGL46S15-0	Square coverglasses, 15 x 15mm, 0 Thickness	100	#N/B
AGL46S15-2	Square coverglasses, 15 x 15mm, 2 Thickness	100	7,10
AGL46S15-3	Square coverglasses, 15 x 15mm, 3 Thickness	100	#N/B
AGL46S18-0	Square coverglasses, 18 x 18mm, 0 Thickness	100	5,50
AGL46S18-1	Square coverglasses, 18 x 18mm, 1 Thickness	100	6,00

Part #	Description	Unit	Price €
AGL46S18-15	Square coverglasses, 18 x 18mm, 1.5 Thickness	100	6,30
AGL46S18-2	Square coverglasses, 18 x 18mm, 2 Thickness	100	15,90
AGL46S18-4	Square coverglasses, 18 x 18mm, 4 Thickness	100	#N/B
AGL46S18-5	Square coverglasses, 18 x 18mm, 5 Thickness	100	#N/B
AGL46S20-1	Square coverglasses, 20 x 20mm, 1 Thickness	100	6,10
AGL46S20-15	Square coverglasses, 20 x 20mm, 1.5 Thickness	100	7,10
AGL46S20-2	Square coverglasses, 20 x 20mm, 2 Thickness	100	#N/B
AGL46S20-5	Square coverglasses, 20 x 20mm, 5 Thickness	100	#N/B
AGL46S22-0	Square coverglasses, 22 x 22mm, 0 Thickness	100	7,50
AGL46S22-1	Square coverglasses, 22 x 22mm, 1 Thickness	100	6,60
AGL46S22-2	Square coverglasses, 22 x 22mm, 2 Thickness	100	20,90
AGL46S24-0	Square coverglasses, 24 x 24mm, 0 Thickness	100	#N/B
AGL46S24-1	Square coverglasses, 24 x 24mm, 1 Thickness	100	9,80
AGL46S24-15	Square coverglasses, 24 x 24mm, 1.5 Thickness	100	8,70
AGL46S24-2	Square coverglasses, 24 x 24mm, 2 Thickness	100	#N/B
AGL46S24-5	Square coverglasses, 24 x 24mm, 5 Thickness	100	#N/B
AGL6101	Coverslip with 10x10 grid , 0.1mm squares, 5 positions. Indexed 1-10 top, A-J c	25	101,70
AGL6102	Coverslip with 10x10 grid , 1mm squares. Each square indexed 0-99.	25	93,70
AGL6103	Coverslip- 20x20Grid 0.5mmSqs Ind:1-20TopA-TSide+Center	25	96,50
AGL6104	Coverslip-10x10 grid 1mmSqs.Index:0-99,18mm Round	25	102,10
AGL6105	Coverslip- 20x20 Grid 0.5mm.1-20Top A-TSide+Cntr 18mm Round	25	99,20
AGLCA301	Rotating Vane Anemometer	each	#N/B
AGLCA501	Vane Anemometer	each	#N/B
AGLT-80	DWS Folding Metal Magnifier 18/04 SGB	each	#N/B
AGM100KIT	Olympus SZ61 Stereo Zoom Microscope with LED illumination st	each	#N/B
AG-M1-100G	M1 100g	each	#N/B
AG-M1-10G	M1 10g	each	#N/B
AG-M1-1G	M1 1g	each	#N/B
AG-M1-1G-100G	M1 1g-100g Weight set	each	#N/B
AG-M1-1G-200G	M1 1g-200g Weight set	each	#N/B
AG-M1-1G-2KG	M1 1g-2kg Weight set	each	#N/B
AG-M1-1G-500G	M1 1g-500g Weight set	500g	#N/B
AG-M1-1KG	M1	1kg	#N/B
AG-M1-200G	M1 200g	each	#N/B
AG-M1-20G	M1 20g	each	#N/B
AG-M1-2G	M1 2g	each	#N/B
AG-M1-2KG	M1 2kg	each	#N/B
AG-M1-500G	M1	500g	#N/B
AG-M1-50G	M1	50g	#N/B
AG-M1-5G	M1 5g	5g	#N/B
AGML6120-SP	Binocular polarising asbestos microscope	each	6.026,40
AGML6520-SP1	Binocular phase contrast microscope	each	4.246,90
AGM-MICFIP	WiFi Digital Microscope - Polarised	each	469,50
AGM-MICFISTD	WiFi Digital Microscope - White Light	each	408,60
AGM-MICFIT6	WiFi Digital Microscope - 600x magnification	each	469,50
AGM-MICFITL	WiFi Digital Microscope - Distance Focus	each	449,80
AGM-MICFIUVW	Wifi Digital Microscope - UV/White Light	each	596,60
AGM-RCKA	Basic Stand for Wifi Digital Microscopes	each	124,00
AGM-RCKBL	Back Light stand for Wifi Digital Microscopes	each	131,00
AGM-RCKX	Vertical Moving Stand for Wifi Digital Microscopes	each	155,00
AGM-RCKY	Vertical & Horizontal Moving Stand for Wifi Digital Microsco	each	167,70
AGMT-25	Power LED Illuminated Magnifier	each	9,50
AGMT4200H	Binocular infinity corrected brightfield microscope with 6V 30W halogen illum	each	2.339,70
AGMT4200L	Binocular infinity corrected brightfield microscope	each	2.339,70
AGMT4210H	Binocular infinity corrected phase contrast microscope with 6V 30W halogen i	each	3.275,60
AGMT4210L	Binocular infinity corrected phase contrast microscope	each	3.275,60
AGMT4300H	Trinocular infinity corrected brightfield microscope with 6V 30W halogen illum	each	2.623,30
AGMT4300L	Trinocular infinity corrected brightfield microscope	each	2.623,30

Part #	Description	Unit	Price €
AGMT4310H	Triocular infinity corrected phase contrast microscope with 6V 30W halogen il	each	3.545,00
AGMT4310L	Triocular infinity corrected phase contrast microscope	each	3.545,00
AGMT5200H	Binocular infinity corrected brightfield microscope with 6V 30W halogen Koeh	each	3.112,50
AGMT5200L	Binocular infinity corrected brightfield microscope	each	3.112,50
AGMT5210H	Phase contrast binocular microscope with 6V 30W halogen Koehler illuminato	each	4.672,30
AGMT5210L	Phase contrast binocular microscope	each	4.672,30
AGMT5300H	Trinocular infinity corrected brightfield microscope with 6V 30W halogen Koel	each	3.403,20
AGMT5300L	Trinocular infinity corrected brightfield microscope	each	3.403,20
AGMT5310H	Phase contrast trinocular microscope with 6V 30W Halogen Koehler illuminatc	each	4.991,30
AGMT5310L	Phase contrast trinocular microscope	each	4.991,30
AGMT6120-SP	Binocular polarising asbestos microscope	each	7.089,90
AGMT6520-SP1	Binocular phase contrast microscope	each	4.246,90
AGMT6820	Binocular asbestos combo microscope	each	8.500,80
AGMT6830	Trinocular asbestos combo microscope	each	9.343,10
AGMT7000H	Binocular metallurgical microscope for reflected light observation with 6V 30V	each	6.522,70
AGMT7000L	Binocular metallurgical microscope for reflected light observation	each	6.522,70
AGMT7100H	Trinocular metallurgical microscope for reflected light observation with 6V 30V	each	6.792,10
AGMT7100L	Trinocular metallurgical microscope for reflected light observation	each	6.792,10
AGMT7520H	Binocular brightfield/darkfield metallurgical microscope	each	7.657,10
AGMT7530H	Trinocular brightfield/darkfield metallurgical microscope	each	7.926,50
AGMT8000	Binocular metallurgical microscope for transmitted and reflective light	each	7.082,80
AGMT8100	Trinocular metallurgical microscope for transmitted and reflective light	each	7.366,40
AGMT8520	Binocular brightfield/darkfield metallurgical microscope	each	9.422,50
AGMT8530	Trinocular brightfield/darkfield metallurgical microscope	each	10.294,80
AGNBL12001E	12000g Nimbus Balance, external cal.	each	#N/B
AGNBL124E	120g Nimbus Balance, external cal.	each	#N/B
AGNBL124I	120g Nimbus Balance, internal cal.	each	#N/B
AGNBL16001E	16000g Nimbus Balance, external cal.	each	#N/B
AGNBL1602E	1600g Nimbus Balance, external cal.	each	#N/B
AGNBL1602I	1620g Nimbus Balance, internal cal.	each	#N/B
AGNBL164E	160g Nimbus Balance, external cal.	each	#N/B
AGNBL164I	160g Nimbus Balance, internal cal.	each	#N/B
AGNBL214E	210g Nimbus Balance, external cal.	each	#N/B
AGNBL214I	210g Nimbus Balance, internal cal.	each	#N/B
AGNBL22001E	22000g Nimbus Balance, external cal.	each	#N/B
AGNBL223E	220g Nimbus Balance, external cal.	each	#N/B
AGNBL223I	220g Nimbus Balance, internal cal.	each	#N/B
AGNBL254E	250g Nimbus Balance, external cal.	250	#N/B
AGNBL254I	250g Nimbus Balance, internal cal.	each	#N/B
AGNBL2602E	2600g Nimbus Balance, external cal.	each	#N/B
AGNBL2602I	2600g Nimbus Balance, internal cal.	each	#N/B
AGNBL3602E	3600g Nimbus Balance, external cal.	each	#N/B
AGNBL3602I	3600g Nimbus Balance, internal cal.	each	#N/B
AGNBL4201E	4200g Nimbus Balance, external cal.	each	#N/B
AGNBL423E	420g Nimbus Balance, external cal.	each	#N/B
AGNBL423I	420g Nimbus Balance, internal cal.	each	#N/B
AGNBL4602E	4600g Nimbus Balance, external cal.	each	#N/B
AGNBL4602I	4600g Nimbus Balance, internal cal.	each	#N/B
AGNBL6201E	6200g Nimbus Balance, external cal.	each	#N/B
AGNBL623E	620g Nimbus Balance, external cal.	each	#N/B
AGNBL623I	620g Nimbus Balance, internal cal.	each	#N/B
AGNBL8201E	8200g Nimbus Balance, external cal.	each	#N/B
AGNBL823E	820g Nimbus Balance, external cal.	each	#N/B
AGNBL823I	820g Nimbus Balance, internal cal.	each	#N/B
AGNBL84E	80g Nimbus Balance, external cal.	each	#N/B
AGNBL84I	80g Nimbus Balance, internal cal.	each	#N/B
AGO6002	Stereoviewer x 4	each	78,20
AGO6091	100mm Glass Scale with Magnifier	each	407,40

Part #	Description	Unit	Price €
AGO613A	Eyeglass X10 plastic	each	#N/B
AGO6170	12W Daylight Table Magnifier Lamp	each	#N/B
AGO6170T	Spare Tube for AGO6170 Magnifier Lamp	each	#N/B
AGO6172	Illuminated magnifiers x1.75, circular lens, standard daylig	each	144,00
AGO6172T	Spare lamp for AGO6172 Magnifier	each	25,20
AGO6173N	Rectangular Long Reach Magnifier Lamp	each	181,80
AGO6173NT	Spare Tube for Rectangular Magnifier Lamp	each	26,20
AGO6173T	Spare Tube for AGO6173 Magnifier Lamp	each	34,50
AGO6175N	Long Reach Slim Line Magnifier Lamp	each	#N/B
AGO6175NT	Spare Tube for Long Reach Magnifier Lamp	each	#N/B
AGO6176N	Classic Round Magnifier Lamp	each	143,80
AGO6176NT	Spare Tube for Round Magnifier Lamp	each	23,00
AGO6177	Bench Bracket for AGO6171 / AGO6172 / AGO6173 Magnifier Lamp	each	45,60
AGO6178	Cordless LED Folding Lamp	each	87,20
AGO6270	Red/Blue Stereo Spectacles Rubber Frame	each	20,20
AGO6271	Red/Blue Clip-On Glasses	each	30,70
AGO6276	Polarised Stereo Spectacles Plastic Frame	each	10,60
AGO6277	Polarised Stereo Spectacles Ophthalmic Frames	each	30,70
AGO6326	Measuring Magnifier x8	each	200,20
AGO6327	Horizontal Graticule for Measuring Magnifier. 25mm, 0.1mm di	each	96,90
AGO6328	Crossed Graticule for AGO6326 Measuring Magnifier	each	95,50
AGO6330	Measuring Magnifier X4	each	398,90
AGO6331	Horizontal Graticule for AGO6330 Measuring Magnifier	each	141,40
AGO6332	Crossed Graticule for AGO6330 Measuring Magnifier	each	167,80
AGO6400	Fibre Optic Illuminator 150w 230/240v	each	1.511,60
AGO6401	Swan-Neck Light Guide - 1 Branch	each	308,90
AGO6402	Swan-Neck Light Guide - 2 Branch	each	614,00
AGO6403	Swan-Neck Light Guide - 3 Branch	each	753,20
AGO642	Red/Green Stereo Spectacles	50	39,70
AGO6421	150W Halogen Light Bulb	each	49,50
AGO6460-1	58mm Six Point Ringlight (900mm Length)	each	o.a.
AGO6460-2	66mm Six Point Ringlight (900mm Length)	each	o.a.
AGO6460-4	40mm Annular Ringlight	each	o.a.
AGO6461	Fibre optic illuminator KL300LED	each	496,30
AGO6500	Viewing Box A4 (234mm x 305mm)	each	#N/B
AGO6501	Viewing Box A3 (305mm x 459mm)	each	#N/B
AGO6502	Viewing Box A2 (459mm x 605mm)	each	#N/B
AGO6503	Viewing Box A1 (605mm x 904mm)	each	#N/B
AGO6504	Viewing Box A0 (904mm x 1226mm)	each	#N/B
AGO652	Glass measuring scale	each	484,20
AGO6521	Swan Neck Light Guide 1 Branch	each	o.a.
AGO6522	Swan Neck Light Guide 2 Branch	each	o.a.
AGO6523	Flexible Light Guide 2 Branch	each	o.a.
AGO6527	Flexible Light Guide 1 Branch	each	o.a.
AGO653	150mm Magscale	each	851,20
AGO654	310mm Single Reader Magscale	each	1.245,50
AGO657	450mm Magscale	each	1.921,30
AGO660	Polarising Spectacles Cardboard Mount	5	12,80
AGO664	Red/Green Stereo Spectacles with Arms	20	23,90
AGO665	Red/Blue Stereo Viewer	50	#N/B
AGO665A	Red/Blue Stereo Viewer	10	14,00
AGO668	Red/Blue Stereo Viewers with Arms	20	19,10
AGO751	Measuring Magnifier, 7x, w/Standard Metric Scale	each	107,10
AGO767	HPRL - High Power Ringlight for Microscopy	each	#N/B
AGO767-1	Polarizing filter set for HPRL	each	#N/B
AGO767-2	Diffuser for HPRL	each	#N/B
AGO767-A1	Ringlight Adapter 66/58 dia	each	#N/B
AGO767-A2	Ringlight Adapter 66/60 dia	each	#N/B

Part #	Description	Unit	Price €
AGO767-A3	Ringlight Adapter 66/62 dia	each	#N/B
AGO767-A4	Ringlight Adapter M48 x 0.75 dia	each	#N/B
AGO767-A5	Ringlight Adapter M55 x 0.75 dia	each	#N/B
AGO768	High Power LED Light Source for Microscopy	each	#N/B
AGO801	Headband Magnifier (x 1.5)	each	116,40
AGO802	Headband Magnifier (x 1.75)	each	#N/B
AGO803	Headband Magnifier (x 2)	each	116,40
AGO804	Headband Magnifier (x 2.5)	each	116,40
AGO805	Headband Magnifier (x 2.75)	each	116,40
AGO806	Headband Magnifier (x 3.5)	each	116,40
AGO807	Optiloupe (Auxiliary Lens)	each	#N/B
AGO808	Lens Plate only (x 1.5)	each	73,20
AGO809	Lens Plate only (x 1.75)	each	73,20
AGO810	Lens Plate only (x 2)	each	76,60
AGO811	Lens Plate only (x 2.5)	each	76,60
AGO812	Lens Plate only (x 2.75)	each	85,10
AGO813	Lens Plate only (x 3.5)	each	85,10
AGO817	LED Rechargeable Task Lamp with Dimmer & USB	each	33,10
AGO818	LED Magnifier Lamp with Warm to Cool light dimmer & 2 lenses	each	206,40
AGP9104	Ilfostop Stop Bath	500ml	20,00
AGP9105	Ilford Wash Aid (1 Litre)	1l	#N/B
AGP9106	* Ilford Phenisol Developer	5l	101,10
AGP9145A	Kodak Photo-Flo Solution	473ml	29,30
AGP9147	Carestream 4489 3.25inchx 4inch Multipack	250	#N/B
AGP9150	Carestream SO163 Film 3.25inch x 4inch	250	#N/B
AGP9152	Ilford Developer PQ	5l	59,20
AGP9154A	Ilford PQ Universal Developer	8 x 500ml	#N/B
AGP9156	Kodak HC-110 Developer	1l	#N/B
AGP9159	Multigrade Developer	1l	30,60
AGP9166A	* Indicator Stop Bath 473ml to Make	1l	#N/B
AGP9169A	Kodak D76 Developer (to Make 3.8 Litre)	each	28,80
AGP9183	Ilford Hypam Fixer	5l	88,30
AGP9190	Ilfotol Wetting Agent	1l	43,20
AGP9252	Ilfospeed 1m Grade 2 12.7cm x 17.8cm	100	61,70
AGP9253	Ilfospeed 1m Grade 3 12.7cm x 17.8cm	100	#N/B
AGP9262	Ilfospeed 1m Grade 2 20.3cm x 25.4cm	100	#N/B
AGP9263	Ilfospeed 1m Grade 3 20.3cm x 25.4cm	100	#N/B
AGP9266	Ilfospeed Multigrade Paper 20.3 x 25.4 Cm	100	#N/B
AGP9281	Ilford Nuclear Emulsion Type K5	50ml	#N/B
AGP9281-1	Ilford Nuclear Emulsion Type K5	100ml	#N/B
AGP9282	Ilford Nuclear Emulsion Type L4	50ml	#N/B
AGP9282-1	Ilford Nuclear Emulsion Type L4	100ml	#N/B
AGP9283	Ilford Nuclear Emulsion Type K2	50ml	#N/B
AGP9283-1	Ilford Nuclear Emulsion Type K2	100ml	#N/B
AGP9284	Kodak NTB emulsion for autoradiography	each	2.651,80
AGP9285	Ilford Nuclear Emulsion K5D, Gel Form	100ml	#N/B
AGP940	Graduated Cylinder Beaker 42 oz	1,2kg	#N/B
AGP941	Mixing Jug 2 Litre	each	#N/B
AGP942	Filter Funnel 5inch	each	7,70
AGP9463	Digital thermometer	each	45,70
AGP960	Carestream SO163 Film 6.5cm x 9cm	100	#N/B
AGP968	Kodak 4489 Film 9 x 12cm	100	265,80
AGP972	Carestream 4489 Film 6.5cm x 9cm	100	#N/B
AGP973	Kodak 4489 Film 3 1/4inch x 3 5/16inch	100	o.a.
AGP976	General Purpose Thermometer	each	17,50
AGP9798	Scotch PMA rolls	each	#N/B
AGP986	* Photomount Aerosol Spray 260g	each	#N/B
AGP987	* Spraymount Aerosol	200ml	#N/B

Part #	Description	Unit	Price €
AGPAN315	Duct Leak Tester (+/-) with NO inst 230V	each	#N/B
AGPAN315-110	Duct Leak Tester (+/-) with NO inst 110V	each	#N/B
AGPAN341	Duct Leak Tester (+/-) with inst 230V	each	#N/B
AGPAN341-110	Duct Leak Tester (+/-) with inst 110V	each	#N/B
AGPH730	Micromanometer Kit	each	#N/B
AGPH731	QAPH731 2ft x 2ft Capture Hood Kit	each	#N/B
AGPH731-STA	QAPH731 with Capture Hood Stand and tabl	each	#N/B
AGPMB202	Moisture Analyzer	200g	#N/B
AGPMB53	Moisture Analyzer	50g	#N/B
AGPVM610	PVM610 Micromanometer	each	#N/B
AGPVM620	PVM620 Micromanometer	each	#N/B
AGPX124	Pioneer Analytical Balance PX124 120g x 0.1mg	each	#N/B
AGPX124-E	Pioneer Analytical Balance PX124/E 120g x 0.1mg	each	#N/B
AGPX224	Pioneer Analytical Balance PX224 220g x 0.1mg	each	#N/B
AGPX224-E	Pioneer Analytical Balance PX224/E 220g x 0.1mg	each	#N/B
AGPX3202	Pioneer Precision Balance PX3202 3200g x 0.01g	each	#N/B
AGPX3202-E	Pioneer Precision Balance PX3202/E 3200g x 0.01g	each	#N/B
AGPX323	Pioneer Precision Balance PX323 320g x 1mg	each	#N/B
AGPX323-E	Pioneer Precision Balance PX323/E 320g x 1mg	each	#N/B
AGPX4201	Pioneer Precision Balance PX4201 4200g x 0.01g	each	#N/B
AGPX4201-E	Pioneer Precision Balance PX4201/E 4200g x 0.01g	each	#N/B
AGPX5202	Pioneer Precision Balance PX5202 3200g x 0.01g	each	#N/B
AGPX5202-E	Pioneer Precision Balance PX5202/E 5200g x 0.01g	each	#N/B
AGPX523	Pioneer Precision Balance PX523 520g x 1mg	each	#N/B
AGPX523-E	Pioneer Precision Balance PX523/E 520g x 1mg	each	#N/B
AGR1000	UA-Zero™ EM Stain	25ml	82,30
AGR1009	Glutaraldehyde 25% Practical Grade	500ml	45,50
AGR1010	Glutaraldehyde 25% EM Grade	500ml	72,00
AGR1011	Glutaraldehyde 25% EM Grade	250ml	41,70
AGR1012	Glutaraldehyde 25% EM Grade	100ml	23,30
AGR1013	* Glutaraldehyde 70% Double Distilled ampoule	5ml	13,00
AGR1014	* Glutaraldehyde 70% Double Distilled ampoule	5ml	9,90
AGR1015	* Osmium Tetroxid ampoules	10 x 0.1g	142,20
AGR1016	* Osmium Tetroxide ampoules	10 x 0.25g	192,10
AGR1017	* Osmium Tetroxide ampoule	1g	85,50
AGR1018	* Paraformaldehyde	500g	18,80
AGR1019	* Osmium Tetroxide 2% Solution	100ml	216,10
AGR1020	*Glutaraldehyde 25% EM Grade ampoule	10 x 10ml	61,70
AGR1021	* Osmium Tetroxide 2% Solution ampoules	5 x 2ml	71,10
AGR1022	* Osmium Tetroxide 2% Solution ampoules	5 x 5ml	78,00
AGR1023	* Osmium Tetroxide 4% Solution ampoules	5 x 2ml	81,80
AGR1024	* Osmium Tetroxide 4% Solution ampoules	5 x 5ml	129,10
AGR1026	* Paraformaldehyde 16% Sol methanol-free	10 x 10ml	70,50
AGR1027	Ruthenium Tetroxide	5 x 10ml	591,30
AGR1030	* Epoxy resin Araldite CY212 Kit	each	146,40
AGR1031	* Agar 100 Resin Kit	each	281,00
AGR1034	* Lowicryl HM20 Kit	each	1.354,50
AGR1035	* Lowicryl K4M Kit	each	1.230,00
AGR1036	* Glycol Methacrylate Resin Kit	each	96,60
AGR1040	Araldite CY212 Resin	500g	49,20
AGR1041	Araldite CY212 Resin	250g	29,70
AGR1042	Araldite CY212 Resin	100g	13,00
AGR1043	Agar 100 Resin	500g	122,60
AGR1044	Agar 100 Resin	250g	66,70
AGR1045	Agar 100 Resin	100g	31,50
AGR1047R	ERL 4221	225ml	64,20
AGR1049	* Methyl Methacrylate	500g	161,50
AGR1050	* N-Butyl Methacrylate	500g	85,90

Part #	Description	Unit	Price €
AGR1051	DDSA (Dodecyl Succinic Anhydride) EM grade redistilled	500g	76,10
AGR1052	DDSA (Dodecyl Succinic Anhydride) EM grade redistilled	250g	40,00
AGR1053	DDSA (Dodecyl Succinic Anhydride) EM grade redistilled	100g	16,80
AGR1054	NSA (Nonenyl Succinic Anhydride) redistilled	500g	89,70
AGR1055	NSA (Nonenyl Succinic Anhydride) redistilled	250g	44,60
AGR1056	NSA (Nonenyl Succinic Anhydride) redistilled	100g	20,70
AGR1058	* Benzoyl Peroxide	100g	95,20
AGR1059	* Benzoyl Peroxide	25g	28,10
AGR1060	* BDMA Benzyl Dimethylamine	500g	73,50
AGR1061	* BDMA Benzyl Dimethylamine	250g	44,50
AGR1062	* BDMA Benzyl Dimethylamine	100g	18,30
AGR1063	* DMP-30 Tri-Dimethylaminomethyl Phenol	500g	35,60
AGR1064	* DMP-30 Tri-Dimethylaminomethyl Phenol	250g	19,80
AGR1065	* DMP-30 Tri-Dimethylaminomethyl Phenol	100g	10,20
AGR1066	* S-1 Dimethylaminoethanol	500g	164,00
AGR1067	* S-1 Dimethylaminoethanol	100g	42,70
AGR1068	Quetol 651	500g	252,10
AGR1070	* Dibutyl Phthalate	500g	41,00
AGR1071	* Dibutyl Phthalate	100g	10,30
AGR1072	DER 736 Diglycidylether of Polypropyleneglycol	500g	57,40
AGR1073	DER 736 Diglycidylether of Polypropyleneglycol	250g	31,50
AGR1074	DER 736 Diglycidylether of Polypropyleneglycol	100g	13,80
AGR1075	* Lowicryl K4M MonoStep Kit 225g	225g	541,40
AGR1076	* Lowicryl HM20 MonoStep Kit 225g	225g	559,10
AGR1078	* Agar Low Viscosity Resin Kit	each	260,00
AGR1080	* Propylene Oxide	500ml	45,20
AGR1081	* MNA Methyl Nadic Anhydride redistilled	500g	47,80
AGR1082	* MNA Methyl Nadic Anhydride redistilled	250g	24,80
AGR1083	* MNA Methyl Nadic Anhydride redistilled	100g	11,90
AGR1084	GMA (2-Hydroxyethyl Methacrylate)	500g	62,60
AGR1085	GMA (2-Hydroxyethyl Methacrylate)	250g	33,80
AGR1086	GMA (2-Hydroxyethyl Methacrylate)	100g	14,20
AGR1087	Carbowax 400 (Polyethylene Glycol 400)	500ml	41,10
AGR1088	* N N-Dimethyl Aniline	100g	23,00
AGR1089	* 2-Butoxyethanol	500g	41,20
AGR1090	Araldite RD-2 Resin	500g	74,80
AGR1091	Araldite RD-2 Resin	250g	39,20
AGR1092	Araldite RD-2 Resin	100g	16,60
AGR1098A	* Xylene	1l	69,40
AGR1099	Carbowax 200 (Polyethylene Glycol 200)	500ml	68,50
AGR1100	* S-Collidine	100ml	146,80
AGR1101	* S-Collidine	25ml	40,80
AGR1102	* Sodium Cacodylate	500g	563,80
AGR1103	* Sodium Cacodylate	250g	336,40
AGR1104	* Sodium Cacodylate	100g	152,30
AGR1105	* Sodium Cacodylate	25g	44,90
AGR1106	Tri-Sodium Citrate EM	500g	94,50
AGR1107	Tri-Sodium Citrate EM	100g	#N/B
AGR1121	Hepes	25g	95,60
AGR1130	* JB-4 Resin Embedding Kit	each	589,50
AGR1131	* Histocryl Resin Kit	each	168,80
AGR1133	* Lowicryl K11M Kit	each	1.098,20
AGR1134	* Lowicryl HM23 Kit	each	982,20
AGR1138	* JB-4 Mini Embedding Kit	each	#N/B
AGR1140	* Agar 100 Premix Kit - Hard	each	164,70
AGR1141	* Agar 100 Premix Kit - Medium	each	164,70
AGR1142	* Agar 100 Premix Kit - Soft	each	155,40
AGR1143	* Araldite CY212 Premix Kit - Hard	each	99,30

Part #	Description	Unit	Price €
AGR1144	* Araldite CY212 Premix Kit - Medium	each	99,30
AGR1145	* Araldite CY212 Premix Kit - Soft	each	93,70
AGR1150	Acridflavin	25g	95,90
AGR1153	Adenosine 5 Monophosphate	1g	17,50
AGR1154	Adenosine 5 Triphosphate	1g	46,80
AGR1155	Ammonium Molybdate	100g	o.a.
AGR1156	Ammonium Molybdate	25g	127,50
AGR1160	Tris Buffer	500g	182,80
AGR1161	Tris Buffer	250g	95,40
AGR1162	Tris Buffer	100g	41,30
AGR1164	* Agar Low Viscosity Premix Kit Hard	each	172,40
AGR1165	* Agar Low Viscosity Premix Kit Medium	each	172,40
AGR1166	* Agar Low Viscosity Premix Kit Soft	each	162,70
AGR1170	Sodium Di-Hydrogen Orthophosphate	500g	156,60
AGR1172	Di-Sodium Hydrogen Orthophosphate	500g	102,50
AGR1180	OCT Compound	118ml	31,00
AGR1190	Cytochrome C	25mg	30,50
AGR1195	Dimethyl Sulphoxide (DMSO)	500ml	93,10
AGR1200	Ferritin, cadmium free	100mg	157,90
AGR1201	Formvar	100g	33,70
AGR1202	Formvar	25g	9,00
AGR1205	Indium Trichloride	10g	234,80
AGR1207	Lanthanum Nitrate	25g	97,50
AGR1209	* Lead Acetate	250g	173,30
AGR1210	* Lead Citrate	50g	72,00
AGR1212	Phosphomolybdic Acid	25g	66,80
AGR1213	* Phosphotungstic Acid	25g	109,30
AGR1214	* Potassium Pyroantimonate	100g	119,30
AGR1215	* Potassium Permanganate	100g	24,60
AGR1217	* Lead Nitrate	500g	111,00
AGR1220	Peroxidase Grade 1	10mg	446,30
AGR1221	Peroxidase Grade 2	100mg	429,70
AGR1222	Peroxidase Grade 2	10mg	59,00
AGR1226	Ruthenium Red	1g	332,70
AGR1227	Ruthenium Red	100mg	45,30
AGR1228	* Hexamethyldisilazane HMDS	100ml	85,80
AGR1230	Sodium Silicotungstate	25g	220,70
AGR1240	* Thiocarbohydrazide	1g	36,40
AGR1241	* Thiosemicarbazide	25g	20,50
AGR1260A	* Uranyl Acetate	25g	389,60
AGR1260C	* Platinum Blue Uranyl Acetate Alternative	1ml	173,40
AGR1260D	EM Stain 336 Uranyl Acetate Alternative	25ml	131,70
AGR1272	* Silver Paint Diluent	100ml	22,30
AGR1275	Pioloform Powder	10g	19,30
AGR1276	Butvar B98 Powder	25g	20,60
AGR1277	Silicon Monoxide	10g	75,80
AGR1278A	Indicating Silica Gel	1kg	164,20
AGR1280	London Resin (Hard)	500g	138,10
AGR1280A	LR White Resin (Hard) 500g Uncatalysed	500g	134,60
AGR1281	London Resin (Medium)	500g	139,20
AGR1281A	LR White Resin (Medium) Uncatalysed	500g	135,90
AGR1281C	Catalyst for London Resin	9,9g	23,60
AGR1282	London Resin (Soft)	500g	146,00
AGR1282A	LR White Resin (Soft) Uncatalysed	500g	146,00
AGR1283	* Accelerator for London Resin	10ml	25,40
AGR1284	LR Gold Resin	500g	173,00
AGR1285	Activator - Benzil	50g	93,30
AGR1286	Polyvinylpyrrolidone	100g	93,30

Part #	Description	Unit	Price €
AGR1287	HAZ Benzoyl Peroxide (stabilised 50%)	50g	47,40
AGR1310	* Glutaraldehyde 25% Vacuum Distilled	500ml	265,60
AGR1311	Glutaraldehyde 25% Vacuum Distilled	100ml	69,80
AGR1312	Glutaraldehyde 25% Vacuum Distilled	10 x 10ml	115,80
AGR1313	* Glutaraldehyde 50% Vacuum Distilled	100ml	120,80
AGR1314	* Glutaraldehyde 50% Vacuum Distilled	10 x 10ml	192,70
AGR1320	Citifluor Glycerol Pbs Solution AF1	25ml	66,20
AGR1320A	Citifluor Glycerol Pbs Solution AF1	100ml	178,20
AGR1321	Citifluor Glycerol Solution AF2	25ml	64,50
AGR1321A	Citifluor Glycerol Solution AF2	100ml	162,70
AGR1322	Citifluor Phosphate Buffered Saline AF3	25ml	63,90
AGR1322A	Citifluor Phosphate Buffered Saline AF3	100ml	181,00
AGR1323	Citifluor non-fluorescent Immersion Oil	10ml	55,20
AGR1324	Citifluor non-fluorescent Immersion Oil	25ml	92,00
AGR1325	PVA Solution 25ml	25ml	53,50
AGR1326	Citifluor Solid Mountant Kit	each	105,80
AGR1327	Citifluor AF4 Mountant Solution	25ml	53,50
AGR1328-1	Citifluor AF100 Antifadent Solution	5ml	12,70
AGR1328-2	Citifluor AF200 Antifadent Solution	10ml	91,60
AGR1328-3	Citifluor AF300 Antifadent Solution	10ml	53,50
AGR1329	Aquatex	50ml	107,10
AGR1330-10	Citifluor CFM-1 Mountant Solution	10ml	43,10
AGR1330-25	Citifluor CFM-1 Mountant Solution	25ml	100,10
AGR1330A-10	Citifluor CFM-1 Mountant Soln. + Antifadent	10ml	53,70
AGR1330A-25	Citifluor CFM-1 Mountant Soln. + Antifadent	25ml	108,20
AGR1331-10	Citifluor CFM-2 Mountant Solution	10ml	40,10
AGR1331-25	Citifluor CFM-2 Mountant Solution	25ml	78,20
AGR1332-10	Citifluor CFM-3 Mountant Solution	10ml	48,10
AGR1332-25	Citifluor CFM-3 Mountant Solution	25ml	93,60
AGR1333-1	Citifluor UVM-1 Mountant Solution	10ml	#N/B
AGR1333-2	Citifluor UVM-2 Mountant Solution	10ml	39,10
AGR1337	* Ethylene Dichloride	500ml	168,00
AGR1339	* Eukitt	100ml	73,40
AGR1340	* DPX contains xylene	500ml	64,40
AGR1343	* Canada Balsam In Xylene	100ml	113,10
AGR1344	* Euparal Mountant	50ml	56,70
AGR1344A	* Euparal Mountant	100ml	98,70
AGR1344B	* Euparal Thinner	100ml	101,10
AGR1346	* Solvent for Polymerised Resin	50ml	37,80
AGR1347	* Solvent for Unpolymerised Resin	500ml	41,30
AGR1348	* Periodic Acid	100g	#N/B
AGR1350	Polyvinyl Alcohol Mol. Wt. 30,000 - 70,000	250g	85,30
AGR1351	* Histomount	100ml	50,40
AGR1352	* Histomount	450ml	96,40
AGR1353	* Histo-Clear II	3,8l	197,20
AGR1357	Omnimount	100ml	63,00
AGR1370	Agar Low Viscosity Resin	500g	115,40
AGR1371	Agar Low Viscosity Resin	250g	67,60
AGR1372	Agar Low Viscosity Resin	100g	32,60
AGR1373	Low Viscosity Hardener VH1	500g	79,90
AGR1374	Low Viscosity Hardener VH1	250g	48,00
AGR1375	Low Viscosity Hardener VH1	100g	20,70
AGR1376	* Low Viscosity Hardener VH2	500g	48,00
AGR1377	* Low Viscosity Hardener VH2	250g	24,90
AGR1378	* Low Viscosity Hardener VH2	100g	12,00
AGR1379	* Low Viscosity Accelerator	100ml	14,80
AGR1380	* Low Viscosity Accelerator	50ml	9,80
AGR1381	* Low Viscosity Accelerator	5 x 2,5ml	21,50

Part #	Description	Unit	Price €
AGR1500	HAZ Sodium Cacodylate Buffer 0.2M, PH 7.4 225ml		37,60
AGR1501	HAZ Sodium Cacodylate Buffer 0.2M, PH 7.4 500ml		69,50
AGR1502	HAZ Sodium Cacodylate Buffer 0.2M, PH 7.4 1ltr		100,90
AGR1503	HAZ Sodium Cacodylate Buffer 0.2M, PH 7.2 500ml		77,80
AGR1504	HAZ Sodium Cacodylate Buffer 0.2M, PH 7.2 1ltr		133,70
AGR1550	HAZGlutaraldehyde 0.625% 0.1M Phosphate Buffer pH 7.4 500ml		170,30
AGR1551	HAZ Glutaraldehyde 0.625% 0.1M Phosphate Buffer pH 7.4 1ltr		298,10
AGR1552	HAZ Glutaraldehyde 2% in 0.1M Phosphate Buffer, pH 7.4 500ml		170,30
AGR1553	HAZ Glutaraldehyde 2% in 0.1M Phosphate Buffer, pH 7.4 1lt		306,80
AGR1554	HAZGlutaraldehyde 2.5% in 0.1M Phosphate Buffer pH 7.4 500ml		177,60
AGR1556	HAZ Glutaraldehyde 2.5% in 0.1M Phosphate Buffer pH 7.4 1ltr		323,30
AGR1700	Acid Fuchsin	25g	78,00
AGR1701	Acridine Orange	25g	307,20
AGR1702	Alcian Blue 8GX	25g	282,30
AGR1703	Aniline Blue (Alcohol Soluble)	25g	69,50
AGR1704	Aniline Blue (Water Soluble)	25g	102,10
AGR1705	Azur II	25g	215,50
AGR1706	Pararosaniline	25g	94,10
AGR1708	Crystal Ponceau	25g	81,90
AGR1709	Eosin (yellowish - Water Soluble)	25g	365,60
AGR1710	Eosin (Spirit Soluble)	25g	133,20
AGR1711	Fast Green FCF	25g	229,70
AGR1712	Fluorescein	25g	91,80
AGR1713	Giemsa's Stain	25g	36,70
AGR1716	Light Green SF	25g	219,70
AGR1717	Luxol Fast Blue MBS	25g	102,00
AGR1718	Methyl Green	25g	170,80
AGR1719	Methylene Blue	25g	74,90
AGR1720	Neutral Red	25g	154,90
AGR1721	Orange G	25g	101,90
AGR1722	Ponceau de Xylidine	25g	116,80
AGR1723	Pyronin Y	10g	543,10
AGR1724	Safranin O	25g	173,40
AGR1725	Sudan Black B	25g	110,60
AGR1726	Tartrazine	25g	29,50
AGR1727	Toluidine Blue	10g	42,50
AGR1728	* Van Gieson's Stain	500ml	72,00
AGR1729	Alizarin Red	25g	76,30
AGR1730	Brilliant Cresyl Blue	25g	92,70
AGR1731	Brilliant Green	25g	227,40
AGR1732	Carmine	25g	227,00
AGR1733	Congo Red	25g	67,20
AGR1734	Methyl Blue	25g	212,20
AGR1735	Nigrosine (Water Soluble)	25g	75,00
AGR1736	Oil Red O	25g	46,00
AGR1737	Scarlet R (Sudan Iv)	25g	32,60
AGR19478	* Mount-Quick, mtg Medium in a Tube	30ccm	22,70
AGS100	Perforated carbon film, 3.05mm grid	each	36,60
AGS104	2160 Lines/mm Diffraction Grating Replica 3mm	each	53,20
AGS106	Cross Grating Replica 3mm	each	62,30
AGS106L	Diffraction grating replica with latex spheres	each	94,60
AGS108	Evaporated Aluminium Film 3mm	each	26,70
AGS110	Evaporated Thallous Chloride 3mm	each	32,40
AGS112	Image rotation, Molybdenum trioxide crystals	each	32,40
AGS114	Evaporated platinum / iridium, 3.05mm grid	each	27,80
AGS116-1	Ultra Smooth Carbon Foils 2.5nm Thick	2	551,90
AGS116-10	Ultra Smooth Carbon Foils 40.0nm Thick	2	513,30
AGS116-11	Ultra Smooth Carbon Foils 50.0nm Thick	2	513,30

Part #	Description	Unit	Price €
AGS116-2	Ultra Smooth Carbon Foils 3.5nm Thick	2	551,90
AGS116-3	Ultra Smooth Carbon Foils 5.0nm Thick	2	551,90
AGS116-4	Ultra Smooth Carbon Foils 7.0nm Thick	2	513,30
AGS116-5	Ultra Smooth Carbon Foils 10.0nm Thick	2	551,90
AGS116-6	Ultra Smooth Carbon Foils 15.0nm Thick	2	513,30
AGS116-7	Ultra Smooth Carbon Foils 20.0nm Thick	2	513,30
AGS116-8	Ultra Smooth Carbon Foils 25.0nm Thick	2	513,30
AGS116-9	Ultra Smooth Carbon Foils 30.0nm Thick	2	572,70
AGS117-1	Quantifoil 7x7um (2um Bar) 200 Mesh Copper	10	131,20
AGS117-2	Quantifoil 7x7um (2um Bar) 300 Mesh Copper	10	131,20
AGS117-3	Quantifoil 7x7um (2um Bar) 400 Mesh Copper	10	134,40
AGS117-4	Quantifoil 7x7um (2um Bar) 200 Mesh Nickel	10	149,60
AGS117-5	Quantifoil 7x7um (2um Bar) 300 Mesh Nickel	10	149,60
AGS117-6	Quantifoil 7x7um (2um Bar) 400 Mesh Nickel	10	149,60
AGS117-7	Quantifoil 7x7um (2um Bar) 200 Mesh Gold	10	149,60
AGS117-8	Quantifoils 7x7um (2um Bar) 300 Mesh Gold	10	143,80
AGS117-9	Quantifoils 7x7um (2um Bar) 400 Mesh Gold	10	149,60
AGS118	Plane spacing 0.56nm, potassium chloroplatinate crystals, 3.05mm grid	each	50,40
AGS122	Plane spacing 0.9nm and 0.45nm, Crocidolite crystals, 3.05mm grid	each	27,80
AGS124	Catalase crystals, 3.05mm grid	each	58,00
AGS126	Ferritin, 3.05mm grid	each	30,30
AGS128A-1	Pd/Pt Shadowed 0.120µm Polystyrene Latex	each	54,30
AGS128A-4	Pd/Pt Shadowed 0.204um Polystyrene Latex	each	54,30
AGS128A-7	Pd/Pt Shadowed 0.945µm Polystyrene Latex	each	50,30
AGS128B	Gold shadowed latex particles	each	44,30
AGS130-02	Polystyrene Latex Particles 0.03µm	10ml	36,30
AGS130-04	Polystyrene Latex Particles 0.08um	10ml	35,10
AGS130-06	Polystyrene Latex Particles 0.09µm	10ml	35,10
AGS130-1	Polystyrene Latex Particles 0.095µm	5ml	30,30
AGS130-4	Polystyrene Latex Particles 0.204um	5ml	30,30
AGS130-7	Polystyrene Latex Particles 0.855µm	5ml	30,30
AGS132	High resolution test specimen gold particles on carbon film	each	75,80
AGS134	Pioloform on 3mm 200 Mesh Cu Grids	25	86,60
AGS134-1	Pioloform on 3mm 100 Mesh Cu Grids	25	139,50
AGS134-3	Pioloform on 3mm 300 Mesh Cu Grids	25	86,60
AGS134-4	Pioloform on 3mm 400 Mesh Cu Grids	25	86,60
AGS134-5	Pioloform on 2mm x 1mm Cu Slot Grids	25	119,20
AGS134-6	Pioloform on 3mm 75 Mesh Cu Grids	25	139,50
AGS134A	Pioloform on 3mm 200 Mesh Au Grids	25	129,20
AGS134A1	Pioloform on 3mm 100 Mesh Au Grids	25	116,50
AGS134A3	Pioloform on 3mm 300 Mesh Au Grids	25	86,60
AGS134A4	Pioloform on 3mm 400 Mesh Au Grids	25	94,30
AGS134A5	Pioloform on 2mm x 1mm Au Slot Grids	25	96,00
AGS134A6	Pioloform on 3mm 75 Mesh Au Grids	25	171,10
AGS134AF1	Pioloform on F1 Au finder grid	25	130,50
AGS134AH7	Pioloform on F7 Au finder grid	25	197,10
AGS134-F1	Pioloform on F1 Cu finder grid	25	104,40
AGS134-H7	Pioloform on F7 Cu finder grid	25	99,90
AGS134N	Pioloform on 3mm 200 Mesh Ni Grids	25	91,60
AGS134N1	Pioloform on 3mm 100 Mesh Ni Grids	25	80,50
AGS134N3	Pioloform on 3mm 300 Mesh Ni Grids	25	161,60
AGS134N4	Pioloform on 3mm 400 Mesh Ni Grids	25	91,60
AGS134N5	Pioloform on 2mm x 1mm Ni Slot Grids	25	176,90
AGS134N6	Pioloform on 3mm 75 Mesh Ni Grids	25	141,80
AGS134NF1	Pioloform on F1 Ni finder grid	25	107,60
AGS134NH7	Pioloform on F7 Ni finder grid	25	119,90
AGS135	Oriented Gold Crystal ,3.05mm grid Plane Spacing 0.204,0.144	each	129,90
AGS136	Plane spacing 1.0nm, Copper phthalocyanine, 3.05mm grid	each	30,30

Part #	Description	Unit	Price €
AGS138	Formvar on 3mm 200 Mesh Cu Grids	25	58,60
AGS138-1	Formvar on 3mm 100 Mesh Cu Grids	25	77,10
AGS138-3	Formvar on 3mm 300 Mesh Cu Grids	25	58,40
AGS138-4	Formvar on 3mm 400 Mesh Cu Grids	25	81,40
AGS138-6	Formvar on 3mm 75 Mesh Cu Grids	25	137,10
AGS138A	Formvar on 3mm 200 Mesh Au Grids	25	o.a.
AGS138A1	Formvar on 3mm 100 Mesh Au Grids	25	o.a.
AGS138A3	Formvar on 3mm 300 Mesh Au Grids	25	o.a.
AGS138A4	Formvar on 3mm 400 Mesh Au Grids	25	93,00
AGS138AF1	Formvar on 3mm F1 Finder Au Grids	25	o.a.
AGS138AH7	Formvar on 3mm H7 Finder Au Grids	25	o.a.
AGS138-F1	Formvar on 3mm F1 Finder Cu Grids	25	o.a.
AGS138-H7	Formvar on 3mm H7 Finder Cu Grids	25	o.a.
AGS138N	Formvar on 3mm 200 Mesh Ni Grids	25	90,10
AGS138N1	Formvar on 3mm 100 Mesh Ni Grids	25	o.a.
AGS138N3	Formvar on 3mm 300 Mesh Ni Grids	25	o.a.
AGS138N4	Formvar on 3mm 400 Mesh Ni Grids	25	o.a.
AGS138N6	Formvar on 3mm 75 Mesh Ni Grids	25	o.a.
AGS138NF1	Formvar on 3mm F1 Finder Ni Grids	25	o.a.
AGS138NH7	Formvar on 3mm H7 Finder Ni Grids	25	o.a.
AGS140	Graphitised Carbon Black , 3.05mm grid Plane Spacing 0.34nm	each	30,30
AGS142	Combined test specimen, 3.05mm grid	each	66,60
AGS143-1	Quantifoil 1.2um dia (1.3µm) 200 Mesh Copper	10	131,20
AGS143-1-100	Quantifoil 1.2um dia (1.3µm) 200 Mesh Copper	100	768,40
AGS143-2	Quantifoil 1.2um dia (1.3µm) 300 Mesh Copper	10	131,20
AGS143-2-100	Quantifoil 1.2um dia (1.3µm) 300 Mesh Copper	100	755,20
AGS143-2-GO	Graphene oxide support film on Quantifoil R1.2/1.3 on 300 mesh	each	203,70
AGS143-3	Quantifoil 1.2um dia (1.3µm) 400 Mesh Copper	10	131,20
AGS143-3-100	Quantifoil 1.2um dia (1.3µm) 400 Mesh Copper	100	830,70
AGS143-4	Quantifoil 1.2um dia (1.3µm) 200 Mesh Nickel	10	149,60
AGS143-5	Quantifoil 1.2um dia (1.3µm) 300 Mesh Nickel	10	149,60
AGS143-6	Quantifoil 1.2um dia (1.3µm) 400 Mesh Nickel	10	149,60
AGS143-7	Quantifoil 1.2um dia (1.3µm) 200 Mesh Gold	10	149,60
AGS143-8	Quantifoil 1.2um dia (1.3µm) 300 Mesh Gold	10	149,60
AGS143-9	Quantifoil 1.2um dia (1.3µm) 400 Mesh Gold	10	149,60
AGS144	Boron carbide standard	each	28,80
AGS145	Medium resolution aluminium-tungsten dendrites	each	277,00
AGS147	Holey Carbon Films on 200 Mesh Copper Grids	50	180,10
AGS147-3	Holey Carbon Films on 300 Mesh Copper Grids	50	180,10
AGS147-3H	Holey Carbon Films on 300 Mesh Copper Grids	25	96,10
AGS147-4	Holey Carbon Films on 400 Mesh Copper Grids	50	180,10
AGS147-4H	Holey Carbon Films on 400 Mesh Copper Grids	25	96,10
AGS147A	Holey Carbon Films on 200 Mesh Gold	50	205,60
AGS147A3	Holey Carbon Films on 300 Mesh Gold	50	205,60
AGS147A3H	Holey Carbon Films on 300 Mesh Gold	25	132,20
AGS147A4	Holey Carbon Films on 400 Mesh Gold	50	205,60
AGS147A4H	Holey Carbon Films on 400 Mesh Gold	25	132,20
AGS147AH	Holey Carbon Films on 200 Mesh Gold	25	132,20
AGS147AH7	Holey Carbon Films on H7 Finder Grids Gold (Pack of 50)		596,00
AGS147AH7H	Holey Carbon Films on H7 Finder Grids Gold	25	284,50
AGS147H	Holey Carbon Films on 200 Mesh Copper Grids	25	96,10
AGS147-H7	Holey Carbon Films on H7 Finder Grids Copper	50	372,50
AGS147-H7H	Holey Carbon Films on H7 Finder Grids Copper	25	220,60
AGS147N	Holey Carbon Films on 200 Mesh Nickel	50	186,80
AGS147N3	Holey Carbon Films on 300 Mesh Nickel	50	186,80
AGS147N3H	Holey Carbon Films on 300 Mesh Nickel	25	120,10
AGS147N4	Holey Carbon Films on 400 Mesh Nickel	50	186,80
AGS147N4H	Holey Carbon Films on 400 Mesh Nickel	25	120,10

Part #	Description	Unit	Price €
AGS147NH	Holey Carbon Films on 200 Mesh Nickel	25	154,00
AGS147NH7	Holey Carbon Films on H7 Finder Grids Nickel	50	360,80
AGS147NH7H	Holey Carbon Films on H7 Finder Grids Nickel	25	193,00
AGS151	1000 Mesh Copper in sandwich grid 3mm	each	30,30
AGS151-5	1500 Mesh Copper in sandwich grid 3mm	each	31,00
AGS155	Standard for HVEM	each	37,90
AGS160	Carbon Films on 200 Mesh Grids Copper	50	126,60
AGS160-3	Carbon Films on 300 Mesh Grids Copper	50	128,90
AGS160-3H	Carbon Films on 300 Mesh Grids Copper	25	71,80
AGS160-4	Carbon Films on 400 Mesh Grids Copper	50	133,50
AGS160-4H	Carbon Films on 400 Mesh Grids Copper	25	77,30
AGS160A	Carbon Films on 200 Mesh Gold Grids	50	222,60
AGS160A3	Carbon Films on 300 Mesh Gold Grids	50	264,00
AGS160A3H	Carbon Films on 300 Mesh Gold Grids	25	158,80
AGS160A4	Carbon Films on 400 Mesh Gold Grids	50	262,60
AGS160A4H	Carbon Films on 400 Mesh Gold Grids	25	137,60
AGS160AF1	Carbon Film on F1 Finder Grid Gold	50	517,70
AGS160AF1H	Carbon Films on F1 Finder grid Gold.	25	276,30
AGS160AH	Carbon Films on 200 Mesh Gold Grids	25	128,70
AGS160AH2	Carbon Films on H2 Finder Grids Gold	50	461,70
AGS160AH2H	Carbon Films on H2 Finder Grids Gold	25	210,00
AGS160AH7	Carbon Film on H7 Finder Grid Gold	50	444,60
AGS160AH7H	Carbon Film on H7 Finder Grid Gold	25	210,00
AGS160-F1	Carbon Films on F1 finder grid Copper	50	342,20
AGS160-F1H	Carbon Films on F1 finder grid Copper.	25	182,50
AGS160H	Carbon Films on 200 Mesh Grids Copper	25	71,80
AGS160-H2	Carbon Films on H2 finder Grids Copper	50	410,30
AGS160-H2H	Carbon Films on H2 finder Grids Copper	25	232,40
AGS160-H7	Carbon Films on H7 finder grid Copper	50	289,70
AGS160-H7H	Carbon Film on H7 finder Grid Copper	25	232,40
AGS160N	Carbon Films on 200 Mesh Nickel Grids	50	141,70
AGS160N3	Carbon Films on 300 Mesh Nickel Grids	50	177,80
AGS160N3H	Carbon Films on 300 Mesh Nickel Grids	25	92,80
AGS160N4	Carbon Films on 400 Mesh Nickel Grids	50	179,80
AGS160N4H	Carbon Films on 400 Mesh Nickel Grids	25	107,70
AGS160NF1	Carbon Film on F1 Finder Grid Nickel	50	170,90
AGS160NF1H	Carbon Films on F1 Finder grid Nickel.	25	232,40
AGS160NH	Carbon Films on 200 Mesh Nickel Grids	25	92,70
AGS160NH2	Carbon Films on H2 Finder Grids Nickel	50	321,10
AGS160NH2H	Carbon Films on H2 Finder Grids Nickel	25	232,40
AGS160NH7	Carbon Film on H7 Finder Grid Nickel	50	170,90
AGS160NH7H	Carbon Film on H7 Finder Grid Nickel	25	232,40
AGS162	Formvar/Carbon on 200 Mesh Copper	50	172,60
AGS162-3	Formvar/Carbon on 300 Mesh Copper	50	172,60
AGS162-3H	Formvar/Carbon on 300 Mesh Copper	25	93,90
AGS162-4	Formvar/Carbon on 400 Mesh Copper	50	195,30
AGS162-4H	Formvar/Carbon on 400 Mesh Copper	25	172,60
AGS162-5	Formvar/Carbon on 2mmx1mm Slot Copper Grids (Pack of 50)	50	317,90
AGS162-6	Formvar/Carbon on 75 Mesh Copper	50	294,10
AGS162-6H	Formvar/Carbon on 75 Mesh Copper	25	194,90
AGS162A	Formvar/Carbon on 200 Mesh Gold	50	266,40
AGS162A1	Formvar/Carbon on 100 Mesh Gold	50	173,50
AGS162A1H	Formvar/Carbon on 100 Mesh Gold	25	96,80
AGS162A3	Formvar/Carbon on 300 Mesh Gold	50	266,40
AGS162A3H	Formvar/Carbon on 300 Mesh Gold	25	213,60
AGS162A4	Formvar/Carbon on 400 Mesh Gold	50	266,40
AGS162A4H	Formvar/Carbon on 400 Mesh Gold	25	213,60
AGS162A5	Formvar/Carbon on 2mmx1mm Gold Slots	50	422,90

Part #	Description	Unit	Price €
AGS162A5H	Formvar/Carbon on 2mmx1mm Gold Slots	25	337,20
AGS162A6H	Formvar/Carbon on 75 Mesh Gold	25	221,10
AGS162AF1	Formvar/Carbon Films on F1 Gold Grids	50	224,30
AGS162AF1H	Formvar/Carbon Films on F1 Gold Grids	25	142,60
AGS162AH	Formvar/Carbon on 200 Mesh Gold	25	133,90
AGS162AH7	Formvar/Carbon Films on H7 Gold Grids	50	383,40
AGS162AH7H	Formvar/Carbon Films on H7 Gold Grids	25	200,90
AGS162-F1	Formvar/Carbon Films on F1 Copper Grids	50	225,50
AGS162-F1H	Formvar/Carbon Films on F1 Copper Grids	25	160,40
AGS162H	Formvar/Carbon on 200 Mesh Copper	25	84,00
AGS162-H7	Formvar/Carbon Films on Copper H7 Grids	50	233,40
AGS162-H7H	Formvar/Carbon Films on Copper H7 Grids	25	109,80
AGS162N	Formvar/Carbon on 200 Mesh Nickel	50	204,20
AGS162N1	Formvar/Carbon on 100 Mesh Nickel	50	234,20
AGS162N1H	Formvar/Carbon on 100 Mesh Nickel	25	105,50
AGS162N3	Formvar/Carbon on 300 Mesh Nickel	50	204,20
AGS162N3H	Formvar/Carbon on 300 Mesh Nickel	25	180,70
AGS162N4	Formvar/Carbon on 400 Mesh Nickel	50	204,20
AGS162N4H	Formvar/Carbon on 400 Mesh Nickel	25	180,70
AGS162N5	Formvar/Carbon on 2mm x 1mm Nickel Slot	50	417,50
AGS162N5H	Formvar/Carbon on 2mmx1mm Nickel Slot	25	284,80
AGS162N6	Formvar/Carbon on 75 Mesh Nickel	50	305,50
AGS162N6H	Formvar/Carbon on 75 Mesh Nickel	25	219,20
AGS162NF1	Formvar/Carbon Films on F1 Nickel Grids	50	181,90
AGS162NF1H	Formvar/Carbon Films on F1 Nickel Grids	25	95,60
AGS162NH	Formvar/Carbon on 200 Mesh Nickel	25	119,80
AGS162NH7	Formvar/Carbon Films on H7 Nickel Grids	50	168,00
AGS162NH7H	Formvar/Carbon Films on H7 Nickel Grids	25	94,20
AGS166	Lacey Carbon Films on 200 Mesh Copper Grids	50	171,00
AGS166-3	Lacey Carbon Films on 300 Mesh Copper Grids	50	171,00
AGS166-3H	Lacey Carbon Films on 300 Mesh Copper Grids	25	97,60
AGS166-4	Lacey Carbon Films on 400 Mesh Copper Grids	50	171,00
AGS166-4H	Lacey Carbon Films on 400 Mesh Copper Grids	25	97,60
AGS166-4H-GO	Graphene oxide support film on Lacey carbon on 400mesh Cu grids	25	217,10
AGS166A	Lacey Carbon Films on 200 Mesh Gold	50	211,90
AGS166A3	Lacey Carbon Films on 300 Mesh Gold	50	211,90
AGS166A3H	Lacey Carbon Films on 300 Mesh Gold	25	136,20
AGS166A4	Lacey Carbon Films on 400 Mesh Gold	50	211,90
AGS166A4H	Lacey Carbon Films on 400 Mesh Gold	25	136,20
AGS166AH	Lacey Carbon Films on 200 Mesh Gold	25	168,00
AGS166AH7	Lacey Carbon Films on H7 finder grid. Gold	50	492,40
AGS166AH7H	Lacey Carbon Films on H7 finder grid. Gold	25	284,50
AGS166H	Lacey Carbon Films on 200 Mesh Copper Grids	25	97,60
AGS166-H2	Lacey Carbon Films on H2 Finder Grid	50	354,70
AGS166-H7	Lacey Carbon Films on H7 finder grid. Copper	50	372,50
AGS166-H7H	Lacey Carbon Films on H7 finder grid. Copper	25	220,60
AGS166N	Lacey Carbon Films on 200 Mesh Nickel	50	209,50
AGS166N3	Lacey Carbon Films on 300 Mesh Nickel	50	193,10
AGS166N3H	Lacey Carbon Films on 300 Mesh Nickel	25	124,20
AGS166N4	Lacey Carbon Films on 400 Mesh Nickel	50	193,10
AGS166N4H	Lacey Carbon Films on 400 Mesh Nickel	25	124,20
AGS166NH	Lacey Carbon Films on 200 Mesh Nickel	25	153,10
AGS166NH7	Lacey Carbon Films on H7 finder grid. Nickel	50	168,00
AGS166NH7H	Lacey Carbon Films on H7 finder grid. Nickel	25	98,70
AGS167	Silicon Dioxide Support Films 40nm	10	511,30
AGS167T	Silicon Dioxide Support Films 8nm	10	736,60
AGS167U	Silicon Dioxide Support Films 18nm	10	606,60
AGS168	Gold-On-Carbon SEM Resolution Test on Pin Stub	each	103,80

Part #	Description	Unit	Price €
AGS168A	Gold-On-Carbon SEM Test on 10mm JEOL Stub	each	103,80
AGS168AZ	Low KV Gold on Carbon Specimen on 10mm JEOL Stub	each	135,00
AGS168B	Gold-On-Carbon SEM Resolution Test on ISI Stub	each	103,80
AGS168BZ	Low KV Gold on Carbon Specimen on ISI/ABT Stub	each	135,00
AGS168C	Gold-On-Carbon SEM Resolution Test on Hitachi Stub	each	103,80
AGS168CZ	Low KV Gold on Carbon Specimen on Hitachi Stub	each	135,00
AGS168D	Gold on Carbon Resolution Test on Customer Stub *	each	103,80
AGS168DZ	Low KV Gold on Carbon Specimen on Customer Stub *	each	135,00
AGS168E	Gold-On-Carbon SEM Test on 12.5mm x 10mm Stub	each	103,80
AGS168EZ	Low KV Gold on Carbon Specimen on 12.5mm JEOL Stub	each	135,00
AGS168S	Gold on Carbon Resolution Test on Short Pin Stub	each	99,80
AGS168SI	Gold on 500um silicon substrate. 5 - 150nm	each	128,90
AGS168T	Gold on Thin Carbon SEM Resolution Test	each	103,80
AGS168TZ	Low KV Gold on Carbon Specimen on Carbon Disc	each	135,00
AGS168U	Gold on Carbon SEM Resolution Test (Unmounted)	each	99,80
AGS168UZ	Low KV Gold on Carbon Specimen Unmounted	each	135,00
AGS168VTSI	Gold on 200µm silicon substrate. 5 - 150nm	each	139,60
AGS168Z	Low KV Gold on Carbon Specimen on 12.5mm Pin Stub	each	135,00
AGS168ZSI	Gold on 500um silicon substrate. 30 - 300nm	each	145,90
AGS168ZVTSI	Gold on 200µm silicon substrate. 30 - 300nm	each	148,60
AGS169-13	Silicon Monoxide/Removable Formvar 300 Mesh, Copper	50	234,20
AGS169-13H	Silicon Monoxide/Removable Formvar 300 Mesh, Copper	25	130,20
AGS171	Chessy test specimen	each	1.496,90
AGS171-1	30nm Si Nitride on 200µm 0.5mm x 0.5mm	10	171,20
AGS171-1A4	30nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	294,10
AGS171-1A5	30nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	o.a.
AGS171-1A7	30nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	752,80
AGS171-1H	30nm Si Nitride on 200µm 1.0mm x 0.25mm	10	171,30
AGS171-1L	30nm Si Nitride on 200µm 1.0mm x 1.0mm	10	166,30
AGS171-1S	30nm Si Nitride on 200µm 0.25mm x 0.25mm	10	166,30
AGS171-1T	30nm Si Nitride on 200µm 0.1mm x 0.1mm	10	166,30
AGS171-2	50nm Si Nitride on 200µm 0.5mm x 0.5mm	10	126,20
AGS171-2A4	50nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	215,60
AGS171-2A5	50nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	o.a.
AGS171-2A7	50nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	701,50
AGS171-2H	50nm Si Nitride on 200µm 1.0mm x 0.25mm	10	126,20
AGS171-2L	50nm Si Nitride on 200µm 1.0mm x 1.0mm	10	130,00
AGS171-2S	50nm Si Nitride on 200µm 0.25mm x 0.25mm	10	134,80
AGS171-2T	50nm Si Nitride on 200µm 0.1mm x 0.1mm	10	134,80
AGS171-3	75nm Si Nitride on 200µm 0.5mm x 0.5mm	10	106,20
AGS171-3A4	75nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	185,10
AGS171-3A5	75nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	o.a.
AGS171-3A7	75nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	584,60
AGS171-3H	75nm Si Nitride on 200µm 1.0mm x 0.25mm	10	105,90
AGS171-3L	75nm Si Nitride on 200µm 1.0mm x 1.0mm	10	109,30
AGS171-3S	75nm Si Nitride on 200µm 0.25mm x 0.25mm	10	105,90
AGS171-3T	75nm Si Nitride on 200µm 0.1mm x 0.1mm	10	109,30
AGS171-4	100nm Si Nitride on 200µm 0.5mm x 0.5mm	10	94,00
AGS171-4A4	100nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	158,60
AGS171-4A5	100nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	257,00
AGS171-4A7	100nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	467,20
AGS171-4H	100nm Si Nitride on 200µm 1.0mm x 0.25mm	10	94,00
AGS171-4L	100nm Si Nitride on 200µm 1.0mm x 1.0mm	10	94,60
AGS171-4S	100nm Si Nitride on 200µm 0.25mm x 0.25mm	10	94,00
AGS171-4T	100nm Si Nitride on 200µm 0.1mm x 0.1mm	10	94,00
AGS171-5	150nm Si Nitride on 200µm 0.5mm x 0.5mm	10	94,60
AGS171-5A4	150nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	149,90
AGS171-5A5	150nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	257,00

Part #	Description	Unit	Price €
AGS171-5A7	150nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	467,20
AGS171-5H	150nm Si Nitride on 200µm 1.0mm x 0.25mm	10	94,60
AGS171-5L	150nm Si Nitride on 200µm 1.0mm x 1.0mm	10	92,80
AGS171-5S	150nm Si Nitride on 200µm 0.25mm x 0.25mm	10	94,00
AGS171-5T	150nm Si Nitride on 200µm 0.1mm x 0.1mm	10	94,00
AGS171-6	200nm Si Nitride on 200µm 0.5mm x 0.5mm	10	94,60
AGS171-6A4	200nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	158,60
AGS171-6A5	200nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	257,00
AGS171-6A7	200nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	467,20
AGS171-6H	200nm Si Nitride on 200µm 1.0mm x 0.25mm	10	94,60
AGS171-6L	200nm Si Nitride on 200µm 1.0mm x 1.0mm	10	94,60
AGS171-6S	200nm Si Nitride on 200µm 0.25mm x 0.25mm	10	94,00
AGS171-6T	200nm Si Nitride on 200µm 0.1mm x 0.1mm	10	91,70
AGS171-7	500nm Si Nitride on 200µm 0.50mm x 0.50mm	10	136,00
AGS171-7A4	500nm Si Nitride 4x4 Array 0.5mm x 0.5mm	each	141,90
AGS171-7A5	500nm Si Nitride 5x5 Array 0.5mm x 0.5mm	each	257,00
AGS171-7A7	500nm Si Nitride 7x7 Array 0.5mm x 0.5mm	each	467,20
AGS171-7H	500nm Si Nitride on 200µm 1.0mm x 0.25mm	10	94,00
AGS171-7L	500nm Si Nitride on 200µm 1.0mm x 1.0mm	10	94,00
AGS171-7S	500nm Si Nitride on 200µm 0.25mm x 0.25mm	10	94,60
AGS171-7T	500nm Si Nitride on 200µm 0.1mm x 0.1mm	10	94,60
AGS171-8	10nm Si Nitride on 200µm 0.5mm x 0.5mm	10	250,80
AGS171-9	20nm Si Nitride on 200µm 0.5mm x 0.5mm	10	198,40
AGS171-9A5	20nm Si Nitride 5x5 Array on 200µm 0.5mm x 0.5mm (Single)	each	535,30
AGS171-9A7	20nm Si Nitride 7x7 Array on 200µm 0.5mm x 0.5mm (Single)	each	1.049,10
AGS171-9H	20nm Si Nitride on 200µm 1.0mm x 0.25mm	10	198,40
AGS171-9S	20nm Si Nitride on 200µm 0.25mm x 0.25mm	10	198,40
AGS171-9T	20nm Si Nitride on 200µm 0.1mm x 0.1mm	10	198,40
AGS172-1	30nm Si Nitride on 100µm 0.50mm x 0.50mm	10	198,40
AGS172-1H	30nm Si Nitride on 100µm 1.0mm x 0.25mm	10	181,50
AGS172-1L	30nm Si Nitride on 100µm 1.0mm x 1.0mm	10	181,50
AGS172-1S	30nm Si Nitride on 100µm 0.25mm x 0.25mm	10	181,50
AGS172-1T	30nm Si Nitride on 100µm 0.1mm x 0.1mm	10	205,20
AGS172-2	50nm Si Nitride on 100µm 0.50mm x 0.50mm	10	166,50
AGS172-2H	50nm Si Nitride on 100µm 1.0mm x 0.25mm	10	150,20
AGS172-2L	50nm Si Nitride on 100µm 1.0mm x 1.0mm	10	156,30
AGS172-2S	50nm Si Nitride on 100µm 0.25mm x 0.25mm	10	150,20
AGS172-2T	50nm Si Nitride on 100µm 0.1mm x 0.1mm	10	166,50
AGS172-3	75nm Si Nitride on 100µm 0.50mm x 0.50mm	10	135,10
AGS172-3H	75nm Si Nitride on 100µm 1.0mm x 0.25mm	10	128,00
AGS172-3L	75nm Si Nitride on 100µm 1.0mm x 1.0mm	10	128,00
AGS172-3S	75nm Si Nitride on 100µm 0.25mm x 0.25mm	10	128,00
AGS172-3T	75nm Si Nitride on 100µm 0.1mm x 0.1mm	10	131,20
AGS172-4	100nm Si Nitride on 100µm 0.50mm x 0.50mm	10	109,30
AGS172-4H	100nm Si Nitride on 100µm 1.0mm x 0.25mm	10	106,20
AGS172-4L	100nm Si Nitride on 100µm 1.0mm x 1.0mm	10	#N/B
AGS172-4S	100nm Si Nitride on 100µm 0.25mm x 0.25mm	10	105,90
AGS172-4T	100nm Si Nitride on 100µm 0.1mm x 0.1mm	10	106,20
AGS172-5	150nm Si Nitride on 100µm 0.50mm x 0.50mm	10	109,30
AGS172-5H	150nm Si Nitride on 100µm 1.0mm x 0.25mm	10	105,90
AGS172-5L	150nm Si Nitride on 100µm 1.0mm x 1.0mm	10	105,90
AGS172-5S	150nm Si Nitride on 100µm 0.25mm x 0.25mm	10	105,90
AGS172-5T	150nm Si Nitride on 100µm 0.1mm x 0.1mm	10	105,90
AGS172-6	200nm Si Nitride on 100µm 0.50mm x 0.50mm	10	109,30
AGS172-6H	200nm Si Nitride on 100µm 1.0mm x 0.25mm	10	105,90
AGS172-6L	200nm Si Nitride on 100µm 1.0mm x 1.0mm	10	105,90
AGS172-6S	200nm Si Nitride on 100µm 0.25mm x 0.25mm	10	105,90
AGS172-6T	200nm Si Nitride on 100µm 0.1mm x 0.1mm	10	105,90

Part #	Description	Unit	Price €
AGS172-7	500nm Si Nitride on 100µm 0.50mm x 0.50mm	10	130,00
AGS172-7H	500nm Si Nitride on 100µm 1.0mm x 0.25mm	10	105,90
AGS172-7L	500nm Si Nitride on 100µm 1.0mm x 1.0mm	10	105,90
AGS172-7S	500nm Si Nitride on 100µm 0.25mm x 0.25mm	10	105,90
AGS172-9	20nm Si Nitride on 100µm 0.5mm x 0.5mm	10	246,50
AGS172-9H	20nm Si Nitride on 100µm 1.0mm x 0.25mm	10	246,50
AGS172-9L	20nm Si Nitride on 100µm 1.0mm x 1.0mm	10	246,50
AGS172-9S	20nm Si Nitride on 100µm 0.25mm x 0.25mm	10	246,50
AGS172-9T	20nm Si Nitride on 100µm 0.1mm x 0.1mm	10	246,50
AGS173-1	Quantifoil 2µm dia (2µm) 200 Mesh Copper	10	145,50
AGS173-1-100	Quantifoils 2µm dia (2µm) 200 Mesh Copper	100	830,70
AGS173-2	Quantifoils 2µm dia (2µm) 300 Mesh Copper	10	145,50
AGS173-2-100	Quantifoils 2µm dia (2µm) 300 Mesh Copper	100	755,20
AGS173-2-GO	Graphene oxide support film on Quantifoil R2/2 on 300mesh Cu	each	181,80
AGS173-3	Quantifoil 2µm dia (2µm) 400 Mesh Copper	10	145,50
AGS173-3-100	Quantifoils 2µm dia (2µm) 400 Mesh Copper	100	830,70
AGS173-4	Quantifoil 2µm dia (2µm) 200 Mesh Nickel	10	149,60
AGS173-5	Quantifoil 2µm dia (2µm) 300 Mesh Nickel	10	149,60
AGS173-6	Quantifoils 2µm dia (2µm) 400 Mesh Nickel	10	149,60
AGS173-7	Quantifoil 2µm dia (2µm) 200 Mesh Gold	10	149,60
AGS173-8	Quantifoil 2µm dia (2µm) 300 Mesh Gold	10	149,60
AGS173-9	Quantifoils 2µm dia (2µm) 400 Mesh Gold	10	149,60
AGS174-1	Quantifoil 2µm dia (1µm) 200 Mesh Copper	10	145,50
AGS174-1-100	Quantifoil 2µm dia (1µm) 200 Mesh Copper	100	842,60
AGS174-2	Quantifoil 2µm dia (1µm) 300 Mesh Copper	10	152,10
AGS174-2-100	Quantifoil 2µm dia (1µm) 300 Mesh Copper	100	842,60
AGS174-3	Quantifoil 2µm dia (1µm) 400 Mesh Copper	10	145,50
AGS174-3-100	Quantifoil 2µm dia (1µm) 400 Mesh Copper	100	919,20
AGS174-4	Quantifoil 2µm dia (1µm) 200 Mesh Nickel	10	163,90
AGS174-5	Quantifoil 2µm dia (1µm) 300 Mesh Nickel	10	163,90
AGS174-6	Quantifoil 2µm dia (1µm) 400 Mesh Nickel	10	153,90
AGS174-7	Quantifoil 2µm dia (1µm) 200 Mesh Gold	10	163,90
AGS174-8	Quantifoil 2µm dia (1µm) 300 Mesh Gold	10	146,90
AGS174-9	Quantifoil 2µm dia (1µm) 400 Mesh Gold	10	163,90
AGS175-1	Quantifoil 2µm dia (4.0µm) 200 Mesh Copper	10	131,20
AGS175-2	Quantifoil 2µm dia (4.0µm) 300 Mesh Copper	10	131,20
AGS175-3	Quantifoil 2µm dia (4.0µm) 400 Mesh Copper	10	131,20
AGS175-4	Quantifoil 2µm dia (4.0µm) 200 Mesh Nickel	10	140,50
AGS175-5	Quantifoil 2µm dia (4.0µm) 300 Mesh Nickel	10	149,60
AGS175-6	Quantifoil 2µm dia (4.0µm) 400 Mesh Nickel	10	140,50
AGS175-7	Quantifoil 2µm dia (4.0µm) 200 Mesh Gold	10	149,60
AGS175-7-100	Quantifoil 2µm dia (4.0µm) 200 Mesh Gold 100 Pk	10	890,30
AGS175-8	Quantifoil 2µm dia (4.0µm) 300 Mesh Gold	10	140,50
AGS175-9	Quantifoil 2µm dia (4.0µm) 400 Mesh Gold	10	149,60
AGS176-1	Quantifoil 1µm dia (4.0µm) 200 Mesh Copper	10	153,60
AGS176-2	Quantifoil 1µm dia (4.0µm) 300 Mesh Copper	10	136,60
AGS176-3	Quantifoil 1µm dia (4.0µm) 400 Mesh Copper	10	145,50
AGS176-4	Quantifoil 1µm dia (4.0µm) 200 Mesh Nickel	10	153,90
AGS176-5	Quantifoil 1µm dia (4.0µm) 300 Mesh Nickel	10	153,90
AGS176-6	Quantifoil 1µm dia (4.0µm) 400 Mesh Nickel	10	163,90
AGS176-7	Quantifoil 1µm dia (4.0µm) 200 Mesh Gold	10	163,90
AGS176-8	Quantifoil 1µm dia (4.0µm) 300 Mesh Gold	10	163,90
AGS176-9	Quantifoil 1µm dia (4.0µm) 400 Mesh Gold	10	163,90
AGS177-1	Quantifoil 3.5µm dia (1.0µm) 200 Mesh Copper	10	145,50
AGS177-2	Quantifoil 3.5µm dia (1.0µm) 300 Mesh Copper	10	145,50
AGS177-3	Quantifoil 3.5µm dia (1.0µm) 400 Mesh Copper	10	145,50
AGS177-4	Quantifoil 3.5µm dia (1.0µm) 200 Mesh Nickel	10	163,90
AGS177-5	Quantifoil 3.5µm dia (1.0µm) 300 Mesh Nickel	10	153,90

Part #	Description	Unit	Price €
AGS177-6	Quantifoil 3.5µm dia (1.0um) 400 Mesh Nickel	10	153,90
AGS177-7	Quantifoil 3.5µm dia (1.0um) 200 Mesh Gold	10	163,90
AGS177-8	Quantifoil 3.5µm dia (1.0um) 300 Mesh Gold	10	165,60
AGS177-9	Quantifoil 3.5µm dia (1.0um) 400 Mesh Gold	10	163,90
AGS179-GO	Suspended Monolayer Graphene on TEM Grids Quantifoil R2/4	each	75,40
AGS179-GO4	Suspended Monolayer Graphene on TEM Grids Quantifoil R2/4	4	183,20
AGS1800-100	100nm Nanospheres	15ml	753,80
AGS1800-125	125nm Nanospheres	15ml	753,80
AGS1800-150	150nm Nanospheres	15ml	783,20
AGS1800-20	20nm Nanospheres	15ml	980,90
AGS1800-200	200nm Nanospheres	15ml	753,80
AGS1800-30	30nm Nanospheres	15ml	735,40
AGS1800-300	300nm Nanospheres	15ml	810,80
AGS1800-40	40nm Nanospheres	15ml	783,20
AGS1800-400	400nm Nanospheres	15ml	864,30
AGS1800-50	50nm Nanospheres	15ml	753,80
AGS1800-500	500nm Nanospheres	15ml	753,80
AGS1800-60	60nm Nanospheres	15ml	810,80
AGS1800-600	600nm Nanospheres	15ml	735,40
AGS1800-70	70nm Nanospheres	15ml	735,40
AGS1800-700	700nm Nanospheres	15ml	735,40
AGS1800-80	80nm Nanospheres	15ml	783,20
AGS1800-800	800nm Nanospheres	15ml	735,40
AGS1800-90	90nm Nanospheres	15ml	735,40
AGS1801-1	1.0um Microspheres	15ml	939,50
AGS1801-10	10µm Microspheres	15ml	939,50
AGS1801-100	100µm Microspheres	15ml	939,50
AGS1801-1000	1000µm Microspheres	1g	939,50
AGS1801-15	15µm Microspheres	15ml	415,30
AGS1801-160	160µm Microspheres	15ml	939,50
AGS1801-2	2.0um Microspheres	15ml	350,60
AGS1801-2.5	2.5µm Microspheres	15ml	1.010,60
AGS1801-20	20µm Microspheres	15ml	1.018,20
AGS1801-200	200µm Microspheres	1g	939,50
AGS1801-25	25µm Microspheres	15ml	916,60
AGS1801-3	3.0um Microspheres	15ml	939,50
AGS1801-30	30µm Microspheres	15ml	415,30
AGS1801-300	300µm Microspheres	1g	916,60
AGS1801-40	40µm Microspheres	15ml	1.018,20
AGS1801-5	5.0um Microspheres	15ml	1.024,10
AGS1801-50	50µm Microspheres	15ml	939,50
AGS1801-500	500um Microspheres	1g	939,50
AGS1801-6	6.0um Microspheres	15ml	916,60
AGS1801-60	60µm Microspheres	15ml	916,60
AGS1801-7	7.0um Microspheres	15ml	916,60
AGS1801-70	70µm Microspheres	15ml	916,60
AGS1801-750	750µm Microspheres	1g	916,60
AGS1801-80	80um Microspheres	15ml	916,60
AGS1802-10	10µm DVB Microspheres	15ml	671,50
AGS1802-15	15µm DVB Microspheres	15ml	671,50
AGS1802-20	20µm DVB Microspheres	15ml	671,50
AGS1802-30	30µm DVB Microspheres	15ml	671,50
AGS1802-40	40µm DVB Microspheres	15ml	671,50
AGS1802-5	5µm DVB Microspheres	15ml	671,50
AGS1803-05	0.5µm Silica Microspheres	15ml	473,90
AGS1803-07	0.7um Silica Microspheres	15ml	492,50
AGS1803-10	1.0um Silica Microspheres	15ml	492,50
AGS1803-16	1.6µm Silica Microspheres	15ml	#N/B

Part #	Description	Unit	Price €
AGS1805	Silicon nitride (Si3N4) ion sputter standard	each	523,60
AGS1807	Nickel / chromium (Ni/Cr) ion sputter standard	each	881,90
AGS1808	38nm Nanogold High Resolution Test on 12.5mm Pin Stub	each	221,80
AGS1808A	38nm Nanogold High Resolution Test on 10mm JEOL Stub	each	o.a.
AGS1808B	38nm Nanogold High Resolution Test on 15mm Topcon/ISI/ABT st	each	o.a.
AGS1808C	38nm Nanogold High Resolution Test on 15mm Hitachi stub	each	221,80
AGS1808D	38nm Nanogold High Resolution Test on customer stub *		o.a.
AGS1808E	38nm Nanogold High Resolution Test on 12.5mm JEOL stub	each	o.a.
AGS1808S	38nm Nanogold High Resolution Test on Short Pin Stub	each	250,60
AGS1808U	38nm Nanogold High Resolution Test Unmounted	each	221,80
AGS1809	27.3nm Nanogold High Resolution Test on 12.5mm Pin Stub	each	221,80
AGS1809A	27.3nm Nanogold High Resolution Test on 10mm JEOL Stub	each	o.a.
AGS1809B	27.3nm Nanogold High Resolution Test on 15mm Topcon/ISI/ABT	each	o.a.
AGS1809C	27.3nm Nanogold High Resolution Test on 15mm Hitachi stub	each	o.a.
AGS1809D	27.3nm Nanogold High Resolution Test on customer stub *		o.a.
AGS1809E	27.3nm Nanogold High Resolution Test on 12.5mm JEOL stub	each	o.a.
AGS1809S	27.3nm Nanogold High Resolution Test on Short Pin Stub	each	250,60
AGS1809U	27.3nm Nanogold High Resolution Test Unmounted	each	221,80
AGS1810	MRS4 Geller Reference Standard Non Traceable	each	o.a.
AGS1810-C	MRS4 Cleaning Only (AGS1810)		o.a.
AGS1810-G	MRS4 GellerRefStd Non Traceable (1" X 3" Glass Slide Mounti		o.a.
AGS1810-O	MRS4 Geller Ref Standard Non Traceable in Optical Retainer	each	o.a.
AGS1810P	MRS Pin ONLY for use with SEM Retainer (MRS3 & MRS4 Range)		o.a.
AGS1810-RXY	MRS4 RXY Recertification & Clean (AGS1810)		o.a.
AGS1810-RXYZ	MRS4 RXYZ Recertification & Clean (AGS1810)		o.a.
AGS1810-S	MRS4 Geller Reference Standard Non Traceable in SEM Retainer	each	o.a.
AGS1810-T	MRS4 GellerRefStd Non Traceable in TEM 3mm x 0/35mm Retaine		o.a.
AGS1811	MRS4 Geller Reference Standard XY Traceable	each	o.a.
AGS1811-C	MRS4 Cleaning Only (AGS1811)		o.a.
AGS1811-G	MRS4 GellerRefStd XY Traceable (1" X 3" Glass Slide Mountin		o.a.
AGS1811-O	MRS4 GellerRefStandard XY Traceable in Optical Retainer	each	o.a.
AGS1811-RXY	MRS4 RXY Recertification & Clean (AGS1811)		o.a.
AGS1811-RXYZ	MRS4 RXYZ Recertification & Clean (AGS1811)		o.a.
AGS1811-S	MRS4 Geller Reference Standard XY Traceable in SEM Retainer	each	o.a.
AGS1811-T	MRS4 GellerRefStd XY Traceable in TEM 3mm x 0/35mm Retainer		o.a.
AGS1812	MRS4 Geller Reference Standard XYZ Traceable	each	o.a.
AGS1812-C	MRS4 Cleaning Only (AGS1812)		o.a.
AGS1812-G	MRS4 GellerRefStd XYZ Traceable (1" X 3" Glass Slide Mounti		o.a.
AGS1812-O	MRS4 GellerRefStandard XYZ Traceable in Optical Retainer	each	o.a.
AGS1812-RXY	MRS4 RXY Recertification & Clean (AGS1812)		o.a.
AGS1812-RXYZ	MRS4 RXYZ Recertification & Clean (AGS1812)		o.a.
AGS1812-S	MRS4 Geller Reference Standard XYZ Traceable in SEM Retainer	each	o.a.
AGS1812-T	MRS4 GellerRefStd XYZ Traceable in TEM 3mm x 0/35mm Retaine		o.a.
AGS1813	MRS-4 Non Certified & TEM Holder	each	o.a.
AGS1814	MRS6 Geller Reference Standard Non Traceable	each	o.a.
AGS1814-C	MRS6 Cleaning Only (AGS1814)		o.a.
AGS1814-P	MRS6 GellerRefStd Non Traceable with Retainer Pin		o.a.
AGS1814-RXY	MRS6 RXY Recertification & Clean (AGS1814)		o.a.
AGS1815	MRS6 Geller Reference Standard XY Traceable	each	o.a.
AGS1815-C	MRS6 Cleaning Only (AGS1815)		o.a.
AGS1815-P	MRS6 GellerRefStd XY Traceable with Retainer Pin		o.a.
AGS1815-RXY	MRS6 RXY Recertification & Clean (AGS1815)		o.a.
AGS1816-002	Glass Microsphere Size Standards - 2um	each	714,10
AGS1816-005	Glass Microsphere Size Standards - 5um	each	714,10
AGS1816-008	Glass Microsphere Size Standards - 8um	each	714,10
AGS1816-010	Glass Microsphere Size Standards - 10um	each	714,10
AGS1816-015	Glass Microsphere Size Standards - 15um	each	714,10
AGS1816-020	Glass Microsphere Size Standards - 20um	each	714,10

Part #	Description	Unit	Price €
AGS1816-030	Glass Microsphere Size Standards - 30µm	each	639,50
AGS1816-040	Glass Microsphere Size Standards - 40µm	each	639,50
AGS1816-050	Glass Microsphere Size Standards - 50µm	each	639,50
AGS1816-060	Glass Microsphere Size Standards - 60µm	each	639,50
AGS1816-070	Glass Microsphere Size Standards - 70µm	each	639,50
AGS1816-080	Glass Microsphere Size Standards - 80µm	each	639,50
AGS1816-090	Glass Microsphere Size Standards - 90µm	each	639,50
AGS1816-100	Glass Microsphere Size Standards - 100µm	each	655,50
AGS1816-1000	Glass Microsphere Size Standards - 1000µm	each	639,50
AGS1816-110	Glass Microsphere Size Standards - 110µm	each	639,50
AGS1816-120	Glass Microsphere Size Standards - 120µm	each	655,50
AGS1816-140	Glass Microsphere Size Standards - 140µm	each	655,50
AGS1816-170	Glass Microsphere Size Standards - 170µm	each	639,50
AGS1816-200	Glass Microsphere Size Standards - 200µm	each	639,50
AGS1816-2000	Glass Microsphere Size Standards - 2000µm	each	639,50
AGS1816-230	Glass Microsphere Size Standards - 230µm	each	639,50
AGS1816-280	Glass Microsphere Size Standards - 280µm	each	639,50
AGS1816-330	Glass Microsphere Size Standards - 330µm	each	639,50
AGS1816-400	Glass Microsphere Size Standards - 400µm	each	639,50
AGS1816-480	Glass Microsphere Size Standards - 480µm	each	639,50
AGS1816-550	Glass Microsphere Size Standards - 550µm	each	639,50
AGS1816-650	Glass Microsphere Size Standards - 650µm	each	639,50
AGS1816-750	Glass Microsphere Size Standards - 750µm	each	639,50
AGS1816-950	Glass Microsphere Size Standards - 950µm	each	639,50
AGS1823	GSR and particle analysis calibration kit	each	o.a.
AGS1830	Universal Set of 20 Standards for XRF	each	o.a.
AGS1831	Set of 10 Rare Earth Standards for XRF	each	o.a.
AGS1832	Set of 5 Precious Metal Standards for XRF	each	o.a.
AGS1833	Standard Set A for XRF	each	o.a.
AGS1834	Standard Set B for XRF	each	o.a.
AGS1835	Single XRF Standard	each	o.a.
AGS185-3	Silicon monoxide on lacey Formvar, 300 Mesh	each	229,00
AGS1861	144nm calibration standard on 12.5 mm pin stub	each	o.a.
AGS1861A	144nm calibration standard on 10 mm JEOL stub	each	o.a.
AGS1861B	144nm calibration standard on 15 mm Topcon/ISI/ABT stub	each	o.a.
AGS1861C	144nm calibration standard on 15 mm Hitachi stub	each	o.a.
AGS1861D	144nm calibration standard on customer stub *		o.a.
AGS1861E	144nm calibration standard on 12.5 mm JEOL stub	each	o.a.
AGS1861S	144nm calibration standard on Short Pin Stub	each	o.a.
AGS1861U	144nm calibration standard, unmounted	each	o.a.
AGS1862	144nm certified calibration standard on 12.5 mm pin stub	each	o.a.
AGS1862A	144nm certified calibration standard on 10 mm JEOL stub	each	o.a.
AGS1862B	144nm certified calibration standard on 15 mm Topcon/ISI/ABT	each	o.a.
AGS1862C	144nm certified calibration standard on 15 mm Hitachi stub	each	o.a.
AGS1862D			o.a.
AGS1862E	144nm certified calibration standard on 12.5 mm JEOL stub	each	o.a.
AGS1862S	144nm certified calibration standard on Short Pin Stub	each	o.a.
AGS1862U	144nm certified calibration standard, unmounted	each	o.a.
AGS1863	300nm calibration standard on 12.5 mm pin stub	each	o.a.
AGS1863A	300nm calibration standard on 10 mm JEOL stub	each	o.a.
AGS1863B	300nm calibration standard on 15 mm Topcon/ISI/ABT stub	each	o.a.
AGS1863C	300nm calibration standard on 15 mm Hitachi stub	each	o.a.
AGS1863D			o.a.
AGS1863E	300nm calibration standard on 12.5 mm JEOL stub	each	o.a.
AGS1863S	300nm calibration standard on Short Pin Stub	each	o.a.
AGS1863U	300nm calibration standard, unmounted	each	o.a.
AGS186-4	3nm Carbon Films on 400 Mesh Copper	50	223,70
AGS1867	292nm resolution standard on 12.5mm pin stub	each	o.a.

Part #	Description	Unit	Price €
AGS1867A	292nm resolution standard on 10mm JEOL stub	each	o.a.
AGS1867B	292nm resolution on 15mm Topcon/ISI/ABT stub	each	o.a.
AGS1867C	292nm resolution standard on 15mm Hitachi stub	each	o.a.
AGS1867D			o.a.
AGS1867E	292nm resolution standard on 12.5mm JEOL stub	each	o.a.
AGS1867S	292nm resolution standard on Short Pin Stub	each	o.a.
AGS1867U	292nm resolution standard, unmounted	each	o.a.
AGS1868	292nm certified resolution standard on 12.5 mm pin stub	each	o.a.
AGS1868A	292nm certified resolution standard on 10 mm JEOL stub	each	o.a.
AGS1868B	292nm certified resolution on 15 mm Topcon/ISI/ABT stub	each	o.a.
AGS1868C	292nm certified resolution standard on 15 mm Hitachi stub	each	o.a.
AGS1868D			o.a.
AGS1868E	292nm certified resolution standard on 12.5 mm JEOL stub	each	o.a.
AGS1868S	292nm certified resolution standard on Short Pin Stub	each	o.a.
AGS186H4	3nm Carbon Films on 400 Mesh Copper	25	120,90
AGS187-3	UltrAuFoil R1.2 / 1.3 on 300 mesh Gold	10	244,80
AGS187-4	Ultra-thin carbon support film, 3nm - on lacey carbon	each	232,40
AGS188-1000A	1000 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-1000C	1000 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-1000N	1000 Mesh Nickel 150mm x 150mm	each	737,60
AGS188-100A	100 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-100C	100 Mesh Copper 150mm x 150mm	each	1.431,60
AGS188-100N	100 Mesh Nickel 150mm x 150mm	each	569,90
AGS188-1500A	1500 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-1500C	1500 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-1500N	1500 Mesh Nickel 150mm x 150mm	each	1.423,10
AGS188-2	UltrAuFoil R2 / 2 on 200 mesh Gold	10	254,50
AGS188-2000N	2000 Mesh Nickel 150mm x 150mm	each	o.a.
AGS188-200A	200 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-200C	200 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-200N	200 Mesh Nickel 150mm x 150mm	each	672,80
AGS188-20A	20 Mesh Gold 165mm x 165mm	each	o.a.
AGS188-20C	20 Mesh Copper 165mm x 165mm	each	1.262,70
AGS188-20N	20 Mesh Nickel 165mm x 165mm	each	569,00
AGS188-300A	300 Mesh Gold 275mm x 275mm	each	o.a.
AGS188-333A	333 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-333C	333 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-333N	333 Mesh Nickel 150mm x 150mm	each	569,50
AGS188-400A	400 Mesh Gold 190mm x 190mm	each	o.a.
AGS188-400C	400 Mesh Copper 190mm x 190mm	each	o.a.
AGS188-400N	400 Mesh Nickel 190mm x 190mm	each	1.021,40
AGS188-500A			o.a.
AGS188-500N	500 Mesh Nickel 150mm x 150mm	each	815,00
AGS188-5A	5 Mesh Gold 140mm x 140mm	each	o.a.
AGS188-5C	5 Mesh Copper 140mm x 140mm	each	1.262,70
AGS188-5N	5 Mesh Nickel 140mm x 140mm	each	569,00
AGS188-670C	670 Mesh Copper 275mm x 275mm	each	o.a.
AGS188-70A	70 Mesh Gold 275mm x 275mm	each	o.a.
AGS188-70C	70 Mesh Copper 275mm x 275mm	each	o.a.
AGS188-70N	70 Mesh Nickel 275mm x 275mm	each	609,00
AGS188-750A			o.a.
AGS188-750C	750 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-750N	750 Mesh Nickel 150mm x 150mm	each	815,00
AGS188-80A	80 Mesh Gold 150mm x 150mm	each	o.a.
AGS188-80C	80 Mesh Copper 150mm x 150mm	each	o.a.
AGS188-80N	80 Mesh Nickel 150mm x 150mm	each	569,50
AGS189-3	UltrAuFoil R0.6 / 1 on 300 mesh Gold	10	304,60
AGS1914	25mm Universal Standards (Set of 45)	each	o.a.

Part #	Description	Unit	Price €
AGS1915	32mm Universal Standards (Set of 45)	each	o.a.
AGS1916	25mm Semiconductor Standards (Set of 21)	each	o.a.
AGS1917	32mm Semiconductor Standards (Set of 21)	each	o.a.
AGS1918	NPL image analysis standard	each	o.a.
AGS1919	SEM Multi-purpose calibration block	each	o.a.
AGS1920	Biological Standard Set - 25mm Block	each	o.a.
AGS1921	Biological Standard Set - 32mm Block	each	o.a.
AGS1928	Basic calibration set	each	o.a.
AGS1929	Faraday cup	each	295,10
AGS1930	Silicon Specimen - Unmounted	each	91,60
AGS1930-CT	Certified Planotec silicon test specimen, unmounted	each	431,30
AGS1931	Silicon Specimens - Unmounted	10	699,50
AGS1932	Silicon Specimen on Pin Stub	each	95,50
AGS1932A	Silicon Specimen on 10 x 10mm JEOL Stub	each	114,90
AGS1932A-CT	Certified Planotec silicon test specimen on 10mm JEOL stub	each	454,90
AGS1932B	Silicon Specimen on 15mm ISI/ABT/Topcon Stub	each	114,90
AGS1932B-CT	Certified Planotec silicon test specimen on 15mm Topcon/ISI/	each	454,90
AGS1932C	Silicon Specimen on Hitachi Stub	each	114,90
AGS1932C-CT	Certified Planotec silicon test specimen on 15mm Hitachi stu	each	454,90
AGS1932-CT	Certified Planotec silicon test specimen on 12.5mm pin stub	each	439,10
AGS1932D	Silicon Specimen on Customer Stub *	each	114,90
AGS1932D-CT	Certified Planotec silicon test specimen on customer stub *	each	454,90
AGS1932E	Silicon Specimen on 12mm JEOL Stub	each	114,90
AGS1932E-CT	Certified Planotec silicon test specimen on 12.5mm JEOL stub	each	454,90
AGS1932S	Silicon Specimen on Short Pin Stub	each	91,90
AGS1934	Silicon Specimen for Incident Light Microscopy	each	128,70
AGS1934-CT	Certified Planotec silicon test specimen for incident light	each	466,10
AGS1936	Mag*I*Cal calibration standard	each	1.539,10
AGS1937	Universal Tin on Carbon Specimen on Pin Stub	each	193,60
AGS1937A	Universal Tin on Carbon Specimen on 10mm JEOL Stub	each	186,20
AGS1937B	Universal Tin on Carbon Specimen on ISI/ABT Stub	each	193,60
AGS1937C	Universal Tin on Carbon Specimen on 15mm Hitachi	each	193,60
AGS1937D	Universal Tin on Carbon Specimen on Customer Stub *	each	193,60
AGS1937E	Universal Tin on Carbon Specimen on 12.5mm JEOL	each	193,60
AGS1937S	Universal Tin on Carbon Specimen on Short Pin Stub	each	186,20
AGS1937SI	Tin on 500um silicon substrate. 5 - 30µm	each	229,30
AGS1937T	Universal Tin on Carbon Specimen (Thin Carbon Disc)	each	193,60
AGS1937TS	Universal Tin on Thin Carbon - Short Pin Stub	each	186,20
AGS1937U	Universal Tin on Carbon Specimen Unmounted	each	190,10
AGS1937VTS	Ultra Thin Tin on Carbon Wafer Test Specimen - Short Pin Stu	each	249,40
AGS1937VTSI	Tin on 200µm silicon substrate. 5 - 30µm	each	239,90
AGS1949	MetroChip calibration standard	each	1.343,40
AGS1950	BSE Reference - Nickel/Copper	each	423,20
AGS1951	BSE Reference - Palladium/Silver	each	409,80
AGS1952	BSE Reference - Platinum/Gold	each	404,60
AGS1953	Duplex reference specimen	each	424,20
AGS1954	BSE Reference - Aluminium/Silicon	each	398,10
AGS1957	Quality control test calibration block	each	o.a.
AGS1959	Chemography test sample, synthetic GSR Paper Process	each	871,40
AGS1959-1	Chemography test sample, synthetic GSR Filter Paper Process	each	871,40
AGS1959R	Chemography test sample Homogeneity & Resolution PhotoPaper	each	871,40
AGS1959R-1	Chemography test sample Homogeneity & Resolution FilterPaper	each	871,40
AGS1960	Re-calibration certificate for 2160 lines/mm SIRA	each	389,30
AGS1961	Re-calibration certificate for 19.7 lines/mm, SIRA	each	389,30
AGS1962	Agar re-calibration certificate	each	389,30
AGS1965	1000 mesh, 3 mm diameter TEM grid for calibration, in vial	each	23,40
AGS1967	Tin on Carbon Test Specimen on Pin Stub	each	150,90
AGS1967A	Tin on Carbon Test Specimen on JEOL Stub	each	150,90

Part #	Description	Unit	Price €
AGS1967B	Tin on Carbon Test Specimen on ISI Stub	each	150,90
AGS1967C	Tin on Carbon Test Specimen on Hitachi Stub	each	145,10
AGS1967D	Tin on Carbon Test Specimen on Special Stub	each	150,90
AGS1967E	Tin on Carbon Test Specimen on 12.5mm JEOL Stub	each	150,90
AGS1967S	Tin on Carbon Test Specimen on Short Pin Stub	each	145,10
AGS1967SI	Tin on 500um silicon substrate. 10 - 100nm	each	175,20
AGS1967T	Tin on Carbon Test Specimen - Unmounted-Thin	each	150,90
AGS1967U	Tin on Carbon Test Specimen - Unmounted	each	147,10
AGS1967VTSI	Tin on 200µm silicon substrate. 10 - 100nm	each	176,40
AGS1968	Prickly Gold on 12.5mm Pin Stub	each	114,90
AGS1968A	Prickly Gold on JEOL Stub	each	114,90
AGS1968B	Prickly Gold on ISI Stub	each	114,90
AGS1968C	Prickly Gold on Hitachi Stub	each	114,90
AGS1968D	Prickly Gold on Customer's Stub	each	267,80
AGS1968E	Prickly Gold on 12.5mm JEOL Stub	each	114,90
AGS1968S	Prickly Gold on Short Pin Stub	each	87,50
AGS1968U	Prickly Gold Unmounted	each	91,70
AGS1969	High Resolution Test (Gold) on 12.5mm Pin Stub	each	169,00
AGS1969A	High Resolution Test (Gold) on JEOL Stub	each	169,00
AGS1969B	High Resolution Test (Gold) on ISI Stub	each	169,00
AGS1969C	High Resolution Test (Gold) on Hitachi Stub	each	169,00
AGS1969D	High Resolution Test (Gold) on Customer Stub *	each	169,00
AGS1969E	High Resolution Test (Gold) on 12.5x10mm	each	169,00
AGS1969S	High Resolution Test (Gold) on Short Pin Stub	each	162,50
AGS1969SI	Gold on 500um silicon substrate. 3 - 50nm	each	179,50
AGS1969T	High Resolution Test (Gold) on Thin Carbon	each	161,00
AGS1969U	High Resolution Test (Gold) Unmounted	each	169,00
AGS1969VTSI	Gold on 200µm silicon substrate. 3 - 50nm	each	181,30
AGS1970	Cobalt Standard 2mm	each	o.a.
AGS1980	Universal TEM X-Ray Standard Set. Set of 25	each	1.181,50
AGS1981	Rare Earth TEM X-Ray Standard Set. Set of 14	each	550,10
AGS1982	Metal Foil Standards. Set of 25	each	o.a.
AGS1983	Metal Foil Standards. Set of 5	each	o.a.
AGS1984	UHV - Electronics Set	each	o.a.
AGS1985	UHV Standards to Customer's Choice	each	o.a.
AGS1987	High Resolution Au-C Test on 12.5mm Pin Stub	each	255,20
AGS1987A	High Resolution Au-C Test on 10mm JEOL Stub	each	245,40
AGS1987B	High Resolution Au-C Test on ISI/ABT/Topcon Stub	each	245,40
AGS1987C	High Resolution Au-C Test on 15mm Hitachi Stub	each	245,40
AGS1987D	High Resolution Au-C Test on Customer Stub *	each	245,40
AGS1987E	High Resolution Au-C Test on 12.5mm JEOL Stub	each	245,40
AGS1987S	High Resolution Au-C Test on Short Pin Stub	each	245,40
AGS1987SI	Gold on 500um silicon substrate. 2 - 30nm	each	338,80
AGS1987T	High Resolution Au-C Test on Carbon Disc	each	245,40
AGS1987U	High Resolution Au-C Test (Unmounted)	each	246,40
AGS1987VTSI	Gold on 200µm silicon substrate. 2 - 30nm	each	332,70
AGS1988	Low KV Tin on Carbon Specimen on 12.5mm Pin Stub	each	181,90
AGS1988A	Low KV Tin on Carbon Specimen on 10mm JEOL Stub	each	181,90
AGS1988B	Low KV Tin on Carbon Specimen on ISI/ABT Stub	each	181,90
AGS1988C	Low KV Tin on Carbon Specimen on 15mm Hitachi Stub	each	181,90
AGS1988D	Low KV Tin on Carbon Specimen on Customer Stub *	each	181,90
AGS1988E	Low KV Tin on Carbon Specimen on 12.5mm JEOL Stub	each	181,90
AGS1988S	Low KV Tin on Carbon Specimen on Short Pin Stub	each	174,90
AGS1988SI	Tin on 500um silicon substrate. 20 - 400nm	each	180,90
AGS1988T	Low KV Tin on Carbon Specimen on Carbon Disc	each	181,90
AGS1988U	Low KV Tin on Carbon Specimen Unmounted	each	181,90
AGS1988VTSI	Tin on 200µm silicon substrate. 20 - 400nm	each	193,50
AGS1990	MRS3 Geller Reference Standard Non Traceable	each	o.a.

Part #	Description	Unit	Price €
AGS1990-C	MRS3 Cleaning Only (AGS1990)		o.a.
AGS1990-G	MRS3 GellerRefStd Non Traceable (1" X 3" Glass Slide Mounti		o.a.
AGS1990-O	MRS3 GellerReferenceStandard NonTraceable in OpticalRetainer	each	o.a.
AGS1990-RXY	MRS3 RXY Recertification & Clean (AGS1990)		o.a.
AGS1990-RXYZ	MRS3 RXYZ Recertification & Clean (AGS1990)		o.a.
AGS1990-S	MRS3 GellerReferenceStandard NonTraceable in SEMRetainer	each	o.a.
AGS1990-T	MRS3 GellerRefStd Non Traceable in TEM 3mm x 0/35mm Retaine		o.a.
AGS1991	MRS3 Geller Reference Standard XY Traceable	each	o.a.
AGS1991-C	MRS3 Cleaning Only (AGS1991)		o.a.
AGS1991-G	MRS3 GellerRefStd XY Traceable (1" X 3" Glass Slide Mountin		o.a.
AGS1991-O	MRS3 GellerRefStandard XY Traceable in OpticalRetainer	each	o.a.
AGS1991-RXY	MRS3 RXY Recertification & Clean (AGS1991)		o.a.
AGS1991-RXYZ	MRS3 RXYZ Recertification & Clean (AGS1991)		o.a.
AGS1991-S	MRS3 GellerRefStandard XY Traceable in SEM Retainer	each	o.a.
AGS1991-T	MRS3 GellerRefStd XY Traceable in TEM 3mm x 0/35mm Retainer		o.a.
AGS1992	MRS3 Geller Reference Standard XYZ Traceable	each	o.a.
AGS1992-C	MRS3 Cleaning Only (AGS1992)		o.a.
AGS1992-G	MRS3 GellerRefStd XYZ Traceable (1" X 3" Glass Slide Mounti		o.a.
AGS1992-O	MRS3 Geller Ref Standard XYZ Traceable in Optical Retainer	each	o.a.
AGS1992-RXY	MRS3 RXY Recertification & Clean (AGS1992)		o.a.
AGS1992-RXYZ	MRS3 RXYZ Recertification & Clean (AGS1992)		o.a.
AGS1992-S	MRS3 Geller Reference Standard XYZ Traceable in SEM Retainer	each	o.a.
AGS1992-T	MRS3 GellerRefStd XYZ Traceable in TEM 3mm x 0/35mm Retaine		o.a.
AGS1995A	CD Calibration Specimen 0.5-10µm Non Certified	each	149,60
AGS1997A	CD Calibration Specimen 0.5-10µm GPTB Certified	each	o.a.
AGS1998	CD Calibration Specimen 500-200-100nm non certified	each	960,00
AGS1999	Agar calibration initial inspection	each	24,40
AGS310	Boron surface, ca. 1mm ² on 1cm ² tria, no adhesive (Standard)	each	494,20
AGS311	Boron surface,ca. 1mm ² on 1cm ² tria, with adhesive (Standard	each	479,90
AGS320	Boron surface, ca. 2mm ² on 1cm ² triangle, no adhesive (Prem)	each	607,30
AGS321	Boron surface ca. 2mm ² on 1cm ² tria, w/ adhesive (Prem)	each	607,30
AGSB100	Grid storage box for 100 grids	each	45,90
AGSEMKIT-FEI	SEM Consumable Kit for FEI with filaments	each	810,00
AGSEMKIT-FEI-2	SEM Consumable Kit for FEI without filaments	each	617,20
AGSEMKIT-FEI-2-H	SEM Consumable Kit for FEI without filaments & without *ar	each	553,50
AGSEMKIT-FEI-H	SEM Consumable Kit for FEI with filaments & without *ardou	each	658,40
AGSEMKIT-HITACHI	* SEM Consumable Kit for Hitachi with filaments	each	1.191,70
AGSEMKIT-HITACHI-2	SEM Consumable Kit for Hitachi without filaments	each	520,90
AGSEMKIT-HITACHI-2	SEM Consumable Kit for Hitachi without filaments & without H	each	461,00
AGSEMKIT-HITACHI-H	SEM Consumable Kit for Hitachi with filaments & without Hazardous items	each	1.148,40
AGSEMKIT-JEOL	SEM Consumable Kit for JEOL with filaments	each	760,40
AGSEMKIT-JEOL-2	SEM Consumable Kit for JEOL without filaments	each	483,10
AGSEMKIT-JEOL-2-H	SEM Consumable Kit for JEOL without filaments & without Hazardous items	each	406,10
AGSEMKIT-JEOL-H	SEM Consumable Kit for JEOL with filaments & without Hazardous items	each	619,30
AGSEMKIT-TESCAN	SEM Consumable Kit for TESCAN with filaments	each	1.191,90
AGSEMKIT-TESCAN-2	SEM Consumable Kit for TESCAN without filaments	each	576,60
AGSEMKIT-TESCAN-2	SEM Consumable Kit for TESCAN without filaments & without Hazardous item:	each	449,80
AGSEMKIT-TESCAN-H	SEM Consumable Kit for TESCAN with filaments & without Hazardous items	each	1.138,40
AGSEMKIT-ZEISS	SEM Consumable Kit for Carl Zeiss with filaments	each	659,50
AGSEMKIT-ZEISS-2	SEM Consumable Kit for Carl Zeiss without filaments	each	600,70
AGSEMKIT-ZEISS-2-H	SEM Consumable Kit for Carl Zeiss without filaments & withou	each	553,50
AGSEMKIT-ZEISS-H	SEM Consumable Kit for Carl Zeiss with filaments & without H	each	613,00
AGT003	Dumont Tweezers Box for 3 Tweezers	each	58,70
AGT005	Dumont Tweezers Box for 5 Tweezers	each	64,40
AGT007	Dumont Tweezers Box for 7 Tweezers	each	69,80
AGT1002	Cryo Gripper, 2mm	each	5,80
AGT1003	Cryo-Gripper Aluminium, 2mm	each	7,60
AGT1003	Cryo-Gripper Aluminium, 2mm	each	7,60

Part #	Description	Unit	Price €
AGT1003	Cryo-Gripper Aluminium, 2mm	each	7,60
AGT1003	Cryo-Gripper Aluminium, 2mm	each	7,60
AGT119	Ultrafine, single Deer Hair with handle	each	38,90
AGT1516	High Precision Tweezers, strong thick flat edged tips	each	34,70
AGT1517	High Precision Tweezers, flat accurate round tips	each	34,70
AGT1518	High Precision Tweezers, curved tips	each	39,90
AGT1519	High Precision Tweezers, very sharp, fine tips	each	35,30
AGT1520	High Precision Tweezers, short fine tips	each	35,30
AGT1521	High Precision Tweezers, extra fine tips	each	37,90
AGT1522	High Precision Tweezers, extra fine tips, superior finish	each	41,30
AGT1523	High Precision Tweezers, extra fine tips, superior finish	each	41,30
AGT1524	High Precision Tweezers, bent tips	each	44,50
AGT1525	High Precision Tweezers, flat angled, fine tips	each	51,60
AGT1526	High Precision Tweezers, very fine, curved tips	each	41,80
AGT1527	High Precision Tweezers, extra fine tips	each	35,30
AGT1529-1-10	ESD safe Sticky swabs, 1mm head	10	25,90
AGT1529-1-25	ESD safe Sticky swabs, 1mm head	25	50,00
AGT1529-2-10	ESD safe Sticky swabs, 2mm head	10	25,90
AGT1529-2-25	ESD safe Sticky swabs, 2mm head	25	50,00
AGT1529-3-10	ESD safe Sticky swabs, 3mm head	10	25,90
AGT1529-3-25	ESD safe Sticky swabs, 3mm head	25	50,00
AGT1529KIT	Mixed Kit of ESD safe sticky swabs, 25x 1/2/3mm (Box of 75)		98,50
AGT1530	High Precision Scissors, Extra Fine, 9cm, Curved -miniature	each	35,40
AGT1531	High Precision Scissors, Extra Fine, 9 cm, Straight -minatur	each	35,40
AGT1532	High Precision Scissors, Extra Fine, 9 cm, Curved - precisio	each	35,40
AGT1533	High Precision Scissors, Extra Fine, 9 cm, Straight - precis	each	35,40
AGT1663-12	SEM Pin Stub Mount Gripper Tweezer	each	37,30
AGT5001	Micro scissors	each	#N/B
AGT5004	Dumont Medical Tweezers serrated, 5 Stainless steel. 0.10 x	each	42,00
AGT5005	Platinum-tipped tweezers	each	774,50
AGT5006	Hexagonal headed nutdriver set, set of 5	each	76,20
AGT5007	Dumont HP crossover tweezers N7 Dumoxel. 0.17 x 0.10mm tip.	each	46,00
AGT5008	Tweezers, opening 1/2 - 3/4 inch	each	53,00
AGT5010	Specimen Pickup Loop with Handle	each	48,20
AGT5011	Spare Platinum Loops	3	78,60
AGT5012	Crosshead screwdrivers, set of 4	each	36,60
AGT5013	Dumont Biology Tweezers 5 Titanium. 0.05 x 0.01mm tip.	each	50,30
AGT5014	Dumont Biology Tweezers 7 Titanium. 0.07 x 0.03mm tip.	each	56,00
AGT5016	Heavy Duty Single Edge Carbon Steel Razor Blades	100	22,50
AGT5017	90mm SEM Gripper, short handled	each	12,30
AGT5018	130mm SEM Gripper, long handled	each	13,90
AGT5019	Blades for Brain Matrices	25	167,40
AGT502	Dumont HP Tweezers 5 Carbon steel. 0.08 x 0.04mm tip.	each	41,70
AGT5022	Long tweezers, 200mm	each	18,20
AGT5023	O-rings.	10	15,80
AGT5024	Disposable Microtome Blades, 76 x 7.9 x 0.3mm	50	227,30
AGT5026	Heavy Duty Cryostat Blades, 76 x 12.6 x 0.52mm	35	359,50
AGT5027	Membrane tweezers	each	44,60
AGT5031	Dumont HP Tweezers 2a Stainless Steel. 1.50 x 0.20mm tip.	each	37,40
AGT5032	Dumont HP Tweezers 3c	each	34,50
AGT5033	Dumont HP Tweezers 4 Stainless Steel. 0.13 x 0.08mm tip.	each	35,90
AGT5034	Dumont HP Tweezers 5 Stainless Steel. 0.10 x 0.06mm tip.	each	37,60
AGT5036	Dumont HP Tweezers 5/45 Stainless Steel. 0.10 x 0.06mm tip.	each	44,40
AGT5038	Dumont HP Tweezers 6 Stainless Steel. 0.17 x 0.10mm tip.	each	52,70
AGT5039	Dumont HP Tweezers 7 Stainless Steel. 0.17 x 0.10mm tip.	each	40,20
AGT5040	Dumont HP Tweezers 7a Stainless Steel. 0.24 x 0.16mm tip.	each	38,80
AGT5042	Dumont HP crossover tweezers N1 Stainless steel. 0.20 x 0.12	each	41,70
AGT5044-SV	Replaceable ESD Safe, PVDF Tip Tweezers, Style 5X, fine tip,	each	72,20

Part #	Description	Unit	Price €
AGT5045	Dumont Tweezers SS - Electronic Tweezers	each	33,10
AGT5046	Spade-end tweezers	each	23,10
AGT5048	Tweezers for 10mm dia stubs, 130mm length	each	30,40
AGT5049	Tweezers for 15mm dia stubs, 140mm length	each	#N/B
AGT5049-25	Tweezers for 25mm dia stubs, 153mm length	each	#N/B
AGT5049-32	Tweezers for 32mm dia stubs, 154mm length	each	36,60
AGT5049-SV	ESD Safe Carbon Wafer Tweezers, Style 4WF for 5.8 - 10.2cm w	each	72,20
AGT5049-SV-1	Replacement tips for ESD Safe Carbon Wafer Tweezers, Style 4	each	29,30
AGT505	Dumont HP Tweezers 4 Dumoxel 0.13 x 0.08mm tip.	each	38,30
AGT5052	Wafer Tweezers Type 2WF	each	39,40
AGT5054	Wafer Tweezers Type 3WF	each	34,00
AGT5055	Wafer Tweezers Type 4WL	each	35,60
AGT5056	Wafer Tweezers Type 4WF	each	39,40
AGT5057	Wafer Tweezers Type 2WFG	each	40,40
AGT5059	Wafer Tweezers Type 4WFG	each	39,40
AGT506	Dumont HP Tweezers 5 Dumoxel 0.10 x 0.06mm tip.	each	39,40
AGT5060	Wafer Tweezers Type 34A	each	48,70
AGT5061	Wafer Tweezers Type 35A	each	33,50
AGT5063	Wafer Tweezers Type 36A	each	51,30
AGT5065	Wafer Tweezers Type 35S	each	#N/B
AGT5066	Component Handling Tweezers Type 571	each	45,60
AGT5070	Component Handling Tweezers Type 578	each	47,40
AGT5071	Dumont HP crossover tweezers NOC Stainless steel. 0.17 x 0.1	each	37,10
AGT5073A	Heat Sink Tweezers	each	23,90
AGT5074	Stainless Steel Scissors with Both Points Blunt	each	28,20
AGT5075	Stainless Steel Scissors with 1 Blunt point and 1 Sharp poin	each	26,80
AGT5076	Electrician's Screwdriver, 75mm blade length, tip width 3.0m	each	5,70
AGT5077	Electrician's Screwdriver, 150mm blade length, tip width 5.0	each	5,00
AGT5078	Stubby Screwdriver, 45mm blade length, tip width 6.3mm	each	#N/B
AGT5079	Adjustable Inspection Mirror	each	#N/B
AGT508	Dumont HP Tweezers 7 Dumoxel 0.17 x 0.10mm tip.	each	40,20
AGT5081	Spatula, Chattaway, Stainless Steel, Micropatten Blades, 100	each	19,20
AGT5082	Spatula with spoon at one end, Nickel, 120mm	each	59,90
AGT5083	Spatula (Palette Knife) with wooden handle and stainless ste	each	114,10
AGT5084	Flat nosed pliers	each	45,70
AGT5085	Top cutting pliers for hard wire	each	15,10
AGT5086	Mini side cutting nipper	each	36,90
AGT509	Dumont HP Tweezers 2a Dumoxel 1.50 x 0.20mm tip.	each	36,00
AGT5090	1/16inch Socket Wrench	each	#N/B
AGT5091	5/64inch Socket Wrench	each	#N/B
AGT5092	3/32inch Socket Wrench	each	#N/B
AGT5093	1/8inch Socket Wrench	each	#N/B
AGT5094	5/32inch Socket Wrench	each	#N/B
AGT5095	3/16inch Socket Wrench	each	#N/B
AGT5096	7/32inch Socket Wrench	each	#N/B
AGT5097	1/4inch Socket Wrench	each	#N/B
AGT5098	1.5mm Socket Wrench	each	1,10
AGT5099	2mm Socket Wrench	each	#N/B
AGT510	Dumont Tweezers AA - Electronic Tweezers	each	28,70
AGT5100	2.5mm Socket Wrench	each	#N/B
AGT5101	3mm Socket Wrench	each	1,30
AGT5102	4mm Socket Wrench	each	#N/B
AGT5103	5mm Socket Wrench	each	#N/B
AGT5104	6mm Socket Wrench	each	#N/B
AGT5106	Safety dispenser for 10 single edge razor blades	each	5,60
AGT5107	Steel rule, 150mm	each	29,00
AGT5112	Perfect Loop with Handle	each	179,00
AGT5112LM	Perfect Loop, LM, Loop and Handle Set	each	185,50

Part #	Description	Unit	Price €
AGT5113	Replacement Loop	each	168,70
AGT5113LM	Perfect Loop, LM, Loop only	each	150,90
AGT5114	WECPREP Blades	each	190,70
AGT5115	FEATHER blades	each	17,20
AGT5116	Paper points	each	35,20
AGT5118	Glass Filled Nylon Pliers Flat Jaws, 14.8cm, white	each	#N/B
AGT512	Dumont Tweezers 24/45 - Electronic tweezers	each	48,80
AGT5128	Dumont Tweezers 5 - with clamping ring	each	71,70
AGT5130	Dumont Biology Tweezers 5 Stainless Steel. 0.05 x 0.01mm tip	each	49,70
AGT5131	Dumont Biology Tweezers 7 Stainless Steel. 0.07 x 0.03mm tip	each	48,80
AGT5142	High Profile Disposable Blades, 76 x 14.1 x 0.3mm	50	227,30
AGT515	Diamond Grips	each	20,40
AGT5150	Ceramic Tweezers - Straight Alloy Shank	each	154,40
AGT5154	Ceramic Tweezers - Curved Polyacetyl Shank	each	94,20
AGT516	Tweezers guards.	20	6,40
AGT5165CP	Fine Tip/Flat Strong Tip 150mm CP Probe	each	9,30
AGT5165SV	Fine Tip/Flat Strong Tip 150mm SV Probe	each	8,40
AGT5166CP	Curved Fine Tip/Flat Strong Tip 148mm CP Probe	each	10,50
AGT5166SV	Curved Fine Tip/Flat Strong Tip 148mm SV Probe	each	8,40
AGT5167CP	Flat Large Fine Tip/Flat Fine Sharp Tip 140mm CP Probe	each	10,50
AGT5167SV	Flat Large Fine Tip/Flat Fine Sharp Tip 140mm SV Probe	each	8,40
AGT5168CP	Set of 3 CP Probes	each	29,80
AGT5168SV	Set of 3 SV Probes	each	23,40
AGT5175	Butterfly tweezers	each	#N/B
AGT5179	Tweezers holders for 5 tweezers	each	42,90
AGT5179-3	Tweezers holders for 3 tweezers	each	36,20
AGT5179-8	Tweezers holders for 8 tweezers	each	53,80
AGT5180	Tweezers holder	each	40,50
AGT5181	Wafer Nylon tips Tweezers, Style 8WNY	each	147,70
AGT5187	All-purpose utility scissors	each	11,90
AGT518A	Plastic Tweezers - Acid Resistant	each	12,70
AGT518B	Plastic Tweezers - Delrin		6,60
AGT518C	Plastic Tweezers - Carbon fiber ESD safe	each	6,60
AGT519	Plastic tweezers, CTFE	each	35,60
AGT520	Dumont tweezers 24 - Electronic tweezers	each	38,80
AGT5201	Gold plated tweezers - Straight, extra fine tips 0.08mm		37,70
AGT5202	Gold plated tweezers - Straight, extra fine tips 0.06mm		#N/B
AGT5203	Gold plated tweezers - Angled, extra fine tips		38,80
AGT5204	Gold plated tweezers - Anti-Acid/Anti-Mag SS - Reverse action - tips: anti capil	each	51,90
AGT5205	Gold plated tweezers - Very fine, curved tips		42,10
AGT5206	Mount gripper for grooved 25mm pin stubs	each	25,60
AGT5207-10	Disposable Scalpel No.10	10	11,90
AGT5207-11	Disposable Scalpel No.11	10	11,90
AGT5207-15	Disposable Scalpel No.15	10	11,90
AGT5207-15A	Disposable Scalpel No.15a	10	11,90
AGT5208	Retractable Blade Scalpel, handle only	each	13,80
AGT5209-61	Scalpel Blades No.61	5	13,80
AGT5209-62	Scalpel Blades No.62	5	13,80
AGT5209-64	Scalpel Blades No.64	5	13,80
AGT5209-65	Scalpel Blades No.65	5	13,80
AGT5209-65A	Scalpel Blades No.65a	5	13,80
AGT5209-67	Scalpel Blades No.67	5	13,80
AGT5209-90	Scalpel Blades No.90	5	78,20
AGT5209-91	Scalpel Blades No.91	5	78,20
AGT521	Dumont Tweezers 3C - Electronic tweezers	each	37,30
AGT5210	Scalpel blade remover, Box of 50	50	30,00
AGT5211	Fine Scalpel Handle 130mm	each	40,80
AGT5212	Bulk non-sterile blade removers	250	135,50

Part #	Description	Unit	Price €
AGT5214	Fine Scalpel Handle 100mm	each	40,80
AGT5215	Fine Scalpel Handle 75mm	each	40,80
AGT5216	Fine Scalpel Handle 50mm	each	40,80
AGT5217-10	Sterile Disposable Scalpels Shape 10	10	9,10
AGT5217-10A	Sterile Disposable Scalpels Shape 10a	10	9,10
AGT5217-11	Sterile Disposable Scalpels Shape 11	10	9,10
AGT5217-12	Sterile Disposable Scalpels Shape 12	10	9,10
AGT5217-15	Sterile Disposable Scalpels Shape 15	10	9,10
AGT5217-15A	Sterile Disposable Scalpels Shape 15a	10	9,10
AGT5217-20	Sterile Disposable Scalpels Shape 20	10	9,80
AGT5217-21	Sterile Disposable Scalpels Shape 21	10	9,80
AGT5217-22	Sterile Disposable Scalpels Shape 22	10	9,80
AGT5217-22A	Sterile Disposable Scalpels Shape 22A	10	6,60
AGT5217-23	Sterile Disposable Scalpels Shape 23	10	9,80
AGT5217-24	Sterile Disposable Scalpels Shape 24	10	9,80
AGT5218	Disposable scalpels.	5	5,90
AGT522	Dumont Tweezers 5 - Electronic tweezers	each	33,10
AGT5220	Curved Vanna's Microscissors	each	181,50
AGT5223	Straight Vanna's Microscissors	each	217,00
AGT5228	Vanna's-Type Micro Scissors Straight 80mm	each	59,80
AGT5229	Vanna's-Type Micro Scissors Curved 80mm	each	59,80
AGT523	Dumont Tweezers 7 - Electronic tweezers	each	37,40
AGT5230	Plastic Tweezers - 6mm Flat Tips	each	#N/B
AGT5231	Plastic Tweezers - 9mm Flat Tips	each	25,60
AGT5232	Plastic Tweezers - 12.5mm Flat Tips	each	25,60
AGT5234	Plastic tweezers, type 2a	each	14,20
AGT5235	Plastic tweezers, angled 45°	each	19,60
AGT5236	Plastic tweezers, heat-resistant	each	25,00
AGT5240P	Dissecting Kit 1 In Canvas Holdall	each	114,70
AGT5242P	Dissecting Kit 2 In Canvas Holdall	each	101,80
AGT5244P	Dissecting Kit 3 In Canvas Holdall	each	123,60
AGT5246P	Dissecting Kit 4 In Canvas Holdall	each	110,50
AGT5248	Holdall for dissecting kits	each	24,60
AGT525	Dumont Tweezers 2A - Electronic tweezers	each	34,50
AGT5251	Dumont HP Tweezers 2 Stainless steel. 0.30 x 0.13mm tip.	each	35,90
AGT5253	Dumont HP Tweezers 3 Stainless Steel. 0.17 x 0.10mm tip.	each	37,40
AGT526	Dumont Tweezers P - Electronic Tweezers	each	#N/B
AGT5262	Dumont HP Tweezers 7a Dumoxel 0.24 x 0.16mm tip.	each	41,70
AGT5263	Dumont HP Tweezers SS Carbon steel. 0.14 x 0.10mm tip.	each	34,70
AGT5264	Dumont HP Tweezers SS Stainless Steel. 0.20 x 0.12mm tip.	each	49,60
AGT5267	Dumont HP crossover tweezers N1 Dumoxel. 0.20 x 0.12mm tip.	each	33,20
AGT5268	Dumont HP crossover tweezers N5 Dumoxel. 0.10 x 0.06mm tip.	each	41,70
AGT5269	Dumont HP crossover tweezers N7 Stainless steel. 0.17 x 0.10	each	51,70
AGT5279	Dumont Medical Tweezers 7a Stainless steel. 0.24 x 0.16mm ti	each	50,30
AGT5288	Dumont Biology Tweezers 4 Stainless Steel. 0.06 x 0.02mm tip	each	40,20
AGT5289	Dumont Biology Tweezers 4 Dumoxel. 0.06 x 0.02mm tip.	each	44,10
AGT5291	Dumont Biology Tweezers 5 Dumoxel. 0.05 x 0.01mm tip.	each	43,10
AGT5293	Dumont Biology Tweezers N5 Stainless Steel. 0.05 x 0.01mm ti	each	53,10
AGT5294	Dumont Biology Tweezers N5 Dumoxel. 0.05 x 0.01mm tip.	each	62,30
AGT5295	Dumont Biology Tweezers N5 Titanium. 0.05 x 0.01mm tip.	each	58,80
AGT5298	Dumont Biology Tweezers N7 Titanium. 0.07 x 0.03mm tip.	each	59,70
AGT5303	Dumont Biology, Medical Tweezers 4 Stainless steel. 0.06 x 0	each	#N/B
AGT5306	Dumont Biology, Medical mirror polished tweezers 4 Stainless	each	#N/B
AGT5314A	Perry tweezers	each	23,60
AGT5315	Tweezers, opening 1/4 - 3/8 inch	each	53,90
AGT5319	Dumont Tweezers N2A - Crossover tweezers	each	41,70
AGT5321	Springbow Dissecting Scissors, Straight blades, 110mm, extra	each	96,80
AGT5322	Springbow Dissecting Scissors, Curved blades, 100mm, extra f	each	83,60

Part #	Description	Unit	Price €
AGT5324	Scalpel SS 38mm Blade	each	15,10
AGT5325	Scalpel SS 45mm Blade	each	13,90
AGT5326	Scalpel SS 50mm Blade	each	13,90
AGT5327	Scalpel Handle No.4	each	6,90
AGT5328	Scalpel Handle No.4L	each	17,70
AGT5329	Scalpel Handle No.3	each	6,90
AGT5330	Scalpel Handle No.3L	each	17,70
AGT5331	Scalpel Handle No.7	each	17,20
AGT5332	Single Edge PTFE coated Stainless Steel Razor Blades	100	30,90
AGT5334	Mini Diamond Grips	each	30,00
AGT5335	Micro Tool - Fork (B)	each	#N/B
AGT5336	Micro Tool - Brush (E)	each	91,70
AGT5337	Micro Tool - Spatula (G)	each	125,00
AGT5338	Micro Tool - 90°ree Hook (H)	each	#N/B
AGT5340	Micro Tool - Straight Needle (N)	each	137,90
AGT5341	Micro Tool - 30°ree Needle (P)	each	160,20
AGT5342	Micro Tool - 90°ree Needle (Q)	each	153,20
AGT5343	Micro Tool - Needle Scribe (R)	each	#N/B
AGT5344	Micro Tool - Handle	each	44,40
AGT5346	Glass scribe	each	65,50
AGT5347	Glass writing diamond	each	101,80
AGT5351	Service engineer's tool kit	each	#N/B
AGT5353	Tool kit for EM lab	each	#N/B
AGT5354	Chinagraph pencils, black	12	42,30
AGT536	Imperial socket screw wrenches 1/16 to 3/8" across flats. Se	each	#N/B
AGT5361	Micro Tool - Needle Long Steel Shank (0.12mm) (A)	each	116,10
AGT5362	Micro Tool - Spade (J)	each	#N/B
AGT5363	Micro Tool - Stainless Steel Mirror (C)	each	#N/B
AGT5364	Micro Tool - Scale 10mm x 0.1mm Divisions (D)	each	#N/B
AGT5366	Micro Tool - Diamond Scribe 60°. Steel Shank (F)	each	99,30
AGT5367	Micro Tool - Saw (L)	each	91,70
AGT5368	Micro Tool - Needle Short Steel Shank (0.12mm) (M)	each	99,00
AGT5369	Micro Tool - Triangular Diamond File (T)	each	#N/B
AGT537	Metric socket wrenches 1.5, 2, 2.5, 3, 4, 5 and 6 mm across	each	#N/B
AGT5372	Springbow dissecting scissors, 95mm	each	27,30
AGT5374	Solvent dispenser 120ml pink	each	23,40
AGT5375	Single Edge, Long Stainless Steel Blades	each	231,00
AGT5377	Solvent dispenser 180ml pink	each	26,00
AGT5384	Dumont HP Tweezers 4 Dumostar 0.13 x 0.08mm tip.	each	52,70
AGT5385	Dumont HP Tweezers 5 Dumostar 0.10 x 0.06mm tip.	each	59,00
AGT5386	Dumont HP Tweezers 7 Dumostar 0.17 x 0.10mm tip.	each	61,90
AGT539	Dumont HP crossover tweezers N5 Stainless steel. 0.10 x 0.06	each	46,00
AGT5390	Dumont Biology Tweezers 5 Dumostar. 0.05 x 0.01mm tip.	each	65,10
AGT5392	Dumont Biology Tweezers 7 Dumostar. 0.07 x 0.03mm tip.	each	79,70
AGT5394NM	Flat Tip 34A.sa	each	45,70
AGT5398	Tweezers for 25mm (1") pin stubs	each	110,30
AGT5398A	Tweezers for 32mm (1 1/4") pin stubs	each	41,20
AGT5400	DumontAnti-magnTweezer w/ ReplaceableTips B2A- Ultrafine		60,10
AGT5400T	Replacement Tip - Ultrafine	each	34,50
AGT5401T	Replacement Tip - Curved	pair	34,50
AGT5402T	Replacement Tip - Straight	pair	34,50
AGT5403R	Dumont Anti-magnetic Tweezer w/ Replaceable Buis Tips- Fine	each	88,00
AGT5403T	Replacement Tip - Fine	each	34,50
AGT5404	Dumont Replaceable Buis Tip Tweezers - Fine Rounded	each	#N/B
AGT5404T	Replacement Tip - Fine Rounded	pair	36,90
AGT541	Diagonal cutters	each	78,30
AGT5410	Clutch Type Scratch Brush 2mm	each	20,20
AGT5410A	Refills for AGT5410	3	5,00

Part #	Description	Unit	Price €
AGT5410B	Refills for AGT5410	12	24,10
AGT5412	Propelling Pencil Scratch Brush 4mm with Brass refills	each	14,80
AGT5412B	Refills for AGT5412 Brass Scratch Brush	10	25,60
AGT5413	Wide pattern scratch brush 8mm	each	28,40
AGT5413A	8mm Refill for AGT5413 8mm Scratch Brush	each	15,30
AGT5415	Set of 6 Fine Probes	each	29,00
AGT5416	Precision oil lubricator	each	14,90
AGT5420	Angled Micro Pick, 190mm, 0.63mm (0.25") tip	each	19,00
AGT5421	Straight Micro Pick, 184mm, 0.63mm (0.25") tip	each	14,80
AGT5433	Superfine Eyelash with handle	each	26,40
AGT5441	9.5mm Disc Punch	each	107,80
AGT5442	11.0mm Disc Punch	each	114,10
AGT5443	12.7mm Disc Punch	each	114,10
AGT5444	Fine Point Diamond Scriber Straight	each	58,00
AGT5445	Fine Point Diamond Scriber 30°ree Bend	each	59,80
AGT5450	Micrometer	each	90,60
AGT5451	Vernier caliper	each	134,60
AGT5453A	Digital Vernier Caliper 0-150mm	each	137,40
AGT5460	Pen light	each	#N/B
AGT5463	Easy-Fit Fine Scalpel Handle 130mm	each	49,90
AGT5464	Scalpel Handle B3	each	40,90
AGT5465	Easy-Fit Fine Scalpel Handle 95mm	each	49,90
AGT5467	Combination Screwdrivers	each	54,20
AGT5470	Loop Handle	each	32,60
AGT5471	1.5mm Nichrome Loops	3	14,10
AGT5472	3.5mm Nichrome Loops	3	14,10
AGT5473	5mm Nichrome Loops	3	14,10
AGT548-1	Mounted Pin In Plastic Handle	each	4,10
AGT5481-A60	Diamond Scriber, bent 30° shaft, tip angle 60°, overall	each	29,60
AGT5481-S	Diamond Scriber, straight shaft, tip angle 90°, overall le	each	42,50
AGT5482	Deluxe Diamond Scribing Pen	each	30,50
AGT548-2	Mounted Pin with Aluminium Handle	each	3,50
AGT5482R	Refill for AGT5482 Deluxe Diamond Scribing Pen	each	24,40
AGT5483	Diamond scriber	each	43,70
AGT5484	Micro-Tools microscopist tool set	each	654,30
AGT549	Spatula, Chattaway, Stainless Steel, 178mm	each	32,20
AGT5490-300	Harris Uni-Core 3.0mm Hole	each	45,30
AGT550	Set of four fine tip probes	each	24,50
AGT5500SV	Fine Tip Type SV, Plastic Tweezers, 0.6mm thick tips	each	7,00
AGT5501SV	Plastic Tweezers - tips: angled, pointed - ESD safe	each	7,00
AGT5502	High precision cutting tweezers, 120mm long	each	49,20
AGT5503	High Precision Miniature Cutting Tweezers, 98mm long	each	38,10
AGT5504	High precision cutting tweezers, 120mm long	each	82,50
AGT5506	High precision super Slim Tweezers, very fine bent tips	each	28,70
AGT5507	High precision super Slim Tweezers, extra fine tips	each	32,70
AGT5509	ESD Soft Grip Tweezers, 115mm, Extra Fine Tips	each	46,90
AGT550A	Stainless Steel Straight Needle Probe	each	8,10
AGT550B	Stainless Steel Hook Curve Needle Probe	each	8,10
AGT550C	Stainless Steel Bent Needle Probe	each	8,10
AGT550D	Stainless Steel Large Curve Needle Probe	each	#N/B
AGT551	Locking Forceps	each	27,30
AGT5510	ESD Soft Grip Tweezers, 120mm, Fine Curved Tips	each	49,40
AGT5511	Carbon fibre replacement tip tweezers, flat, round tips	each	30,80
AGT5511-T	Tips for plastic tipped tweezers - flat round tips	each	9,30
AGT5512	Carbon fibre replacement tip tweezers, very fine tips, rever	each	34,90
AGT5512-T	Tips for plastic tipped tweezers - flat fine tips	each	9,30
AGT5513	Carbon fibre replacement tip tweezers, very fine, curved tip	each	30,80
AGT5513-T	Tips for AGT5513 Curved Tip Tweezers	each	9,30

Part #	Description	Unit	Price €
AGT5514	Ceramic Tipped Tweezers - Straight, round tips	each	64,30
AGT5514-T	Tips for Ceramic Tipped Tweezers - Straight	each	40,80
AGT5515	Ceramic Tipped Tweezers - Curved, fine curved tips	each	64,30
AGT5515-T	Tips for Ceramic Tipped Tweezers - Curved	each	46,00
AGT5516	Ceramic Tipped Tweezers - Fine Straight, strong tips	each	64,30
AGT5516-T	Tips for Ceramic Tipped Tweezers - Fine Straight	each	40,80
AGT5520	Magnetic pick-up mirror	each	20,30
AGT5521	Micro powder spatulas, set of 4	each	60,10
AGT5522	Blue micro spatula, plastic handle, fine pointed tip 0.3mm	each	3,90
AGT5523	Black micro spatula, plastic handle, extra fine tip 0.2mm	each	3,90
AGT5524	Red micro spatula, plastic handle, larger tip 0.8mm	each	4,00
AGT5525	Green micro spatula, plastic handle, medium tip 0.5mm	each	4,00
AGT5526	Micro spatulas - kit of four	each	14,00
AGT5527	Vampire Vacuum Pick-Up	each	71,70
AGT5528	Pocket Vacuum Syringe with 3 suction cups	each	16,10
AGT5529	Rubber suction cups kit: 3.5mm, 6.5mm, 9.5mm dia	each	9,60
AGT553	Surgical Scissors	each	39,70
AGT5530	Tweezers kit	each	141,50
AGT5530D	Tweezers kit - Diamond coated		240,90
AGT5530T	Tweezers kit - Titanium		179,20
AGT5531	Wafer Tweezers Type 8WF	each	39,60
AGT5532	Ceramic Tipped Tweezers - Straight, fine tips	each	64,30
AGT5532-T	Tips for Ceramic Tipped Tweezers - Straight, fine	each	40,80
AGT5533	Ceramic Tipped Tweezers - Straight, fine tips	each	64,30
AGT5533-T	Tips for Ceramic Tipped Tweezers - Straight, strong, fine	each	40,40
AGT5534	Ceramic Tipped Tweezers - Straight, fine tips	each	63,90
AGT5534-T	Tips for Ceramic Tipped Tweezers - Straight, fine	each	40,80
AGT5536	Fine ESD Epoxy coated Tweezers SS - Anti-Magnetic 110mm	each	29,70
AGT5537	Dumont Tweezers Inox 08, Tip 0.06 thick x 0.10wide	each	84,80
AGT5538	ESD epoxy coated tweezers, Anti-Acid/Anti-Mag SS, Extra Fine	each	29,70
AGT554	Dissecting Scissors	each	36,60
AGT5540	Kit of four stainless steel probes	each	14,50
AGT5541	Kit of four stainless steel spatulas	each	15,30
AGT5542	Kit of six stainless steel probes	each	22,10
AGT5543	Stainless steel probe with straight needle tip	each	4,60
AGT5544	Stainless steel probe with angled needle tip	each	4,60
AGT5545	Stainless steel probe with single bend tip	each	4,60
AGT5546	Stainless steel probe with hook tip	each	4,60
AGT5547	Stainless steel probe with triple bend tip	each	4,60
AGT5548	Stainless steel probe with flat tip	each	4,60
AGT5549	Stainless steel spatula with short pointed tip and long curv	each	4,60
AGT555	Scalpel Handle No 4, plus 5 Blades, shape 24	each	10,00
AGT5550	Stainless steel spatula with short rounded tip and large fla	each	4,60
AGT5551	Stainless steel spatula with short flat squared tip and long	each	4,60
AGT5552	Stainless steel spatula with long curved pointed tip and lon	each	4,60
AGT5553	ESD Soft Grip Tweezers, 125mm, Fine Tips		45,20
AGT5559	Retractable Disposable Scalpel (Sterile), No.10	25	73,20
AGT5560	Retractable Disposable Scalpel (Sterile), No.11P	25	73,20
AGT5561	Retractable Disposable Scalpel (Sterile), No.15	25	73,20
AGT5562	Retractable Disposable Scalpel (Sterile), No.20	25	102,30
AGT556-20	Scalpel blades, shape 20.	5	2,10
AGT556-21	Scalpel blades, shape 21.	5	2,10
AGT556-22	Scalpel blades, shape 22.	5	2,10
AGT556-22A	Scalpel blades, shape 22A.	5	2,10
AGT556-22H	Scalpel blades, shape 22.	100	33,30
AGT556-23	Scalpel blades, shape 23.	5	2,10
AGT556-24	Scalpel blades, shape 24.	5	2,10
AGT556-24H	Scalpel blades shape 24.	100	33,30

Part #	Description	Unit	Price €
AGT556-25	Scalpel blades, shape 25.	5	2,10
AGT556-26	Scalpel blades, shape 26.	5	2,10
AGT556-27	Scalpel blades, shape 27.	5	2,10
AGT5563	Retractable Disposable Scalpel (Sterile), No.21	25	102,30
AGT5564	Retractable Disposable Scalpel (Sterile), No.22	25	102,30
AGT5565	Retractable Disposable Scalpel (Sterile), No.23	25	102,30
AGT5566	Retractable Disposable Scalpel (Sterile), No.24	25	102,30
AGT5567	Retractable Disposable Scalpel (Sterile), No.25A	25	102,30
AGT5568	Retractable Disposable Scalpel (Sterile), StitchCutter	25	73,20
AGT557	Scalpel Handle No 3, plus 1 blade of each shape: 10, 10A, 11	each	11,30
AGT5580	KLEEN BMS, Fits No.3 Handle, No.10	50	135,40
AGT5581	KLEEN BMS, Fits No.3 Handle, No.110	50	135,40
AGT558-10	Blades No.10 for No.3 Handle	5	1,90
AGT558-10A	Blades No.10A for No.3 Handle	5	1,90
AGT558-11	Blades No.11 for No.3 Handle	5	1,90
AGT558-11P	Blades No.11P for No.3 Handle	5	1,90
AGT558-12	Blades No.12 for No.3 Handle	5	1,90
AGT558-15	Blades No.15 for No.3 Handle	5	1,90
AGT558-15A	Blades No.15A for No.3 Handle	5	1,90
AGT558-15T	Blades No.15T for No.3 Handles	5	1,90
AGT558-16	Blades No.16 for No.3 Handle	5	1,90
AGT5582	KLEEN BMS, Fits No.3 Handle, No.15	50	135,40
AGT5583	KLEEN BMS, Fits No.3 Handle, No.11P	50	135,40
AGT5584	KLEEN BMS, Fits No.3 Handle, No.15T	50	135,40
AGT5585	KLEEN BMS, Fits No.4 Handle, No.20	50	139,20
AGT5586	KLEEN BMS, Fits No.4 Handle, No.21	50	139,20
AGT5587	KLEEN BMS, Fits No.4 Handle, No.22 0	50	139,20
AGT5588	KLEEN BMS, Fits No.4 Handle, No.23 0	50	139,20
AGT5589	KLEEN BMS, Fits No.4 Handle, No.24 0	50	139,20
AGT560	Pin vice, 75mm long	each	#N/B
AGT561	Piercing Saw	each	52,50
AGT5610	Carbon fibre replacement tip tweezers, flat round tips	each	34,90
AGT5611	Carbon fibre replacement tip tweezers, fine tips	each	30,80
AGT5611-T	Tips for AGT5611 Fine Tip Tweezers	each	9,30
AGT5612	Carbon fibre replacement tip tweezers, flat, very fine tips	each	30,80
AGT5613	Carbon fibre replacement tip tweezers, flat, small round tip	each	30,80
AGT5613-T	Tips for AGT5513 small round Tip Tweezers	each	9,30
AGT5614	Carbon fibre replacement tip tweezers, flat, round tips	each	30,80
AGT5614-T	Tips for AGT5614 flat round tip tweezers	each	9,30
AGT5615	Carbon fibre replacement tip tweezers, fine tips	each	30,80
AGT5615-T	Tips for AGT5615 fine tip tweezers	each	9,30
AGT5621	Dumont Type GG Tweezers, strong tips SS (Cryo TEM)	each	42,80
AGT562-A	Blades for Piercing Saw	10	14,60
AGT567	Razor Blade Holder	each	42,70
AGT568	Replacement blades, carbon steel, Dispenser of 20	each	16,20
AGT569	Replacement blades, carbon steel, for AGT567	250	52,60
AGT569A	Replacement blades, stainless steel	each	76,20
AGT569T	Replacement blades, PTFE coated stainless steel, Dispenser	each	12,20
AGT570	Demagnetiser	each	782,50
AGT570-1	Demagnetiser	each	551,20
AGT571	Arkansas Stone, triangular, 75 x 6mm	each	34,90
AGT5710	Diamond sharpening stone, Coarse Mini Stone	each	24,70
AGT5711	Diamond sharpening stone, Fine Mini Stone	each	#N/B
AGT5712	Diamond sharpening stone, Extra Fine Mini Stone	each	#N/B
AGT5713	Set of 3 Mini Stones	each	#N/B
AGT572	Arkansas Stone, square, nominal size 100 x 10 x 10mm	each	33,90
AGT5720	Pen-Vac tweezers, complete kit	each	182,40
AGT5720-1A	Vacuum Cup Kit: Vacuum cup, Black Buna ESD of each size (8 total)	each	58,20

Part #	Description	Unit	Price €
AGT5721	Pen-Vac vacuum pick-up systems, 220V	each	252,20
AGT5721-11	Spare tip for T5721 pickup System 25 Gauge (0.2mm)	each	23,60
AGT5721-9	Spare Tip for AGT5721 Pickup System 18 Gauge (0.8mm)	each	11,40
AGT5722	Pen-Vac Vacuum Pump only	each	129,90
AGT5722-1	Pen-Vac standard with 6 probes and cups, 146mm long, 12.7mm	each	182,40
AGT5724	Dumoxel anti-capillary N4 reverse action, straight tweezers		84,80
AGT5725	Dumoxel anti-capillary N7 precision tipped tweezers		70,70
AGT572A	Arkansas Stone, flat, minimum size 75 x 20 x 6mm	each	49,80
AGT577	Iris Scissors 90mm	each	42,30
AGT577A	Angled Iris Scissors 115mm	each	83,60
AGT579	Carbon rod sharpener	each	77,90
AGT581	Miniature Saw	each	43,30
AGT5811	Diamond Coated Tweezer 00 -Tips: Straight	each	#N/B
AGT5814	Diamond Coated Tweezer 3C -Tips: Straight Sharp Fine	each	#N/B
AGT5815	Diamond Coated Tweezer 2A -Tips: Straight Flat	each	50,30
AGT5819	Diamond Coated Tweezer -Tips: Straight extra fine Sharp	each	54,80
AGT582	Spare Blades for Miniature Saw AGT581	10	25,70
AGT5821	Diamond Coated Tweezer 7 - Tips: Very fine, Curved	each	57,00
AGT583	Single Bristle Holder	each	27,50
AGT584	Spare Bristles for AGT583	10	2,60
AGT585	Single Edge Carbon Steel Razor Blades	100	15,90
AGT586	Single Edge Stainless Steel Razor Blades	100	29,90
AGT588	Dumont Tweezers N5AC	each	76,10
AGT592	Hand tool for carbon rod shaping	each	180,60
AGT597	Ballpoint Tip Socket Screwdriver (Metric), Set of 7	each	147,00