

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: AF4 Antifadent Solution**
- **Article number: AGR1327**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
A glycerol- based antioxidant solution for use in fluorescence microscopy
- **Application of the substance / the preparation:** No further relevant information.
- **1.3 Details of the supplier of the safety data sheet**
- **Supplier.**  
Agar Scientific  
Elektron Technology UK Ltd  
Parsonage Lane  
Stansted CM24 8GF  
United Kingdom  
sales@agarscientific.com  
Tel: +44 (0) 1279 813 519
- **Further information obtainable from:** Technical Support
- **1.4 Emergency telephone number:** 24 hours: +44 (0)1856 407333

## SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



GHS07

Skin Sens. 1 H317 May cause an allergic skin reaction.

- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**



Xi; Sensitising

R43: May cause sensitisation by skin contact.

- **Information concerning particular hazards for human and environment:**

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

- **Classification system:**

The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

- **Hazard pictograms**



GHS07

- **Signal word** Warning
- **Hazard-determining components of labelling:**  
propyl gallate
- **Hazard statements**  
H317 May cause an allergic skin reaction.

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

**Trade name: AF4 Antifadent Solution**

(Contd. of page 1)

### • Precautionary statements

- P101 If medical advice is needed, have product container or label at hand.  
 P102 Keep out of reach of children.  
 P103 Read label before use.  
 P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.  
 P321 Specific treatment (see on this label).  
 P363 Wash contaminated clothing before reuse.  
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
 P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

### • 2.3 Other hazards

### • Results of PBT and vPvB assessment

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

## SECTION 3: Composition/information on ingredients

### • 3.2 Chemical characterisation: Mixtures

- **Description:** Mixture of substances listed below with nonhazardous additions.

#### • Dangerous components:

CAS: 56-81-5 EINECS: 200-289-5	glycerol substance with a Community workplace exposure limit	95-99%
CAS: 121-79-9 EINECS: 204-498-2	propyl gallate ☒ Xn R22; ☒ Xi R43 ☠ Acute Tox. 4, H302; Skin Sens. 1, H317	1-5%

- **Additional information:** For the wording of the listed risk phrases refer to section 16.

## SECTION 4: First aid measures

### • 4.1 Description of first aid measures

#### • After inhalation:

If unconscious, place the patient on their side in the recovery position and ensure they are breathing. Clean nose and mouth with water. Artificial respiration may be administered by suitably qualified first-aiders. Get medical attention if any discomfort continues.

#### • After skin contact:

Remove affected person from source of contamination. Remove contaminated clothing and shoes. Wash the skin immediately with soap and water. Take especial care to clean folds, crevices, crease and groin. Get medical attention immediately.

#### • After eye contact:

Check for contact lenses which must be removed from the eyes before rinsing. Promptly rinse eyes with plenty of clean water while lifting the eyelids. Continue to rinse for at least 15 minutes. Continue until the eyes are free of all traces of contamination. Get medical attention if any discomfort or irritation persists.

#### • After swallowing:

If swallowed do NOT induce vomiting. Never give anything by mouth to an unconscious person. If person vomits, keep head low to prevent vomit entering the lungs. Rinse mouth thoroughly. Get medical attention if symptoms persist. Keep patient at rest

### • 4.2 Most important symptoms and effects, both acute and delayed

No further relevant information available.

### • 4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

GB

(Contd. on page 3)

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

Trade name: **AF4 Antifadent Solution**

(Contd. of page 2)

## SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
- **Suitable extinguishing agents:**  
Extinguish with alcohol-resistant foam, carbon dioxide, or dry powder.
- **For safety reasons unsuitable extinguishing agents:**  
Do not use a water jet as this may spread the fire
- **5.2 Special hazards arising from the substance or mixture**  
If heated, vapours/gases hazardous to health (e.g. CO, NOx,) may be formed.
- **5.3 Advice for firefighters**
- **Protective equipment:**  
Wear self contained breathing apparatus and full protective clothing. Keep all unnecessary people away. fire water run-off must not be allowed to contaminate ground or enter drains, sewers or water courses. Provide bunding against fire water run-off.

## SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**  
Wear protective gloves
- **6.2 Environmental precautions:** Prevent discharge into drains, water courses or onto the ground
- **6.3 Methods and material for containment and cleaning up:**  
Stop leak if possible without risk. Absorb in vermiculite, dry sand or earth and place into containers. Flush with plenty of water to clean spillage area. Do not contaminate sewers or water courses.
- **6.4 Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**  
Avoid spilling, skin and eye contact. Wash thoroughly after handling.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**  
Store in tightly closed original container in a dry, cool and well-ventilated place.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **7.3 Specific end use(s)** No further relevant information available.

## SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:**  
Provide adequate ventilation. Observe occupational exposure limits and minimise the risk of inhalation of vapours.
- **8.1 Control parameters**

• <b>Ingredients with limit values that require monitoring at the workplace:</b>
----------------------------------------------------------------------------------

<b>56-81-5 glycerol</b>
-------------------------

WEL	Long-term value: 10 mg/m <sup>3</sup>
-----	---------------------------------------

(Contd. on page 4)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

**Trade name: AF4 Antifadent Solution**

(Contd. of page 3)

- **Additional information:** The lists valid during the making were used as basis.
- **8.2 Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
 DO NOT SMOKE IN WORK AREA! Wash hands at the end of each work shift and before eating, smoking and using the toilet. Wash promptly if skin becomes wet or contaminated. Promptly remove any clothing that becomes contaminated. Use appropriate skin cream to prevent drying of skin. When using do not eat, drink or smoke.
- **Respiratory protection:**  
 In case of inadequate ventilation or when the product is heated use suitable respiratory equipment with gas filter (type A2).
- **Protection of hands:**



Protective gloves

Wear suitable protective gloves conforming to EN 374. Seek recommendation from manufacturer or supplier. After using gloves the hands should be washed and dried thoroughly and a suitable moisturiser applied.

- **Material of gloves**  
 The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**  
 The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Use approved safety goggles or face shield.
- **Body protection:** Wear appropriate clothing to prevent any possibility of skin contact.

## SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**

Form:	Liquid
Colour:	Colourless
- **Odour:** Characteristic
- **Odour threshold:** Not determined.
- **pH-value:** Not determined.
- **Change in condition**

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gaseous):** Not applicable.
- **Ignition temperature:**
- **Decomposition temperature:** Not determined.
- **Self-igniting:** Product is not selfigniting.
- **Danger of explosion:** Product does not present an explosion hazard.

(Contd. on page 5)

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

Trade name: **AF4 Antifadent Solution**

(Contd. of page 4)

· <b>Explosion limits:</b>	
<b>Lower:</b>	Not determined.
<b>Upper:</b>	Not determined.
· <b>Vapour pressure:</b>	
Not determined.	
· <b>Density:</b>	
Not determined.	
· <b>Relative density</b>	
Not determined.	
· <b>Vapour density</b>	
Not determined.	
· <b>Evaporation rate</b>	
Not determined.	
· <b>Solubility in / Miscibility with water:</b>	
Fully miscible.	
· <b>Partition coefficient (n-octanol/water):</b>	
Not determined.	
· <b>Viscosity:</b>	
<b>Dynamic:</b>	Not determined.
<b>Kinematic:</b>	Not determined.
· <b>9.2 Other information</b>	
No further relevant information available.	

## SECTION 10: Stability and reactivity

- **10.1 Reactivity**
- **10.2 Chemical stability** Stable under normal conditions of storage and use.
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** Avoid exposure to heat in any form including sunlight
- **10.5 Incompatible materials:** Oxidising agents, strong acids.
- **10.6 Hazardous decomposition products:**  
Decomposition can lead to vapours/gases hazardous to health (e.g. CO, NO<sub>x</sub>,) being formed.

## SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
- **Acute toxicity**
- **Skin corrosion/irritation** May cause skin sensitisation by skin contact
- **Serious eye damage/irritation** May cause temporary irritation
- **Respiratory or skin sensitisation** Sensitisation possible through skin contact.
- **Other information (about experimental toxicology):**  
Inhalation: In high concentrations or on heating, vapours may irritate the throat and respiratory system causing coughing.  
Ingestion: May cause discomfort if swallowed.
- **Additional toxicological information:**  
The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:  
Irritant

## SECTION 12: Ecological information

- **12.1 Toxicity**  
The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.

(Contd. on page 6)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

**Trade name: AF4 Antifadent Solution**

(Contd. of page 5)

- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** This product is expected to be biodegradable
- **12.3 Bioaccumulative potential**  
This product does not contain any substances expected to be bioaccumulating
- **12.4 Mobility in soil**  
This product is completely water soluble and therefore may spread into water systems
- **Additional ecological information:**
- **General notes:**  
Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water  
Do not allow product to reach ground water, water course or sewage system.  
Danger to drinking water if even small quantities leak into the ground.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects** No further relevant information available.

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**  
Product is hazardous waste. Do not allow into drains, sewers or water courses. Disposal must be by means of a waste contractor
- **Recommendation:**  
Small amounts may be flushed with water to sewer. Larger volumes must be sent to approved plant for destruction.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

### SECTION 14: Transport information

- |                                                                                        |                 |
|----------------------------------------------------------------------------------------|-----------------|
| · <b>14.1 UN-Number</b><br>· <b>ADR, IMDG, IATA</b>                                    | Void            |
| · <b>14.2 UN proper shipping name</b><br>· <b>ADR, IMDG, IATA</b>                      | Void            |
| · <b>14.3 Transport hazard class(es)</b><br>· <b>ADR, IMDG, IATA</b><br>· <b>Class</b> | Void            |
| · <b>14.4 Packing group</b><br>· <b>ADR, IMDG, IATA</b>                                | Void            |
| · <b>14.5 Environmental hazards:</b><br>· <b>Marine pollutant:</b>                     | No              |
| · <b>14.6 Special precautions for user</b>                                             | Not applicable. |
| · <b>14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b>  | Not applicable. |
| · <b>UN "Model Regulation":</b>                                                        | -               |

GB

(Contd. on page 7)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 12.10.2015

Revision: 12.10.2015

**Trade name: AF4 Antifadent Solution**

(Contd. of page 6)

### SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

This information is based on our present knowledge and should assist the user with the safe handling of this material when properly applied. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**  
H302 Harmful if swallowed.  
H317 May cause an allergic skin reaction.  
.....  
R22 Harmful if swallowed.  
R43 May cause sensitisation by skin contact.
- **Department issuing MSDS:** Sales department
- **Contact:**  
sales@agarscientific.com  
Tel: +44 (0) 1279 813 519
- **Abbreviations and acronyms:**  
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonised System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)  
VOC: Volatile Organic Compounds (USA, EU)  
Acute Tox. 4: Acute toxicity, Hazard Category 4  
Skin Sens. 1: Sensitisation - Skin, Hazard Category 1

GB